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

#### Copurification of miRNA and total RNA

#### Separate purification of miRNA-enriched fraction and total RNA\*

\* Separate purification using the miRNeasy Mini Kit requires additionally the RNeasy MinElute Cleanup Kit. For separate purification using the miRNeasy 96 Kit, an additional RNeasy 96 plate is required. For economical purchase of an additional RNeasy 96 plate, we recommend ordering an RNeasy 96 Kit (订购信息见第10页)。

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Total RNA was purified from A) 10<sup>2</sup>-10<sup>7</sup> Jurkat cells using the miRNeasy 96 Kit or B) a dilution series of rat lung tissue homogenate from 20 mg to 200 pg tissue, using the miRNeasy Mini Kit. miRNA-enriched fractions (<200 nt) were also isolated from the same samples. Purified RNA was used as a template in quantitative, real-time RT-PCR assays for the miRNA miR-16.

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Total RNA including miRNA was purified from 1x10<sup>6</sup> Jurkat cells using the miRNeasy Mini Kit or precipitation from TRIzol Reagent. OD260 was lower when using TRIzol, indicating phenol carryover. In addition, the OD230 measurement was higher in the TRIzol prepared RNA, indicating salt carryover from the TRIzol Reagent.


图3. miRNeasy Kit纯化的结果优于TRIzol。

Real-time RT-PCR assays for the miRNA miR-16. Results showed that CT values were lower after purification using the miRNeasy Kit, indicating that higher amounts of miRNA were purified than when using TRIzol. miRNA was effectively purified from as little as 1x10<sup>2</sup> cells using the miRNeasy Kit. In contrast, no miRNA was detected after 40 PCR cycles from 1x10<sup>2</sup> cells when TRIzol was used for purification.

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

Total RNA including miRNA was purified from 25 mg of a range of RNAlater® stabilized rat tissues using the miRNeasy 96 Kit. Purified RNA was used as a template in quantitative, real-time RT-PCR assays for the miRNA miR-16 and for the larger mRNA of the PGK1 gene. Results showed successful detection of both PGK1 mRNA (large RNA) and miR-16 (small RNA) from the same eluates.

## miScript Reverse Transcription Kit

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
**图5. 单次cDNA合成反应可用于检测多种miRNAs和mRNAs**  
Total RNA was prepared from HeLa S3 cells using the miRNeasy Mini Kit. The miScript System was used for real-time PCR analysis of 4 miRNAs(Let-7a, miR-16, miR-21, and miR-25). QuantiTect Primer Assays were used for real-time PCR analysis of 4 mRNAs (GAPDH, PPIA, CDC2, and Erk2).

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## miScript System — miScript SYBR Green PCR Kit

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Synthetic miRNAs of each Let-7 isoform were used in cDNA synthesis reactions performed with the miScript Reverse Transcription Kit. An aliquot of the resultant cDNA was used as a template in real-time PCR reactions with a miScript Primer Assay for each isoform and the miScript SYBR Green PCR Kit. The % relative detection was calculated using the differences between the CT values achieved from the mismatching miScript Primer Assays and those from the perfectly matching miScript Primer Assays (% relative detection =  $2^{-\Delta CT} \times 100$ ).

miRNA sequence	
Let-7a	UGAGGUAGUAGGUUGUAUAGUU
Let-7b	UGAGGUAGUAGGUUGUGUGUU
Let-7c	UGAGGUAGUAGGUUGUAUGUU
Let-7d	<u>A</u> GAGGUAGUAGGUUG <u>C</u> AUAGU•
Let-7e	UGAGGUAG <u>G</u> AGGUUGUAUAGU•
Let-7f	UGAGGUAGUAG <u>A</u> UUGUAUAGUU
Let-7g	UGAGGUAGUAG <u>U</u> UUGUA <u>C</u> AGU•
Let-7i	UGAGGUAGUAG <u>U</u> UUGU <u>G</u> CUGU•

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These sequences show the Let-7 isoforms. Base changes are red and underlined. Changes in length are indicated by a red dot.

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RNA was purified from HeLa S3 cells using the miRNeasy Mini Kit. A range of amounts of RNA from 10 pg to 1 μg were used in cDNA synthesis reactions using the miScript Reverse Transcription Kit. CDNA was used as a template in quantitative, real-time PCR assays for 5 miRNAs (miR-16, miR-20a, miR-25, miR-28, and miR-26a).

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## miScript miRNA Mimics and miScript miRNA Inhibitors

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
HDAC4 is predicted to be a target of miR-1 and has 2 miR-1 binding sites in the 3' UTR (see Ref. 2). HeLa S3 cells ( $8 \times 10^4$  cells/well) were transfected using HiPerFect Transfection Reagent and a range of concentrations of miR-1 Mimic. After 24 and 48 hours, **A** HDAC4 protein level was analyzed by western blotting (indicated with arrows). Tubulin was also analyzed as an internal control (not shown). No effect was observed after 24 hours. After 48 hours, protein level was significantly decreased. **B** HDAC4 transcript level was analyzed by real-time PCR using the miScript System after transfection using 50 nM miR-1 Mimic. The negative control was AllStars Negative Control siRNA. miR-1 #1, miR-1 #2, and miR-1 #3 are different batches of miR-1 Mimic. Results were normalized to GAPDH and to the negative control and are presented as  $2^{-\Delta\Delta CT}$ . Transcript levels did not change significantly at either the 24-hour or 48-hour time point.

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miR-1 is not endogenously expressed in HeLa S3 cells. miR-16 is endogenously expressed at high levels in HeLa S3 cells (see Ref. 1). Using HiPerFect Transfection Reagent, HeLa S3 cells ( $6 \times 10^4$  cells/well) were transfected with a luciferase reporter construct with a binding site for miR-1 in the 3' UTR; cotransfected with a luciferase reporter construct with a binding site for miR-1 in the 3' UTR and with an miR-1 Mimic; or transfected with a luciferase reporter construct with a binding site for miR-16 in the 3' UTR. Twenty-four hours later, luciferase activity was measured. miR-1 Mimic resulted in a similar decrease in luciferase level as endogenous miR-16.

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

**图11. 转染miScript miRNA Inhibitor后下调了内生miRNA的表达水平。** HeLa S3 cells ( $6 \times 10^4$  cells/well) were transfected using HiPerFect Transfection Reagent with miScript miR-16 Inhibitor or a negative control. Endogenous miR-16 levels were quantified by real-time PCR using the miScript System. U6B was used as the reference RNA for data normalization and miR-16 levels were expressed as  $2^{-\Delta\Delta CT}$  values relative to the miR-16 levels in untransfected cells. **A** Transfection of miScript miR-16 Inhibitor strongly reduced the amount of detectable endogenous miR-16. **B** The amplification plots indicated a significant decrease in the level of detectable miR-16 as a result of miR-16 Inhibitor transfection. Dissociation curve analysis (inset) reflects the specificity of the real-time PCR detection using the miScript System.

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HeLa S3 cells express miR-16 at high levels and do not express miR-1. In a cotransfection experiment using HiPerFect Transfection Reagent, HeLa S3 cells ( $6 \times 10^4$  cells/well) were cotransfected with a luciferase reporter construct with an miR-16 binding site in the 3' UTR together with miR-16 Inhibitor. miR-16 Inhibitor was used at varying concentrations in the experiment to evaluate the optimal inhibitor concentration required to see the inhibitory effect. Alternatively, cells were transfected with miR-1 Inhibitor alone as a control. A luciferase construct without an miRNA binding site (WT) was transfected as a positive control. An increase in luciferase expression following transfection of the miR-16 Inhibitor indicated that it prevented endogenous miR-16 from downregulating the reporter gene.

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**表3. 优化24孔板miRNA Mimics/miRNA Inhibitors转染**

Amount (conc.) of miRNA mimic	(10 nM)75 ng	(10 nM)75 ng	(10 nM)75 ng
Amount (conc.) of miRNA inhibitor	(100 nM)750 ng	(100 nM)750 ng	(100 nM)750 ng
Volume of HiPerFect Reagent	1.5 $\mu$ l	3 $\mu$ l	4.5 $\mu$ l
Amount (conc.) of miRNA mimic	(5 nM)37.5 ng	<b>(5 nM)37.5 ng</b>	(5 nM)37.5 ng
Amount (conc.) of miRNA inhibitor	(50 nM)375 ng	<b>(50 nM)375 ng</b>	(50 nM)375 ng
Volume of HiPerFect Reagent	1.5 $\mu$ l	<b>3 <math>\mu</math>l</b>	4.5 $\mu$ l
Amount (conc.) of miRNA mimic	(1 nM)7.5 ng	(1 nM)7.5 ng	(1 nM)7.5 ng
Amount (conc.) of miRNA inhibitor	(25 nM)187.5 ng	(25 nM)187.5 ng	(25 nM)187.5 ng
Volume of HiPerFect Reagent	1.5 $\mu$ l	3 $\mu$ l	4.5 $\mu$ l

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## Attractene Transfection Reagent **New**

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产品	规格	货号
miRNeasy Mini Kit (50)	50 RNeasy Mini Spin Columns, Collection Tubes (1.5 ml and 2 ml), QIAzol Lysis Reagent, RNase-Free Reagents and Buffers	217004
miRNeasy FFPE Kit (50)	50 RNeasy MinElute Spin Columns, 50 gDNA Eliminator Spin Columns, Collection Tubes, Proteinase K, RNase-Free Reagents and Buffers	217404
miRNeasy 96 Kit (4)	4 RNeasy 96 Plates, Collection Microtubes (racked), Elution Microtubes CL, Caps, S-Blocks, AirPore Tape Sheets, QIAzol Lysis Reagent, RNase-Free Reagents and Buffers	217061
miScript Reverse Transcription Kit (10)	For 10 reactions: miScript Reverse Transcriptase Mix, miScript RT Buffer, RNase-Free Water	218060
miScript Reverse Transcription Kit (50)	For 50 reactions: miScript Reverse Transcriptase Mix, miScript RT Buffer, RNase-Free Water	218061
miScript Primer Assay (100)	10x miScript Primer Assay (contains one miRNA-specific primer)	Varies*
miScript SYBR Green PCR Kit (200)	For 200 reactions: QuantiTect SYBR Green PCR Master Mix, miScript Universal Primer	218073
miScript SYBR Green PCR Kit (1000)	For 1000 reactions: QuantiTect SYBR Green PCR Master Mix, miScript Universal Primer	218075
miScript miRNA Mimic (1 nmol)	1 nmol (for approx. 400 transfections in 24-well plates)* cell-culture-grade miRNA mimic	Varies*
miScript miRNA Mimic (5 nmol)	5 nmol (for approx. 2000 transfections in 24-well plates)* cell-culture-grade miRNA mimic	Varies*
miScript miRNA Inhibitor (1 nmol)	1 nmol (for approx. 40 transfections in 24-well plates)* cell-culture-grade miRNA inhibitor	Varies*
miScript miRNA Inhibitor (5 nmol)	5 nmol (for approx. 200 transfections in 24-well plates)* cell-culture-grade miRNA inhibitor	Varies*
<b>单独纯化miRNA组分所额外需要的试剂</b> RNeasy® MinElute® Cleanup Kit (50)	50 RNeasy MinElute Spin Columns, Collection Tubes (1.5 ml and 2 ml), RNase-Free Reagents and Buffers	74204
RNeasy 96 Kit (4) <sup>†</sup>	4 RNeasy 96 Plates, Elution Microtubes CL, Caps, S-Blocks, AirPore Tape Sheets, RNase-Free Reagents and Buffers	74181
<b>经优化的用于miRNA Mimics和miRNA Inhibitors的转染试剂</b> HiPerFect Transfection Reagent (1 ml)	HiPerFect Transfection Reagent for up to 333 transfections in 24-well plates or up to 1333 transfections in 96-well plates	301705
Attractene Transfection Reagent (1 ml)	Attractene Transfection Reagent for up to 660 transfections in 24-well plates	301005
<b>相关产品</b> QuantiTect Primer Assay(200)	10x QuantiTect Primer Assay(contains a mix of forward and reverse primers for a specific target)	Varies*
QuantiTect SYBR Green PCR Kit(200) <sup>†</sup>	For 200 X 50 µl reactions: 3 x 1.7ml 2x QuantiTect SYBR Green PCR Master Mix, 2 x 2 ml RNase-Free Water	204143

\*请登陆[www.qiagen.com/GeneGlobe](http://www.qiagen.com/GeneGlobe)查询所需miRNA引物。†更大包装规格及‡实验操作手册，请咨询QIAGEN

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