

Labware List — QIASymphony® DNA Investigator® Kit

Sample and eluate tubes/racks that can be used with the QIASymphony DNA Investigator Kit and the QIASymphony SP (software version 4.0; labware package SOW-516-6)

Revision management

This document is the QIASymphony DNA Investigator Kit Labware List for labware package SOW-516-6, Revision 1.

Important: Before using this Labware List, confirm that it is consistent with the labware package installed on your QIASymphony system.

Legend

- Recommended: These tubes are for use in combination with this protocol.
- User defined: User is responsible for defining and/or validating as the user deems appropriate.
- Not recommended: These tubes are not for use in combination with this protocol.
- x/■ Minimum sample volume (μ l) required per sample per protocol (denoted by x)/clot detection possible.
- x/□ Minimum sample volume (μ l) required per sample per protocol (denoted by x)/clot detection not possible.

Note: Be aware that other tubes are available that can be used in combination with other QIASymphony Kits, but are not recommended for any protocols that can be used with this kit.

“Sample” drawer, tube carrier

Note: Ensure that you remove swabs before using tubes on the QIASymphony SP.

Supplier	Material	Example cat. no.	Name in touchscreen	Insert	Protocol	
					Investigator*	Investigator HE (High Efficiency) [†]
Eppendorf	1.5 ml Eppendorf Safe-Lock Tube	0030 120.086	EP#0030120.086 T1.5 Snap Cap	Insert 5A		
Eppendorf	2.0 ml Eppendorf Safe-Lock Tube	0030 120.094	EP#0030120.094 T2.0 Snap Cap	Insert 5A		
QIAGEN	2 ml QIAGEN Collection Tube	19201	QIA#19201 Collection	Insert 03/3B		

Table continued on next page

Labware List — QIASymphony DNA Investigator Kit

“Sample” drawer, tube carrier (continued)

Note: Ensure that you remove swabs before using tubes on the QIASymphony SP.

Supplier	Material	Example cat. no.	Name in touchscreen	Insert	Protocol	
					Investigator*	Investigator HE (High Efficiency)†
QIAGEN	2 ml QIAGEN Collection Tube	19201	QIA#19201 Collection	Insert 03/3B		
Sarstedt	1.5 ml Snap Cap Micro tube	72.690	SAR#72.690 T1.5 Snap Cap	Insert 5A		
Sarstedt	2 ml Snap Cap Micro tube	72.691	SAR#72.691 T2.0 Snap Cap	Insert 5A		
Sarstedt	2 ml Micro tube, PP, NON-SKIRTED	72.693	SAR#72.693 T2.0 Screw	Insert 03/3B		
Sarstedt	2 ml Micro tube, PP, SKIRTED	72.694	SAR#72.694 T2.0 ScrewSkirt	Insert 3B		
Sarstedt	2 ml Snap Cap Safe Seal tube	72.695	SAR#72.695 T2.0 Snap Cap	Insert 5A		
Sarstedt	1.5 ml Snap Cap Safe Seal tube	72.706	SAR#72.706 T1.5 Snap Cap	Insert 5A		
Treff	ClickFit 2.0 ml PP tube	96.09329.9.01	TR#96.09329.9.01 T2.0 SnapCap	Insert 5A		

* Reference, Casework, and Casework ADV protocols.

† Casework HE and Casework ADV HE protocols.

Labware List — QIASymphony DNA Investigator Kit

"Sample" drawer, plate carrier

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Sample adapter	Protocol	
						Investigator*	Investigator HE (High Efficiency) [†]
ABgene	ABgene 2.2 ml Storage Plate, Mark II	0932	Deep Well	AB#0932 2.2mlSPMarkII	No adapter needed		
QIAGEN	96-well S-block with 2.2 ml wells	19585	Deep Well	QIA#19585 S-Block96	No adapter needed		
Sarstedt	1.5 ml Micro tube, PP, NON-SKIRTED	72.607	Tube 1.5 ml	SAR#72.607 T1.5 Screw	Adapter, tubes, 2 ml		
Sarstedt	2 ml Micro tube, PP, NON-SKIRTED	72.693	Tube 2.0 ml	SAR#72.693 T2.0 Screw	Adapter, tubes, 2 ml		
Sarstedt	2 ml Micro tube, PP, SKIRTED	72.694	Tube 2.0 ml	SAR#72.694 T2.0 ScrewSkirt	Adapter, tubes, 2 ml		

* Reference, Casework, and Casework ADV protocols.

[†] Casework HE and Casework ADV HE protocols.

Labware List — QIASymphony DNA Investigator Kit

"Eluate" drawer

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots [†]	Protocol	
						Investigator [§]	Investigator HE (High Efficiency) [¶]
ABgene	Thermo-Fast 96 non-skirted	AB-0600	PCR Plate	AB#0600 *PCR96	PCR Plate 96 QS (cooling slot 1) Elution slots 2–4 are not used in this procedure		
ABgene	Thermo-Fast 96 non-skirted, low-profile	AB-0700	PCR Plate	AB#0700 *PCR96 LowPro	PCR Plate 96 QS (cooling slot 1) Elution slots 2–4 are not used in this procedure		
ABgene	Thermo-Fast 96 non-skirted, low-profile	AB-0800	PCR Plate	AB#0800 *PCR96 SkirtLowPro	PCR Plate 96 QS (cooling slot 1)		
				AB#0800 PCR96 SkirtLowPro	No adapter needed (non-cooling slots 2–3)		
ABgene	ABgene 2D Storage Plate 1.2 ml	AB-1185	Deep Well	AB#AB-1185 *2.0ml2DPlate	SBS Universal QS V2 (cooling slot 1)		
				AB#AB-1185 2.0ml2DPlate	No adapter needed (non-cooling slots 2–3)		
ABgene	ABgene deep well plate 0.8 ml	AB-0765	Deep Well	AB#0765 *0.8ml StoragePlate	SBS Universal QS V2 (cooling slot 1) ^{††}		
				AB#0765 0.8ml Storage Plate	No adapter needed (non-cooling slots 2–3)		
Abbott	Abbott 96 Deep Well Plate	4J7130	Deep Well	ABT#4J7130 DeepWell96 RB	No adapter needed (non-cooling slots 2–3) Elution slot 1 is not used in this procedure		
Corning	Polystyrene Microplate, 96 well, flat bottom	3370	Micro Plate	CO#3370 *MTP96 FB	SBS Universal QS V2 (cooling slot 1)	†	†
				CO#3370 MTP96 FB	No adapter needed (non-cooling slots 2–3)		
Corning	Polystyrene Microplate, 96 well, round bottom	3795	Micro Plate	CO#3795 *MTP96 RB	SBS Universal QS V2 (cooling slot 1)	†	†
				CO#3795 MTP96 RB	No adapter needed (non-cooling slots 2–3)		

Table continued on next page

Labware List — QIASymphony DNA Investigator Kit

"Eluate" drawer (continued)

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots [†]	Protocol	
						Investigator [§]	Investigator HE (High Efficiency) [¶]
Corning	Polystyrene Microplate, 96 well, v-bottom	3897	Micro Plate	CO#3897 *MTP96 VB	SBS Universal QS V2 (cooling slot 1)	‡	‡
				CO#3897 MTP96 VB	No adapter needed (non-cooling slots 2–3)		
Eppendorf	Eppendorf Deepwell Block 96	951033405	Deep Well	EP#951033405 *Block96	SBS Universal QS V2 (cooling slot 1)	§	§
				EP#951033405 Block96	No adapter needed (non-cooling slots 2–3)		
Eppendorf	Eppendorf twin.tec PCR Plate 96, semi-skirted	0030 128.575	PCR Plate	EP#128.575 *PCR96 SemiSkirt	PCR Plate 96 QS (cooling slot 1) Elution slots 2–4 are not used in this procedure	§	¶
Eppendorf	Eppendorf twin.tec PCR Plate 96, skirted	128648	PCR Plate	EP#128.648 *PCR96 Skirt	PCR Plate 96 QS (cooling slot 1) Elution slots 2–4 are not used in this procedure	§	¶
Eppendorf	1.5 ml Eppendorf LoBind Snap Cap Safe-Lock tube	0030 108.051	Tube 1.5 ml	EP#0030108.051 *T1.5 Snap Cap	Micro Tube Snap Cap (cooling slot 1) ^{††}	§	§
				EP# 0030108.051 T1.5 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Eppendorf	1.5 ml Eppendorf Snap Cap tube	0030 125.150	Tube 1.5 ml	EP#0030125.150 *T1.5 Snap Cap	Micro Tube Snap Cap (cooling slot 1) ^{††}	§	§
				EP#0030125.150 T1.5 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Eppendorf	1.5 ml Eppendorf Snap Cap Safe-Lock tube	0030 120.086	Tube 1.5 ml	EP#0030120.086 *T1.5 Snap Cap	Micro Tube Snap Cap (cooling slot 1) ^{††}	§	§
				EP#0030120.086 T1.5 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Eppendorf	2.0 ml Eppendorf Snap Cap Safe-Lock tube	0030 120.094	Tube 2.0 ml	EP#0030120.094 *T2.0 Snap Cap	Micro Tube Snap Cap (cooling slot 1) ^{††}	§	§
				EP#0030120.094 T2.0 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		

Table continued on next page

Labware List — QIASymphony DNA Investigator Kit

"Eluate" drawer (continued)

Supplier	Material	Example no.	cat.	Category	Name in touchscreen	Adapter on Elution slots [†]	Protocol	
							Investigator [§]	Investigator HE (High Efficiency) [¶]
Greiner Bio-One	Polystyrene Microplate, 96 well, u-bottom	650161		Micro Plate	GR#650161 *MTP96 RB	96-Well Round Bottom QS (cooling slot 1)	‡	‡
					GR#650161 MTP96 RB	No adapter needed (non-cooling slots 2–3)		
Micronic	Micronic 1.4 TraXis tube with box	M52621		Deep Well	MIC#M52621 *1.4TraxisTubes	SBS Universal QS V2 (cooling slot 1)	§	¶
					MIC#M52621 1.4TraxisTubes	No adapter needed (non-cooling slots 2–3)		
NUNC	Nunc U96 DeepWell Plate1	260252		Deep Well	NU#260252 *DeepWell96 RB	SBS Universal QS V2 (cooling slot 1)	§	¶
					NU#260252 DeepWell96 RB	No adapter needed (non-cooling slots 2–3)		
QIAGEN	Elution Microtubes CL 96	19588		Deep Well	QIA#19588 *EMTR	Elution Microtube Rack QS (cooling slot 1)	‡	‡
					QIA#19588 EMTR	No adapter needed (non-cooling slots 2–3)		
QIAGEN	96-well S-block with 2.2 ml wells	19585		Deep Well	QIA#19585 *S-Block96	96-Well Round Bottom QS (cooling slot 1) ^{††}	‡	‡
					QIA#19585 S-Block96	No adapter needed (non-cooling slots 2–3)		
QIAGEN	MinElute 96 UF PCR Purification Plates	28053		Micro Plate	QIA#28053 *MiniElute 96PCR	SBS Universal QS V2 (cooling slot 1) ^{††}	§	¶
					QIA#28053 MiniElute 96PCR	No adapter needed (non-cooling slots 2–3)		
REMP Sample Management	REMP Plate 96 Tube 300	STBR96-300		Deep Well	RE#STBR96-300 *TubeRack300	SBS Universal QS V2 (cooling slot 1)	§	¶
					RE#STBR96-300 TubeRack300	No adapter needed (non-cooling slots 2–3)		
Sarstedt	1.5 ml Snap Cap tube	72.690		Tube 1.5 ml	SAR#72.690 *T1.5 Snap Cap	Micro Tube Snap Cap (cooling slot 1) ^{††}	§	¶
					SAR#72.690 T1.5 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		

Table continued on next page

Labware List — QIASymphony DNA Investigator Kit

"Eluate" drawer (continued)

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots†	Protocol	
						Investigator§	Investigator HE (High Efficiency)¶
Sarstedt	1.5 ml Snap Cap Safe-Lock tube	72.706	Tube 1.5 ml	SAR#72.706 *T1.5 Snap Cap	Micro Tube Snap Cap (cooling slot 1)††		
				SAR#72.706 T1.5 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Sarstedt	2.0 ml Snap Cap tube	72.691	Tube 2.0 ml	SAR#72.691 *T2.0 Snap Cap	Micro Tube Snap Cap (cooling slot 1)††		
				SAR#72.691 T2.0 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Sarstedt	2.0 ml Snap Cap Safe-Lock tube	72.695	Tube 2.0 ml	SAR#72.695 *T2.0 Snap Cap	Micro Tube Snap Cap (cooling slot 1)††		
				SAR#72.695 T2.0 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		
Sarstedt	2 ml Micro tube, PP, NON-SKIRTED	72693	Tube 2.0 ml / Tube_2.0ml AdapterV1 (no BC)	SAR#72.693 **T2.0 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.693 *T2.0 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.693 T2.0 Screw	24-Tube 1.5/2.0 ml QS (non-cooling slots 2–4)		
Sarstedt	1.5 ml Micro tube, PP, NON-SKIRTED	72607	Tube, 1.5 ml / Tube_1.5ml AdapterV1 (no BC)	SAR#72.607 **T1.5 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.607 *T1.5 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.607 T1.5 Screw	24-Tube 1.5/2.0 ml QS (non-cooling slots 2–4)		
Sarstedt	2 ml Micro tube, PP, SKIRTED	72694	Tube 2.0 ml / Tube_2.0ml AdapterV1 (no BC)	SAR#72.694 **T2.0 ScrewSkirt	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.694 *T2.0 ScrewSkirt	Micro Tube Screw Cap QS (cooling slot 1)		
				SAR#72.694 T2.0 ScrewSkirt	24-Tube 1.5/2.0 ml QS (non-cooling slots 2–4)		

Table continued on next page

Labware List — QIASymphony DNA Investigator Kit

"Eluate" drawer (continued)

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots [†]	Protocol	
						Investigator [§]	Investigator HE (High Efficiency) [¶]
Starlab	1.5 ml Graduated Conical Tube, Natural Standard Screw Cap	E1415-2231	Tube, 1.5 ml / Tube_1.5ml AdapterV1 (no BC)	SL#E1415-2231 **T1.5 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SL#E1415-2231 *T1.5 Screw	Micro Tube Screw Cap QS (cooling slot 1)		
				SL#E1415-2231 T1.5 Screw	24-Tube 1.5/2.0 ml QS (non-cooling slots 2–4)		
Thermo Scientific	Matrix 2D Barcoded ScrewTop Storage Tubes 500 µl	3744	Deep Well	TS#3744 *2DStorageTubes	SBS Universal QS V2 (cooling slot 1)		
				TS#3744 2DStorageTubes	No adapter needed (non-cooling slots 2–3)		
Thermo Scientific	Matrix 2D Barcoded Storage Tubes 1.4 ml, Flat Bottom Tubes	3711	Deep Well	TS#3711 *1.4ml 2DStorageT	SBS Universal QS V2 (cooling slot 1)		
				TS#3711 1.4ml 2DStorageT	No adapter needed (non-cooling slots 2–3)		
Thermo Scientific	Matrix Storage Tubes 1.4 ml, Round Bottom Tubes	4247	Deep Well	TS#4247 *1.4ml RBStorageT	SBS Universal QS V2 (cooling slot 1)		
				TS#4247 1.4ml RBStorageT	No adapter needed (non-cooling slots 2–3)		
Thermo Scientific	Matrix 2D Barcoded Storage Tubes 500 µl 2D Barcoded, V Bottom Tubes	3734	Deep Well	TS#3734 *0.5ml 2D VBStorageT	SBS Universal QS V2 (cooling slot 1) ^{††}		
				TS#3734 0.5ml 2D VBStorageT	No adapter needed (non-cooling slots 2–3)		
Treff	ClickFit 2.0 ml PP tube	96.09329.9.01	Tube 2.0 ml	TR#9609329901 T2.0 Snap Cap	Micro Tube Snap Cap (non-cooling slots 2–4)		

* Indicates labware that can be cooled using a cooling adapter with bar code.

** Indicates labware that can be cooled using a cooling adapter without bar code (non-transferable and not usable on QIASymphony AS).

† Do not use 96-well plates on "Elution slot 4" because the robotic arm cannot access all positions.

‡ A maximum elution volume of 200 µl can be used with this tube/plate format. If larger elution volumes are required, please select another type of tube/plate.

§ Reference, Casework, and Casework ADV protocols.

¶ Casework HE and Casework ADV HE protocols.

†† Labware that cannot be transferred to or used on QIASymphony AS.

Labware List — QIASymphony DNA Investigator Kit

For ordering information, visit www.qiagen.com/goto/QIASymphony.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIASymphony®, Investigator®, MinElute® (QIAGEN Group); Abbott (Abbott® Laboratories); Corning® (Corning, Inc.); Eppendorf®, eppendorf twin.tec® (Eppendorf AG); Bio-One® (Greiner Bio-One GmbH); Micronic™, TraXis® (Micronic Holding B.V.); Remp® (Remp AG); Sarstedt® (Sarstedt AG and Co.); Starlab® (Starlab Group); 2D Matrix®, ABgene®, DeepWell™, Nunc®, Thermo-Fast®, Thermo Scientific®, (Thermo Fisher Scientific or its subsidiaries). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Jan-16 HB-0986-L03-001 © 2012–2016 QIAGEN, all rights reserved.



Sample & Assay Technologies