September 2016

BLU-V[®] System User Manual

For neutralization of DNA of dead organisms

Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Contents

1. Introduction	5
1.1 Intended Use of the BLU-V system	5
1.2 General Information	5
2. Safety Information	7
2.1 Proper use	9
2.2 Electrical safety	10
2.3 Environment	13
2.4 Biological safety	14
2.5 Chemicals	17
2.6 Safety maintenance	19
2.7 Waste disposal	20
2.8 Symbols on the BLU-V System	21
3. General Description	23
3.1 Features of the BLU-V System	23
4. Installation	27
4.1 Unpack the instrument	27
4.2 Site requirements	27
4.3 Power cable connection	27
4.4 Grounding requirements	28
4.5 LED status lights	28
5. Operating Procedures	29
5.1 Preparation	29
5.2 Performing the run	31
6. Maintenance	33
6.1 Cleaning procedure	33
6.2 Replacing the main fuse	33
Appendix A	35
Technical Data	35
Waste Electrical and Electronic Equipment	36

Warranty	37
Appendix B: Declaration of Conformity	38
Appendix C: BLU-V System Accessories	38
Appendix D: Related Products	39

1. Introduction

Thank you for choosing the BLU-V System. We are confident that this instrument will become an integral part of your laboratory.

Before using the BLU-V System, it is essential that you read this user manual carefully. Following the instructions and safety information in this user manual will ensure safe operation and maintain the system in a safe condition.

1.1 Intended Use of the BLU-V system

The BLU-V System is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

The BLU-V System allows the neutralization of the DNA of dead organisms through the use of photo-activation on intercalating agents. Neutralizing the DNA of dead organisms allows the downstream detection of DNA of living cells by conventional PCR or real-time PCR.

The BLU-V System is intended for use by professional users trained in molecular biology techniques and the operation of the BLU-V System instrument.

1.2 General Information

1.2.1 Scope of delivery

The delivery includes the following items.

- BLU-V System
- BLU-V Sample Holder
- Power supply cable*

1.2.2 Technical assistance

At QIAGEN, we pride ourselves on the quality and availability of our technical support. Our Technical Services Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN products. If you have any questions or experience any

* For details, see Appendix C.

difficulties regarding the BLU-V System or QIAGEN products in general, do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists, as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance, contact one of the QIAGEN Technical Services Departments or local distributors (see back cover).

1.2.3 Policy statement

It is the policy of QIAGEN to improve products as new techniques and components become available. QIAGEN reserves the right to change the specifications of products at any time.

In an effort to produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact QIAGEN Technical Services with any feedback.

1.2.4 Requirements for BLU-V System users

Table 1 covers the general level of competence for the use and servicing of the BLU-V System.

Table 1. Level of competence required for the BLU-V System

Task	Personnel	Training and experience
Routine use	Laboratory technicians or equivalent	Training in molecular biology techniques for laboratory instrument operation
Servicing	QIAGEN Service Specialists only	Trained, certified, and authorized by QIAGEN

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<p>WARNING</p> 	<p>Electrical hazard [W2]</p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument unsafe.</p> <p>Intentional interruption is prohibited.</p> <p>Lethal voltages inside the instrument</p> <p>When the instrument is connected to line power, its terminals may be live. Opening covers or removing parts is likely to expose live parts.</p>
---	--

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<p>WARNING</p> 	<p style="text-align: right;">[W5]</p> <p>Samples containing infectious agents</p> <p>Some samples used with the BLU-V System may contain infectious agents. Handle such samples in accordance with the required safety regulations.</p> <p>The responsible person(s) (e.g., laboratory manager) must take the necessary precautions to ensure that the workplace is safe and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents, as defined in the applicable Safety Data Sheets (SDSs) or OSHA*, ACGIH†, or COSHH‡ documents.</p> <p>Venting of fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	--

*OSHA: Occupational Safety and Health Administration (United States of America).

† ACGIH: American Conference of Government Industrial Hygienists (United States of America).

‡ COSHH: Control of Substances Hazardous to Health (United Kingdom).

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*OSHA: Occupational Safety and Health Administration (United States of America).

† ACGIH: American Conference of Government Industrial Hygienists (United States of America).

‡ COSHH: Control of Substances Hazardous to Health (United Kingdom).

2.5 Chemicals

<p>WARNING</p> 	<p style="text-align: right;">[W6]</p> <p>Hazardous chemicals</p> <p>Some chemicals used with the BLU-V System may be hazardous or may become hazardous after completion of the protocol run.</p> <p>Always wear safety glasses, gloves, and a lab coat.</p> <p>The responsible person(s) (e.g., laboratory manager) must take the necessary precautions to ensure that the workplace is safe and that the instrument operators are suitably trained and not exposed to hazardous levels of toxic substances (chemical or biological), as defined in the applicable Safety Data Sheets (SDSs) or OSHA*, ACGIH†, or COSHH‡ documents.</p> <p>Venting of fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	--

*OSHA: Occupational Safety and Health Administration (United States of America).

† ACGIH: American Conference of Government Industrial Hygienists (United States of America).

‡ COSHH: Control of Substances Hazardous to Health (United Kingdom).

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
3. General Description

The BLU-V System is a thermally stable instrument that generates a constant and uniform light dosage to allow simultaneous photo-activation of 12 samples in a simple and efficient manner. Its design minimizes intra-and inter-laboratory variability.


3.1 Features of the BLU-V System

- | | |
|---|------------------------|
| 1 | Socket for power cable |
| 2 | On/off switch |
| 3 | Control panel |
| 4 | Protective shield |
| 5 | BLU-V Sample Holder |
| 6 | Air inlet fan |
| 7 | Air outlet slits |

3.1.1 Control panel

- 1 Numerical display
- 2 Time set button and LED
- 3 Light intensity button and LED
- 4 Start button and operation LED

3.1.2 Protective shield

The protective shield in the operating position.

Note: You may only perform a run when the protection shield is closed. A sensor detects whether the shield is closed or not.


3.1.3 Air inlet fan

The air inlet fan is located on the right side of the instrument.

3.1.4 Power supply

- 1 Socket for power cable
- 2 Fuse box
- 3 ON/OFF switch

Note: Details on changing the fuse are given in the “Maintenance” section.

3.1.5 Power cable

The BLU-V System is equipped with a power cable with a plug suitable for the destination country.

3.1.6 BLU-V Sample Holder


The BLU-V Sample Holder enables the insertion of 12 sample tubes, each with a tube volume of up to 2 ml into the BLU-V System instrument.

CAUTION 	<p>Damage to the instrument [C1]</p> <p>Never put the sample tubes directly into the holes in the BLU-V instrument. Always use the BLU-V Sample Holder to put sample tubes into the instrument.</p>
DE	<p>Beschädigung des Gerätes</p> <p>Setzen Sie nie die Probengefäße direkt in das BLU-V Gerät. Benutzen Sie immer den Probenhalter (BLU-V sample holder), um die Proben in das Gerät zu setzen.</p>
FR	<p>Détérioration de l'appareil</p> <p>Ne pas placer les échantillons directement dans emplacements pour tubes du BLU-V system. Toujours utiliser le support d'échantillon du BLU-V system prévu à cet effet pour placer les échantillons dans l'instrument.</p>

3.1.7 Using the BLU-V Incubation Box (optional)

The BLU-V Incubation Box is intended for light-protected incubation of reaction solutions containing photo-reactive reagents. The BLU-V Incubation Box and a second BLU-V Sample Holder for use with it must be ordered separately. For details, see Appendix B: BLU-V System Accessories.

Tubes in BLU-V Sample Holder within BLU-V Incubation Box

BLU-V Incubation Box with closed protective shield

The procedure for use of the BLU-V Incubation Box is as follows:

1. **Set up the reaction solutions in tubes with volumes up to 2 ml.**
2. **Open the BLU-V Incubation Box and insert the BLU-V Sample Holder.**
3. **Place the tubes into the BLU-V Sample Holder.**
4. **Close the protective shield and keep it closed during incubation.**
5. **When incubation is complete, the BLU-V Sample Holder containing the tubes can be transferred into the BLU-V System.**

4. Installation

4.1 Unpack the instrument

The packaging of the BLU-V System can be stored for reuse.

4.2 Site requirements

Place the BLU-V System instrument on a level, stable surface, ideally a firmly secured bench top, near an earthed/grounded electrical outlet. Unobstructed access to the electrical outlet to which the instrument is connected must be ensured at all times. Allow at least 10 cm (4 in) of clearance on all sides of the unit for proper ventilation. While operating, make sure that the instrument does not touch other objects. Locate it away from vibration-sensitive instruments, such as an analytical balance.

4.3 Power cable connection

The socket for connecting the power cable is on the left of the BLU-V System instrument.

1. **Connect the supplied power cord to the socket on the instrument.**
2. **Plug the instrument into the electrical outlet.**
3. **Switch the instrument ON.**
4. **The green LED beside the start button on the control panel flashes green when the instrument is ON.**

When the BLU-V System is not in use for a longer period of time, we recommend disconnecting the power cable.

4.4 Grounding requirements


To protect operating personnel, the instrument is equipped with a 3-conductor AC power cord. To preserve this protection feature, do not operate the instrument from an AC power outlet that has no ground (earth) connection.

For replacement of the 3-conductor AC power cord, contact QIAGEN Technical Services for an authorized spare part.

4.5 LED status lights

There are three light emitting diode (LED) status lights on the control panel.

- 1 This LED is on and red when the display shows the selected time, and this setting can be changed by pushing the time set button.
- 2 This LED is on and red when the display shows the selected light intensity and this setting can be changed by pushing the light intensity button
- 3 This LED flashes green when the instrument is switched ON. It changes to red and stops flashing when the start button is pressed

5. Operating Procedures

5.1 Preparation

Before switching on the instrument:

1. **Prepare the reaction solutions.**
2. **Insert the reaction tubes into the BLU-V Sample Holder.**
3. **Place the BLU-V Sample Holder with the sample tubes into the BLU-V System.**

Note: If the BLU-V Incubation Box with the BLU-V Sample Holder has been used, this holder can be directly transferred into the instrument.

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5.1.1 Setting the time


- The display shows the selected time in minutes.
- The default setting of the timer when the instrument is first used is 15 minutes.
- The time setting can be changed by pushing the time set button until the required value is reached. It can be pushed step by step or continuously. The time increments are 1 minute each.
- The maximum possible time is 60 minutes. The time setting is a cyclic push button — after 60, it begins again with 1.
- Adjust the time settings based on the recommendations given in the corresponding BLU-V Kit handbooks.

- When the start button is pressed, the entered time disappears and the count starts at zero.

Note: The setting of the last completed run is saved as the new default.

5.1.2 Setting the light intensity

- Press the light intensity button to set the internal LED light intensity.
- The default light intensity is 100%.
- The light intensity setting can be changed by pushing the light intensity button until the required value is reached. The increments are 10% each.
- The light intensity button is a cyclic push button — after 100%, it begins again with 10%.
- Adjust the light intensity settings based on the recommendations given in the corresponding BLU-V Kit handbooks. For standard illumination settings keep the light power at 100%.

Note: The setting of the last completed run is saved as the new default.

5.1.3 Using the protective shield

- Close the protective shield before starting the run.
- Keep the protective shield closed until the run is completed.

If the protective shield is opened during a run, the instrument will interrupt automatically and the photo-activation LEDs will switch off. The time will also stop. The run and the time will only continue when the protective shield is closed again.

5.2 Performing the run

1. To start the run, press the start button. The color of the indicator LED will change from green to red. The time display will show the elapsed time from the beginning of the operation in minutes.

- 2. Throughout the run, the indicator LED will be red to show that the internal LEDs are working.**
- 3. When the run time is complete, the internal LEDs are shut off, but the indicator LED remains red and the fan stays on to cool the instrument.**
- 4. After 30 seconds the indicator LED will change from red to flashing green and the set time is displayed again, open the protective shield and remove the BLU-V Sample Holder with the sample tubes.**
- 5. Keep the instrument on until the fan stops. This helps to cool the instrument. It takes around 30 seconds.**
- 6. The instrument is ready for a new set of samples when the fan stops.**

6. Maintenance

6.1 Cleaning procedure

Disconnect the BLU-V System from the power supply for cleaning. Use a lightly moistened cloth or an antistatic cloth. Do not use cleaning agents or abrasive solvent products.

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
6.2 Replacing the main fuse

To replace the main fuse, follow this procedure.

A

B

1. Switch the instrument OFF.

- 2. Remove the power cord from the socket.**
- 3. Remove the fuse holder (indicated by the red arrow A).**
- 4. Remove the fuse (indicated by the red arrow B) from the fuse holder.**
- 5. Insert the new fuse and insert the fuse holder into the instrument.**

Specifications of the fuse are given in Appendix A "Technical data".

Appendix A

Technical Data

QIAGEN reserves the right to change specifications at any time.

Environmental conditions

Operating conditions

Power range	100–240 V AC Mains supply voltage fluctuations must not exceed +/- 10% of nominal supply voltages
Frequency range	50/60 Hz
Fuse	T2A/250V Time-Delay
Maximum power	60 W
Overvoltage category	II
Air temperature	15–32°C (59–89.6°F)
Relative humidity	15–75% (noncondensing) Maximum of 75% relative humidity for temperatures up to 31°C (88°F), decreasing linearly to 50% humidity at 32°C (89.6°F)
Altitude	Up to 2000 m (6500 ft.)
Place of operation	For indoor use only
Pollution degree	2

Transportation conditions

Air temperature	–25°C to 60°C (–13°F to 140°F) in manufacturer's packaging
Relative humidity	Maximum of 75% (noncondensing)

Storage conditions

Air temperature	5°C to 40°C (41°F to 104°F) in manufacturer's packaging
Relative humidity	Maximum of 85% (noncondensing)

Dimensions and Weight

Dimensions	Width: 355 mm (14 in.)
	Height: 83 mm (3.3 in.)
	Depth: 90 mm (3.5 in.)
Weight	1.2 kg (2.65 lb.)

Waste Electrical and Electronic Equipment

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

Recycling can be provided by QIAGEN upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by QIAGEN, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact your local QIAGEN sales office for the required return form. Once the form is submitted, you will be contacted by QIAGEN either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

FCC declaration

The "United States Federal Communications Commission" (USFCC) (in 47 CFR 15. 105) declared that the users of this product must be informed of the following facts and circumstances.

"This device complies with part 15 of the FCC:

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation."

"This Class B digital apparatus complies with Canadian ICES-0003."

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the accompanying documentation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES 003 for digital apparatus. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/T.V. technician for help.

QIAGEN GmbH Germany is not responsible for any radio television interference caused by unauthorized modifications of this equipment or the substitution or attachment of connection cables and equipment other than those specified by QIAGEN GmbH, Germany. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

Warranty

The BLU-V System instrument has a 24-month warranty period which starts when the final customer purchases it.

The warranty is void when misuse of the equipment can be proved.

Damage or manufacturing faults caused by impacts, chemical or corrosive products, liquids, damp, or other external factors, such as radiation, fire, or inadequate transport are not included.

In addition, the warranty will not apply if the equipment has been handled, repaired, or modified by not qualified or specifically designated personnel.

Appendix B: Declaration of Conformity

Name and address of the legal manufacturer

QIAGEN GmbH
QIAGEN Strasse 1
40724 Hilden
Germany

An up-to-date Declaration of Conformity can be requested from QIAGEN Technical Services.

Appendix C: BLU-V System Accessories

The following accessories are not included in the delivery. State the order number (cat. no.) when ordering the accessory.

Product	Contents	Cat. no.
BLU-V Incubation Box	Box for up to 12 tubes (tube volumes up to 2 ml) for light protection during incubation	9022908
Fuse T2A, 250V	Set of fuses	9240628
Power cord	Power cord for use in Australia	9232548
Power cord	Power cord for use in Switzerland	9232546
Power cord	Power cord for use in India	9243305

The following accessory is included in the delivery and can be ordered separately. State the order number (cat. no.) when ordering the accessory.

Product	Contents	Cat. no.
BLU-V Sample Holder	Sample holder for up to 12 tubes (tube volumes up to 2 ml)	9022909

Appendix D: Related Products

The following related products are not included in the delivery. State the order number (cat. no.) when ordering the accessory.

Product	Contents	Cat. no.
BLU-V Viability PMA Kit	PMA Reagent 2 x 0.7 mg, EB Buffer, RNase-Free water	296015

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services.

Trademarks: QIAGEN[®], Sample to Insight[®], BLU-V[®] (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-1761-006 © 2014–2016 QIAGEN, all rights reserved.

Notes

Notes

Notes

www.qiagen.com

Australia ■ techservice-au@qiagen.com

Austria ■ techservice-at@qiagen.com

Belgium ■ techservice-bnl@qiagen.com

Brazil ■ suportetecnico.brasil@qiagen.com

Canada ■ techservice-ca@qiagen.com

China ■ techservice-cn@qiagen.com

Denmark ■ techservice-nordic@qiagen.com

Finland ■ techservice-nordic@qiagen.com

France ■ techservice-fr@qiagen.com

Germany ■ techservice-de@qiagen.com

Hong Kong ■ techservice-hk@qiagen.com

India ■ techservice-india@qiagen.com

Ireland ■ techservice-uk@qiagen.com

Italy ■ techservice-it@qiagen.com

Japan ■ techservice-jp@qiagen.com

Korea (South) ■ techservice-kr@qiagen.com

Luxembourg ■ techservice-bnl@qiagen.com

Mexico ■ techservice-mx@qiagen.com

The Netherlands ■ techservice-bnl@qiagen.com

Norway ■ techservice-nordic@qiagen.com

Singapore ■ techservice-sg@qiagen.com

Sweden ■ techservice-nordic@qiagen.com

Switzerland ■ techservice-ch@qiagen.com

UK ■ techservice-uk@qiagen.com

USA ■ techservice-us@qiagen.com