

QuantiNova[®] Multiplex RT-PCR Kit

The QuantiNova Multiplex RT-PCR Kit (cat. nos. 208552, 208554, 208556) should be stored immediately upon receipt at -15 to -30°C in a constant-temperature freezer and protected from light. QuantiNova Multiplex RT-PCR Master Mix can also be stored protected from light at 2 – 8°C for up to 12 months, depending on the expiration date.

Further information

- *QuantiNova Multiplex RT-PCR Kit Handbook*: www.qiagen.com/HB-2313
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- This protocol is optimized for quantification of RNA targets in a multiplex format (up to 5-plex), using TaqMan[®] probes with any real-time cyclers and condition for fluorescence normalization. ROX[™] dye is required for various cyclers at the following concentrations:
No requirement for ROX dye: Rotor-Gene[®], Bio-Rad[®] CFX, Roche[®] LightCycler[®] 480 and Agilent[®] Technologies Mx instruments.
Low concentration of ROX dye: Applied Biosystems[®] 7500, ViiA[®]7 and QuantStudio[™] Real-Time PCR Systems.
High concentration of ROX dye: ABI PRISM[®] 7000, Applied Biosystems 7300, 7900 and StepOne[™] Real-Time PCR Systems.
- QuantiNova ROX Reference Dye should be used as a 20x concentrated solution for a 1x reaction when using an instrument requiring a high ROX dye concentration. For instruments requiring a low ROX dye concentration, use the dye as a 200x concentrate.
- The reference dye in QuantiNova Yellow Template Dilution Buffer allows tracking of pipetted samples in the qRT-PCR. When template is added to the blue QuantiNova Multiplex RT-PCR Master Mix, the color changes from blue to green.

It is provided as a 100x concentrate and should be diluted (using water) to obtain a 1x final concentration within the sample*. The buffer does not affect sample stability or qPCR.

- For the highest efficiency in real-time RT-PCR using TaqMan probes, amplicons should ideally be 60–150 bp in length.
- Before performing multiplex analysis, choose suitable combinations of reporter dyes and quenchers that are compatible with multiplex analysis using the detection optics of your real-time cycler. We strongly recommend using dual-labeled probes with non-fluorescent quenchers.
- Always start with the cycling conditions and primer concentrations specified in this protocol.
- The RT-PCR must start with an initial reverse transcription step of 10 min at 50°C to transcribe the RNA into cDNA, followed by an incubation step of 2 min at 95°C to activate the QuantiNova DNA Polymerase and to inactivate the HotStaRT-Script Reverse Transcriptase.
- For ease of use, we recommend preparing a 20x primer–probe mix consisting of 16 μM forward primer, 16 μM reverse primer and 5 μM probe in TE buffer for each target.
- For information on suitable combinations of dyes for multiplex PCR using various cyclers, please refer to the *QuantiNova Multiplex RT-PCR Handbook*.
- The QuantiNova Internal Control RNA (QN IC RNA) is optionally used to test successful reverse transcription/amplification and to report instrument or chemistry failures, or the presence of inhibitors. It is detected as a 200 bp internal control (IC) in the yellow (Rotor-Gene Q) or in the VIC®/HEX dye channel (other real-time PCR instruments) when using the QuantiNova IC Probe Assay (cat. no. 205813). Alternatively, it is detected in the red (Rotor-Gene Q) or in the Cy®5 dye channel (other real-time PCR instruments) when using the QuantiNova IC Probe Assay Red 650 (cat. no. 205824). Before use, add 180 μl (or 900 μl) of RNase-free water to 20 μl (or 100 μl) of QN IC RNA provided in the kit and mix thoroughly by vortexing.

* Example: Add 0.5 μl Yellow Template Dilution Buffer to a 50 μl sample, which can be used as template in various PCR runs, regardless of the volume added to each reaction. Yellow Template Dilution Buffer can be pre-diluted using water. In this example, add 5 μl of 1:10 pre-diluted Yellow Template Dilution Buffer.

1. Thaw 4x QuantiNova Multiplex RT-PCR Master Mix, QuantiNova Yellow Template Dilution Buffer, template RNA, primers, probes, QN ROX Reference Dye (if required), and RNase-free water. Mix the individual solutions. Put QuantiNova Multiplex Reverse Transcription Mix on ice.
2. Prepare a reaction mix according to Table 1. Due to the hot start, it is not necessary to keep samples on ice during reaction setup or while programming the real-time cycler.

Table 1. Reaction mix setup

| Component | 96-well block, Rotor-Gene | 384-well block | Final concentration |
|---|------------------------------|-----------------------------|---|
| 4x Multiplex RT-PCR Master Mix | 5 μ l | 2.5 μ l | 1x |
| QN ROX Reference Dye (AB instruments only) | 1 μ l/0.1 μ l* | 0.5 μ l/0.05 μ l* | 1x |
| 100x Multiplex Reverse Transcription Mix | 0.2 μ l | 0.1 μ l | 1x |
| 20x primer–probe mix [†] (for each of up to 5 targets) | 1 μ l | 0.5 μ l | 0.8 μ M forward primer 0.8 μ M reverse primer 0.25 μ M TaqMan probe |
| QuantiNova IC Probe Assay 10x (optional) | 2 μ l | 1 μ l | 1x |
| QN IC RNA (optional) | 1 μ l | 1 μ l | 1x |
| RNase-Free Water | Variable | Variable | – |
| Template RNA (added at step 4) | Variable | Variable | \leq 800 ng/reaction |
| Total reaction volume | 20 μl | 10 μl | – |

*Results in a 1:20 dilution for high ROX dye cyclers (i.e., ABI PRISM 7000, Applied Biosystems 7300, 7900 and StepOne Real-Time PCR Systems) and a 1:200 dilution for low-ROX dye cyclers (i.e., Applied Biosystems 7500 and ViiA7 Real-Time PCR Systems) in the final 1x reaction.

[†] A 20x primer–probe mix consists of 16 μ M forward primer, 16 μ M reverse primer, and 5 μ M probe in TE buffer for each target. Primers and probes can either be pre-mixed and added simultaneously, or primer-probe mixes for each target can be added separately. If concentration of primer-probe mix(es) differ, the respective added volume needs to be adjusted, to achieve a final concentration of 0.8 μ M for each primer and 0.25 μ M for each probe.

3. Mix the reaction thoroughly and dispense appropriate volumes into PCR tubes or wells of a PCR plate.
4. Add template RNA (\leq 800 ng/reaction) to the individual PCR tubes or wells containing the reaction mix.
5. Program the real-time cycler according to Table 2.

Note: Data acquisition should be performed during the combined annealing/extension step.

Table 2. Cycling conditions

| Step | Time | Temperature | Ramp rate |
|------------------------------|------------|-------------|-------------------|
| Reverse transcription step | 10 min | 50°C | Maximal/fast mode |
| PCR initial heat activation | 2 min | 95°C | Maximal/fast mode |
| 2-step cycling | | | |
| Denaturation | 5 s | 95°C | Maximal/fast mode |
| Combined annealing/extension | 30 s | 60°C | Maximal/fast mode |
| Number of cycles | 40* | | |

* The number of cycles depends on the amount of template RNA.

6. Place the PCR tubes or plates in the real-time cycler and start the cycling program.
7. Perform data analysis.

Note: Before performing data analysis, select the analysis settings for each probe (i.e., baseline settings and threshold values). Optimal analysis settings are a prerequisite for accurate quantification data.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, Sample to Insight®, QuantiNova®, Rotor-Gene® (QIAGEN Group); Agilent® (Agilent Technologies, Inc.); ABI PRISM®, Applied Biosystems®, QuantStudio™, StepOne™, VIC®, ViiA® (Thermo Fisher Scientific or its subsidiaries); Bio-Rad® (Bio-Rad Laboratories, Inc.); LightCycler®, Roche®, TaqMan® (Roche Group); SmartCycler® (Cepheid); Cy® (GE Healthcare). 1101759 12/2016 HB-2208-001 © 2016 QIAGEN, all rights reserved.