

# HotStarTaq<sup>®</sup> DNA Polymerase

HotStarTaq DNA Polymerase (cat. nos. 203203, 203205, 203207 and 203209), including buffers and reagents, should be stored immediately upon receipt at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

## Further information

- *HotStarTaq PCR Handbook*: [www.qiagen.com/HB-0452](http://www.qiagen.com/HB-0452)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- HotStarTaq DNA Polymerase requires a heat-activation step of 15 min at  $95^{\circ}\text{C}$  (see step 5).
- The PCR Buffer provides a final concentration of 1.5 mM  $\text{MgCl}_2$  in the reaction mix, which will give satisfactory results in most cases. However, in some cases, reactions may be improved by increasing the final  $\text{Mg}^{2+}$  concentration. If a higher  $\text{Mg}^{2+}$  concentration is required, add the appropriate volume of 25 mM  $\text{MgCl}_2$  to the reaction mix as described in the *HotStarTaq PCR Handbook*.
- If required, prepare a dNTP mix containing 10 mM of each dNTP. Store this mix in aliquots at  $-20^{\circ}\text{C}$ . High-quality, PCR-grade dNTP mix (10 mM) is available from QIAGEN (cat. no. 201900).
- HotStarTaq DNA Polymerase is provided with Q-Solution<sup>®</sup>, which facilitates amplification of templates that have a high degree of secondary structure or that are GC-rich by modifying the melting behavior of DNA. When using Q-Solution for the first time for a particular primer–template pair, always perform parallel reactions with and without Q-Solution.

- It is not necessary to keep PCR tubes on ice as nonspecific DNA synthesis cannot occur at room temperature due to the inactive state of HotStarTaq DNA Polymerase.
  - A No Template Control (NTC) should always be included.
1. Thaw 10x PCR Buffer, dNTP mix, primer solutions, Q-Solution (if required) and 25 mM MgCl<sub>2</sub> (if required). Mix thoroughly before use to avoid localized differences in salt concentration.
  2. Prepare a reaction mix according to Table 1. The reaction mix typically contains all the components needed for PCR except the template DNA. Prepare a volume of reaction mix 10% greater than that required for the total number of PCR assays to be performed.

**Table 1. Reaction setup using HotStarTaq DNA Polymerase**

Component	Volume/reaction	Final concentration
<b>Reaction mix</b>		
10x PCR Buffer*	10 µl	1x
<b>Optional:</b> 5x Q-Solution <sup>†</sup>	20 µl	1x
dNTP mix (10 mM of each)	2 µl	200 µM of each dNTP
Primer A	Variable	0.1–0.5 µM
Primer B	Variable	0.1–0.5 µM
HotStarTaq DNA Polymerase	0.5 µl	2.5 units/reaction
Distilled water	Variable	–
<b>Template DNA</b> (added at step 4)	Variable	≤1 µg/reaction
<b>Total reaction volume</b>	100 µl <sup>‡</sup>	

\* Contains 15 mM MgCl<sub>2</sub>.

<sup>†</sup> For templates with GC-rich regions or complex secondary structure.

<sup>‡</sup> If using different reaction volumes, adjust the volume of each component accordingly.

3. Mix the reaction mix gently but thoroughly, for example, by pipetting up and down a few times. Dispense appropriate volumes into PCR tubes.
4. Add template DNA ( $\leq 1 \mu\text{g}/100 \mu\text{l}$  reaction) to the individual PCR tubes containing the reaction mix. For RT-PCR, add an aliquot from the reverse transcriptase reaction. This should not exceed 10% of the final PCR volume.
5. Program the thermal cycler according to the manufacturer's instructions.

**Note:** Each PCR program must start with an initial heat-activation step at 95°C for 15 min. A typical PCR cycling program is outlined in Table 2. For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.

6. Place the PCR tubes in the thermal cycler and start the cycling program.

**Note:** After amplification, samples can be stored overnight at 2–8°C, or at –20°C for longer storage.

**Table 2. Optimized cycling conditions**

Step	Time	Temperature	Comment
<b>Initial heat activation</b>	15 min	95°C	Activates HotStarTaq DNA Polymerase.
<b>3-step cycling:</b>			
Denaturation	0.5–1 min	94°C	
Annealing	0.5–1 min	50–68°C	Approximately 5°C below $T_m$ of primers.
Extension	1 min	72°C	For PCR products longer than 1 kb, use an extension time of approximately 1 min per kb DNA.
Number of cycles	25–35		
<b>Final extension</b>	10 min	72°C	



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, HotStarTaq®, Q-Solution® (QIAGEN Group). 1101226 03/2016 HB-0654-002 © 2016 QIAGEN, all rights reserved.