

AdnaTest BreastCancerSelect

AdnaTest BreastCancerSelect must be stored at 2–8°C. The components must not be used beyond the expiration date.

Further information

- *AdnaTest BreastCancerSelect, BreastCancerDetect and ER/PR-Detect Handbook:*
www.qiagen.com/HB-2453
- Safety Data Sheets: **www.qiagen.com/safety**
- Technical assistance: toll-free 00800-22-44-6000 or **support.qiagen.com**

Protocol

1. Resuspend BreastSelect Beads thoroughly and transfer 100 µl for each blood sample into a 1.5 ml reaction tube.
2. Wash BreastSelect Beads with 3 x 1 ml PBS.
3. Resuspend BreastSelect Beads in 100 µl PBS per blood sample.
4. When using standard EDTA tubes, transfer 5 ml whole blood into a 15 ml Collection Tube (provided).
When using ACD-A blood in a BD Vacutainer® ACD-A Tube, transfer 5 ml of blood into an AdnaTube.
5. Add 100 µl of washed BreastSelect Beads to each blood sample.
6. Incubate for 30 min at room temperature under tilting and rotation at approximately 5 rpm.
7. Place tube for 3 min in the AdnaMag-L to separate the beads.
Release any blood drops captured in the cap by swinging the AdnaMag-L downwards.

-
8. Remove blood supernatant.
 9. Wash beads with 3 x 5 ml PBS.
 10. Resuspend beads in 1 ml PBS and transfer into a new 1.5 ml reaction tube (not provided).
 11. Separate beads in the AdnaMag-S and remove supernatant.
 12. Resuspend beads in 200 µl AdnaTest Lysis/Binding Buffer by pipetting at least five times.
 13. Place reaction tubes into the AdnaMag-S and transfer supernatant into a new 1.5 ml reaction tube (provided).
 14. Continue immediately with the AdnaTest BreastCancerDetect or store at -20°C for a maximum of 2 weeks.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN[®], Sample to Insight[®] (QIAGEN Group); Vacutainer[®] (Becton Dickinson and Company). 1106856 03/2017 HB-2454-001
© 2017 QIAGEN, all rights reserved.