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# QuantiFERON<sup>®</sup>-TB Gold Plus (QFT<sup>®</sup>-Plus) 仿單

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1. Andersen, P. et al. (2000) Specific immune-based diagnosis of tuberculosis. *Lancet* 356, 1099.
2. Balcells, M.E. et al. (2008) A comparative study of two different methods for the detection of latent tuberculosis in HIV-positive individuals in Chile. *Int. J. Infect. Dis.* 12, 645.
3. Bartalesi, F. et al. (2009) QuantiFERON-TB Gold and TST are both useful for latent TB screening in autoimmune diseases. *Eur. Respir. J.* 33, 586.
4. Bocchino, M. et al. (2008) Performance of two commercial blood IFN-gamma release assays for the detection of *Mycobacterium tuberculosis* infection in patient candidates for anti-TNF-alpha treatment. *Eur. J. Clin. Microbiol. Infect. Dis.* 27, 907.
5. Brock, I. et al. (2006) Latent tuberculosis in HIV positive, diagnosed by the *M. tuberculosis* specific interferon-gamma test. *Respir. Res.* 7, 56.
6. Chun, J.K. et al. (2008) The role of a whole blood interferon gamma assay for the detection of latent tuberculosis infection in Bacille Calmette-Guerin vaccinated children. *Diagn. Microbiol. Infect. Dis.* 62, 389.
7. Connell, T.G. et al. (2008) A three-way comparison of tuberculin skin testing, QuantiFERON-TB Gold and T-SPOT.TB in children. *PLoS ONE* 3, e2624.

8. Detjen, A.K. et al. (2007) Interferon-gamma release assays improve the diagnosis of tuberculosis and nontuberculous mycobacterial disease in children in a country with a low incidence of tuberculosis. *Clin. Infect. Dis.* 45, 322.
9. Diel, R. et al. (2009) Comparative performance of tuberculin skin test, QuantiFERON-TB Gold In Tube assay, and T-Spot.TB test in contact investigations for tuberculosis. *Chest* 135, 1010.
10. Diel, R. et al. (2008) Predictive value of a whole-blood IFN- $\gamma$  assay for the development of active TB disease. *Am. J. Respir. Crit. Care Med.* 177, 1164.
11. Diel, R. et al. (2006) Tuberculosis contact investigation with a new, specific blood test in a low-incidence population containing a high proportion of BCG-vaccinated persons. *Respir. Res.* 7, 77.
12. Dogra, S. et al. (2007) Comparison of a whole blood interferon-gamma assay with tuberculin skin testing for the detection of tuberculosis infection in hospitalized children in rural India. *J. Infect.* 54, 267.
13. Drobniowski, F. et al. (2007) Rates of latent tuberculosis in health care staff in Russia. *PLoS Med.* 4, e55.
14. Gerogianni, I. et al. (2008) Whole-blood interferon-gamma assay for the diagnosis of tuberculosis infection in an unselected Greek population. *Respirology* 13, 270.
15. Harada, N. et al. (2008) Comparison of the sensitivity and specificity of two whole blood interferon-gamma assays for *M. tuberculosis* infection. *J. Infect.* 56, 348.
16. Higuchi, K. et al. (2009) Comparison of performance in two diagnostic methods for tuberculosis infection. *Med. Microbiol. Immunol.* 198, 33.

17. Kang, Y.A. et al. (2005) Discrepancy between the tuberculin skin test and the whole-blood interferon gamma assay for the diagnosis of latent tuberculosis infection in an intermediate tuberculosis-burden country. *JAMA* 293, 2756.
18. Katiyar, S. K. et al. (2008) Use of the QuantiFERON-TB Gold In-Tube test to monitor treatment efficacy in active pulmonary tuberculosis. *Int. J. Tuberc. Lung Dis.* 12, 1146.
19. Kipfer, B. et al. (2008) Tuberculosis in a Swiss army training camp: contact investigation using an Interferon gamma release assay. *Swiss Med. Wkly.* 138, 267.
20. Luetkemeyer, A. et al. (2007) Comparison of an interferon-gamma release assay with tuberculin skin testing in HIV-infected individuals. *Am. J. Respir. Crit. Care Med.* 175, 737.
21. Mackensen, F. et al. (2008) QuantiFERON-TB Gold – a new test strengthening long-suspected tuberculous involvement in serpiginous-like choroiditis. *Am. J. Ophthalmol.* 146, 761.
22. Manuel, O. et al. (2007) Comparison of QuantiFERON-TB Gold with tuberculin skin test for detecting latent tuberculosis infection prior to liver transplantation. *Am. J. Transplant.* 7, 2797.
23. Matulis, G. et al. (2007) Detection of latent tuberculosis in immunosuppressed patients with autoimmune diseases performance of a Mycobacterium tuberculosis antigen specific IFN-gamma assay. *Ann. Rheum. Dis.* 67, 84.
24. Mirtskhulava, V. et al. (2008) Prevalence and risk factors for latent tuberculosis infection among health care workers in Georgia. *Int. J. Tuberc. Lung Dis.* 12, 513.
25. Nakaoka, H. et al. (2006) Risk for tuberculosis among children. *Emerg. Infect. Dis.* 12, 1383.

26. Pai, M. et al. (2005) *Mycobacterium tuberculosis* infection in health care workers in rural India: comparison of a whole-blood, interferon-gamma assay with tuberculin skin testing. *JAMA* 293, 2746.
27. Ponce de Leon, D. et al. (2008) Comparison of an interferon-gamma assay with tuberculin skin testing for detection of tuberculosis (TB) infection in patients with rheumatoid arthritis in a TB-endemic population. *J. Rheumatol.* 35, 776.
28. Richeldi, L. et al. (2008) Prior tuberculin skin testing does not boost QuantiFERON-TB results in paediatric contacts. *Eur. Respir. J.* 32, 524.
29. Rothel, J.S., Andersen, P. (2005) Diagnosis of latent *Mycobacterium tuberculosis* infection: Is the demise of the Mantoux test imminent? *Expert Rev. Anti Infect. Ther.* 3, 981.
30. Schoepfer, A.M. et al. (2008) Comparison of interferon-gamma release assay versus tuberculin skin test for tuberculosis screening in inflammatory bowel disease. *Am. J. Gastroenterol.* 103, 2799.
31. Silverman, M.S. et al. (2007) Use of an interferon-gamma based assay to assess bladder cancer patients treated with intravesical BCG and exposed to tuberculosis. *Clin. Biochem.* 40, 913.
32. Stebler, A. et al. (2008) Whole-blood interferon-gamma release assay for baseline tuberculosis screening of healthcare workers at a Swiss university hospital. *Infect. Control Hosp. Epidemiol.* 29, 681.
33. Turner, J. et al. (1996) Stimulation of human peripheral blood mononuclear cells with live *Mycobacterium bovis* BCG activates cytolytic CD8<sup>+</sup> T cells in vitro. *Immunology* 87, 339.

- 
34. Brookes, R.H. et al. (2003) CD8+ T cell-mediated suppression of intracellular *Mycobacterium tuberculosis* growth in activated human macrophages. *Eur. J. Immunol.* 33, 3293.
  35. Stenger, S. et al. (1998) An antimicrobial activity of cytolytic T cells mediated by granulysin. *Science* 282, 121.
  36. Lalvani, A. et al. (1998) Human cytolytic and interferon gamma-secreting CD8+ T lymphocytes specific for *Mycobacterium tuberculosis*. *Proc. Natl. Acad. Sci. U.S.A.* 95, 270.
  37. Lewinsohn, D.M. et al. (2001) Classically restricted human CD8+ T lymphocytes derived from *Mycobacterium tuberculosis*-infected cells: definition of antigenic specificity. *J. Immunol.* 166, 439.
  38. Lewinsohn, D.A. et al. (2007) Immunodominant tuberculosis CD8 antigens preferentially restricted by HLA-B. *PLoS Pathol.* 3, 1240.
  39. Barcellini, L. et al. (2016) First independent evaluation of QuantiFERON-TB Plus performance. *Eur. Respir. J.* 47, 1587.
  40. Day, C.L. et al. (2011) Functional capacity of *Mycobacterium tuberculosis*-specific T cell responses in humans is associated with mycobacterial load. *J. Immunol.* 187, 2222.
  41. Rozot, V. et al. (2013) *Mycobacterium tuberculosis*-specific CD8+ T cells are functionally and phenotypically different between latent infection and active disease. *Eur. J. Immunol.* 43, 1568.

42. Nikolova, M. et al. (2013) Antigen-specific CD4- and CD8-positive signatures in different phases of *Mycobacterium tuberculosis* infection. *Diagn. Microbiol. Infect. Dis.* 75, 277.
43. Chicchio, T. et al. (2014) Polyfunctional T-cells and effector memory phenotype are associated with active TB in HIV-infected patients. *J. Infect.* 69, 533.
44. Ongaya, A. et al. (2013) *Mycobacterium tuberculosis*-specific CD8+ T cell recall in convalescing TB subjects with HIV co-infection. *Tuberculosis* 93, S60.
45. Lanicioni, C. et al. (2012) CD8+ T cells provide an immunologic signature of tuberculosis in young children. *Am. J. Respir. Crit. Care Med.* 185, 206.
46. Mazurek, G.H., et al.; IGRA Expert Committee; Centers for Disease Control and Prevention (CDC). (2010) Updated guidelines for using Interferon Gamma Release Assays to detect *Mycobacterium tuberculosis* infection – United States, 2010. *MMWR Recomm. Rep.* 59, 1.
47. Chiappini, E. et al. (2014). Interferon-gamma release assay sensitivity in children younger than 5 years is insufficient to replace the use of tuberculin skin test in western countries. *Pediatric Infectious Disease Journal*, 33(12), 1291–1293.
48. Garazzino, S. et al. (2014). Performance of interferon- $\gamma$  release assay for the diagnosis of active or latent tuberculosis in children in the first 2 years of age: a multicenter study of the Italian Society of Pediatric Infectious Diseases. *Pediatric Infectious Disease Journal*, 33(9).
49. Rose, W. et al. (2014). Quantiferon gold-in-tube assay for TB screening in HIV infected children: Influence of Quantitative Values. *BMC Infectious Diseases*, 14:516.

- 
50. Pavić, I. et al. (2015). Discordance between tuberculin skin test and interferon- $\gamma$  release assay in children younger than 5 years who have been vaccinated with bacillus calmette-guérin. *Laboratory Medicine*, 46(3), 200–206.
  51. Howley, M. M. et al. (2015). Evaluation of Quantiferon-TB gold in-tube and tuberculin skin tests among immigrant children being screened for latent tuberculosis infection. *Pediatric Infectious Disease Journal*, 34(1), 35–39.
  52. Rose, W. et al. (2014). Relating tuberculosis (TB) contact characteristics to quantiferon-tb-gold and tuberculin skin test results in the Toronto Pediatric TB Clinic. *Journal of the Pediatric Infectious Diseases Society*, 4(2), 96–103.
  53. Grinsdale, J. A. et al. (2014). Interferon-gamma release assays and pediatric public health tuberculosis screening: The San Francisco Program Experience 2005 to 2008. *Journal of the Pediatric Infectious Diseases Society*, 5(2), 122–130.
  54. Guðjónsdóttir, M. J. et al. (2016). Relation between BCG vaccine scar and an interferon-gamma release assay in immigrant children with “positive” tuberculin skin test ( $\geq 10$  mm). *BMC Infectious Diseases*, 16(1):540.
  55. Nejat, S. & Bennet, R. (2016). Interferon-gamma release assays can effectively screen migrants for the tuberculosis infection, but urgent, active cases need clinical recognition. *Acta Paediatrica*, 105(6), 671–675.
  56. Dehority, W. et al. (2017). Comparison of the Quantiferon TB gold in-tube assay with tuberculin skin test for the diagnosis of latent tuberculosis infection among HIV-infected and uninfected children. *Pediatric Infectious Disease Journal*, 36(12): e317-e321.

57. Marino, A. et al. (2017). Prebiologic therapy tuberculosis screening experience in a pediatric rheumatology center: TST and IGRA are both necessary. *Pediatric Infectious Disease Journal*, 36(4), 440–441.
58. Petrucci, R. et al. (2017). Quantiferon-TB gold in-tube improves tuberculosis diagnosis in children. *Pediatric Infectious Disease Journal*, 36(1), 44–49.
59. Hormi, M. et al. (2018). Performance of the quantiferon-TB gold assay among HIV-infected children with active tuberculosis in France. *Pediatric Infectious Disease Journal*, 37(4), 339–344.
60. Kay, A. W. et al. (2018). Interferon- $\gamma$  release assay performance for tuberculosis in childhood. *Pediatrics*, 141(6):e20173918.
61. Stout, J.E. et al. (2018). Tuberculosis Epidemiologic Studies Consortium. Evaluating latent tuberculosis infection diagnostics using latent class analysis. *Thorax*, 73(11):1062-1070.
62. Chiappini, E. et al. (2019). QuantiFERON-TB Gold In-Tube test performance in a large pediatric population investigated for suspected tuberculosis infection. *Paediatr Respir Rev*. 32:36-47.
63. Debulpaep, S. et al. (2019). Contribution of QuantiFERON-TB Gold-in-Tube to the Diagnosis of Mycobacterium tuberculosis Infection in Young Children in a Low TB Prevalence Country. *Front Pediatr*. 18;7:291.
64. Lombardi, G. et al. (2019). QuantiFERON-TB Performs Better in Children, Including Infants, than in Adults with Active Tuberculosis: a Multicenter Study. *J Clin Microbiol*. 24;57(10):e01048-19.

65. Mastrolia, M.V. et al. (2019). Utility of tuberculin skin test and IGRA for tuberculosis screening in internationally adopted children: Retrospective analysis from a single center in Florence, Italy. *Travel Med Infect Dis.* 28:64-67.
66. Venkatappa, T.K. et al. (2019). Comparing QuantiFERON-TB Gold Plus with Other Tests To Diagnose Mycobacterium tuberculosis Infection. *J Clin Microbiol.* 23;57(11):e00985-19.
67. Ahmed, A. et al. (2020). Tuberculosis Epidemiologic Studies Consortium. Interferon- $\gamma$  Release Assays in Children <15 Years of Age. *Pediatrics.* 145(1):e20191930.
68. Velasco-Arnaiz, E. et al. (2020). Impact of Baseline Tuberculin Skin Test and Isoniazid Chemoprophylaxis on Subsequent Quantiferon-TB Gold In-Tube Performance in Young Children Assessed After Tuberculosis Contact in Catalonia. *Pediatr Infect Dis J.* 39(2):e22-e25.
69. Buonsenso, D. et al. (2020). Accuracy of QuantiFERON-TB Gold Plus Test for Diagnosis of Mycobacterium tuberculosis Infection in Children. *J Clin Microbiol.* 26;58(6):e00272-20.
70. Quintana-Ortega, C. et al. (2020). False-positive Results of Quantiferon-Tb-Gold Assay in Children. *Pediatr Infect Dis J.* 39(7):620-623.
71. Worjolah, A. et al. (2011). Interferon gamma release assay compared with the tuberculin skin test for latent tuberculosis detection in pregnancy. *Obstet Gynecol.* 118(6):1363-70.
72. Lighter-Fisher, J. & Surette, A.M. (2012). Performance of an interferon-gamma release assay to diagnose latent tuberculosis infection during pregnancy. *Obstet Gynecol.* 119(6):1088-95.

73. Aichelburg, MC. et al. (2014). Reversion and conversion of interferon- $\gamma$  release assay results in HIV-1-infected individuals. *J Infect Dis.* 209(5):729-33.
74. Aichelburg, MC. et al. (2014). The association of smoking with IGRA and TST results in HIV-1-infected subjects. *Int J Tuberc Lung Dis.* 18(6):709-16.
75. Pullar, ND. et al. (2014). Low prevalence of positive interferon-gamma tests in HIV-positive long-term immigrants in Norway. *Int J Tuberc Lung Dis.* 18(2):180-7.
76. Pullar, ND. et al. (2014). HIV patients with latent tuberculosis living in a low-endemic country do not develop active disease during a 2 year follow-up; a Norwegian prospective multicenter study. *BMC Infect Dis.* 14:667.
77. Sauzullo, I. et al. (2014). Interferon- $\gamma$  release assay in HIV-infected patients with active tuberculosis: impact of antituberculous drugs on host immune response. *New Microbiol.* 37(2):153-61.
78. Sester, M. et al. (2014). Risk assessment of tuberculosis in immunocompromised patients. A TBNET study. *Am J Respir Crit Care Med.* 190(10):1168-76.
79. Bourgarit, A. et al. (2015). IGRAVIH Study Group. Latent Tuberculosis Infection Screening and 2-Year Outcome in Antiretroviral-Naive HIV-Infected Patients in a Low-Prevalence Country. *Ann Am Thorac Soc.* 12(8):1138-45.
80. Parrella, R. et al. (2015). Interferon gamma release assays and tuberculin skin test performance in different settings of HIV immunodeficiency. *In Vivo.* 29(1):137-40.
81. Kurti, Z. et al. (2015). Tuberculin Skin Test and Quantiferon in BCG Vaccinated, Immunosuppressed Patients with Moderate-to-Severe Inflammatory Bowel Disease. *J Gastrointestin Liver Dis.* 24(4):467-72.

- 
82. Perifanou, D. et al. (2018). Screening for Latent Tuberculosis Infection in Patients with Autoimmune Diseases Before Initiating TNF- $\alpha$  Inhibitors Therapy. *Mater Sociomed.* 30(1):32-37.
  83. Chiacchio, T. et al. (2019). Characterization of QuantiFERON-TB-Plus results in latent tuberculosis infected patients with or without immune-mediated inflammatory diseases. *J Infect.* 79(1):15-23.
  84. Diamantopoulos, P.T. et al. (2017). Evaluation of the performance of tuberculin skin test and Quantiferon-TB gold in tube test in patients with hematologic malignancies. *Infect Dis (Lond).* 49(7):545-548.

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© 2023 QIAGEN，保留所有權利