

# QIAGEN Supplementary Protocol

---

## Purification of DNA amplified using REPLI-g® Kits

This protocol is optimized for the purification of 2–60 µg REPLI-g amplified DNA using a simple centrifugation step. In general, recovery of at least 80% of REPLI-g amplified DNA can be achieved using this protocol. Purification with silica-column-based cleanup methods are not recommended due to recovery loss of up to 50%, which is caused by the large fragment size (2000–70000 bp) of the amplified DNA.

Purification of REPLI-g amplified DNA is necessary if residual primers or nucleotides interfere with downstream analysis (e.g., direct labeling of REPLI-g amplified DNA using labeled nucleotides).

**IMPORTANT:** Please refer to the handbooks supplied with the respective REPLI-g Kits for general information on the handling and storage of kit components.

## Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- 1.5 or 2 ml microcentrifuge tubes
- Microcentrifuge
- Pipet tips (pipet tips with aerosol barriers for preventing cross-contamination are recommended)
- TE buffer (10 mM Tris·Cl; 1 mM EDTA, pH 8.0)
- 96–100% ethanol
- 70% ethanol

## Important points before starting

- All centrifugation steps should be performed at room temperature (15–25°C).
- Equilibrate REPLI-g amplified DNA to room temperature (15–25°C).
- Wear gloves throughout the entire procedure. In case of contact between gloves and sample, change gloves immediately.



# QIAGEN Supplementary Protocol

---

## Procedure

1. **After the REPLI-g reaction, equilibrate REPLI-g amplified DNA to room temperature (15–25°C).**

This protocol is suitable for use immediately following amplification using REPLI-g Kits, or with REPLI-g amplified DNA that has been stored at –20°C.

**Note:** If necessary, add the reaction mixture to a new 1.5 ml microcentrifuge tube.

2. **Add 150 µl ethanol (96–100%). Mix carefully by flicking the tube.**

**Note:** Do not mix by pipetting up and down.

3. **Centrifuge at maximum speed for 1–2 min.**

**Note:** Do not centrifuge for longer than 2 min. Centrifugation for longer periods would condense the precipitate, making it more difficult to dissolve.

4. **Aspirate the supernatant carefully using a pipet and discard the supernatant.**

5. **Add 100 µl of 70% ethanol to the pellet.**

6. **Centrifuge at maximum speed for 1–2 min.**

**Note:** Do not centrifuge for longer than 2 min. Centrifugation for longer periods will cause the precipitate to condense, making it more difficult to dissolve.

7. **Aspirate the supernatant carefully using a pipet and discard the supernatant.**

8. **Centrifuge briefly to collect residual supernatant at the bottom of the tube.**

9. **Aspirate the residual supernatant carefully using a pipet and discard the supernatant.**

10. **Incubate the microcentrifuge tube containing the precipitate of REPLI-g amplified DNA for 5 min at room temperature (15–25°C).**

11. **Add 50 µl 1x TE buffer (pH 8.0) to dissolve REPLI-g amplified DNA.**

**Note:** Mix carefully to avoid shearing of REPLI-g amplified DNA.

12. **Purified DNA can be directly used for all downstream applications.**

Optical density (OD) measurements can be performed to accurately determine the concentration of purified REPLI-g amplified DNA.

# QIAGEN Supplementary Protocol

---

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Selected handbooks can be downloaded from [www.qiagen.com/literature](http://www.qiagen.com/literature). Safety data sheets (SDS) for any QIAGEN product can be downloaded from [www.qiagen.com/safety](http://www.qiagen.com/safety).

Trademarks: QIAGEN®, REPLI-g® (QIAGEN Group).  
RG21 Nov-12 © 2012 QIAGEN, all rights reserved.

**www.qiagen.com**

**Australia** ■ 1-800-243-800

**Austria** ■ 0800/281010

**Belgium** ■ 0800-79612

**Canada** ■ 800-572-9613

**China** ■ 021-51345678

**Denmark** ■ 80-885945

**Finland** ■ 0800-914416

**France** ■ 01-60-920-930

**Germany** ■ 02103-29-12000

**Hong Kong** ■ 800 933 965

**Ireland** ■ 1800 555 049

**Italy** ■ 800 787980

**Japan** ■ 03-6890-7300

**Korea (South)** ■ 1544 7145

**Luxembourg** ■ 8002 2076

**The Netherlands** ■ 0800 0229592

**Norway** ■ 800-18859

**Singapore** ■ 65-67775366

**Spain** ■ 91-630-7050

**Sweden** ■ 020-790282

**Switzerland** ■ 055-254-22-11

**UK** ■ 0808-2343665

**USA** ■ 800-426-8157



---

Sample & Assay Technologies