

***pigtype*[®] Trichinella Ab Handbook**



1 (catalog no. 273501)



5 (catalog no. 273503)



20 (catalog no. 273505)*

Multi-species ELISA kit for the detection of antibodies to *Trichinella* spp.

REF

273501, 273503, 273505*



QIAGEN Leipzig GmbH, Deutscher Platz 5b, 04103
Leipzig, Germany



* Available only on request.

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

In addition, QIAGEN provides high-quality, easy-to-use, and sensitive molecular solutions to enable veterinary pathogen detection and animal pathogen research. The QIAGEN veterinary portfolio includes a broad range of pathogen-specific PCR-assays and an extensive and growing ELISA portfolio. For more information, visit www.qiagen.com/Animal-and-Veterinary-Testing.

Contents

| | |
|--|-----------|
| Kit Contents | 4 |
| Intended Use | 5 |
| Symbols | 5 |
| Storage | 6 |
| Safety Information | 6 |
| Quality Control | 7 |
| Introduction | 8 |
| Principle | 9 |
| Equipment and Reagents to Be Supplied by User | 10 |
| General precautions | 11 |
| Protocol: | |
| ■ ELISA | 13 |
| Data Interpretation | 15 |
| Validation Criteria | 15 |
| Troubleshooting Guide | 17 |
| Ordering Information | 18 |
| Quick guide | 20 |
| Data interpretation | 20 |

Kit Contents











| <i>pigtype</i> Trichinella Ab | | | |
|---|---------------|---------------|----------------|
| Catalog no. | 273501 | 273503 | 273505* |
| Number of plates | 1 | 5 | 20 |
| Test Plate: microtiter plate with 96 wells, coated with non-infectious <i>Trichinella</i> E/S-antigen | 1 | 5 | 20 |
| Sample diluent, ready-to-use | 1 x 60 ml | 2 x 125 ml | 2 x 500 ml |
| Negative Control, ready-to-use | 1 x 1.5 ml | 1 x 3.5 ml | 2 x 3.5 ml |
| Positive Control, ready-to-use | 1 x 1.5 ml | 1 x 3.5 ml | 2 x 3.5 ml |
| Wash Buffer (10x) | 1 x 125 ml | 2 x 125 ml | 2 x 500 ml |
| Conjugate, ready-to-use | 1 x 12 ml | 1 x 60 ml | 1 x 240 ml |
| TMB Substrate, ready-to-use | 1 x 12 ml | 1 x 60 ml | 1 x 240 ml |
| Stop solution, ready-to-use | 1 x 12 ml | 1 x 60 ml | 1 x 240 ml |
| Handbook | 1 | 1 | 1 |

* Available only on request.

Intended Use

The pigtype *Trichinella* Ab is an indirect enzyme immunoassay (ELISA). It is intended for the detection of antibodies to *Trichinella* spp. in pig and wild boar serum, plasma, and meat juice samples. Pig and wild boar serum and plasma samples can be tested as pools of up to 10 individual samples. The test kit may also be used for other mammalian species e.g., horse, fox, and rodents. For veterinary use only.

Symbols

| | |
|---|-------------------------------------|
|  | Contains reagents for <N> plates |
|  | Legal manufacturer |
|  | Lot number |
|  | Use by date |
|  | Temperature limitations for storage |
|  | Handbook |
|  | Catalog number |
|  | Material number |
|  | Protect from light |
|  | For pig and wild boar samples |

Storage

The components of the *pigtype* Trichinella Ab ELISA should be stored at 2–8°C and are stable until the expiration date stated on the label. Wash Buffer (10x) and Stop Solution may be stored at room temperature (18–25°C) to avoid salt crystallization. Store the remaining unused test strips in the resealed foil pouch with desiccant at 2–8°C until next use. The test strips can be stored for at least 6 weeks after opening the plate pouch.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.



**CAUTION: The Stop Solution contains
0.5 M sulphuric acid.**

All sample residues and objects which have come into contact with samples must be decontaminated or disposed as potentially infectious material.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of *pigtype* Trichinella Ab is tested against predetermined specifications to ensure consistent product quality.

Introduction

The *pigtype* *Trichinella* Ab is a highly sensitive solution for the detection of antibodies to *Trichinella* spp. Trichinellosis is a zoonosis which is caused by the nematode *Trichinella*. In humans the infection can cause mild to lethal illness. Humans are infected by consumption of pork containing *Trichinella*-larvae. Most infected pigs do not show clinical symptoms. To prevent human infections, all pigs slaughtered for the food chain are tested by artificial digestion.

The European Commission Regulation (EC) No. 2075/2005 allows holdings to be certified as *Trichinella*-free under certain conditions. One of these conditions is a monitoring program. Serological testing can provide a more cost-effective method, to the commonly used digestion-method, for the monitoring of *Trichinella* antibodies in pork from integrated production systems. The *pigtype* *Trichinella* Ab can be used for monitoring *Trichinella* infections as part of such programs.

Principle

pigtype Trichinella Ab is an indirect ELISA. The microtiter plate is coated with inactivated *Trichinella*-antigen (E/S-antigen). During sample incubation *Trichinella*-specific antibodies bind to the immobilized antigen. Unbound material is removed by rinsing. Serum antibodies bound to the antigen are detected by a multi-species horseradish peroxidase (HRP) conjugate. Unbound conjugate is removed by rinsing. A colorimetric reaction is initiated by adding Substrate Solution and stopped after 10 minutes. In the presence of *Trichinella*-specific antibodies, within the sample, HRP catalyzes a blue color development, which turns yellow after adding the Stop Solution. The optical density (OD) is measured using a spectrophotometer at 450 nm. The OD values correlate with the concentration of *Trichinella*-specific antibodies in the sample.

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- Beakers
- Measuring cylinders
- Pipets (adjustable)
- Multichannel pipets (adjustable)
- Aluminum or adhesive foil for covering the Test Plate
- Device for delivery and aspiration of wash solution (optional)
- Microtiter plate absorbance reader
- Tubes or plates for diluting the samples
- Distilled water

General precautions

The user should always pay attention to the following:

- Do not expose the TMB Substrate Solution to intense light or to sunlight during the performance of the test.
- Components of the test kit should not be contaminated or mixed with components from other batches.
- Do not use the components of the test kit past expiration date.
- Water from ion-exchange systems used for diluting the Wash Solution (10x) may interfere with the assay if not pure enough. Water quality of double distilled water or highly purified water (Milli-Q) is suitable.
- The use of clean glass devices, careful pipetting and rinsing during the test, and strict adherence to the indicated incubation times is essential for precise test results.

Things to do before starting

- Bring reagents to room temperature (18–25°C) immediately before use. In case of precipitated salt crystals in the Wash Buffer (10x), dissolve by gentle swirling and warming.

Wash Buffer: Dilute Wash Buffer (10x) 1:10 in distilled water, e.g., for one test plate dilute 25 ml Wash Buffer (10x) in 225 ml distilled water and mix.

Serum/plasma: Prior to sample analysis, with serum/plasma samples, dilute **1:100** in Sample Diluent (e.g., dilute 5 μ l sample in 495 μ l Sample Diluent) and mix well. Use plastic tubes or uncoated microtiter plates for dilution. Change pipet tips for each sample.

Pool samples: Serum or plasma samples can also be tested as pools of 10 single samples (e.g., pool 10 μl of each of 10 samples). Dilute the pool samples **1:20** with Sample Diluent (e.g., dilute 10 μl sample in 190 μl Sample Diluent) and mix.

Meat juice: Prior to sample analysis, with meat juice samples, dilute **1:10** in Sample Diluent (e.g., dilute 25 μl sample in 225 μl Sample Diluent) and mix well.

Alternatively, meat juice samples can be diluted directly in the Test Plate (see Procedure 1a).

To extract meat juice samples, freeze approximately 10 g blood- and fat-free muscle meat (e.g., from the muscular pillars of lumbar diaphragm) in a meat juice sampling device and then thaw. Take the meat juice released from the thawed samples and store at 2–8°C. Samples stored at 2–8°C should be analyzed within 24 hours (alternatively, meat juice samples can be stored at –20°C for several months until analysis).

- **Controls** are ready-to-use and do not require dilution.

Protocol: ELISA

Please read “Things to do before starting”, page 11.

Procedure

1. Pipet 100 μ l of each of the ready-to-use Negative Control (in duplicates) and Positive Control (in duplicates) and the 1:10 diluted meat juice samples and/or 1:100 diluted serum or plasma samples into the Test Plate wells.

1a. Alternatively, pipet 90 μ l of Sample Diluent in each sample well and add 10 μ l of the undiluted meat juice sample and/or of the 1:10 pre-diluted serum or plasma sample. Mix well.

Record the positions of the controls and samples in a test protocol. The use of a multichannel pipet is recommended for the transfer of samples. Cover the Test Plate.

2. Incubate for 60 min at room temperature (18–25°C) or overnight (12–18 hours) at 2–8°C.
3. Remove solution from the wells by aspiration or tapping.
4. Rinse each well 3 x with 300 μ l of prepared Wash Buffer. Remove the buffer after each rinse.
5. Pipet 100 μ l ready-to-use Conjugate to each well and incubate for 30 min at room temperature.
6. Remove solution from wells by aspiration or tapping.
7. Rinse each well 3 x with 300 μ l of prepared Wash Buffer. Remove the buffer after each rinse.
8. Pipet 100 μ l TMB Substrate Solution to each well.
9. Incubate for 10 min at room temperature in the dark. Begin timing after the first well is filled.

- 10. Stop the reaction by adding 100 μ l Stop Solution per well. Add the Stop Solution in the same order as the Substrate Solution was added.**
- 11. Measure the OD in the plate reader at 450 nm within 20 min after stopping the reaction.**
Measuring at a reference wavelength (620–650 nm) is optional.

Data Interpretation

Validation Criteria

The results are valid if the following criteria are met:

- The mean OD of the Positive Control must be ≥ 0.7
- The mean OD of the Negative Control must be ≤ 0.2

In case of invalid assays the test should be repeated after a thorough review of the instructions for use.

Calculation

Calculate the mean values (MV) of the measured OD for the Negative Control (NC) and the Positive Control (PC).

The ratio (S/P) of sample OD to mean OD of the Positive Control is calculated according to the following equation:

$$S/P = \frac{OD_{\text{sample}} - MV OD_{\text{NC}}}{MV OD_{\text{PC}} - MV OD_{\text{NC}}}$$

Interpretation short protocol (60 min sample incubation)

Samples with S/P-ratio ≥ 0.3 are positive.

Specific antibodies to *Trichinella* spp. were detected.

Samples with S/P-ratio < 0.3 are negative.

Specific antibodies to *Trichinella* spp. could not be detected.

Interpretation overnight protocol (O/N sample incubation)

Samples with S/P-ratio ≥ 0.45 are positive.

Specific antibodies to *Trichinella* spp. were detected.

Samples with S/P-ratio < 0.45 are negative.

Specific antibodies to *Trichinella* spp. could not be detected

Troubleshooting Guide

The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit www.qiagen.com).

Ordering Information

| Product | Contents | Cat. no. |
|---|--|-----------------|
| <i>pigtype</i> Trichinella Ab (1) | For 96 reactions: 1 Test Plate (strips), Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273501 |
| <i>pigtype</i> Trichinella Ab (5) | For 480 reactions: 5 Test Plates (strips), Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273503 |
| <i>pigtype</i> Trichinella Ab (20)* | For 1920 reactions: 20 Test Plates, Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273505 |
| Related products | | |
| <i>pigtype</i> Salmonella Ab (5) [†] | For 480 reactions: 5 Test Plates (strips), Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273003 |
| <i>pigtype</i> Yersinia Ab (1) | For 96 reactions: 1 Test Plate (strips), Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273801 |

* Available only on request.

[†] Other kit sizes are available; see www.qiagen.com.

| Product | Contents | Cat. no. |
|---|--|-----------------|
| <i>pigtype</i> Toxoplasma Ab (5)* | For 480 reactions: 5 Test Plates (strips), Wash Buffer, Sample Diluent, Positive Control, Negative Control, Conjugate, TMB Substrate Solution, Stop Solution | 273403 |

* Other kit sizes are available; see www.qiagen.com.

QIAGEN offer a range of ELISA kits and real-time PCR and real-time RT-PCR kits for the detection of animal pathogens. Visit www.qiagen.com/Animal-and-Veterinary-Testing for more information about the *bactotype*[®], *cador*[®], *cattletype*[®], *flocktype*[®], *pigtype*, and *virotype*[®] products.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Quick guide

Sample dilution:

Serum/plasma 1:100

Meat juice 1:10

| Step | Short protocol | Overnight protocol |
|---------------|----------------|--------------------|
| 1. Sample | | 100 µl/well |
| 2. Incubation | 60 min. RT | Overnight 2–8°C |
| 3. Wash | | 3 x 300 µl |
| 4. Conjugate | | 100 µl/well |
| 5. Incubation | | 30 min. RT |
| 6. Wash | | 3 x 300 µl |
| 7. TMB | | 100 µl/well |
| 8. Incubation | | 10 min. RT |
| 9. Stop | | 100 µl/well |
| 10. Read | | 450 nm |

Data interpretation

| | Negative | Positive |
|--------------------|------------|------------|
| Short protocol | S/P < 0.3 | S/P ≥ 0.3 |
| Overnight protocol | S/P < 0.45 | S/P ≥ 0.45 |

Notes

Notes

Trademarks: QIAGEN®, *bactotype*®, *cador*®, *cattletype*®, *flocktype*®, *pigtype*®, *virotype*® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

© 2015 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ techservice-au@qiagen.com

Austria ■ techservice-at@qiagen.com

Belgium ■ techservice-bnl@qiagen.com

Brazil ■ suportetecnico.brasil@qiagen.com

China ■ techservice-cn@qiagen.com

Denmark ■ techservice-nordic@qiagen.com

Finland ■ techservice-nordic@qiagen.com

France ■ techservice-fr@qiagen.com

Germany ■ techservice-de@qiagen.com

Hong Kong ■ techservice-hk@qiagen.com

India ■ techservice-india@qiagen.com

Ireland ■ techservice-uk@qiagen.com

Italy ■ techservice-it@qiagen.com

Japan ■ techservice-jp@qiagen.com

Korea (South) ■ techservice-kr@qiagen.com

Luxembourg ■ techservice-bnl@qiagen.com

Mexico ■ techservice-mx@qiagen.com

The Netherlands ■ techservice-bnl@qiagen.com

Norway ■ techservice-nordic@qiagen.com

Singapore ■ techservice-sg@qiagen.com

Sweden ■ techservice-nordic@qiagen.com

Switzerland ■ techservice-ch@qiagen.com

UK ■ techservice-uk@qiagen.com

