

# DNA quality measurement and somatic mutation profiling in PAXgene tissue samples with qBiomarker Somatic Mutation PCR Arrays



Samuel Long<sup>1</sup>, Daniel Groelz<sup>1</sup>, Lili Yang<sup>2</sup>, Sally Yousif<sup>1</sup>, Karen Rodland<sup>1</sup>, Sunita Badola<sup>2</sup> and Eric Lader<sup>1</sup>  
<sup>1</sup>QIAGEN Sciences Inc., 6951 Executive Way, Frederick, MD 27103; <sup>2</sup>Millennium Pharmaceuticals, Inc., 40 Landsdowne St., Cambridge, MA 02139

## Abstract

Current formaldehyde-based tissue fixation methods are of limited use for molecular analysis due to degradation and chemical modification of nucleic acids that can cause inhibition in PCR or real-time PCR. The PAXgene Tissue Kit is a novel, powerful system that enables both molecular and traditional pathology analysis from the same specimen. The system consists of a tissue collection device (tissue container for collection, stabilization, storage, and transportation of tissue specimens) and kits for purification of nucleic acids. The fixation and stabilization method used in this system preserves tissue morphology and the integrity of nucleic acids for an extended period without the destructive cross-linking and degradation found in formalin-fixed tissues. In addition, samples stabilized using the PAXgene Tissue Kit can be embedded in paraffin for histological studies. Nucleic acids can be isolated from the stabilized samples either before or after embedding in paraffin.



In the present study, the specialized qBiomarker Somatic Mutation PCR Array for Human DNA QC assessment, which consists of a panel of 7 copy number PCR control assays, was employed to compare the quality of DNA extracted from fresh frozen, formaldehyde-fixed, and PAXgene tissues that were prepared at the same time, from the same set of mouse xenograft samples of human melanoma. The DNA quality of PAXgene tissue samples was found to be similar to that of fresh frozen samples (i.e. intact DNA), while more than 80% of the DNA fraction in formaldehyde-fixed samples exists in non-amplifiable fractions under typical real-time PCR conditions.

The Human EGFR Pathway qBiomarker Somatic Mutation PCR Array was used to profile the mutation status in these melanoma xenograft samples for the following genes: EGFR, BRAF, KRAS, NRAS, HRAS, c-MET, PIK3CA, AKT and PTEN. While the identical set of BRAF and NRAS mutations was identified in DNAs from samples of all three preparation methods, at all DNA input doses tested (200 ng per sample as the lowest dose for screening 87 mutations, and up to 1200 ng per sample), the mutation call was at least 4-fold more sensitive in DNA extracted from PAXgene samples and fresh frozen samples than that from the FFPE DNA. In addition, when duplicate samples were compared, among the three sample preparation methods, PAXgene sample DNA quality showed the least variability between duplicate samples. This consistency in sample quality also translated into the highest reproducibility in mutation calls among the three sample preparation methods.


The PAXgene Tissue System, with its excellent tissue morphology preservation and close-to-intact level of nucleic acid preservation, is an ideal system for molecular and histological testing of research samples. The combination of the PAXgene Tissue Kit with disease- or pathway-focused qBiomarker Somatic Mutation PCR Arrays allows specific, highly sensitive, and accurate somatic mutation detection in tumor samples (especially for mutations occurring at a low frequency), and has distinct advantages over other systems that are used for sample stratification.

The applications presented here are for research use only. Not for use in diagnostic procedures.

## PAXgene Tissue System workflow and experimental design




## Performance of PAXgene Tissue Kit in preserving morphology of surgically resected tissues




**Figure 2. Tissue morphology preservation by the PAXgene Tissue Kit.** Morphology of H&E (hematoxylin and eosin) stained, mirrored human spleen, liver, pancreas and large intestine tissue sections (4 μm), preserved either by PFPE (PAXgene Tissue fixed, paraffin embedded) or FFPE (formalin-fixed). Surgically resected tissues were collected by a commercial tissue provider (Cureline) with prior written informed consent from respective patients. Spleen and liver: original magnifications x 10; pancreas and large intestine: original magnifications x 20.

## qBiomarker Somatic Mutation PCR Assay principle and performance




**Figure 3. qBiomarker Somatic Mutation PCR Assay principle and performance.** [A] qBiomarker Somatic Mutation PCR Assays employ ARMS primers to differentiate between mutant and wild-type alleles using template matches and mismatches at a primer's 3' end. An ARMS primer that has a match to the mutant template at the 3' end, which is mismatched with the wild-type, will extend on the mutant instead of the wild-type template. Additional mismatch(es) can be included in the ARMS primer to increase allele discrimination. Amplicons containing the mutation are subsequently detected by a probe. [B] Assay sensitivity test for p53 R280K qBiomarker Somatic Mutation PCR Assay. A series of 10 ng genomic DNA samples, which contain 4, 21, 106, 528, 2640 and 13200 copies of mutant DNA template respectively, from MDA-MB-231 cell line (mutant harboring the p53 R280K mutation) mixed with genomic DNA from the Coriell GM00131 cell line (wild-type), were tested with the p53 R280K mutation assay. Amplification plots for duplicate reactions are shown. Mutation detection limit for this assay is determined to be 0.03%.

## qBiomarker Somatic Mutation PCR Array workflow




**Figure 4. qBiomarker Somatic Mutation PCR Array workflow.** The procedure involves DNA extraction, qPCR detection on qBiomarker Somatic Mutation PCR Arrays, and data analysis (using the qBiomarker Somatic Mutation Data Analysis Template or web portal). The arrays can be used on all real-time cyclers except the Cepheid® SmartCycler® and the Roche® LightCycler® 2.0.

## DNA integrity of PAXgene-preserved tissues is far superior to that of FFPE samples



**Figure 5. DNA integrity comparison among fresh frozen, PAXgene tissue, and FFPE samples.** [A] Analysis of 300 ng DNA in agarose gel (0.8%) electrophoresis. DNA was isolated from mirrored human tumor tissues (left: colorectal cancer; right: breast cancer) that were fresh frozen (liquid nitrogen snap frozen; LN2), PAXgene Tissue fixed and embedded in paraffin (PFPE), or FFPE. [B] qBiomarker Somatic Mutation PCR Array DNA QC plate layout. 7 probe-based gene copy number reference assays are used to interrogate sample DNA quality, while the SMPC assay monitors the presence of PCR inhibitors in the sample. Each 96-well plate can be used to QC 12 DNA samples, and each 384-well plate can QC 48 DNA samples. [C] Real-time PCR DNA quality analysis of DNAs from 4 mouse xenografts of human melanoma (color-coded gray, light blue, pink and light green) that were preserved by liquid nitrogen snap frozen (fresh frozen), PAXgene tissue kit followed by paraffin embedding (PAXgene), or formalin fixation (FFPE) at the same time. For each preservation method, each sample was prepared in duplicates. About 8 weeks after preservation processing, DNA was isolated from the samples using QIAamp DNA Mini Kit (fresh frozen), PAXgene Tissue DNA Kit, and QIAamp DNA FFPE Tissue Kit, respectively. Purified DNA samples were then analyzed on the qBiomarker Somatic Mutation PCR Array DNA QC plate at an input of 2 ng DNA per reaction. Sample QC C<sub>T</sub> average: the average of the C<sub>T</sub> values of the 7 gene copy number reference assays for the sample. Set QC C<sub>T</sub> average: the average of the "sample QC C<sub>T</sub> average" for the 8 samples with the same preservation method. On average, more than 80% of the DNA fraction in FFPE samples exists in non-amplifiable fractions under typical real-time PCR conditions, as compared to fresh frozen and PAXgene sample.

## Somatic mutation detection in PAXgene tissue samples using qBiomarker Somatic Mutation PCR Array



**Figure 6. Somatic mutation profiling in mirrored fresh frozen, PAXgene and FFPE samples using the EGFR pathway qBiomarker Somatic Mutation PCR Array.** [A] Assay layout of the EGFR pathway qBiomarker Somatic Mutation PCR Array. Three categories of probe-based assays are included on this array. Firstly, mutation detection assays are included for the most frequently occurring somatic mutations for genes in the EGFR pathway. Gene copy number reference assays are also incorporated for all genes covered on this array. These assays provide DNA quality and gene dosage information during data analysis. Finally, SMPC assays are included for monitoring PCR conditions, such as the presence of PCR inhibitors. [B] More robust performance of PAXgene samples compared to FFPE samples in real time PCR-based somatic mutation analysis. The samples in Figure 5(C) were profiled on the EGFR pathway qBiomarker Somatic Mutation PCR Array (4ng DNA per reaction) and C<sub>T</sub> values for selected loci are shown. On average, mutation call was at least 4 fold more sensitive (i.e. C<sub>T</sub>s ≥ 2 Ct lower) in DNA extracted from PAXgene and fresh frozen samples than FFPE DNA. C<sub>T</sub> values in black: wildtype calls; C<sub>T</sub> values in green: mutant.

## Conclusions

- Tissue morphology preservation in PAXgene tissue samples is comparable to or better than FFPE samples of the same age.
- DNA isolated from PAXgene tissues is mostly in 10 to 23 kb molecular weight range, which translates to significantly higher amplifiable fractions in real-time PCR.
- qBiomarker Somatic Mutation PCR Arrays consist of pathway- or disease-focused collections of mutation detection assays with sensitivity as high as 0.03%.
- Mutation detection in PAXgene samples is at least 4-fold more sensitive than in FFPE samples; the combination of the PAXgene Tissue Kit with qBiomarker Somatic Mutation PCR Arrays allows specific, highly sensitive, and accurate somatic mutation detection in tumor samples, and has distinct advantages over other systems that are used for sample stratification.