

Protocol Sheet

Applied Biosystems® StepOnePlus (for software version 2.0) instrument setup instructions for EpiTect® ChIP PCR Arrays

Important points before starting

- Please read the handbook supplied with the EpiTect ChIP PCR Array, paying careful attention to the “Safety Information” and “Important Notes” sections, before beginning this procedure.

Procedure

1. Open the ABI StepOnePlus software on the desktop of the computer that is connected to the ABI StepOnePlus system.
2. Select “New Experiment” on the upper toolbar.
3. “Define: Experiment Properties”
 - Label experiment.
 - Type in “Experiment Name.”
 - Type in “Barcode,” “User Name,” “Comments” (optional).
 - Select instrument.
 - “StepOnePlus Instrument (96 wells)”
 - Select experiment type.
 - “Quantitation”
 - Click “Next” on bottom of the screen.
4. “Define: Methods & Materials”
 - “Quantitation Method”
 - “Standard Curve”
 - “Reagents to Detect Target Sequence”
 - “SYBR® Green Reagents”
 - Keep “Melt Curve” checked.
 - “Ramp Speed”
 - “Standard (~2 hours to complete a run)”
 - “Template Type”
 - “cDNA” (complementary DNA)
 - Click “Next” on bottom of the screen.
5. “Set Up: Targets”



- **“How Many Targets Do You Want to Quantify?”**
 - **“1”**
- **Uncheck “Setup Standards”**
 - **“Target Name”:**
 - **“Target 1”**
 - **“Reporter”**
 - **“SYBR”**
 - **“Quencher”**
 - **“None”**
- **Click “Next” on bottom of the screen.**
- **Ignore the warning — click “OK.”**

6. **“Set Up: Standards”**

- **“How Many Points?”**
 - **“2”**
- **“How Many Replicates?”**
 - **“1”**
- **Click “Next” on bottom of the screen.**

7. **“Set Up: Samples”**

- **“How Many Samples?”**
 - **“96”**

Note: If the instrument is not recognizing all 96 wells, please see additional instructions on the last page.

- **“How Many Replicates?”**
 - **“1”**
- **“How Many Negative Controls?”**
 - **“0”**
- **“Which Sample/Target Reactions Do You Want To Set Up?”**
 - **Select: “ALL Sample/Target Reactions”**
- **Verify all wells in “Plate Layout” view have the “U” symbol (“U” = unknown).**
- **Click “Next” on bottom of the screen.**

8. **“Set Up: Run Method”**

- **This setting should default to the run protocol with melting curve.**
 - **“Verify Data Capture” icon is present at:**
 - **“Cycling Stage”:** 60°C (1 minute step)
 - **“Melting Curve Stage”:** During ramp from 60°C to 95°C

- Set "Reaction Volume" to 25 µl.
 - Verify "Number of Cycles" is set to 40.
9. Click "Finish Designing Experiment."
 10. Ignore warning.
 11. Click "OK" when prompted "You did not set up standards on the plate."
 12. Load your plate into the instrument.
 13. Start the run for this experiment.
 14. Save your experiment before starting the run.

Note: For those customers whose instruments do not recognize all 96 wells of the PCR arrays, please use the following instructions.

ABI StepOnePlus — modified setup

1. Open the ABI StepOnePlus software on the desktop of the computer that is connected to the ABI StepOnePlus system.
2. Select "Advanced Setup."
3. "Define: Experiment Properties"
 - Label the experiment.
 - Type in "Experiment Name."
 - Type in "Barcode", "User Name", "Comments" (optional).
 - Select instrument.
 - "StepOnePlus Instrument (96 wells)"
 - Select experiment type.
 - "Quantitation — Standard Curve"
 - Select reagents.
 - "SYBR Green"
 - Select ramp speed.
 - "Standard (~2 hours to complete)"
4. Click "Plate Setup" (on left).
 - Click "Assign Targets and Samples" tab.
 - Highlight the entire plate.
 - Check the box next to "Target 1" under "Assign Targets to the Selected Wells."
 - Verify all wells in "Plate Layout" view have the "U" symbol ("U" = unknown).
5. Click "Run Method" (this setting should default to run protocol with melting curve.)

- Verify “Data Capture” icon is present at:
 - “Cycling Stage”: 60°C (1 minute step)
 - “Melting Curve Stage”: during ramp from 60°C to 95°C
- Set “Reaction Volume” to 25 µl.
- Verify “Number of Cycles” is set to 40.

6. Click “Start Run.”

The EpiTect ChIP PCR Arrays are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected handbooks can be downloaded from www.qiagen.com/literature. Safety data sheets (SDS) for any QIAGEN product can be downloaded from www.qiagen.com/safety.

Trademarks: EpiTect®, QIAGEN® (QIAGEN Group); Applied Biosystems®, SYBR® (Life Technologies Corporation).

Oct-12 © 2012 QIAGEN, all rights reserved.

www.qiagen.com
Australia ■ 1-800-243-800
Austria ■ 0800/281010
Belgium ■ 0800-79612
Canada ■ 800-572-9613
China ■ 021-51345678
Denmark ■ 80-885945
Finland ■ 0800-914416

France ■ 01-60-920-930
Germany ■ 02103-29-12000
Hong Kong ■ 800 933 965
Ireland ■ 1800 555 049
Italy ■ 800 787980
Japan ■ 03-6890-7300
Korea (South) ■ 1544 7145
Luxembourg ■ 8002 2076

The Netherlands ■ 0800 0229592
Norway ■ 800-18859
Singapore ■ 65-67775366
Spain ■ 91-630-7050
Sweden ■ 020-790282
Switzerland ■ 055-254-22-11
UK ■ 01293-422-911
USA ■ 800-426-8157

