

Advancing cancer research: A novel digital PCR tool for simultaneous detection of multiple hallmark mutations in *BRAF* and *EGFR*



Claudia Kappmeier¹, Sherina Edward¹, Corinna Hochstein¹, Jette Ellehauge², Stina Christensen², Tae Jung Kim³, Dávid Kis⁴, Orsolya Biró⁴, Ellen Bruske¹, Francesca Di Pasquale¹, Ronny Kellner¹

¹QIAGEN GmbH, QIAGEN Strasse 1, 40724, Hilden, Germany; ²Klinisk Biokemisk Afdeling, Ringstedgade 57, 4700 Næstved, Denmark;

³Seoul St. Mary's Hospital, Department of Pathology, 222 Banpo-Daero, Seocho-Gu, Seoul, Korea; ⁴Clinomics Europe Ltd., 1094 Budapest, Mihálkovichs utca 10, Hungary

Faster and easier monitoring of hallmark cancer mutations

Precision oncology requires the identification of the specific genetic mutations driving cancer. Some hallmark mutations, such as those found in *BRAF* and *EGFR* genes, are important in the analysis of multiple types of cancer. The capability to concurrently and accurately detect multiple hallmark mutations would greatly benefit cancer research.

Digital PCR (dPCR) technology provides an opportunity for highly sensitive, quantitative mutation analysis. Here, we introduce the new and innovative dPCR PanCancer Kits (RUO), which are designed to detect multiple hallmark mutations in *BRAF* or *EGFR*.

Each kit is tailored to target a spectrum of mutations associated with these genes, facilitating comprehensive mutation analysis. The dPCR PanCancer Kit BRAF V600 FAM (200) targets 8 BRAF V600 mutations in a single assay and the dPCR PanCancer Kit EGFRex19del FAM (200) targets 23 EGFR exon 19 deletions. Each PanCancer assay contains a reference assay for the human single-copy gene *AP3B1* to determine the number of genome copies in the sample and to control for dPCR efficiency.

Here, we present our initial data for various sample types, including blood, plasma, stool and FFPE samples. Through a meticulously optimized dPCR setup, we have achieved exceptional sensitivity and specificity, enabling the detection of multiple mutations in a single channel at allelic frequencies below 1%.


dPCR PanCancer Kit (RUO) assay design and function

dPCR PanCancer Kits (RUO) offer duplex assays to simultaneously detect multiple hallmark mutations in cancer-specific genes and a reference gene. Mutations are detected in the green channel and the reference gene is detected in the yellow channel. Combined with an optimized dPCR Master Mix, these assays reduce the cost and time required for the comprehensive analysis of hallmark mutations.

To comprehensively assess the performance of the dPCR PanCancer Kits (RUO) across varied sample types and methodologies, three reputable institutes conducted rigorous testing. The findings from all three independent evaluations validate the efficacy of the assay design, demonstrating robust compatibility with diverse samples and concordance with testing results from established industry standard methods.

List of mutations targeted by the PanCancer assays

Assay	Mutation aa	Mutation nucleotide	COSMIC ID
dPCR PanCancer Kit BRAF V600 FAM (200)	p.V600K	c.1799_1799delinsAA	COSV56057713
	p.V600R	c.1798_1799delinsAG	COSV56058419
	p.V600E	c.1799_1800delinsAA	COSV56059110
	p.V600E	c.1799T>A	COSV56056643
	p.V600D	c.1799_1800delinsAT	COSV56059623
	p.V600G	c.1799T>G	COSV56080151
	p.V600M	c.1798G>A	COSV56075762
	p.V600R	c.1798_1799delinsCG	COSV56288820
dPCR PanCancer Kit EGFRex19del FAM (200)	p.E745_E749del	c.2231_2247del	COSV51769442
	p.E746_A750delinsP	c.2235_2248delinsAATTC	COSV51817953
	p.E746_A750del	c.2235_2249del	COSV51765119
	p.E746_T751delinsP	c.2235_2251delinsAATTC	COSV51782151
	p.E746_T751delinsI	c.2235_2252delinsAAT	COSV51850034
	p.E746_A750del	c.2236_2250del	COSV51765066
	p.E746_T751delinsA	c.2237_2251del	COSV51769364
	p.E746_T751delinsV	c.2237_2252delinsT	COSV51769364
	p.E746_T751delinsVA	c.2237_2253delinsTTGCT	COSV51771891
	p.E746_S752delinsV	c.2237_2253delinsT	COSV51765862
	p.L747_A750delinsP	c.2238_2248delinsGCA	COSV51782279
	p.L747_T751delinsQ	c.2238_2253delinsGCA	COSV51863059
	p.E746_S752delinsD	c.2238_2255del	COSV51772418
	p.L747_E749del	c.2239_2247del	COSV51780076
	p.L747_A750delinsP	c.2239_2248delinsC	COSV51765099
	p.L747_T751delinsP	c.2239_2251delinsP	COSV51765856
	p.L747_S752del	c.2239_2256del	COSV51767308
p.L747_S752delinsQ	c.2239_2256delinsCAA	COSV51778874	
p.L747_T753delinsG	c.2239_2258delinsCA	COSV51785746	
p.L747_T753delinsS	c.2240_2251del	COSV51768180	
p.L747_T753del	c.2240_2254del	COSV51766247	
p.L747_A750delinsS	c.2240_2248del	COSV51810296	
p.L747_T753delinsS	c.2240_2257del	COSV51767961	



Field test results

These results from the Klinisk Biokemisk Afdeling are for experiments using the dPCR PanCancer Kits (RUO) with DNA extracted from stool samples and compared to NGS screening of DNA extracted from tumor tissue. DNA extraction for dPCR was performed on collected stool samples from donors with known and unknown status of BRAF V600 mutations and EGFR exon 19 deletions using extraction kits optimized for stool samples. Donors with known mutation status were previously analyzed for the presence of BRAF and EGFR mutations by applying NGS to DNA extracted from tumor tissue. In the dPCR reaction, 1–50 ng of extracted stool DNA was used.

Screening results for DNA from stool samples.

ID	Expected mutation	PanCancer Kit result
Stool 1	Unknown	BRAF V600 Mut
Stool 2	Unknown	BRAF V600 Mut
Stool 4	Unknown	BRAF V600 WT
Stool 5	BRAF V600 WT	BRAF V600 WT
Stool 6	Unknown	BRAF V600 WT
Stool 8	BRAF V600 Mut	BRAF V600 Mut
Stool 11	BRAF V600 WT	BRAF V600 Mut
Stool 12	Unknown	BRAF V600 Mut
Stool 1	EGFR WT	n.a.
Stool 3	EGFR WT	EGFR WT
Stool 5	EGFR WT	EGFR WT
Stool 6	EGFR WT	n.a.
Stool 7	EGFR WT	EGFR WT
Stool 8	EGFR WT	EGFR WT
Stool 9	EGFR WT	EGFR WT
Stool 10	EGFR WT	EGFR WT
Stool 11	EGFR WT	EGFR WT
Stool 12	EGFR WT	n.a.

Evaluation of positive predictive validity (PPV) and negative PV (NPV) for dPCR PanCancer Kits (RUO).

PanCancer Kit (dPCR from stool samples)	Standard method (NGS from tumor tissue)		PPV = 100%	NPV = 100%
	Pos	Neg		
Pos	1	1	5	3
Neg	0	8	3	2
	1	9		

Screening results for DNA from FFPE samples.

ID	Expected mutation	PanCancer Kit result
Sample 1	BRAF V600 Mut	BRAF V600 Mut
Sample 2	BRAF V600 WT	BRAF V600 WT
Sample 3	BRAF V600 Mut	BRAF V600 Mut
Sample 4	BRAF V600 WT	BRAF V600 WT
Sample 5	BRAF V600 Mut	BRAF V600 Mut
Sample 6	BRAF V600 WT	BRAF V600 WT
Sample 7	EGFR exon 19 deletion	EGFR exon 19 deletion
Sample 8	Unknown	EGFR exon 19 deletion

Evaluation of PPV and NPV for dPCR PanCancer Kits (RUO).


PanCancer Kit (dPCR from stool samples)	Standard method (NGS from tumor tissue)		PPV = 100%	NPV = 100%
	Pos	Neg		
Pos	5	0	3	2
Neg	0	3	3	2
	5	3		

Screening results for DNA from human reference cell lines.


Cell line	dPCR template input	Expected mutation	PanCancer Kit result
SK-MEL28	70 ng	BRAF V600E	BRAF V600 Mut
SK-MEL28	14 ng	BRAF V600E	BRAF V600 Mut
SK-MEL28	2.8 ng	BRAF V600E	BRAF V600 Mut
A-431	70 ng	BRAF V600 WT	BRAF V600 WT
A-431	14 ng	BRAF V600 WT	BRAF V600 WT
A-431	2.8 ng	BRAF V600 WT	BRAF V600 WT
C5693 P1 mix	70 ng	EGFR exon 19 deletion	EGFR exon 19 deletion
C5693 P1 mix	14 ng	EGFR exon 19 deletion	EGFR exon 19 deletion
C5693 P1 mix	2.8 ng	EGFR exon 19 deletion	EGFR exon 19 deletion
OVCAR-3	70 ng	EGFR exon 19 WT	EGFR 19 WT
OVCAR-3	14 ng	EGFR exon 19 WT	EGFR 19 WT
OVCAR-3	2.8 ng	EGFR exon 19 WT	EGFR 19 WT

Evaluation of PPV and NPV for dPCR PanCancer Kits (RUO).


PanCancer Kit (dPCR from cell line DNA)	Standard method (qPCR from cell line DNA)		PPV = 100%	NPV = 100%
	Pos	Neg		
Pos	2	0	2	2
Neg	0	2	2	2
	2	2		



Examples of 1D scatterplots for individual samples analyzed with the dPCR PanCancer Kit BRAF V600 FAM (200).




Examples of 1D scatterplots for individual samples analyzed with the dPCR PanCancer Kit BRAF V600 FAM (200).



Barplots of concentrations (measured copies/μL) for human DNA samples extracted from reference cell lines. Analyses were performed with the dPCR PanCancer Kits for BRAF (top) and EGFR (bottom). DNA input was titrated, and synthetic DNA (Gblocks) were used as the positive control (PC). Mean cp./μL values from replicates are given above the bars.

Complete workflow with blood samples

Testing the compatibility of the dPCR PanCancer Kits (RUO) with various DNA extraction protocols and workflows is essential prior to applying it in cancer mutation research. To demonstrate compatibility with one manual and two automated extraction protocols (QIAAsymphony® SP, EZ2®), human gDNA was extracted from blood samples from five healthy donors and subsequently analyzed using the dPCR PanCancer Kits (RUO) for BRAF (below left) and EGFR (below right). The results showed stable DNA recovery across the three extraction methods as measured with the AP3B1 reference genes and consistent amplification and quantification of spiked-in mutation templates for both dPCR PanCancer Kits (RUO). This experiment was performed at the QIAGEN GmbH facilities in Germany.



Detection of somatic mutations in genomic DNA extracted from blood and analyzed using the dPCR PanCancer Kits on the QIAcuity® platform. DNA extraction was performed on collected blood procured from five healthy donors (1–5) employing three distinct methodologies. Blood collection for donors 1 and 2 utilized PAXgene® Blood DNA Tubes, while donors 3–5 underwent blood collection using EDTA blood collection tubes. Manual extraction was executed using the QIAamp® DNA Blood Mini Kit, while automated extractions were conducted using an EZ2 with the EZ1.2® Blood 350 μL Kit and a QIAAsymphony SP with the QIAAsymphony DSP DNA Midi Kit. To simulate the presence of mutated DNA content, 2000 copies of synthetic DNA (Gblocks) were spiked into the dPCR reaction, resulting in an expected 2000 copies/μL. In the subsequent dPCR reactions, 1 μL of the eluate served as the input template. This quantity represented 1 μL of blood from the manual preparations, 3.5 μL of blood from the EZ2 preparations, and 2.5 μL of blood from the QIAAsymphony preparations. Error bars depict the standard deviation from 2 replicates.

Conclusions

Both dPCR PanCancer Kits (RUO) have potential applications in research prescreening samples, e.g., prior to next-generation sequencing (NGS), or research into monitoring cancer cells. The assays simultaneously assess multiple mutations, reducing time and costs and saving sample material. Additionally, this novel technology is adaptable for other cancer-associated genes, so similar assays can potentially be developed.

Overall, we have demonstrated that our dPCR PanCancer kits (RUO) provide a robust, fast and efficient technology to identify critical mutations, ultimately enhancing our understanding of BRAF- and EGFR-driven cancers. There is a strong concordance between dPCR PanCancer Kit (RUO) assay results and results from industry standard methods for cancer mutation detection. The assay design is compatible with various sample types.

The dPCR PanCancer Kit is for research use only. Not for the diagnosis, prevention, or treatment of a disease. For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit instructions for use or user operator manual. QIAGEN instructions for use and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services (or your local distributor).

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAcuity®, QIAAsymphony®, EZ1.2®, EZ2® (QIAGEN GmbH); PAXgene® (PreAnalytix GmbH). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. QPRO 6708 04/2024 © 2024 QIAGEN, all rights reserved.