

## Quick-Start Protocol

# QIAwave DNA Blood & Tissue Kit

The QIAwave DNA Blood & Tissue Kit (cat.no. 69554 and cat. no 69556) can be stored at room temperature (15–25°C) for up to 1 year after delivery.

### Further information

- QIAwave DNA Blood & Tissue Handbook: [www.qiagen.com/HB-2987](http://www.qiagen.com/HB-2987)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)
- The QIAwave DNA Blood & Tissue Kit can be automated on the QIAcube Connect using the DNeasy Blood & Tissue Kit protocols that be downloaded at [www.qiagen.com/qiacubeprotocols](http://www.qiagen.com/qiacubeprotocols).

### Notes before starting

- Perform all centrifugation steps at room temperature (15–25°C).
- Redissolve any precipitates in Buffer AL and Buffer ATL.
- Equilibrate frozen tissue or cell pellets to room temperature.
- Preheat an incubator to 56°C.
- Refer to the handbook for pretreatment of fixed tissue, insect, bacterial, or other materials.
- Preassemble DNeasy® Mini Spin Columns with Waste Tubes.
- **Preparation of final buffers from concentrates:** Transfer the entire volume of buffer concentrates from the 2 mL tube or 15 mL bottle into a glass bottle appropriate for the final volume (Table 1), either by using a pipette or by pouring. Add ultrapure water

and/or ethanol (96–100%) according to Table 1. To label the glass bottle, use the enclosed label and transfer it onto the glass bottle.

**Table 1. Preparation of Buffer concentrates**

Kit (cat. no.)	Final buffer	Buffer*	Ultrapure water	Ethanol (96-100%)	Final volume
69554	AW1	AW1/C	–	20 mL	35 mL
	AW2	AW2/C	15 mL	40 mL	56.5 mL
	AE	AE/C	22 mL	–	24 mL
69556	AW1	AW1/C	–	130 mL	228 mL
	AW2	AW2/C	60 mL	160 mL	226 mL
	AE	AE/C	110 mL	–	120 mL

## Procedure

1. **Tissue:** Cut tissue ( $\leq 10$  mg spleen or  $\leq 25$  mg other tissue) into small pieces, and place in a 1.5 ml microcentrifuge tube (not provided). For rodent tails, use 1 (rat) or 2 (mouse) 0.4–0.6 cm lengths of tail. Add 180  $\mu$ l Buffer ATL. Add 20  $\mu$ l proteinase K, mix by vortexing and incubate at 56°C until completely lysed. Vortex occasionally during incubation. Vortex 15 s directly before proceeding to step 2.

**Nonnucleated blood:** Pipet 20  $\mu$ l proteinase K into a 1.5 ml or 2 ml microcentrifuge tube (not provided). Add 50–100  $\mu$ l anticoagulant-treated blood. Adjust volume to 220  $\mu$ l with PBS. Proceed to step 2.

**Nucleated blood:** Pipet 20  $\mu$ l proteinase K into a 1.5 ml or 2 ml microcentrifuge tube (not provided). Add 5–10  $\mu$ l anticoagulant-treated blood. Adjust volume to 220  $\mu$ l with PBS. Proceed to step 2.

**Cultured cells:** Centrifuge a maximum of  $5 \times 10^6$  cells for 5 min at  $300 \times g$  (190 rpm). Resuspend in 200  $\mu$ l PBS. Add 20  $\mu$ l proteinase K. Proceed to step 2.

2. Add 200  $\mu$ l Buffer AL. Mix thoroughly by vortexing. Incubate blood samples at 56°C for 10 min.

3. Add 200  $\mu$ l ethanol (96–100%). Mix thoroughly by vortexing.
  4. Pipet the mixture into a DNeasy Mini spin column placed in a 2 ml Waste Tube (supplied). Centrifuge at  $\geq 6000 \times g$  (8000 rpm) for 1 min. Discard the flow-through and reuse the Waste Tube.
  5. Add 500  $\mu$ l Buffer AW1, and centrifuge for 1 min at  $\geq 6000 \times g$ . Discard the flow-through and reuse the Waste Tube.
  6. Add 500  $\mu$ l Buffer AW2 and centrifuge for 3 min at  $20,000 \times g$  (14,000 rpm). Discard the flow-through and Waste Tube.
  7. Transfer the spin column to a new 1.5 ml or 2 ml microcentrifuge tube (not provided).
  8. Elute the DNA by adding 200  $\mu$ l Buffer AE to the center of the spin column membrane. Incubate for 1 min at room temperature. Centrifuge for 1 min at  $\geq 6000 \times g$ .
- Optional:** Repeat step 8 for increased DNA yield.

## Document Revision History

Date	Changes
01/2022	Initial release
06/2023	Addition of cat. no. 69554 and necessary procedures Updated QSP according to new brand template



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, DNeasy® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

06/2023 HB-2988-002 © 2023 QIAGEN, all rights reserved.