PartoSure® Test

Instructions For Use

Assess the Risk of Preterm Birth

EN English	3
BG Български	4
CS Český	5
DA Dansk	6
DE Deutsch	7
ES Español	8
FI Finish	9
FR Français	10
HR Hrvatski	11
HU Magyar	12
ID Indonesian	13
IT Italiano	14
KK Kazakh	15
LT Lietuvių	16
LV Latviešu	17
MK Macedonian	18
NL Nederlands	19
NO Norsk	20
PL Polish	21
PT Portuguese (Basil berian)l	22
PT Português	23
RO Română	24
RU Русский	25
SK Slovenčina	26
SR Srpski	27
SV Svenska	28
TR Türkçe	29

QIAGEN
19300 Germantown Road
Germantown, MD 20874 USA
+1-800-426-8157
www.QIAGEN.com

QIAGEN GmbH
QIAGEN Strasse 1
40724 Hilden, GERMANY

TTDT-1-5-OUS
TTDT-1-20-OUS

1108422 Rev. 04

PATENT AND TRADEMARK INFORMATION

©2018 QIAGEN

The PartoSure Test and its use are covered by one or more of the following patents granted or licensed to QIAGEN and/or its subsidiaries; US Patent No. 7,709,272; 8,114,610; 9,891,233, 9,429,580; corresponding foreign patents and other patents pending.

K-Y is a registered trademark of Reckitt Benckiser.

Surgilube is a registered trademark of Fougera.

Monistat is a registered trademark of Insight Pharmaceuticals.

Betadine is a registered trademark of Purdue Products.

PartoSure and associated logos are trademarks and/or registered trademarks of QIAGEN and/or its subsidiaries in the United States and/or other countries.

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INTENDED USE

The PartoSure Test is a rapid, non-instrumented, qualitative immunochromatographic test for the *in vitro* detection of placental alpha microglobulin-1 (PAMG-1) in vaginal secretions of pregnant women. The device is designed as an aid to rapidly assess the risk of preterm delivery in ≤ 7 or ≤ 14 days from the time of cervicovaginal sample collection in pregnant women with signs and symptoms of early preterm labor, intact amniotic membranes, and minimal cervical dilatation (≤ 3 cm), sampled between 20 weeks, 0 days and 36 weeks, 6 days gestation.

SUMMARY AND EXPLANATION OF THE TEST

An accurate risk assessment of preterm birth is clinically important among pregnancies with threatened preterm labor. This is particularly true with respect to both the administration of corticosteroids, which have an optimal benefit within 7 days of administration,¹ as well as the transfer of patients to a tertiary care center capable of caring for the birth of a premature infant. Clinical evaluation alone, including the measurement of cervical length and dilatation, is not sufficiently predictive of imminent delivery.² Similarly, currently available biomarker tests, such as that based on the detection of fetal fibronectin, have extremely poor predictive values for imminent delivery.³ Thus, the accurate prediction of time to delivery in patients with threatened preterm labor remains an important clinical concern. The high negative and positive predictive values of the PartoSure Test may allow for the timely identification of patients that would benefit from available interventions such as antenatal corticosteroids. Additionally, a more accurate identification of these patients will likely reduce the costs associated with unnecessary admissions caused by the poor positive predictive values of existing methods. The test kit is a self-contained system that provides rapid, qualitative results without outside instrumentation or an invasive speculum examination.

PRINCIPLE OF THE TEST

The PartoSure Test is a lateral flow, immunochromatographic assay designed to identify the presence of human placental alpha microglobulin-1 (PAMG-1). The test employs monoclonal antibodies sufficiently sensitive to detect 1 ng/ml of PAMG-1. For the analysis, a sample of cervicovaginal discharge collected by vaginal swab is extracted into a solvent. The presence of PAMG-1 antigen is then detected by inserting a lateral-flow test strip into the vial. The sample flows from an absorbent pad to a nitrocellulose membrane, passing through a reactive area containing monoclonal anti-PAMG-1 antibodies conjugated to a gold particle. The antigen-antibody complex flows to the test region where it is immobilized by a second anti-PAMG-1 antibody. This event leads to the appearance of the test line. Unbound antigen-antibody complexes continue to flow along the test strip and are immobilized by a second antibody. This leads to the appearance of the internal control line.

REAGENTS AND COMPONENTS

The PartoSure Test kit includes the following components: the PartoSure Test strip in a foil pouch with desiccant, a sterile flocked vaginal swab, and a plastic vial with solvent solution (0.9% NaCl, 0.05% Na₂S₂O₅, 0.01% Triton X100). No additional materials are required to perform the PartoSure Test.

PRECAUTIONS AND WARNINGS

- Use of any swab or solvent solution other than the one provided with the test kit is prohibited.
- Specimens should be collected prior to collection of culture specimens. Collection of vaginal specimens for microbiologic culture frequently requires aggressive collection techniques that may abrade the cervical or vaginal mucosa and may potentially interfere with sample preparation.
- Care must be taken not to contaminate the swab or cervicovaginal secretions with lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser). These substances may interfere with absorption of the specimen by the swab or with the antibody-antigen reaction of the PartoSure Test and lead to invalid test results.
- If it is suspected that the patient has applied a topical disinfectant (e.g. Monistat®, miconazole nitrate cream) to the vaginal area within 24 hours, delay specimen collection until 24 hours from application of the topical disinfectant have passed as these products can lead to false-negative test results.
- The PartoSure Test is not intended for use in women with moderate or gross vaginal bleeding. The presence of vaginal bleeding may contribute to difficulty in interpreting the PartoSure Test result. Testing a moderately to grossly bloody sample may lead to false-positive results. If upon visual examination you are concerned about the presence of moderate or gross vaginal blood, it is recommended that the sample be collected following the cessation of active vaginal bleeding.
- Specimens should not be obtained from patients with suspected or known placental abruption or placenta previa.
- Do not use the PartoSure Test after the expiration date, which is printed on the product packaging.
- PartoSure Test kit components are for single-use only.
- Specimens not tested within 24 hours of collection must be stored refrigerated at 2° to 8°C and tested within five days of collection.
- The PartoSure Test is for *in vitro* diagnostic use only and no component of the test should be taken internally.
- Do not use the kit if the swab or test strip package integrity is compromised or if the solvent vial has leaked.
- Do not bend or fold the test strip or the foil pouch with the test strip in it; doing so may damage the strip and lead to inaccurate results.
- The PartoSure Test kit includes a plastic vial with solvent solution (0.9% NaCl, 0.05% Na₂S₂O₅, 0.01% Triton X100). Sodium azide may react with plumbing to form potentially explosive metal azides. Avoid contact with skin, eyes, and clothing. In case of contact with any of these reagents, wash area thoroughly with water. If disposing of this reagent, always flush the drain with large volumes of water to prevent azide build-up.
- Safety precautions should be observed when collecting, handling, and disposing of test samples.
- Used PartoSure Test kits are biohazardous; dispose of all components using necessary precautions.


STORAGE AND STABILITY

Store the PartoSure Test kit in a dry place at 15 to 25°C (59 to 77°F). The test should not be frozen. When stored in the foil pouch at the recommended temperature, the test is stable until the expiration date printed on the pouch. The PartoSure Test strip should be used within six (6) hours after removing the test strip from the foil pouch.


QUALITY CONTROL

The PartoSure Test strip contains an internal procedural control mechanism that ensures analytical functionality. The appearance of one or two lines in the results region of the test strip verifies the integrity of the test procedure and components.

TEST PROCEDURE

1. Take the solvent vial by its cap and ensure all liquid in the vial has dropped to the bottom. Open the solvent vial and place it in a vertical position.
2. To collect a sample from the vagina, use only the sterile flocked swab provided with the PartoSure Test kit. Remove the swab from its package following the instructions on the packaging. The tip of the swab should not touch anything prior to insertion into the vagina. Hold the swab by the middle of its shaft and, while the patient is lying on her back, carefully insert the tip of the swab into the vagina until the fingers contact the skin (no more than 5-7 cm deep). Withdraw the swab from the vagina **after 30 seconds**.
3. After the swab has been removed from the vagina, immediately place the tip into the provided solvent vial and rinse by rotating for **30 seconds**.
4. Remove the swab from the vial and dispose of it.
5. Tear open the foil pouch at the tear notches and remove the PartoSure Test strip.
6. Insert the white end of the test strip (marked with arrows facing downward) into the vial with solvent.
7. **Remove the test strip from the vial if two lines are clearly visible in the test region or after 5 minutes sharp.** Read the results by placing the test strip on a clean, dry, and flat surface in a well-lit environment via either natural or fluorescent lighting. A positive result is indicated by two lines in the test region, while a negative result is indicated by a single line in the test region. Do not read or interpret the results after 10 minutes have passed since inserting the test strip into the vial.

Positive: Two Lines	Negative: One Control Line	Invalid: No Lines or Test Line Only
Imminent delivery within ≤ 7 or ≤ 14 days is highly likely	Imminent delivery within ≤ 7 or ≤ 14 days is highly unlikely	Results not valid; retest

The intensity of the lines may vary; the test result is valid even if the lines are faint or uneven. Do not interpret the test result based on the intensity of the lines.

LIMITATIONS OF THE TEST

- The PartoSure Test result should not be interpreted as absolute evidence for the presence or absence of a process that will result in delivery ≤ 7 or ≤ 14 days from specimen collection.
- The PartoSure Test result should always be used in conjunction with information available from the clinical evaluation of the patient and other diagnostic procedures such as cervical examination, assessment of uterine activity, and evaluation of other risk factors.
- These instructions for use must be followed exactly; failure to do so may lead to inaccurate results.
- PartoSure Test performance has been characterized from specimens taken from the vaginal cavity. Samples obtained from other locations should not be used. A speculum examination is not required.
- Results should be interpreted with caution when a specimen is obtained from a patient with unconfirmed gestational age.
- PartoSure Test results are qualitative and not quantitative. No quantitative interpretation should be made based on the strength of the test or control lines.
- The PartoSure Test should only be used in patients with signs and symptoms of preterm labor.
- Invalid test results may occur if lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser) have been used by the patient.

EXPECTED VALUES

PartoSure detects trace amounts of human PAMG-1. PAMG-1 concentrations greater than 1 ng/ml indicate an elevated risk of delivering within the next 7 days. The cutoff of 1 ng/ml was validated in an international multicenter clinical trial evaluating the association of PAMG-1 present in cervicovaginal secretions and the test-to-delivery time interval.

PERFORMANCE CHARACTERISTICS

Precision and reproducibility were determined using three lots of the PartoSure Test at three intended-use sites by three different intended users at each site. Five replicates of seven different PAMG-1 concentrations above and below the limit of detection of the PartoSure Test were used, including an absolute zero (0.0 ng/ml), a low-negative (0.2 ng/ml), a high-negative (0.5 ng/ml), the limit of detection (1.0 ng/ml), a low-positive (2.0 ng/ml), and two concentrations within the C₅-C₉s interval for the device (0.7 ng/ml and 0.9 ng/ml). The results of this study demonstrated 100% negative readings at ≤ 0.5 ng/ml and 100% positive readings at ≥ 1 ng/ml PAMG-1 levels by intended users.

The clinical performance of the PartoSure Test was determined in an international multi-site study. Patients who (i) reported signs, symptoms, or complaints suggestive of preterm labor, (ii) were between 20 and 36⁺ weeks of gestation, and (iii) had clinically intact membranes were invited to participate in the trial. In total, 203 patients were included in the final analysis. The following performance estimates with corresponding 95% confidence intervals were calculated:

Performance Measure	Delivery Within	
	≤ 7 Days	≤ 14 Days
Sensitivity (SN) (95% CI)*	80% (63.1% – 91.6%)	63% (47.4% – 76.1%)
Specificity (SP) (95% CI)*	95% (90.1% – 97.5%)	96% (90.9% – 98.2%)
Negative Predictive Value (NPV) (95% CI)*	96% (91.5% – 98.3%)	89% (83.4% – 93.5%)
Positive Predictive Value (PPV) (95% CI)*	76% (58.8% – 88.2%)	81% (64.8% – 92.0%)

*95% confidence intervals (CI) computed via the Clopper-Pearson procedure
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

CROSS REACTIVITY

The PartoSure Test was evaluated using a panel of potentially cross-reactive protein substances likely to be found in vaginal secretions, including human chorionic gonadotropin, trophoblastic beta-2 glycoprotein, human placental lactogen, alpha-fetoprotein, IGFBP-3, and human serum albumin. Four (4) replicates of each sample containing the potentially cross-reactive substance were tested using a high-negative PAMG-1 sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml). Each potentially cross-reactive substance was tested at the highest concentration of substance that was considered clinically relevant. None of the potentially cross-reactive substances tested demonstrated cross-reactivity with the PartoSure Test.

INTERFERENCE STUDIES

Lubricants, Disinfectants, Antiseptics, Soaps & Creams
Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against lubricant jelly, miconazole nitrate disinfectant cream, antiseptic cleanser, body wash soap, and vaginal cream. While all high-negative replicates were correctly identified in the presence of miconazole nitrate, low-positive replicates were not. The presence of lubricating jelly led to several invalid results for both high-negative and low-positive replicates, though in cases where a negative or positive result was reported, the result was correct. The presence of antiseptic cleanser led to invalid test results in all replicates. Therefore, care must be taken not to contaminate the swab or cervicovaginal secretions with lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser). If it is suspected that the patient has applied a topical disinfectant (e.g. Monistat®, miconazole nitrate cream) to the vaginal area within 24 hours, delay specimen collection until 24 hours from application have passed as these products can lead to false-negative test results. The presence of soap or cream did not interfere with the PartoSure Test.

Pharmacological Agents

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against 10 non-biologically based interfering substances including: 17-OH progesterone (50 µg/ml), ampicillin (152 µmol/L), cephalaxin (337 µg/ml), erythromycin (81.6 µmol/L), gentamycin (21 µmol/L), dexamethasone (1.53 µmol/L), magnesium sulfate (50 µg/ml), oxytocin (58 µU/ml), terbutaline (1 mg/ml), and ritodrine (100 µg/ml). None of these non-biological substances interfered with the PartoSure Test.

Vaginal Bacterial Pathogens

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against 3 vaginal bacterial pathogens, including Gardnerella vaginalis, Candida albicans, and Trichomonas vaginalis. None of these vaginal infection pathogens interfered with the PartoSure Test.

Maternal Bleeding

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against ten (10) individual maternal blood samples at the three (3) lowest admixture levels determined to represent "trace," "moderate," or "gross" levels of maternal bleeding on the vaginal collection swab. Moderate or gross vaginal bleeding may therefore contribute to difficulty in interpreting the PartoSure Test result and may lead to false-positive results.

Semen & Urine

A high-negative (0.2 ng/ml) sample and a low-positive sample (2.0 ng/ml) were tested against 10 individual samples of semen and 10 individual samples of maternal urine. Neither semen nor maternal urine demonstrated interference with the PartoSure Test.

BIBLIOGRAPHY

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. *Clinical management guidelines for obstetrician-gynecologists.* Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement. NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. <http://www.cochrane.org>

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*95% intervaly spolehlivosti (CI) byly vypočteny pomocí postupu Clopper-Pearson Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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LITERATURA

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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LITTERATURLISTE

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin nr. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, 4. udgave. Art. nr.: CD006843. [Systematic review]

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*95%-ige Konfidenzintervalle (KI) berechnet nach Clopper-Pearson
Nikolova T, Bayev Q, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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LITERATUR

1. Antenatal Corticosteroids Revisited: Repeat Courses. *NIH Consensus Statement*: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin* no. 127: Management of Preterm Labor. *Clinical management guidelines for obstetrician-gynecologists*. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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LÄHTEET

1. Antenatal Corticosteroids Revisited: Report Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308–1316.
3. Berghella V, Hayes E, Vasthine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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*Intervalles de confiance (IC) 95 % calculés selon la procédure de Clopper-Pearson
 Nikolaova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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BIBLIOGRAPHIE

1. Antenatal Corticosteroids Revisited: Report Consensus. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintini J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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
*Intervali pouzdanosti (CI) od 95 % izračunati su Clopper-Pearsonovim postupkom Nikolova T, Bayev O, Nikolaeva N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. lipanj 2015;43(4):395-402.

KRIŽNA REAKTIVNOST

Test PartoSure procijenjen je pomoću panela potencijalno križno reaktivnih proteinskih tvari koje se mogu pronaći u vaginalnim uzorcima, uključujući humani korionski gonadotropin, trofoblastični beta-2 glikoprotein, humani placentalni laktogen, alfa-fetoprotein, IGFBP-3 i humani serumski albumin. Testirano je četiri (4) ponavljanja svakog uzorka koji sadrži potencijalno križno reaktivnu tvar primjenom uzorka s visoko negativnom razinom spoja PAMG-1 (0,2 ng/ml) i uzorka s nisko pozitivnom razinom spoja (2,0 ng/ml). Svaka potencijalno križno reaktivna tvar testirana je pri najvećoj koncentraciji tvari koja se smatrala klinički relevantnom. Nijedna od testiranih potencijalno reaktivnih tvari nije pokazala križnu reaktivnost s testom PartoSure.

ISPITIVANJA INTERFERENCIJE

Lubrikanti, sredstva za dezinfekciju, antiseptici, sapuni i kreme Četiri (4) ponavljanja visoko negativnog uzorka (0,2 ng/ml) i nisko pozitivnog uzorka (2,0 ng/ml) testirano je kako bi se odredila interferencija s lubrikantom u gelu, dezinfekcijskom kremom s mikonazolom nitratom, antiseptičkim sredstvom za ispiranje, sapunom za pranje tijela i vaginalnom kremom. Iako su sva visoko negativna ponavljanja točno identificirana u prisutnosti mikonazol nitrata, nisko pozitivna ponavljanja nisu. Prisutnost lubrikanta u gelu uzrokovala je nekoliko nevažeci rezultata za visoko negativna i nisko pozitivna ponavljanja, iako je u slučajevima u kojem je dobiven negativan ili pozitivan rezultat, taj rezultat bio točan. Prisutnost antiseptičkog sredstva za ispiranje dovela je do nevažeci rezultata testa za sva ponavljanja. Stoga, treba voditi brigu o tome da se štapić za uzimanje brisa ili cervovaginalni isjettak ne kontaminiraju lubrikantima ili antisepticima (npr. lubrikanti u gelu K-Y® ili Surgilube®, sredstvo za ispiranje Betadine®). Ako postoji sumnja da je pacijentica lokalno primijenila sredstvo za dezinfekciju (npr. Monistat®, kremu s mikonazolom nitratom) na vaginalno područje unutar posljednja 24 sata, odgovodite prikupljanje uzorka dok ne prođe 24 sata od primjene jer takvi proizvodi mogu uzrokovati lažno negativne rezultate testa. Prisutnost sapuna ili kreme nije utjecala na test PartoSure.

Farmakološka sredstva

Četiri (4) ponavljanja visoko negativnog uzorka (0,2 ng/ml) i nisko pozitivnog uzorka (2,0 ng/ml) testirano je kako bi se odredila interferencija s 10 interferirajućih tvari koje nisu na biološkoj bazi uključujući: 17-OH progesteron (50 µg/ml), ampicilin (152 µmol/L), cefaleksin (337 µmol/L), eritromicin (81,6 µmol/L), gentamicin (21 µmol/L), deksametazon (1,53 µmol/L), magnezijev sulfat (50 µg/ml), oksticion (58 µg/ml), terbutalin (1 mg/ml) i ritodrin (100 µg/ml). Nijedna od tih nebioloških tvari nije interferirala s testom PartoSure.

Vaginalni bakterijski patogeni

Četiri (4) ponavljanja visoko negativnog uzorka (0,2 ng/ml) i nisko pozitivnog uzorka (2,0 ng/ml) testirano je kako bi se odredila interferencija 3 vaginalna bakterijska patogena, konkretno Gardnerella vaginalis, Candida albicans i Trichomonas. Nijedna od tih vaginalnih bakterijskih patogena nije interferirala s testom PartoSure.

Maternalno krvarenje

Četiri (4) ponavljanja visoko negativnog uzorka (0,2 ng/ml) i nisko pozitivnog uzorka (2,0 ng/ml) testirano je kako bi se odredila interferencija s deset (10) individualnih uzorka maternalne krvi na tri (3) najniže razine mješavine za koje je utvrđeno da predstavljaju razinu maternalnog krvarenja „u tragovima“, umjerenju razinu odnosno „obilnu“ razinu na štapiću za uzimanje vaginalnog brisa. Umjerenju ili obilno vaginalno krvarenje stoga može doprinijeti poteškoćama s tumačenjem rezultata testa PartoSure i uzrokovati lažno pozitivne rezultate.

Sperma i urin

Visoko negativni (0,2 ng/ml) i nisko pozitivni uzorak (2,0 ng/ml) testirani su kako bi se odredila interferencija s 10 individualnih uzorka sperme i 10 individualnih uzorka urina trudnice. Nije se pokazala interferencija niti sperme niti urina trudnice s testom PartoSure.

BIBLIOGRAFIJA

1. Antenatal Corticosteroids Revised: Report Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Vastine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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BIBLIOGRÁFIA

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308–1316.
3. Berghella V, Hoyer E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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
*95% interval kepercayaan (CI/Confidence Intervals) terkomputasi melalui prosedur Clopper-Pearson
Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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*Intervalli di confidenza (IC) del 95% calcolati con la procedura di Clopper-Pearson. Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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BIBLIOGRAFIA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrical-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Berghele V, Hayes E, Vistinne J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119: 1308-1316.

- Berghella V, Hayes E, Vitinsine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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* betrouwbaarheidsinterval (CI) van 95% berekend met de Clopper-Pearson-methode Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art.Nr.: CD006843. (Systematic review)

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

*95 % konfidensintervaller (CI) beregnet via Clopper-Perison-prosedyren Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* Juli 2015;43(4):395-402.

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BIBLIOGRAFI

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hobas E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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BIBLIOGRAFIA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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* Intervalos de confiança de 95% (CI) foram calculados através do procedimento de Clopper-Pearson
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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BIBLIOGRAFIA

1. Antenatal Corticosteroids Revisited: Repeat Courses. *NIH Consensus Statement*; NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin* no. 127: Management of Preterm Labor. *Clinical management guidelines for obstetrician-gynecologists*. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visamine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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*Interval de încredere (II) 95% calculat prin procedura Clopper-Pearson
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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BIBLIOGRAFIE

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins—Obstetrics. ACOG Practice Bulletin no. 122: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012;119(1):1308-1316.
- Borghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Annual Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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1. Antenatal Corticosteroids Revisited: Repeat Courses. *NIH Consensus Statement*. NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin* no. 127: Management of Preterm Labor. *Clinical Management guidelines for obstetrician-gynecologists*. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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

*95% interval pouzdanosti (IP) izračunat pomoću Klopper-Pisonovog (Clopper-Pearson) postupka
 Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

UNAKRSNA REAKTIVNOST

Procena PartoSure testa je izvršena pomoću panela potencijalno unakrsno reaktivnih proteina koji su česti u vaginalnom uzorcima, uključujući humani horionski gonadotropin, beta-2 glikoprotein trofoblasta, humani placentalni laktogen, alfafetoprotein, IGFBP-3 i humani serumski albumin. Četiri (4) replikata svakog uzorka sa potencijalno unakrsno reaktivnom supstancom su testirana pomoću uzorka PAMG-1 visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml). Svaka potencijalno unakrsno reaktivna supstanca je testirana pri najvećoj koncentraciji supstance koja je ocenjena kao klinički relevantna. Nijedna testirana potencijalno unakrsno reaktivna supstanca nije iskazala unakrsnu reaktivnost sa PartoSure testom.

STUDIJE SMETNJI

Lubrikanti, sredstva za dezinfekciju, antiseptici, sapuni i kreme
 Četiri (4) replikata uzorka visoko negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na lubrikant, krem sa mikonazol-nitratom, antiseptik, sapun za telo i vaginalnu kremu. Dok je prisustvo svih replikata visoke negativne vrednosti utvrđeno pri korišćenju mikonazol-nitrate, to se nije desilo u slučaju replikata niske negativne vrednosti. Prisustvo lubrikanta proizvelo je nekoliko nevažećih rezultata za replikate sa visokom negativnom i niskom pozitivnom vrednošću, ali je rezultat bio tačan u slučaju prijavljivanja negativnog ili pozitivnog rezultata. Prisustvo antiseptika je proizvelo nevažeće rezultate testa u svim replikatima. Dakle, vodite računa da ne kontaminirate štapić za uzorkovanje ili cervikovaginalni iscedak lubrikantima ili antisepticima (npr. lubrikant K-Y® ili Surgilube®, antiseptik Betadine®). Ako sumnjate da je pacijentkinja nanela lokalno sredstvo za dezinfekciju (npr. Monistat®, krem sa mikonazol-nitratom) u područje vagine tokom prethodna 24 sata, odložite prikupljanje uzorka 24 sata nakon primene jer ovi proizvodi mogu da proizvedu lažno negativne rezultate testa. Prisustvo sapuna ili kreme nije uticalo na PartoSure test.

Farmakološka sredstva

Četiri (4) replikata uzorka visoko negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na 10 nebioloških interferentnih supstanci, uključujući: 17-OH progesteron (50 µg/ml), ampicilin (152 µmol/l), cefaleksin (337 µmol/l), eritromicin (81,6 µmol/l), gentamicin (21 µmol/l), deksametazon (1,53 µmol/l), magnezijum sulfat (50 µg/ml), oksitocin (58 µU/ml), terbutalin (1 mg/ml) i ritodrin (100 µg/ml). Nijedna od ovih nebioloških supstanci nije uticala na test.

Vaginalni bakterijski patogeni

Četiri (4) replikata uzorka visoko negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na 3 vaginalna bakterijska patogena, uključujući Gardnerella vaginalis, Candida albicans i Trichomonas vaginalis. Nijedan od ovih vaginalnih infektivnih patogena nije uticao na test.

Krvarenje trudnice

Četiri (4) replikata uzorka visoko negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na deset (10) pojedinačnih uzoraka krvi trudnice u tri (3) najniža nivoa mešavine koji su predstavljali „mestimično“, „umereno“ ili „obilno“ krvarenje trudnice na štapiću za uzorkovanje iz vagine. Umereno ili obilno vaginalno krvarenje može dovesti do poteškoća pri tumačenju rezultata PartoSure testa i posledično do lažno pozitivnih rezultata.

Semena tečnost i urin

Uzorak visoko negativne vrednosti (0,2 ng/ml) i uzorak niske pozitivne vrednosti (2,0 ng/ml) su testirani sa 10 pojedinačnih uzoraka semene tečnosti i 10 pojedinačnih uzoraka urina. Ni semena tečnost ni urin nisu uticali na rezultate PartoSure testa.

LITERATURA

- Antenatal Corticosteroids Revisited: Report Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrical-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Borghello V, Hayes E, Vasthine J, Baxter JK. Fetal Fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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*95 % konfidensintervall (CI) beräknat via Clopper-Pearson-proceduren
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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REFERENSER

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Artnr.: CD006843. <http://dx.doi.org/10.1002/rev.1>

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