

# AdnaTest ProstateCancerSelect

AdnaTest ProstateCancerSelect must be stored at 2–8°C. The components must not be used beyond the expiration date.

## Further information

- *AdnaTest ProstateCancerSelect and ProstateCancerDetect Handbook:*  
**[www.qiagen.com/HB-2410](http://www.qiagen.com/HB-2410)**
- Safety Data Sheets: **[www.qiagen.com/safety](http://www.qiagen.com/safety)**
- Technical assistance: toll-free 00800-22-44-6000 or **[support.qiagen.com](mailto:support.qiagen.com)**

## Protocol

1. Resuspend ProstateSelect Beads thoroughly and transfer 100 µl for each blood sample into a 1.5 ml reaction tube.
2. Wash ProstateSelect Beads with 3 x 1 ml PBS.
3. Resuspend ProstateSelect Beads in 100 µl PBS per blood sample.
4. When using standard EDTA tubes, transfer 5 ml whole blood into a 15 ml Collection Tube (provided).  
When using ACD-A blood in a BD Vacutainer® ACD-A Tube, transfer 5 ml of blood into an AdnaTube.
5. Add 100 µl of washed ProstateSelect Beads to each blood sample.
6. Incubate for 30 min at room temperature under tilting and rotation at approximately 5 rpm.
7. Place tube for 3 min in the AdnaMag-L to separate the beads.  
Release any blood drops captured in the cap by swinging the AdnaMag-L downwards.

- 
8. Remove blood supernatant.
  9. Wash beads with 3 x 5 ml PBS.
  10. Resuspend beads in 1 ml PBS and transfer into a new 1.5 ml reaction tube (not provided).
  11. Separate beads in the AdnaMag-S and remove supernatant.
  12. Resuspend beads in 200 µl AdnaTest Lysis/Binding Buffer by pipetting at least five times.
  13. Place reaction tubes into the AdnaMag-S and transfer supernatant into a new 1.5 ml reaction tube (provided).
  14. Continue immediately with the AdnaTest ProstateCancerDetect or store at  $-20^{\circ}\text{C}$  for a maximum of 2 weeks.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup> (QIAGEN Group); Vacutainer<sup>®</sup> (Becton Dickinson and Company). 1106858 03/2017 HB-2411-001  
© 2017 QIAGEN, all rights reserved.