

# EpiTect<sup>®</sup> Plus DNA Bisulfite Kit – Protocol 1

MinElute<sup>®</sup> DNA spin columns, DNA Protect Buffer and Buffer BD from the EpiTect Plus DNA Bisulfite Kit (cat. no. 59124) should be stored at 2–8°C. All other buffers and Bisulfite Mix should be stored at room temperature (15–25°C) for up to 6 months if not otherwise stated on label. Dissolved Bisulfite Mix can be stored at –20°C for up to 4 weeks.

## Further information

- *EpiTect Plus Bisulfite Conversion Handbook*: [www.qiagen.com/HB-0388](http://www.qiagen.com/HB-0388)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](mailto:support.qiagen.com)

## Notes before starting

- Add 30 ml ethanol (96–100%) to Buffer BW and store at room temperature (15–25°C).
- Add 27 ml ethanol (96–100%) to Buffer BD and store at 2–8°C.
- Add 310 µl RNase-free water to carrier RNA and store in aliquots at –20°C.
- Carrier RNA (step 1) is not necessary if >100 ng DNA are used.

## Bisulfite conversion of DNA

1. Add 800 µl RNase-free water to each aliquot of Bisulfite Mix needed, and vortex until Bisulfite Mix is completely dissolved. This may take up to 5 min. Dissolving the Bisulfite Mix may require heating the solution to 60°C.
2. Set up the bisulfite reactions in 200 µl PCR tubes according to Table 1. Add each component in the order listed.
3. Close the PCR tubes and mix the bisulfite reactions thoroughly. DNA Protect Buffer should turn blue indicating sufficient mixing and correct pH.

4. Program the thermal cycler according to Table 2. Use a cycler with a heated lid. If using a thermal cycler that does not allow you to enter the reaction volume (140  $\mu$ l), set the instrument to the largest volume setting available.
5. Place the PCR tubes in the thermal cycler and start the incubation. Converted DNA can be left in the thermal cycler overnight without loss of performance.
6. Proceed to Protocol 2 “Cleanup of converted DNA” included in the EpiTect Plus DNA Bisulfite Kit.

**Table 1. Bisulfite reaction setup**

Component	High concentration samples (1 ng – 2 $\mu$ g) Volume per reaction ( $\mu$ l)	Low concentration samples (1 ng – 500 ng) Volume per reaction ( $\mu$ l)
DNA solution	Variable* (maximum 20 $\mu$ l)	Variable† (maximum 40 $\mu$ l)
RNase-free water	Variable*	Variable†
Bisulfite Mix	85	85
DNA Protect Buffer	35	15
<b>Total volume</b>	<b>140</b>	<b>140</b>

\* The combined volume of DNA solution and RNase-free water must total 20  $\mu$ l.

† The combined volume of DNA solution and RNase-free water must total 40  $\mu$ l.

**Table 2. Bisulfite conversion thermal cycler conditions**

Step	Time	Temperature
Denaturation	5 min	95°C
Incubation	25 min	60°C
Denaturation	5 min	95°C
Incubation	85 min (1 h 25 min)	60°C
Denaturation	5 min	95°C
Incubation	175 min (2 h 55 min)	60°C
Hold	Indefinite*	20°C

\* Converted DNA can be left in the thermal cycler overnight without any loss of performance.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, Sample to Insight®, EpiTect®, MinElute® (QIAGEN Group). 1101284 03/2016 HB-0645-002 © 2016 QIAGEN, all rights reserved.