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# QuantStudio™ 5 Real-Time PCR System 用 Investigator® Quantiplex® Pro ハンドブック

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The screenshot displays the 'Assign Targets and Samples' window in the software. It is divided into two main sections: 'Targets' and 'Samples', and a 96-well plate grid on the right.

**Targets Table:**

	Name	Reporter	Quencher	Comp	Task	Quantity
<input checked="" type="checkbox"/>	Human	QPP_FAM	None		N	
<input checked="" type="checkbox"/>	Male	QPP_ATT0647N	None		N	
<input checked="" type="checkbox"/>	Degradat...	QPP_ATT0550	None		N	
<input checked="" type="checkbox"/>	IC	QPP_JOE	None		U	

**Samples Table:**

	Sample Name	Comments
<input checked="" type="checkbox"/>	NTC	
<input type="checkbox"/>	Standard -1	
<input type="checkbox"/>	Standard -2	
<input type="checkbox"/>	Standard -3	

The 96-well plate grid shows columns 1-12 and rows A-H. In row A, wells 9 and 10 have a grey 'N' button highlighted with a red box. In row B, wells 9 and 10 also have a grey 'N' button highlighted with a red box. All other wells have a blue 'U' button.

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The screenshot shows the 'Assign Targets and Samples' window in the software. The 'Advanced Setup' tab is selected. On the left, there are two tables: 'Targets' and 'Samples'. The 'Targets' table has columns for Name, Reporter, Quencher, Cor, Task, and Quantity<sup>1</sup>. The 'Task' column is highlighted with a red box, and the 'S' button is selected for 'Degradat...', 'Human', and 'Male', while the 'U' button is selected for 'IC'. The 'Samples' table lists NTC and Standard -1, -2, -3. On the right, a well plate grid is shown with wells A1-A12. A red dashed box highlights wells A1-A8, which contain 'S' buttons, and wells A9-A12, which contain 'U' buttons. The status bar at the bottom indicates 'Wells: 96 8 2' and '0 Empty'.

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The screenshot shows the 'Assign Targets and Samples' window in the QuantStudio 5 software. The interface is divided into several sections:


- Targets Table:** A table with columns: Name, Reporter, Quencher, Com..., Tas..., and Quantity. The 'Quantity' column is highlighted with a red box. The data is as follows:
 

Name	Reporter	Quencher	Com...	Tas...	Quantity
Human	QPP_FAM	None		S	50.0
Male	QPP_ATTO647N	None		S	50.0
Degradat...	QPP_ATTO550	None		S	50.0
IC	QPP_JOE	None		U	
- Samples Table:** A table with columns: Sample Name and Comments. The 'Standard -1' row is highlighted with a red box.
 

Sample Name	Comments
NTC	
Standard -1	
Standard -2	
Standard -3	
- Assay Plate Grid:** A 12x8 grid of wells. A tooltip is visible over well B3, showing target and sample information:
  - Target: Degradation
  - Task: Standard
  - Dye: QPP\_ATTO550-None
  - Quantity: 50.000
  - Target: Human
  - Task: Standard
  - Dye: QPP\_FAM-None
  - Quantity: 50.000
  - Target: IC
  - Task: Unknown
  - Dye: QPP\_JOE-None
  - Quantity: 50.000
  - Target: Male
  - Task: Standard
  - Dye: QPP\_ATTO647N-None
  - Quantity: 50.000
  - Sample: Standard -1

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

**Default Ct Settings**  
 Default Ct settings are used to calculate the Ct for targets without custom settings. To edit the default settings, click **Edit Default Settings**.

Threshold: AUTO Baseline Start Cycle: AUTO Baseline End Cycle: AUTO **Edit Default Settings**

Target	Threshold	Baseline Start	Baseline End
Degradation	0.125	3	15
Human	0.2	3	15
IC	0.05	3	15
Male	0.2	3	15

**Ct Settings for Male**  
 Ct Settings to Use:  Default Settings  
 Automatic Threshold  
 Threshold: 0.2  
 Automatic Baseline  
 Baseline Start Cycle: 3 End Cycle: 15

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Instructions

## Quantification Assay Data Handling Configuration

### Import QIASymphony Rack File Data for Quantification Setup

Define the root directory where the QIASymphony Rack File is stored

Root directory

Defined root directory

To import QIASymphony Rack files

Browse...

The screenshot shows the 'Instructions' section of the 'Configuring root directories, Quantiplex Pro specific settings, Diluent naming & Archive options' window. The window has a breadcrumb trail: Contents | Configuration | Quantification Setup | Quantification Reagent Volumes | Importing Quantification Data | Virtual Quantification Standards | Standard Assessment | PCR Setup | PCR Reagent Volumes | CE Setup on QIAGility. The 'Instructions' section contains the following text:

**Description**  
The Quant Assay Data Handling tool (QDHT) reads, reformats and exports data from defined locations, these locations are specified in the "Configuration" sheet and must be configured prior to using the associated function. The tool also performs analysis of exported data from the 7500, QuantStudio and RQ2 for Quantiplex Pro and supports CE set-up on the QIAGility. Archive options can be set for the reported Quantification data and the PCR Setup sheets.

**Instructions**  
The configuration sheet is presented in 5 main sections as described below. To specify the directory for the import and export functions click the browse button and navigate to the appropriate folder.

- Quantification Data Processing**  
The QDHT will produce plate record files for use with the RQ2 and Life Technologies 7500 (Q25/HD) and QuantStudio. These files contain data relating to standards, samples and locations on the plate. Specify here where these files should be saved.
- Quant Result Import**  
Specify here where in the directory the quantification result files are located.  
Following the import of quantification run results, for Investigator Quantiplex and Investigator Quantiplex HTm the QDHT will calculate PCR Setup volumes. For Investigator Quantiplex Pro the QDHT will perform a quality analysis assessment. Investigator Quantiplex Pro data may then be used for calculating PCR Setup volumes.
- Quantiplex Pro Plate Setup and Analysis Criteria**  
For Quantiplex Pro it is possible to specify the target name for each of the targets in the result file, this option is only required if the targets have been changed from the default values in the Quantiplex Pro handbook. The Quality Assessment thresholds may also be modified, see the Quantiplex Pro handbook for more details.  
Note: Default handbook values may be entered for the Target Names and Threshold Quality Criteria by pressing the "Enter Defaults" button located in each table as required.  
For Quantiplex Pro it is possible to define a Quantification control sample(s) in the run, entering the name of the control here will allow the control to be filtered, if required, from the imported sample data. If filtering is not required delete the entry from this field.  
To configure Virtual Standards, please refer to the "Virtual Quantification Standards" tab at the top of this window.
- Normalization Diluent Naming**  
The diluent used at PCR stage may be defined here, the PCR Setup sheet will then present the name entered into this field.
- Archive Settings**  
The archive format, archive location and a filename prefix (which will be appended with a timestamp) may be defined using the options presented in this section.
- Export Quantification Result to Instrument for PCR Setup**

Save the QDHT to retain configuration settings for future use, to do this, select "File/Save".

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## QIAGEN Quantification Assay Data Handling and STR Setup Tool

Release Date: 16.12.2022 Version: 4.0.1

**This tool will generate CSV (comma-separated value) text files for use on Real-Time PCR Instruments  
For compatibility the regional settings of the PC running this tool must be set to US or UK English**

This tool enables the creation of sample records for use with the following assays as follows:

Investigator Quantiplex® – 7500 SDS, 7500 HID, and RGQ  
Investigator Quantiplex HYres – 7500 SDS, 7500 HID, and RGQ  
Investigator Quantiplex Pro – 7500 SDS, 7500 HID, and QuantStudio™  
Investigator Quantiplex Pro RGQ – Rotor-Gene® Q Instruments

Result data may be exported from the Rotor-Gene Q, AB 7500, QuantStudio, and Bio-Rad® CFX instruments and imported using this tool.

The import process will format the data removing standards and NTC data.

In the case of AB 7500 Quantiplex HYres, the tool will prompt which result to import (human or male).

From the sample quantification data, information is provided for the setup of QIAGEN STR reactions including normalization if required.

Quantiplex Pro users on ABI 7500 & QuantStudio cyclers – Click to perform a One Time update

**Note: This document requires Excel macros to be enabled in order to function.**

**This tool has been tested with Excel 2010 and Excel 2016.**

**Excel 2007 has not been tested but may be compatible; versions prior to 2007 are not compatible.**

Instructions for each function can be found by pressing the "Instructions" button located on each page.

**QIAGEN shall not be held liable in connection with the use of this Excel tool by any user.**

... Opening Page Configuration Quantification Setup Quant Component Volumes Quant Result Import Quant Standard Plots PCR Setup

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Instructions

### Quantification Assay Data Handling Configuration

#### Import QIASymphony Rack File Data for Quantification Setup

Define the root directory where the QIASymphony Rack File is stored  
Root directory  Defined root directory  
To import QIASymphony Rack files

#### Quantification Data Processing

Define the root directory where the quantification RT-PCR plate record will be saved  
Root directory  Defined root directory  
To save Quant batch files

#### Quant Result Import

Define the root directory where concentration data result file from the RGQ/AB7500/QuantStudio will be saved  
Root directory  Defined root directory  
To import Quant result files   
Use Virtual Standard for Data Analysis?   Yes  No

#### Quantiplex Pro Plate Setup and Analysis Criteria Options

Result File Target Name Specification (Quantiplex Pro only)

Target Description	Name in Results File	<input type="button" value="Enter Defaults"/>
Human Target	Human	
Male Target	Male	
Human Degradation Target	Degradation	
Internal Positive Control (IPC)	IC	
Male Degradation Target (RGQ Only)	Male Degradation	

#### Threshold Specification for Quality Assessment in Quantiplex Pro only

Quality Assessment	Threshold	<input type="button" value="Enter Defaults"/>
Mixture Index (Human/Male)	2	
Human Degradation Index (Human/Human Degradation)	10	
Inhibition Index (IC Shift)	1	
Male Degradation Index (Male/Male Degradation)	10	

#### Quantification NTC & QC Control Specification, to be removed during data import (Quantiplex Pro only)

Control Name in Result File

Opening Page **Configuration** Quantification Setup Quant Component Volumes Quant Result Import Quant Stand

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**Instructions** Quantification Assay Data Handling Configuration

**Import QIASymphony Rack File Data for Quantification Setup**  
 Define the root directory where the QIASymphony Rack File is stored  
 Root directory Defined root directory  
 To import QIASymphony Rack files **Browse...**

**Quantification Data Processing**  
 Define the root directory where the quantification RT-PCR plate record will be saved  
 Root directory Defined root directory  
 To save Quant batch files **Browse...**

**Quant Result Import**  
 Define the root directory where concentration data result file from the RGQ/AB7500/QuantStudio will be saved  
 Root directory Defined root directory  
 To import Quant result files **Browse**  
 Use Virtual Standard for Data Analysis? [What's This?](#)  Yes  No **Edit Standard Data** **Delete Saved Standard**

**Quantiplex Pro Plate Setup and Analysis Criteria Options**  
 Result File Target Name Specification (Quantiplex Pro only)

Target Description	Name in Results File	<b>Enter Defaults</b>
Human Target	Human	
Male Target	Male	
Human Degradation Target	Degradation	
Internal Positive Control (IPC)	IC	
Male Degradation Target (RGQ Only)	Male Degradation	

**Threshold Specification for Quality Assessment in Quantiplex Pro only**

Quality Assessment	Threshold	<b>Enter Defaults</b>
Mixture Index (Human/Male)	2	
Human Degradation Index (Human/Human Degradation)	10	
Inhibition Index (IC/Strain)	1	
Male Degradation Index (Male/Male Degradation)	10	

**Quantification NTC & QC Control Specification, to be removed during data import (Quantiplex Pro only)**  
 Control Name in Result File 3948

Opening Page **Configuration** Quantification Setup Quant Component Volumes Quant Result Import Quant Stand

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Instructions Quantification Assay Data Handling Configuration

Import QIASymphony Rack File Data for Quantification Setup

Define the root directory where the QIASymphony Rack File is stored

Root directory

To import QIASymphony Rack files

Quantification Data Processing

Define the root directory where the quantification RT-PCR plate record will be saved

Root directory

To save Quant batch files

Quant Result Import

Define the root directory where concentration data result file from the RIGQIAB7500QuantStudio will be saved

Root directory

To import Quant result files

Use Virtual Standard for Data Analysis? [What's This?](#)  Yes  No

Quantplex Pro Plate Setup and Analysis Criteria Options

Result File Target Name Specification (Quantplex Pro only)

Target Description	Name in Results File	<input type="button" value="Enter Defaults"/>
Human Target	Human	
Male Target	Male	
Human Degradation Target	Degradation	
Internal Positive Control (IPC)	IC	
Male Degradation Target (RIGQ Only)	Male Degradation	

Threshold Specification for Quality Assessment in Quantplex Pro only

Quality Assessment	Threshold	<input type="button" value="Enter Defaults"/>
Mixture Index (Human/Male)	2	
Human Degradation Index (Human/Human Degradation)	10	
Inhibition Index (IC Shift)	1	
Male Degradation Index (Male/Male Degradation)	10	

Quantification NTC & QC Control Specification, to be removed during data import (Quantplex Pro only)

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Instructions

## Quantification Assay Data Handling Configuration

### Import QIASymphony Rack File Data for Quantification Setup

Define the root directory where the QIASymphony Rack File is stored

Root directory

Defined root directory

To import QIASymphony Rack files

Browse...

### Quantification Data Processing

Define the root directory where the quantification RT-PCR plate record will be saved

Root directory

Defined root directory

To save Quant batch files

Browse...

### Quant Result Import

Define the root directory where concentration data result file from the PGQA87500/QuantStudio will be saved

Root directory

Defined root directory

To import Quant result files

Browse...

Use Virtual Standard for Data Analysis?

What's This?

Yes

No

Edit Standard Data

Delete Saved Standard

### Quantiplex Pro Plate Setup and Analysis Criteria Options

Result File Target Name Specification (Quantiplex Pro only)

#### Target Description

#### Name in Results File

Enter Defaults

Human Target

Human

Male Target

Male

Human Degradation Target

Degradation

Internal Positive Control (IPC)

IC

Male Degradation Target (PGQ Only)

Male Degradation

#### Threshold Specification for Quality Assessment in Quantiplex Pro only

#### Quality Assessment

#### Threshold

Enter Defaults

Mixture Index (Human/Male)

2

Human Degradation Index (Human/Human Degradation)

10

Inhibition Index (IC Shift)

1

Male Degradation Index (Male/Male Degradation)

10

#### Quantification NTC & QC Control Specification, to be removed during data import (Quantiplex Pro only)

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
Result Summary	Human			Human Degradation			IPC			Male			Quality Assessment			
	WT	Quantity	Quantity Mean	CT	Quantity	Quantity Mean	CT	Quantity	Quantity Mean	Mixture Index	Mixture Threshold	Degradation Index	Degradation Threshold	Inhibition Index	Inhibition Threshold	
B1 HE 10xJ1	30.132	0.359	0.361	29.664	0.106	0.116	23.422	27.025	0.320	0.343	0.62	Below Threshold	0.96	Below Threshold	-1.68	Possible Inhibition
C1 HE 10xJ1	29.516	0.334	0.346	27.357	0.2970	0.2934	22.318	26.509	0.9953	0.940	1.76	Below Threshold	1.16	Below Threshold	-0.66	Below Threshold
D1 Mixture 3	21.197	51.182	51.267	20.247	51.6025	52.879	22.599	33.145	0.0014	0.0012	37666.79	Possible Mixture	0.99	Below Threshold	-0.66	Below Threshold
E2 HE 10xJ1	30.124	0.364	0.361	28.941	0.105	0.116	23.106	27.125	0.197	0.243	0.89	Below Threshold	0.96	Below Threshold	-1.17	Possible Inhibition
C2 HE 10xJ1	28.577	0.378	0.346	27.282	0.279	0.294	22.231	26.507	0.198	0.240	1.78	Below Threshold	1.14	Below Threshold	-0.29	Below Threshold
D2 Mixture 3	21.195	51.5732	51.267	20.772	54.2694	52.879	22.278	33.405	0.0011	0.0012	40256.34	Possible Mixture	0.96	Below Threshold	-0.34	Below Threshold
C3 HE 20xJ1	30.000	0.191	0.195	29.809	0.0480	0.058	24.634	28.804	0.920	0.930	0.76	Below Threshold	2.42	Below Threshold	-2.70	Possible Inhibition
D3 Male DNA 500bp	28.673	0.2971	0.309	28.731	0.1053	0.107	22.422	28.662	0.964	0.930	1.78	Below Threshold	2.42	Below Threshold	-0.49	Below Threshold
D3 Male 100ngJ1	19.362	109.246	109.279	18.789	149.099	150.949	21.940	17.751	120.242	120.275	1.09	Below Threshold	0.99	Below Threshold	0.93	Below Threshold
D4 HE 10xJ1	29.921	0.329	0.196	29.298	0.0971	0.096	23.914	28.662	0.900	0.910	0.82	Below Threshold	1.76	Below Threshold	-1.88	Possible Inhibition
C4 Male DNA 500bp	28.628	0.307	0.309	28.621	0.191	0.197	22.260	28.642	0.175	0.190	1.79	Below Threshold	2.69	Below Threshold	-0.32	Below Threshold
D4 Male 100ngJ1	18.986	138.8052	138.979	18.753	152.8842	150.949	21.547	17.709	132.389	130.275	1.05	Below Threshold	0.91	Below Threshold	0.39	Below Threshold
E5 HS 25ngJ1	29.378	0.303	0.321	28.720	0.1062	0.103	23.135	27.257	0.3095	0.132	1.66	Below Threshold	1.70	Below Threshold	-1.20	Possible Inhibition
E5 Male DNA 300bp	28.395	0.2902	0.290	28.930	0.1429	0.1486	21.989	28.748	0.896	0.974	1.57	Below Threshold	5.70	Below Threshold	-0.05	Below Threshold
D6 Female 100ngJ1	29.943	87.505	88.251	19.297	102.8036	104.7978	22.060	Undetermined	0.0000	0.0000	0.00	Below Threshold	0.66	Below Threshold	-0.12	Below Threshold
E6 HS 25ngJ1	29.350	0.339	0.321	29.619	0.1143	0.103	23.001	27.146	0.1978	0.132	1.58	Below Threshold	1.61	Below Threshold	-1.08	Possible Inhibition
C6 Male DNA 300bp	28.862	0.258	0.2550	29.769	0.0434	0.0466	22.003	28.778	0.953	0.974	1.67	Below Threshold	5.26	Below Threshold	-0.06	Below Threshold
D6 Female 100ngJ1	29.568	91.938	89.251	19.245	109.7120	104.7978	21.902	Undetermined	0.0000	0.0000	0.00	Below Threshold	0.86	Below Threshold	0.04	Below Threshold
D7 HS 10ngJ1	29.428	0.140	0.154	29.395	0.0663	0.0635	24.105	27.460	0.164	0.151	1.47	Below Threshold	2.82	Below Threshold	-2.17	Possible Inhibition
D7 Male DNA 100bp	29.401	0.175	0.172	28.154	0.0095	0.0095	21.951	27.385	0.940	0.925	1.71	Below Threshold	373.63	Below Threshold	-0.01	Below Threshold
D7 Mixture 1400.000	19.485	97.8239	200.0430	19.404	245.4578	245.3752	21.741	16.442	0.0002	0.0003	690289.60	Possible Mixture	0.81	Below Threshold	0.20	Below Threshold
D7 HS 33ngJ1	29.406	0.178	0.154	29.489	0.0606	0.0635	24.208	27.266	0.978	0.131	1.64	Below Threshold	2.82	Below Threshold	-2.35	Possible Inhibition

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04-2023 HB-3355-001 © 2023 QIAGEN、無断複写・転載を禁じます。