# RT<sup>2</sup> Profiler PCR Array

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Housekeeping genes      Genomic DNA control      Reverse transcription controls      Positive PCR controls

图 1. RT<sup>2</sup> Profiler PCR Array 中的每孔测定了与通路或疾病状态相关的基因的表达。A typical 96-well format is shown. This is also available in a 384-well plate, 100-well disc, and 96 × 96 chip format.


图 2. RT<sup>2</sup> SYBR Green vs. TaqMan 可比的生物学结果。\* Gene expression analysis was compared between RT<sup>2</sup> Profiler PCR Arrays (SYBR<sup>®</sup> Green-based) and the TaqMan<sup>®</sup> platform. Regression analysis of fold differences, with data normalized against POLR2A, demonstrate that both platforms yield similar biological results.


图 3. RT<sup>2</sup> SYBR Green vs. TaqMan 试剂的灵敏度。\* PCR amplicons detected using the same primer pair with or without TaqMan probes in either SYBR Green or TaqMan chemistry. SYBR green chemistry yields earlier C<sub>s</sub> for each dilution, demonstrating better sensitivity than TaqMan chemistry.


图 4. Stress and Toxicity PathwayFinder RT<sup>2</sup> Profiler PCR Array 揭示了不同的基因表达谱。RNA from HepG2 cells treated with three different glitazone PPAR $\gamma$  agonists for type 2 diabetes mellitus was characterized, and the results were compared to that of a vehicle (DMSO) control. The drug Rezulin induces very different changes in the expression of stress-related genes than Avandia and Actos.



图 5. Common Cytokine RT<sup>2</sup> Profiler PCR Array 鉴定出外周血单核细胞经刺激后的上调和下调基因。Triplicate total RNA samples from human peripheral blood mononuclear cells (either untreated or stimulated with 50 ng/ml PMA and 1 mg/ml ionomycin for 6 hours) were characterized with the human Common Cytokine RT<sup>2</sup> Profiler PCR Array. Twenty-three cytokine genes are up-regulated (> 5-fold, p < 0.0005), while six interleukin and TNF ligand genes are down-regulated.

\* Arikawa, E., et al. (2008) Cross-platform comparison of SYBR Green real-time PCR with TaqMan PCR, microarrays and other gene expression measurement technologies evaluated in the MicroArray Quality Control (MAQC) study. BMC Genomics 9, 378.

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
**图 6. A. 高效去除基因组 DNA 污染。** RNA from HEK 293T cells, mouse spinal tissue, mouse brain tissue, or rat brain tissue was characterized on RT<sup>2</sup> Profiler PCR Arrays before (light blue bars) and after (dark blue bars) treatment with gDNA Elimination Buffer from the RT<sup>2</sup> First Strand Kit. **B. 监控逆转录中的抑制。** Human universal RNA was added with magnesium salt to simulate RNA degradation or added with TRIzol reagent to simulate contamination that inhibits enzyme activity. RT<sup>2</sup> First Strand Kit was used for cDNA synthesis.

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

图 7. MECP2 小鼠模型中下丘脑的基因表达变化。Validation of expression changes for 66 genes by qPCR analysis. Gene expression levels from microarray analyses were validated in four MECP2-Tg males and four MeCP2-null males.\*

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
图 9. RT<sup>2</sup> Profiler PCR Array 可检测低至 1 ng RNA。Different amounts of universal total RNA were characterized using the Human Inflammatory Cytokines and Receptors RT<sup>2</sup> Profiler PCR Array (PAHS-011) with or without PreAMP.

\* Chahrour, M., et al. (2008) MeCP2, A key contributor to neurological disease, activates and represses transcription. Science 320, 1224.

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
图 10. RT<sup>2</sup> Profiler PCR Array 可检测宽动态范围内的 RNA。Ten-fold serial dilutions of human CHRNA5 were characterized with the respective RT<sup>2</sup> qPCR Primer Assay.

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
图 11. 火山图显示了基因表达变化的统计显著性。The x-axis plots the log<sub>2</sub> of the fold-differences, while the y-axis plots their p-values based on student's t-test of your replicate raw C<sub>T</sub> data. The blue and red symbols outside the gray area conveniently have the same meaning as the scatter plot. Symbols in the volcano plot above the dashed line readily identify fold-differences at least as statistically significant as a threshold that can be defined.

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## RT<sup>2</sup> lncRNA PCR Array product list

Pathway VC*	Name
LAHS-001Z	Human lncRNA Finder
LAHS-002Z	Human lncRNA Cancer PathwayFinder
LAHS-003Z	Human lncRNA Cell Differentiation & Development
LAHS-004Z	Human lncRNA Inflammatory Response & Autoimmunity

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Pathway VC*	Name
LAMM-001Z	Mouse lncRNA Finder
LAMM-002Z	Mouse lncRNA Cancer PathwayFinder
LAMM-003Z	Mouse lncRNA Cell Differentiation & Development
LAMM-004Z	Mouse lncRNA Inflammatory Response & Autoimmunity

\* **Array Format/Plates per order:** 96 well -2, 96 well -6, 96 well -12, 96 well -24, 384 4x96 (16), Fluidigm (5).

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产品	内容	货号
RT <sup>2</sup> lncRNA PCR Arrays	Pathway, disease, or custom panels of gene assays	330721
RT <sup>2</sup> lncRNA qPCR Assay	Individual laboratory-verified qPCR assay for each lncRNA	337001
RT <sup>2</sup> SYBR Green qPCR Mastermixes	Reagents for real-time PCR reactions (available with ROX, fluorescein, or no internal reference dye)	varies
RT <sup>2</sup> First Strand Kit (12) or (50)	Reagents for cDNA synthesis reactions for 12 or 50 samples	330401/330404

可兼容用于下列制造商生产的 real-time PCR 仪器: QIAGEN、Agilent/Stratagene、Applied Biosystems®、Bio-Rad®、Eppendorf®、Fluidigm、Roche®、TaKaRa。

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产品	规格	货号
RNeasy Mini Kit (50)*	For purification of RNA from cells & easy-to-lyse tissues	74104
RNase-Free DNase Set (50)	For DNase digestion RNA purification	79254
RNeasy Fibrous Tissue Mini Kit (50)*	For purification of RNA from fiber-rich tissues	74704
RNeasy Plus Universal Mini Kit (50)*	For purification of RNA from all tissue types	73404
RNeasy 96 Kit (12)	For purification of RNA from cells in 96-well format	74182
QIAzol® Lysis Reagent (200 ml)	For lysis of fatty and standard tissues before RNA isolation	79306

\* Automatable on the QIAcube®. Find out more at [www.qiagen.com/goto/QIAcube](http://www.qiagen.com/goto/QIAcube).

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
产品	规格	货号
Allprotect Tissue Reagent (100 ml)	For stabilization of RNA, DNA, & protein in tissues	76405
AllPrep <sup>®</sup> DNA/RNA Mini Kit (50)	For simultaneous purification of DNA and RNA	80204
RNA <sub>later</sub> RNA Stabilization Reagent (250 ml)	For stabilization of RNA in tissues	76106
RNeasy Protect Mini Kit (50)	For stabilization of RNA in tissues and RNA purification	74124
RNAprotect <sup>®</sup> Cell Reagent (250 ml)	For stabilization of RNA in cells	76526
RNAprotect Animal Blood Tubes (50 x 100 µl)	For collection of 100 µl animal blood with RNA stabilization	76544
RNeasy Protect Animal Blood Kit (50)	For purification of RNA from blood collected in RNAprotect Animal Blood Tubes	73224
miRNeasy Protect Animal Blood Kit (50)	For purification of RNA, including miRNA, from blood collected in RNAprotect Animal Blood Tubes	217304

**图 15. 高效的 RNA 稳定。** Rat tissues were stored at 25°C for 2–24 hours in Allprotect Reagent or PBS prior to real-time RT-PCR analysis. **A.** Rat lung tissue was analyzed for c-fos expression. **B.** Rat intestine tissue was analyzed for Madh7 expression. Transcript levels relative to those in liquid nitrogen stabilized tissues were calculated. Changes in transcript levels were prevented by Allprotect Reagent.

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"RNAlater<sup>®</sup>" is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, QIAxcel<sup>®</sup>, QIAzol<sup>®</sup>, AllPrep<sup>®</sup>, Allprotect<sup>®</sup>, RNAprotect<sup>®</sup>, RNeasy<sup>®</sup> (QIAGEN Group); Applied Biosystems<sup>®</sup>, SYBR<sup>®</sup> (Life Technologies Corporation); Bio-Rad<sup>®</sup> (Bio-Rad Laboratories, Inc.); Eppendorf<sup>®</sup> (Eppendorf AG); Excel<sup>®</sup> (Microsoft Corporation); Roche<sup>®</sup>, TaqMan<sup>®</sup> (Roche Group).

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

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© 2016 QIAGEN, all rights reserved.