

ForenSeq DNA Signature Prep Kit with ForenSeq Enhanced PCR1 Buffer System

Reference Guide

VEROGEN PROPRIETARY

Document # VD2021046 Rev. A

January 2022

Revision History

| Document # | Date | Description of Change |
|------------------|--------------|-----------------------|
| VD2021046 Rev. A | January 2022 | Initial release |

Table of Contents

| | |
|---|-----------|
| Revision History | 2 |
| Table of Contents | 3 |
| Chapter 1. Overview | 5 |
| Introduction | 5 |
| Protocol Steps | 6 |
| DNA Input Recommendations | 7 |
| Acronyms | 7 |
| Additional Resources | 8 |
| Chapter 2. Protocol | 9 |
| Introduction | 9 |
| Tips and Techniques | 10 |
| Amplify and Tag Targets | 11 |
| Enrich Targets | 13 |
| Purify Libraries | 15 |
| Normalize Libraries | 17 |
| Pool Libraries | 19 |
| Denature and Dilute Libraries | 20 |
| Chapter 3. Materials | 22 |
| Kit Contents and Storage | 22 |
| Index Adapter Sequences | 25 |
| Consumables and Equipment | 26 |
| Appendix A. Amplicon Information | 29 |
| Loci Detected with DPMA and DPMB | 29 |
| Interpreting Loci D22S1045 | 37 |
| Interpreting Loci DYS392 | 39 |
| Technical Support | 41 |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

This document and its contents are proprietary to Verogen, Inc. and its affiliates, and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Verogen. Verogen does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document. Patents: US Patent No. US10422002B2, Chinese Patent No. CN106164298B, and pending patents.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

VEROGEN DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2021 Verogen, Inc. All rights reserved.

All other names, logos, and other trademarks are the property of their respective owners.

For Research, Forensic, or Paternity Use Only– not for any clinical or therapeutic use in humans or animals.

This product includes 2800M Control DNA manufactured by Promega Corporation for distribution by Verogen, Inc.



Notice to Purchaser: Limited License

The purchase price of the product includes a limited, nontransferable license under US and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this Product.

The Product is sold under license from Affibody AB, Sweden.

Phusion DNA Polymerase is manufactured by Thermo Fisher Scientific. Phusion is a trademark or registered trademark of Thermo Fisher Scientific, or its subsidiaries.


Overview

Introduction 5
 Protocol Steps 6
 DNA Input Recommendations 7
 Acronyms 7
 Additional Resources 8

Introduction

The ForenSeq® DNA Signature Prep Kit generates 384 dual-indexed libraries for sequencing. Each sample is combined with a primer mix that contains a pair of tagged oligos for each target sequence. PCR cycles link the tags to copies of each target, forming DNA templates consisting of the regions of interest flanked by universal primer binding sequences. The tags are then used to attach index adapters and the resulting library is amplified, purified, and pooled for sequencing. The ForenSeq DNA Signature Prep Kit can be used with the standard PCR1 buffer that is available with the kit or an enhanced buffer system (ePCR1) that is available as an add-on kit. This Reference Guide summarizes the protocol for using the ForenSeq DNA Signature Prep Kit with the ForenSeq Enhanced PCR1 Buffer System.

Figure 1 Assay overview



Kit Features

The ForenSeq DNA Signature Prep Kit offers the following features:

- Simultaneous preparation of up to 96 libraries. Each library is a collection of amplified DNA fragments from one sample.
- Amplify short tandem repeat (STR) and single-nucleotide polymorphism (SNP) amplicons from challenging samples in one reaction.
- Choose between two targeted primer mixes, DNA Primer Mix A (DPMA) or DNA Primer Mix B (DPMB), to prepare samples for databasing or casework.

ForenSeq Enhanced PCR1 Buffer System offers the following features:

- An optimized buffer system that support the simultaneously amplification of 96 challenging samples in the presence of high concentrations of inhibitors, such as humic acid.

Protocol Steps

The following diagram lists the steps to prepare libraries with challenging examples, using the enhanced PCR1 buffer, along with hands-on times, total times, and reagents. Safe stopping points are marked between steps.

Figure 2 Overview of the ForenSeq DNA Signature Prep protocol with the ForenSeq Enhanced PCR1 Reaction Mix



DNA Input Recommendations

Use 1 ng purified human genomic DNA (gDNA) as input. Before starting the protocol, quantify the input using a fluorometric-based method or qPCR and assess quality.

Controls

Each preparation must include at least one positive amplification control and at least one negative amplification control. If these controls are not included, troubleshooting support is limited.

The kit includes Control DNA 2800M (2800M) for use as the positive template control and the negative amplification control is nuclease-free water. The protocol includes instructions to prepare each control.

Acronyms

| Acronym | Definition |
|---------|--|
| DNL | Diluted Normalized Libraries |
| DPMA | DNA Primer Mix A |
| DPMB | DNA Primer Mix B |
| ePCR1 | Enhanced PCR1 Reaction Mix |
| gDNA | Genomic DNA |
| FEM | Enzyme Mix |
| FSP | ForenSeq Sample Plate |
| HP3 | 2 N NaOH |
| HSC | Human Sequencing Control |
| HT1 | Hybridization Buffer |
| LNA1 | Library Normalization Additives 1 |
| LNB1 | Library Normalization Beads 1 |
| LNS2 | Library Normalization Storage Buffer 2 |
| LNW1 | Library Normalization Wash 1 |
| 2800M | Control DNA 2800M |
| NLP | Normalized Library Plate |
| NWP | Normalization Working Plate |
| PBP | Purification Bead Plate |
| PCR2 | PCR2 Reaction Mix |
| PLP | Purified Library Plate |
| PNL | Pooled Normalized Libraries |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Acronym | Definition |
|---------|---------------------------|
| ProK | Proteinase K |
| RSB | Resuspension Buffer |
| SPB | Sample Purification Beads |

Additional Resources

This guide provides comprehensive information on the use of ForenSeq DNA Signature Prep Kit with detailed protocol instructions for use of the ForenSeq Enhanced PCR1 Buffer System. Visit the [Documentation page](#) on the Verogen website to download additional kit documentation and access the latest versions.

| Resource | Description |
|---|---|
| <i>ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Checklist (document # VD2021047)</i> | Provides concise protocol instructions for the experienced user. |
| <i>ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Materials List (document # VD2021048)</i> | Lists the consumables and equipment needed to perform the protocol. |

Protocol

| | |
|-------------------------------------|----|
| Introduction | 9 |
| Tips and Techniques | 10 |
| Amplify and Tag Targets | 11 |
| Enrich Targets | 13 |
| Purify Libraries | 15 |
| Normalize Libraries | 17 |
| Pool Libraries | 19 |
| Denature and Dilute Libraries | 20 |

Introduction

This chapter describes the ForenSeq DNA Signature Prep protocol with step-by-step instructions to prepare libraries for sequencing using the ForenSeq Enhanced PCR1 buffer system. For an overview of the protocol with reagents and durations for each step, see [Protocol Steps \(page 6\)](#).

Before starting, confirm kit contents and make sure that you have the necessary reagents, consumables, and equipment. For a list of items, see [Materials \(page 22\)](#).

Number of Samples

Process at least eight samples at a time, including positive and negative amplification controls. Preparing master mixes for fewer than eight samples can introduce pipetting inaccuracies due to small volumes.

Use the following table to determine the maximum number of libraries to pool for a run, depending on primer mix, sample type, and MiSeq FGx[®] reagent kit. Casework recommendations are intended for samples where DNA mixtures are possible or challenging samples with < 1 ng gDNA available and partial degradation.

Table 1 Maximum number of libraries

| Primer Mix | Sample Type | MiSeq FGx Reagent Micro Kit | MiSeq FGx Reagent Kit |
|------------|-----------------------|-----------------------------|-----------------------|
| DPMA | Database or reference | 36 | 96 |
| | Casework | 12 | 32 |
| DPMB | Database or reference | 12 | 32 |
| | Casework | 12 | 32 |

Primer Mixes

The kit includes two primer mixes: DPMA and DPMB. Both primer mixes detect identity informative SNPs (iiSNPs), autosomal STRs (aSTRs), and X- and Y-STRs. In addition to these targets, DPMB detects ancestry-informative SNPs (aiSNPs) and phenotypic-informative SNP (piSNPs).

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

The ForenSeq DNA Signature Prep Kit supports analysis of these SNPs and STRs from gDNA ranging from high-quality, single-source samples to difficult. One reaction with integrated indexing enables sequencing of up to 96 database samples using DPMA or 32 casework samples using DPMB in one run with the MiSeq FGx Reagent Kit.

Tips and Techniques

Protocol Continuity

- Follow the steps in the order indicated using the specified volumes and incubation parameters.
- Unless a safe stopping point is specified, proceed immediately to the next step.

Plate Setup

- Create a sample sheet to record the position of each sample, control, and index adapter.
- Reference the sample sheet throughout the protocol to ensure proper plate setup.

The *ForenSeq Universal Analysis Software v1.3 Reference Guide (document # VD2018007)* provides detailed information on sample sheets and input of sample information.

Preventing Cross-Contamination

- Set up the [Amplify and Tag Targets \(facing page\)](#) process in a pre-PCR environment. Perform all other processes in a post-PCR environment.
- When adding or transferring samples, change tips between **each sample**.
- When adding adapters or primers, change tips between **each well**.
- Do not reuse index adapter caps.
- Remove unused index adapter tubes from the working area.

Sealing the Plate

- Apply a microseal to cover the plate and seal with a rubber roller. After each use, discard seals from plates.
- Use Microseal 'A' pressure film for thermal cycling. When using fewer than 96 wells, you can cut the film to size.
- Use Microseal 'B' adhesive film for shaking, centrifuging, and long-term storage. These seals are effective at -40°C to 110°C.

Handling Beads

- For optimal performance and yield, confirm that beads are at room temperature before use.
- Aspirate and dispense beads slowly due to viscosity.
- Do not centrifuge plates and tubes containing beads, except when indicated.
- Vortex beads before use and frequently throughout the protocol to resuspend. Resuspended beads are evenly distributed and homogenous in color.
- If beads aspirate into pipette tips during supernatant removal, dispense back to the plate on the magnetic stand and wait until the liquid is clear (~2 minutes).

Amplify and Tag Targets

This process uses an oligonucleotide primer mix with regions specific to the DNA sequences upstream and downstream of STRs and SNPs to tag and amplify the input gDNA.

Consumables

- 2800M (Control DNA 2800M)
- DPMA (DNA Primer Mix A) or DPMB (DNA Primer Mix B)
- FEM (Enzyme Mix)
- ePCR1 (Enhanced PCR1 Reaction Mix)
- 1.7 ml microcentrifuge tubes (2)
- 96-well PCR plate, semiskirted
- Input gDNA
- Microseal 'A' film
- Nuclease-free water
- [Optional] RNase/DNase-free 8-tube strip and caps

About Reagents

- Use PCR1 for standard samples and ePCR1 for inhibited samples.
- ePCR1 should not be used with crude lysates or FTA card punches.
- For information on DPMA and DPMB, see [Loci Detected with DPMA and DPMB \(page 29\)](#).

Preparation

1. Prepare the following consumables:

| Item | Storage | Instructions |
|--------------|----------------|--|
| 2800M | 2°C to 8°C | Let stand for 30 minutes to bring to room temperature. Invert three times to mix, and then centrifuge briefly. |
| DPMA or DPMB | -25°C to -15°C | Thaw at room temperature. Invert three times to mix, and then centrifuge briefly. |
| FEM | -25°C to -15°C | Thaw at room temperature, and then centrifuge briefly. Return to storage immediately after use. |
| ePCR1 | -25°C to -15°C | Thaw at room temperature for 30 minutes. Vortex to mix, and then centrifuge briefly. |

2. Save the following PCR1 program on the thermal cycler in the post-amplification area. See [Table 2](#) for ramp modes.
 - Choose the preheat lid option. See [Table 2](#) for lid temperatures.
 - 98°C for 3 minutes
 - 8 cycles of:
 - 96°C for 45 seconds

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

- 80°C for 30 seconds
- 54°C for 2 minutes, with applicable ramp mode
- 68°C for 2 minutes, with applicable ramp mode
- 10 cycles of:
 - 96°C for 30 seconds
 - 68°C for 3 minutes, with applicable ramp mode
- 68°C for 10 minutes
- Hold at 10°C

Table 2 Thermal Cycler lid temperature and ramp modes

| Thermal Cycler | Temperature Mode | Lid Temperature | Ramp Mode | Vessel Type |
|-------------------------------|----------------------------|---------------------------|------------------|--------------------------------|
| ABI LTI thermal cycler 9700 | 9600 emulation | Heated | 8% | Polypropylene plates and tubes |
| Bio-Rad | Calculated | Heated, constant at 100°C | 0.2°C per second | Polypropylene plates and tubes |
| Eppendorf Mastercycler Pro S | Gradient S, Simulated Tube | Heated | 2% | Plate |
| Proflex 96-well PCR System | Not applicable | Heated, constant at 105°C | 0.2°C per second | Polypropylene plates and tubes |
| Veriti 96-well thermal cycler | Standard | Heated, constant at 105°C | 4% | Polypropylene plates and tubes |

The PCR1 program takes ~3.5 hours and can be run overnight.

3. Label a new PCR plate FSP for ForenSeq Sample Plate.
4. Label a new 1.7 ml tube per your input type:

| Input Type | Label |
|--------------|------------|
| Purified DNA | Master Mix |

Procedure

Purified DNA

1. Using nuclease-free water, dilute 1 ng purified DNA input to 0.2 ng/μl.
2. In the Master Mix tube, combine the following volumes. Multiply each volume by the number of samples and add 10% for overage.
 - ePCR1 (4.7 μl)
 - FEM (0.3 μl)
 - DPMA or DPMB (5 μl)

For example, for eight samples prepare 88 μl master mix: 41.4 μl ePCR1, 2.6 μl FEM, and 44 μl DPMA or DPMB.
3. Pipette to mix, and then cap and centrifuge briefly.
4. [Optional] Evenly distribute the master mix among each well of an 8-tube strip. Use a multichannel pipette to dispense.

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

5. Add 10 μ l master mix to each well of the FSP.
6. In a new 1.7 ml tube, combine the following volumes to dilute 2800M:
 - 2800M (2 μ l)
 - Nuclease-free water (98 μ l)
7. Cap and gently invert three times to mix, and then centrifuge briefly.
8. Add 5 μ l diluted 2800M to at least one well of the FSP as a positive amplification control.
9. Pipette to mix.
10. Add 5 μ l nuclease-free water to at least one well of the FSP as a negative amplification control.
11. Pipette to mix.
12. Add 5 μ l 0.2 ng/ μ l DNA to each well of the FSP.
13. Seal and centrifuge at 1000 \times g for 30 seconds.
14. Transport to the post-PCR area.
15. Place on the preprogrammed thermal cycler and run the PCR1 program.
16. Unless you are stopping, proceed to [Enrich Targets \(below\)](#).

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

Enrich Targets

This process amplifies the DNA and adds Index 1 (i7) Adapters, Index 2 (i5) Adapters, and the sequences needed for cluster generation. The index adapters tag the DNA with a unique combination of sequences that identify each sample for analysis.

When preparing eight samples, you can perform this process using an 8-tube strip instead of the 96-well PCR plate.

Consumables

- Index 1 Adapters (R7XX) and orange caps
- Index 2 Adapters (A50X) and white caps
- PCR2 (PCR2 Reaction Mix)
- ForenSeq Index Plate Fixture
- Microseal 'A' film
- [Optional] 1.7 ml microcentrifuge tubes (1 per index adapter tube)

About Reagents

- Dispense PCR2 slowly to prevent bubbles.
- Centrifuge index adapter tubes in the 1.7 ml tubes, if necessary.

Preparation

1. Prepare the following consumables:

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Item | Storage | Instructions |
|----------------|----------------|--|
| Index adapters | -25°C to -15°C | Remove only the index adapters being used. Thaw at room temperature for 20 minutes. Vortex each tube to mix, and then centrifuge briefly. |
| PCR2 | -25°C to -15°C | Thaw at room temperature for 20 minutes, and then invert to mix. |

2. Save the following PCR2 program on the thermal cycler:

- Choose the preheat lid option and set to 100°C
- 98°C for 30 seconds
- 15 cycles of:
 - 98°C for 20 seconds
 - 66°C for 30 seconds
 - 68°C for 90 seconds
- 68°C for 10 minutes
- Hold at 10°C


Total program time is ~46 minutes.

Procedure

1. Centrifuge the sealed FSP at 1000 × g for 30 seconds.
2. Place the Index 1 Adapter tubes (orange caps) in each column of the ForenSeq Index Plate Fixture.
3. Place the Index 2 Adapter tubes (white caps) in each row of the ForenSeq Index Plate Fixture.
4. Place the FSP on the ForenSeq Index Plate Fixture.

Figure 3 shows the fixture setup with index adapter tubes and the FSP.

Figure 3 Fixture setup for 96 libraries



5. Using a multichannel pipette, add index adapters to the FSP:

- a. Add 4 µl R7XX down each column.

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

- b. Replace the cap on each Index 1 Adapter tube with a new orange cap.
 - c. Add 4 μ l A50X across each row.
 - d. Replace the cap on each Index 2 Adapter tube with a new white cap.
6. Vortex PCR2, and then centrifuge briefly.
7. [Optional] Evenly distribute PCR2 among each tube of an 8-tube strip. Use a multichannel pipette to dispense.
8. Add 27 μ l PCR2 to each well.
9. Seal and centrifuge at 1000 \times g for 30 seconds.
10. Place on the preprogrammed thermal cycler and run the PCR2 program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 7 days. Alternatively, leave on the thermal cycler overnight.

Purify Libraries

This process uses purification beads to purify the amplified libraries from other reaction components.

Consumables

- RSB (Resuspension Buffer)
- SPB (Sample Purification Beads)
- ProK (Proteinase K)
- Freshly prepared 80% ethanol (EtOH)
- 96-well midi plate
- 96-well PCR plate, skirted or semiskirted
- PVC reagent reservoir
- Microseal 'B' film

Preparation

1. Prepare the following consumables:

| Item | Storage | Instructions |
|------|----------------|--|
| ProK | -25°C to -15°C | Thaw at room temperature, and then centrifuge briefly. Return to storage immediately after use. |
| RSB | 2°C to 8°C | Let stand for 30 minutes to bring to room temperature. Vortex and invert to mix. |
| SPB | 2°C to 8°C | Let stand for 30 minutes to bring to room temperature. Vortex for ≥ 1 minute and invert to mix. |

2. Label plates as follows.

| Plate Type | Label |
|------------|---------------------------------|
| Midi | PBP for Purification Bead Plate |
| PCR | PLP for Purified Library Plate |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

3. Calculate the appropriate volume for the SPB Master Mix:

| SPB Beads Master Mix | SPB | ProK |
|----------------------|---------|---------|
| for 96 samples | 5.50 ml | 25.0 ul |
| for 64 samples | 3.50 ml | 16.4 ul |
| for 48 samples | 2.70 ml | 12.5 ul |
| for 32 samples | 1.75 ml | 8.2 ul |
| for 24 samples | 1.30 ml | 6.2 ul |

Procedure

1. Add 45 µl SPB Master Mix to each well of the PBP.
2. Centrifuge the sealed FSP at 1000 × g for 30 seconds.
3. Transfer 45 µl reaction from each well of the FSP to the corresponding well of the PBP.
4. Discard the FSP plate.
5. Seal the PBP and shake at 1800 rpm for 2 minutes.
6. Incubate at room temperature for 5 minutes.
7. Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
8. Remove and discard all supernatant.
9. Keep on the magnetic stand and wash as follows.
 - a. Add 200 µl fresh 80% EtOH to each well.
 - b. Incubate for 30 seconds.
 - c. Remove and discard all supernatant.
10. Wash a **second** time.
11. Seal and centrifuge at 1000 × g for 30 seconds.
12. Place on the magnetic stand.
13. With a 20 µl pipette, remove residual EtOH from each well.
14. Remove from the magnetic stand.
15. Add 52.5 µl RSB to each well.
16. Seal and shake at 1800 rpm for 2 minutes.
17. If the beads are not fully resuspended, pipette to mix or reshake at 1800 rpm for 2 minutes.
18. Incubate at room temperature for 2 minutes.
19. Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
20. Transfer 50 µl supernatant from each well of the PBP to the corresponding well of the PLP.
21. Seal and centrifuge at 1000 × g for 30 seconds.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 1 year.

Normalize Libraries

This process normalizes the concentration of each library for even representation without post-PCR quantification and individual normalization. Samples of varying types and input amounts achieve consistent cluster density, optimizing the resolution of each library in the pool.

Consumables

- HP3 (2 N NaOH)
- LNA1 (Library Normalization Additives 1)
- LNB1 (Library Normalization Beads 1)
- LNS2 (Library Normalization Storage Buffer 2)
- LNW1 (Library Normalization Wash 1)
- Nuclease-free water
- One of the following tubes:
 - 1.7 ml microcentrifuge tube
 - 15 ml conical tube
- 1.7 ml microcentrifuge tube
- 96-well midi plate
- 96-well PCR plate, skirted or semiskirted
- PVC reagent reservoir
- Microseal 'B' film

About Reagents

- The volumes combined in the LNA1/LNB1 Master Mix tube and the 0.1 N HP3 tube include overage, so calculating additional overage is not necessary.

Warning: This set of reagents contains potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in accordance with applicable regional, national, and local laws and regulations. **For complete environmental, health, and safety information, see the safety data sheets (SDS) at verogen.com/product-documentation.**

Preparation

1. Prepare the following consumables:

| Item | Storage | Instructions |
|------|----------------|--|
| HP3 | -25°C to -15°C | Thaw at room temperature for ≥ 30 minutes. Vortex to mix, and then centrifuge briefly. |
| LNA1 | -25°C to -15°C | Thaw at room temperature for ≥ 30 minutes. Vortex with intermittent inversion |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Item | Storage | Instructions |
|------|--------------|---|
| LNB1 | 2°C to 8°C | Let stand for 30 minutes to bring to room temperature. Vortex for at least 1 minute, inverting 5 times every 15 seconds. Pipette to mix until the bead pellet at the bottom is resuspended. |
| LNW1 | 2°C to 8°C | Let stand for 30 minutes to bring to room temperature. |
| LNS2 | 15°C to 30°C | Remove from storage. |

2. Label vessels as follows.

| Vessel | Label |
|-----------------------------------|-------------------------------------|
| 1.7 ml tube or 15 ml conical tube | LNA1/LNB1 Master Mix |
| Midi plate | NWP for Normalization Working Plate |
| PCR plate | NLP for Normalization Library Plate |

3. Dedicate separate hazardous waste disposal containers for liquids and solids.

Procedure

- In the LNA1/LNB1 Master Mix tube, combine the following volumes. Multiply each volume by the number of samples, but do not add overage.
 - LNA1 (46.8 µl)
 - LNB1 (8.5 µl)For example, for eight samples, combine 374.4 µl LNA1 and 68 µl LNB1.
- Vortex, and then invert several times to mix.
- Transfer the entire volume to a reagent reservoir.
- Add 45 µl LNA1/LNB1 Master Mix to each sample well of the NWP.
- To clear any beads that might have aspirated, place the PLP on the magnetic stand and wait until the liquid is clear (~2 minutes).
- Transfer 20 µl supernatant from each well of the PLP to the corresponding well of the NWP.
- Seal the NWP and shake at 1800 rpm for 30 minutes.
- While the plate is shaking, perform steps 9–11 to save time later in the process.
- In the 0.1 N HP3 tube, combine the following volumes. Multiply each volume by the number of samples, but do not add overage.
 - Nuclease-free water (33.3 µl)
 - HP3 (1.8 µl)For example, eight samples require 266.4 µl nuclease-free water and 14.4 µl HP3.
- Invert several times to mix, and then set aside.
- Add 30 µl LNS2 to each sample well of the NLP.
- Immediately after shaking, place the NWP on the magnetic stand and wait until the liquid is clear (~2 minutes).
- Remove and discard all supernatant.
- Remove from the magnetic stand.
- Wash as follows.

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

- a. Add 45 μ l LNW1 to each well.
 - b. Seal and shake at 1800 rpm for 5 minutes.
 - c. Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
 - d. Remove and discard all supernatant.
 - e. Remove from the magnetic stand.
16. Wash a **second** time.
 17. Seal and centrifuge at 1000 \times g for 30 seconds.
 18. Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
 19. With a 20 μ l pipette, remove residual LNW1 from each well.
 20. Remove from the magnetic stand.
 21. Add 32 μ l freshly prepared 0.1 N HP3 to each well.
 22. Seal and shake at 1800 rpm for 5 minutes.
 23. If the beads are not fully resuspended, pipette to mix or reshake at 1800 rpm for 5 minutes.
 24. Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
 25. Transfer 30 μ l supernatant from the NWP to the corresponding well of the NLP.
 26. Pipette to mix.
 27. Seal and centrifuge at 1000 \times g for 30 seconds.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 30 days.

Pool Libraries

This process combines equal volumes of each normalized library to create a pool of libraries that are sequenced together on the same flow cell.

Consumables

- 1.7 ml microcentrifuge tube
- RNase/DNase-free 8-tube strip and caps
- Microseal 'B' film

Preparation

1. Select libraries to pool for sequencing.
For recommendations, see [Number of Samples \(page 9\)](#).
2. Label the 1.7 ml tube PNL for Pooled Normalized Libraries.

Procedure

1. Using a multichannel pipette, transfer 5 μ l of each library to a new 8-tube strip.
2. Seal the NLP and store in the post-PCR area at -25°C to -15°C for \leq 30 days.
3. Transfer libraries from each well of the 8-tube strip to the PNL tube.

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

4. Cap and vortex to mix, and then centrifuge briefly.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 30 days.

Denature and Dilute Libraries

This process dilutes libraries to the loading concentration, adds a sequencing control, and uses a heat-based method to denature the libraries for sequencing.

Start this process when you are ready to prepare sequencing reagents and set up the run. Delays can impact template loading.

Consumables

- HP3 (2 N NaOH)
- HSC (Human Sequencing Control)
- MiSeq FGx Reagent Kit contents:
 - HT1 (Hybridization Buffer)
 - Reagent cartridge
- Nuclease-free water
- Pooled libraries
- 1.7 ml microcentrifuge tubes (2)

Preparation

1. Prepare the following consumables:

| Item | Storage | Instructions |
|-------------------|----------------|--|
| HP3 | -25°C to -15°C | Thaw at room temperature for ≥ 30 minutes, and then centrifuge briefly. |
| HSC | -25°C to -15°C | Let stand for 30 minutes to bring to room temperature. Invert to mix, and then centrifuge. |
| HT1 | -25°C to -15°C | Thaw at room temperature, and then vortex to mix. |
| Reagent cartridge | -25°C to -15°C | Thaw in a water bath at room temperature. |

2. Preheat the microheating system to 96°C.
3. Label two new 1.7 ml tubes:
 - Denatured HSC
 - DNL for Denatured Normalized Libraries

Procedure

1. In the Denatured HSC tube, combine the following volumes:
 - HSC (2 µl)
 - HP3 (2 µl)
 - Nuclease-free water (36 µl)

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

2. Pipette gently to mix. Cap and centrifuge briefly to mix.
3. Incubate at room temperature for 5 minutes.
4. Add 591 μ l HT1 to the DNL tube.
5. Transfer 7 μ l library from the PNL tube to the DNL tube.
6. Pipette to mix.
7. Cap the PNL tube and store at -25°C to -15°C for \leq 30 days.
Exceeding 30 days significantly reduces cluster density.
8. Add 4 μ l denatured HSC to the DNL tube.
9. Pipette to mix.
You can store the denatured HSC at room temperature for \leq 1 day.
10. Cap and vortex to mix, and then centrifuge briefly.
11. Place on the preheated microheating system and incubate for 2 minutes.
12. Invert several times to mix.
13. Immediately cool on the benchtop cooler or in the ice-water bath for 5 minutes.
14. Immediately transfer the entire volume to the reagent cartridge per instructions in the *MiSeq FGx Sequencing System Reference Guide (document # VD2018006)*.

Materials

| | |
|---------------------------------|----|
| Kit Contents and Storage | 22 |
| Index Adapter Sequences | 25 |
| Consumables and Equipment | 26 |

Kit Contents and Storage

Make sure that you have the reagents identified in this section before starting the protocol. When you receive the kit, promptly store reagents at the indicated temperatures.

| Kit Name | Part # |
|---|-----------|
| ForenSeq DNA Signature Prep Kit (96 Reactions) | V16000023 |
| ForenSeq DNA Signature Prep Kit (384 Reactions) | 15066151 |
| ForenSeq Enhanced PCR1 Buffer System | V16000137 |

All reagents in a box are shipped at the same temperature. When a reagent has a different storage temperature than most other reagents in the box, you can initially store the reagent at the same temperature as the other reagents. After first use, store the reagent at the indicated temperature.

ForenSeq DNA Signature Prep Kit (96 Reactions) (V16000023)

Pre-PCR Box 1

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|-------------------|--------|----------------|
| 1 | 2800M | Control DNA 2800M | Black | 2°C to 8°C* |
| 2 | DPMA | DNA Primer Mix A | Blue | -25°C to -15°C |
| 2 | DPMB | DNA Primer Mix B | Red | -25°C to -15°C |
| 2 | FEM | Enzyme Mix | Yellow | -25°C to -15°C |
| 2 | PCR1 | PCR1 Reaction Mix | Green | -25°C to -15°C |

* Shipped at -25°C to -15°C

Post-PCR Box 2

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|--|--------|-------------------|
| 1 | HP3 | 2 N NaOH | Orange | -25°C to -15°C |
| 1 | HSC | Human Sequencing Control | Pink | -25°C to -15°C |
| 1 | LNA1 | Library Normalization Additives 1 | Clear | -25°C to -15°C |
| 1 | LNS2 | Library Normalization Storage Buffer 2 | Clear | Room temperature* |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|------------------------------|--------|----------------|
| 2 | LNW1 | Library Normalization Wash 1 | Clear | 2°C to 8°C* |
| 2 | PCR2 | PCR2 Reaction Mix | Purple | -25°C to -15°C |
| 1 | A501 | A501 Index Adapter | White | -25°C to -15°C |
| 1 | A502 | A502 Index Adapter | White | -25°C to -15°C |
| 1 | A503 | A503 Index Adapter | White | -25°C to -15°C |
| 1 | A504 | A504 Index Adapter | White | -25°C to -15°C |
| 1 | A505 | A505 Index Adapter | White | -25°C to -15°C |
| 1 | A506 | A506 Index Adapter | White | -25°C to -15°C |
| 1 | A507 | A507 Index Adapter | White | -25°C to -15°C |
| 1 | A508 | A508 Index Adapter | White | -25°C to -15°C |
| 1 | R701 | R701 Index Adapter | Orange | -25°C to -15°C |
| 1 | R702 | R702 Index Adapter | Orange | -25°C to -15°C |
| 1 | R703 | R703 Index Adapter | Orange | -25°C to -15°C |
| 1 | R704 | R704 Index Adapter | Orange | -25°C to -15°C |
| 1 | R705 | R705 Index Adapter | Orange | -25°C to -15°C |
| 1 | R706 | R706 Index Adapter | Orange | -25°C to -15°C |
| 1 | R707 | R707 Index Adapter | Orange | -25°C to -15°C |
| 1 | R708 | R708 Index Adapter | Orange | -25°C to -15°C |
| 1 | R709 | R709 Index Adapter | Orange | -25°C to -15°C |
| 1 | R710 | R710 Index Adapter | Orange | -25°C to -15°C |
| 1 | R711 | R711 Index Adapter | Orange | -25°C to -15°C |
| 1 | R712 | R712 Index Adapter | Orange | -25°C to -15°C |

* Shipped at -25°C to -15°C

Post-PCR Box 3

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|-------------------------------|--------|------------|
| 1 | LNB1 | Library Normalization Beads 1 | White | 2°C to 8°C |
| 1 | RSB | Resuspension Buffer | Purple | 2°C to 8°C |
| 1 | SPB | Sample Purification Beads | Red | 2°C to 8°C |

ForenSeq DNA Signature Prep Kit (384 Reactions) (15066151)

Pre-PCR Box 1

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|-------------------|--------|----------------|
| 2 | 2800M | Control DNA 2800M | Black | 2°C to 8°C* |
| 8 | DPMA | DNA Primer Mix A | Blue | -25°C to -15°C |
| 8 | DPMB | DNA Primer Mix B | Red | -25°C to -15°C |
| 8 | FEM | Enzyme Mix | Yellow | -25°C to -15°C |
| 8 | PCR1 | PCR1 Reaction Mix | Green | -25°C to -15°C |

* Shipped at -25°C to -15°C

Post-PCR Box 2

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|--|--------|-------------------|
| 3 | HP3 | 2 N NaOH | Orange | -25°C to -15°C |
| 1 | HSC | Human Sequencing Control | Pink | -25°C to -15°C |
| 4 | LNA1 | Library Normalization Additives 1 | Clear | -25°C to -15°C |
| 4 | LNS2 | Library Normalization Storage Buffer 2 | Clear | Room temperature* |
| 8 | LNW1 | Library Normalization Wash 1 | Clear | 2°C to 8°C* |
| 8 | PCR2 | PCR2 Reaction Mix | Purple | -25°C to -15°C |
| 1 | A501 | A501 Index Adapter | White | -25°C to -15°C |
| 1 | A502 | A502 Index Adapter | White | -25°C to -15°C |
| 1 | A503 | A503 Index Adapter | White | -25°C to -15°C |
| 1 | A504 | A504 Index Adapter | White | -25°C to -15°C |
| 1 | A505 | A505 Index Adapter | White | -25°C to -15°C |
| 1 | A506 | A506 Index Adapter | White | -25°C to -15°C |
| 1 | A507 | A507 Index Adapter | White | -25°C to -15°C |
| 1 | A508 | A508 Index Adapter | White | -25°C to -15°C |
| 1 | R701 | R701 Index Adapter | Orange | -25°C to -15°C |
| 1 | R702 | R702 Index Adapter | Orange | -25°C to -15°C |
| 1 | R703 | R703 Index Adapter | Orange | -25°C to -15°C |
| 1 | R704 | R704 Index Adapter | Orange | -25°C to -15°C |
| 1 | R705 | R705 Index Adapter | Orange | -25°C to -15°C |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|--------------------|--------|----------------|
| 1 | R706 | R706 Index Adapter | Orange | -25°C to -15°C |
| 1 | R707 | R707 Index Adapter | Orange | -25°C to -15°C |
| 1 | R708 | R708 Index Adapter | Orange | -25°C to -15°C |
| 1 | R709 | R709 Index Adapter | Orange | -25°C to -15°C |
| 1 | R710 | R710 Index Adapter | Orange | -25°C to -15°C |
| 1 | R711 | R711 Index Adapter | Orange | -25°C to -15°C |
| 1 | R712 | R712 Index Adapter | Orange | -25°C to -15°C |

* Shipped at -25°C to -15°C

Post-PCR Box 3

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|-------------------------------|--------|------------|
| 4 | LNB1 | Library Normalization Beads 1 | White | 2°C to 8°C |
| 1 | RSB | Resuspension Buffer | Purple | 2°C to 8°C |
| 2 | SPB | Sample Purification Beads | Red | 2°C to 8°C |

ForenSeq Enhanced PCR1 Buffer System (96 Reactions) (V16000137)

| Quantity | Reagent | Description | Cap | Storage |
|----------|---------|-------------------------------------|--------|----------------|
| 1 | ePCR1 | ForenSeq Enhanced PCR1 Reaction Mix | Orange | -25°C to -15°C |
| 1 | ProK | Proteinase K | Clear | -25°C to -15°C |

Index Adapter Sequences

The following tables list the 8 bp sequences for the index adapters included in the kit.

Index 1 (i7)

| Index Name | Sequence |
|------------|----------|
| R701 | ATCACGAT |
| R702 | CGATGTAT |
| R703 | TTAGGCAT |
| R704 | TGACCAAT |
| R705 | ACAGTGAT |
| R706 | GCCAATAT |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Index Name | Sequence |
|------------|----------|
| R707 | CAGATCAT |
| R708 | ACTTGAAT |
| R709 | GATCAGAT |
| R710 | TAGCTTAT |
| R711 | GGCTACAT |
| R712 | CTTGAAT |

Index 2 (i5)

| Index Name | Sequence |
|------------|----------|
| A501 | TGAACCTT |
| A502 | TGCTAAGT |
| A503 | TGTTCTCT |
| A504 | TAAGACAC |
| A505 | CTAATCGA |
| A506 | CTAGAACA |
| A507 | TAAGTTCC |
| A508 | TAGACCTA |

Consumables and Equipment

Make sure that you have the following user-supplied consumables and equipment before starting the protocol. These items supplement the library prep reagents and index adapters provided in the kit.

The protocol is optimized and validated using the items listed. Comparable performance is not guaranteed when using alternate consumables and equipment.

Consumables

| Consumable | Supplier |
|--|-----------------------------------|
| 1.7 ml microcentrifuge tubes | General lab supplier |
| 15 ml conical tube | General lab supplier |
| 20 µl barrier pipette tips | General lab supplier |
| 200 µl barrier pipette tips | General lab supplier |
| 96-well deep well storage plates (midi plates) | Fisher Scientific, part # AB-0859 |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Consumable | Supplier |
|--|---|
| 96-well twin.tec PCR plates, semiskirted | One of the following suppliers: <ul style="list-style-type: none"> • Eppendorf, catalog # 951020303 • VWR, catalog # 89136-706 |
| 96-well twin.tec PCR plate, skirted, 150 µl | Eppendorf, catalog # 951020401 |
| Ethyl alcohol, pure | Sigma-Aldrich, catalog # E7023 |
| Microseal 'A' sealing film | Bio-Rad, catalog # MSA5001 |
| Microseal 'B' sealing film, adhesive, optical | Bio-Rad, catalog # MSB1001 |
| One of the following kits: <ul style="list-style-type: none"> • MiSeq FGx Reagent Kit • MiSeq FGx Reagent Micro Kit | Verogen part #: <ul style="list-style-type: none"> • 15066817 • 20021681 |
| Multichannel reagent reservoirs, PVC, disposable | Labcor, part # 730-001 |
| Nuclease-free water | General lab supplier |
| RNase/DNase-free 8-tube strips and caps | General lab supplier |

Equipment

| Equipment | Supplier | Pre-PCR | Post-PCR |
|---|---|---------|----------|
| 20 µl multichannel pipettes (8-channel) | General lab supplier | | X |
| 200 µl multichannel pipettes (8-channel) | General lab supplier | | X |
| Benchtop microcentrifuge | General lab supplier | X | X |
| Heating system, 96-well, 1.5 ml | General lab supplier | | X |
| Magnetic stand-96 | Life Technologies, part # AM10027 | | X |
| Microplate centrifuge | General lab supplier | X | X |
| Thermal cycler, 96-well with heated lid | See Thermal Cyclers (below) | | X |
| Thermoshaker, one of the following: <ul style="list-style-type: none"> • BioShake iQ • BioShake XP | QInstruments, item #: <ul style="list-style-type: none"> • 1808-0506 • 1808-0505 | | X |
| Rubber roller | General lab supplier | X | X |
| Vortexer | General lab supplier | X | X |
| [Optional] Benchtop cooler | VWR, catalog # 414004-286 | | X |

Thermal Cyclers

The following table lists supported thermal cyclers with recommended settings. If your laboratory has an unlisted thermal cycler, evaluate the thermal cycler before performing the protocol.

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Thermal Cycler | Temperature Mode | Lid Temperature | Vessel Type |
|--|----------------------------|---------------------------|--------------------------------|
| ABI LTI thermal cycler 9700 ¹ | 9600 emulation | Heated | Polypropylene plates and tubes |
| Bio-Rad | Calculated | Heated, constant at 100°C | Polypropylene plates and tubes |
| Eppendorf Mastercycler Pro S | Gradient S, Simulated Tube | Heated | Plate |
| Proflex 96-well PCR System ² | Not applicable | Heated, constant at 105°C | Polypropylene plates and tubes |
| Veriti 96-well thermal cycler ² | Standard | Heated, constant at 105°C | Polypropylene plates and tubes |

¹ Only gold heat blocks are supported.

² Settings were verified after developmental validation of the ForenSeq DNA Signature Prep Kit.

Amplicon Information

| | |
|--|----|
| Loci Detected with DPMA and DPMB | 29 |
| Interpreting Loci D22S1045 | 37 |
| Interpreting Loci DYS392 | 39 |

Loci Detected with DPMA and DPMB

The following tables list loci detected with DPMA or DPMB. **Loci in the piSNPs and aiSNPs are exclusive to DPMB.** All other loci are detected with both primer mixes.

- Amplicon lengths exclude 120 bp for adapter sequences.
- Amplicon start and end positions are the one-base endpoints of the entire amplicon, including the sequence that matches primers on the hg19 human reference genome.
- Amelogenin is a genetic marker that confirms the gender of the biological sample donor. The size range is 106–112 bp and the control DNA is male.

Autosomal STRs

| Locus | Repeats | Amplicon Length (bp) | Chromosome | 2800M Control Alleles |
|---------------------|---------|----------------------|------------|-----------------------|
| D1S1656 | 7–21.3 | 133–192 | 1 | 12,13 |
| TPOX | 4–16 | 61–109 | 2 | 11,11 |
| D2S441 | 7–17 | 137–177 | 2 | 10,14 |
| D2S1338 | 10–33.1 | 110–203 | 3 | 22,25 |
| D3S1358 | 8–22 | 138–194 | 3 | 17,18 |
| D4S2408 | 8–13 | 98–118 | 4 | 9,9 |
| FGA | 12.2–53 | 150–312 | 4 | 20,23 |
| D5S818 | 4–20 | 98–162 | 5 | 12,12 |
| CSF1PO | 5–17 | 72–120 | 5 | 12,12 |
| D6S1043 | 8–26 | 154–226 | 6 | 12,20 |
| D7S820 ¹ | 5–21.1 | 118–183 | 7 | 8,11 |
| D8S1179 | 6–20 | 82–138 | 8 | 14,15 |
| D9S1122 | 8–15 | 104–132 | 9 | 12,12 |
| D10S1248 | 7–20 | 124–176 | 10 | 13,15 |
| TH01 | 3–14 | 96–140 | 11 | 6,9,3 |
| vWA | 11–26 | 135–195 | 12 | 16,19 |
| D12S391 | 13–28 | 229–289 | 12 | 18,23 |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Repeats | Amplicon Length (bp) | Chromosome | 2800M Control Alleles |
|-----------------------|---------|----------------------|------------|-----------------------|
| D13S317 | 5-17 | 138-186 | 13 | 9,11 |
| PentaE | 5-28.4 | 362-481 | 15 | 7,14 |
| D16S539 | 4-17 | 132-184 | 16 | 9,13 |
| D17S1301 | 9-15 | 130-154 | 17 | 11,12 |
| D18S51 | 6-40 | 136-272 | 18 | 16,18 |
| D19S433 | 4-27 | 148-240 | 19 | 13,14 |
| D20S482 | 9-17 | 125-157 | 20 | 14,15 |
| D21S11 | 12-41.2 | 147-265 | 21 | 29,31.2 |
| PentaD | 1.1-19 | 209-298 | 21 | 12,13 |
| D22S1045 ² | 8-19 | 201-245 | 22 | 16,16 |

¹ Might include a low-level plus 0.1 base pair artifact with one T addition at the end of the STR repeat sequence of the parent allele. For example, 8,8.1 or 11,11.1.

² Interpret with caution. See [Interpreting Loci D22S1045 \(page 37\)](#) for more information.

X-STRs

| Locus | Repeats | Amplicon Length (bp) | Chromosome | 2800M Control Alleles |
|----------|---------|----------------------|------------|-----------------------|
| DXS10074 | 7-22 | 184-244 | X | 21 |
| DXS10103 | 14-21 | 157-185 | X | 18 |
| DXS10135 | 15.3-34 | 239-312 | X | 28 |
| DXS7132 | 11-20 | 175-211 | X | 13 |
| DXS7423 | 10-18 | 188-220 | X | 15 |
| DXS8378 | 8-14 | 434-458 | X | 12 |
| HPRTB | 8-17 | 193-229 | X | 12 |

Y-STRs

| Locus | Repeats | Amplicon Length (bp) | Chromosome | 2800M Control Alleles |
|-----------|---------|----------------------|------------|-----------------------|
| DYF387S1 | 30-44 | 207-263 | Y | 37,38 |
| DYS19 | 9-19 | 269-309 | Y | 14 |
| DYS385a-b | 7-28 | 232-316 | Y | 13,16 |
| DYS389I | 9-17 | 236-268 | Y | 14 |
| DYS389II | 24-34 | 283-323 | Y | 31 |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Repeats | Amplicon Length (bp) | Chromosome | 2800M Control Alleles |
|-----------|---------|----------------------|------------|-----------------------|
| DYS390 | 17-28 | 290-334 | Y | 24 |
| DYS391 | 5-16 | 119-163 | Y | 10 |
| DYS392* | 6-17 | 318-362 | Y | 13 |
| DYS437 | 10-18 | 194-226 | Y | 14 |
| DYS438 | 6-16 | 129-179 | Y | 9 |
| DYS439 | 6-17 | 167-211 | Y | 12 |
| DYS448 | 14-26 | 330-402 | Y | 19 |
| DYS460 | 7-14 | 348-376 | Y | 11 |
| DYS481 | 17-32 | 129-174 | Y | 22 |
| DYS505 | 9-15 | 162-186 | Y | 11 |
| DYS522 | 8-17 | 298-334 | Y | 12 |
| DYS533 | 7-17 | 186-226 | Y | 12 |
| DYS549 | 10-14 | 210-226 | Y | 13 |
| DYS570 | 10-26 | 142-206 | Y | 17 |
| DYS576 | 10-25 | 163-223 | Y | 18 |
| DYS612 | 26-33 | 275-296 | Y | 29 |
| DYS635 | 15-30 | 242-302 | Y | 21 |
| DYS643 | 7-15 | 141-181 | Y | 10 |
| Y-GATA-H4 | 8-15 | 159-187 | Y | 11 |

* Interpret with caution. See [Interpreting Loci *DYS392* \(page 39\)](#) for more information.

Identity Informative SNPs

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Control Alleles |
|------------|----------------------|------------|-------------------------|-----------------------|-----------------------|
| rs10495407 | 109 | 1 | 238439234 | 238439342 | GG |
| rs1294331 | 85 | 1 | 233448359 | 233448443 | GA |
| rs1413212 | 64 | 1 | 242806767 | 242806830 | GG |
| rs1490413 | 98 | 1 | 4367256 | 4367353 | AA |
| rs560681 | 90 | 1 | 160786641 | 160786730 | AG |
| rs891700 | 115 | 1 | 239881850 | 239881964 | AG |
| rs1109037 | 118 | 2 | 10085691 | 10085808 | GG |
| rs12997453 | 100 | 2 | 182413195 | 182413294 | AA |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Control Alleles |
|------------|----------------------|------------|-------------------------|-----------------------|-----------------------|
| rs876724 | 119 | 2 | 114945 | 115063 | CC |
| rs907100 | 115 | 2 | 239563542 | 239563656 | CG |
| rs993934 | 120 | 2 | 124109120 | 124109239 | CC |
| rs1355366 | 119 | 3 | 190806041 | 190806159 | AG |
| rs1357617 | 120 | 3 | 961696 | 961815 | AT |
| rs2399332 | 157 | 3 | 110300999 | 110301155 | AC |
| rs4364205 | 98 | 3 | 32417576 | 32417673 | GG |
| rs6444724 | 120 | 3 | 193207306 | 193207425 | TT |
| rs1979255 | 102 | 4 | 190318007 | 190318108 | GG |
| rs2046361 | 120 | 4 | 10968994 | 10969113 | AA |
| rs279844 | 167 | 4 | 46329584 | 46329750 | AT |
| rs6811238 | 120 | 4 | 169663541 | 169663660 | GG |
| rs13182883 | 169 | 5 | 136633252 | 136633420 | AG |
| rs159606 | 104 | 5 | 17374845 | 17374948 | AA |
| rs251934 | 97 | 5 | 174778619 | 174778715 | TT |
| rs338882 | 157 | 5 | 178690599 | 178690755 | CC |
| rs717302 | 110 | 5 | 2879333 | 2879442 | GG |
| rs13218440 | 170 | 6 | 12059928 | 12060097 | AG |
| rs1336071 | 120 | 6 | 94537182 | 94537301 | GG |
| rs214955 | 120 | 6 | 152697629 | 152697748 | GG |
| rs727811 | 115 | 6 | 165045254 | 165045368 | AA |
| rs321198 | 165 | 7 | 137029715 | 137029879 | TT |
| rs6955448 | 120 | 7 | 4310285 | 4310404 | CT |
| rs737681 | 120 | 7 | 155990742 | 155990861 | TT |
| rs917118 | 109 | 7 | 4456953 | 4457061 | CC |
| rs10092491 | 116 | 8 | 28411037 | 28411152 | CT |
| rs2056277 | 104 | 8 | 139399038 | 139399141 | CC |
| rs4606077 | 151 | 8 | 144656710 | 144656860 | CT |
| rs763869 | 85 | 8 | 1375576 | 1375660 | CT |
| rs1015250 | 117 | 9 | 1823702 | 1823818 | GG |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Control Alleles |
|------------|----------------------|------------|-------------------------|-----------------------|-----------------------|
| rs10776839 | 103 | 9 | 137417271 | 137417373 | GG |
| rs1360288 | 119 | 9 | 128967994 | 128968112 | CC |
| rs1463729 | 99 | 9 | 126881396 | 126881494 | GA |
| rs7041158 | 115 | 9 | 27985907 | 27986021 | CC |
| rs3780962 | 94 | 10 | 17193284 | 17193377 | CC |
| rs735155 | 170 | 10 | 3374133 | 3374302 | AA |
| rs740598 | 120 | 10 | 118506839 | 118506958 | AG |
| rs826472 | 153 | 10 | 2406511 | 2406663 | TT |
| rs964681 | 105 | 10 | 132698394 | 132698498 | CT |
| rs10488710 | 118 | 11 | 115207134 | 115207251 | CG |
| rs1498553 | 111 | 11 | 5708981 | 5709091 | CT |
| rs2076848 | 118 | 11 | 134667502 | 134667619 | AT |
| rs901398 | 90 | 11 | 11096173 | 11096262 | TT |
| rs10773760 | 99 | 12 | 130761623 | 130761721 | AG |
| rs2107612 | 103 | 12 | 888262 | 888364 | AG |
| rs2111980 | 94 | 12 | 106328186 | 106328279 | GG |
| rs2269355 | 65 | 12 | 6945881 | 6945945 | CC |
| rs2920816 | 157 | 12 | 40862976 | 40863132 | TT |
| rs1058083 | 76 | 13 | 100038193 | 100038268 | AG |
| rs1335873 | 109 | 13 | 20901665 | 20901773 | TT |
| rs1886510 | 116 | 13 | 22374646 | 22374761 | CT |
| rs354439 | 170 | 13 | 106938320 | 106938489 | TT |
| rs1454361 | 118 | 14 | 25850765 | 25850882 | AT |
| rs4530059 | 170 | 14 | 104769099 | 104769268 | GG |
| rs722290 | 101 | 14 | 53216686 | 53216786 | GG |
| rs873196 | 114 | 14 | 98845506 | 98845619 | CT |
| rs1528460 | 115 | 15 | 55210664 | 55210778 | TT |
| rs1821380 | 118 | 15 | 39313343 | 39313460 | GG |
| rs8037429 | 63 | 15 | 53616876 | 53616938 | TT |
| rs1382387 | 89 | 16 | 80106318 | 80106406 | GT |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Control Alleles |
|-----------|----------------------|------------|-------------------------|-----------------------|-----------------------|
| rs2342747 | 104 | 16 | 5868645 | 5868748 | AG |
| rs430046 | 119 | 16 | 78016980 | 78017098 | CC |
| rs729172 | 104 | 16 | 5606153 | 5606256 | CC |
| rs740910 | 113 | 17 | 5706552 | 5706664 | AA |
| rs8078417 | 143 | 17 | 80461847 | 80461989 | CT |
| rs938283 | 98 | 17 | 77468433 | 77468530 | TT |
| rs9905977 | 170 | 17 | 2919324 | 2919493 | GG |
| rs1024116 | 98 | 18 | 75432317 | 75432414 | AA |
| rs1493232 | 75 | 18 | 1127945 | 1128019 | AA |
| rs1736442 | 153 | 18 | 55225698 | 55225850 | GG |
| rs9951171 | 119 | 18 | 9749789 | 9749907 | GG |
| rs576261 | 76 | 19 | 39559780 | 39559855 | AC |
| rs719366 | 170 | 19 | 28463281 | 28463450 | TT |
| rs1005533 | 158 | 20 | 39487066 | 39487223 | AA |
| rs1031825 | 126 | 20 | 4447416 | 4447541 | CC |
| rs1523537 | 117 | 20 | 51296076 | 51296192 | CC |
| rs445251 | 119 | 20 | 15124865 | 15124983 | CG |
| rs221956 | 97 | 21 | 43606933 | 43607029 | CC |
| rs2830795 | 114 | 21 | 28608089 | 28608202 | AA |
| rs2831700 | 79 | 21 | 29679639 | 29679717 | AA |
| rs722098 | 101 | 21 | 16685561 | 16685661 | AG |
| rs914165 | 156 | 21 | 42415865 | 42416020 | AG |
| rs1028528 | 78 | 22 | 48362256 | 48362333 | AG |
| rs2040411 | 68 | 22 | 47836378 | 47836445 | AA |
| rs733164 | 120 | 22 | 27816711 | 27816830 | AG |
| rs987640 | 120 | 22 | 33559450 | 33559569 | AT |

Phenotypic Informative SNPs

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Control Alleles |
|---------------------|----------------------|------------|-------------------------|-----------------------|-----------------------|
| rs28777 | 92 | 5 | 33958916 | 33959007 | AA |
| rs12203592 | 110 | 6 | 396273 | 396382 | CC |
| rs4959270 | 161 | 6 | 457655 | 457815 | AC |
| rs683 | 120 | 9 | 12709246 | 12709365 | AC |
| rs1042602 | 113 | 11 | 88911659 | 88911771 | AC |
| rs1393350 | 99 | 11 | 89010977 | 89011075 | GG |
| rs12821256 | 119 | 12 | 89328278 | 89328396 | CT |
| rs12896399 | 73 | 14 | 92773627 | 92773699 | GG |
| rs2402130 | 120 | 14 | 92801169 | 92801288 | AA |
| rs1800407 | 119 | 15 | 28230246 | 28230364 | GG |
| N29insA | 112 | 16 | 89985688 | 89985799 | CC |
| rs1110400 | 173 | 16 | 89986044 | 89986216 | TT |
| rs11547464 | 173 | 16 | 89986044 | 89986216 | GG |
| rs1805005 | 213 | 16 | 89985774 | 89985986 | GG |
| rs1805006 | 213 | 16 | 89985774 | 89985986 | CC |
| rs1805007 | 173 | 16 | 89986044 | 89986216 | CC |
| rs1805008 | 173 | 16 | 89986044 | 89986216 | CC |
| rs1805009 | 227 | 16 | 89986484 | 89986710 | GG |
| rs201326893_Y152OCH | 173 | 16 | 89986044 | 89986216 | CC |
| rs2228479 | 213 | 16 | 89985774 | 89985986 | GG |
| rs885479 | 173 | 16 | 89986044 | 89986216 | GG |
| rs2378249 | 118 | 20 | 33218028 | 33218145 | AA |

Ancestry Informative SNPs

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Allele Control |
|-----------|----------------------|------------|-------------------------|-----------------------|----------------------|
| rs2814778 | 120 | 1 | 159174650 | 159174769 | AA |
| rs3737576 | 98 | 1 | 101709521 | 101709618 | AA |
| rs7554936 | 106 | 1 | 151122413 | 151122518 | CT |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Allele Control |
|-------------|----------------------|------------|-------------------------|-----------------------|----------------------|
| rs10497191 | 101 | 2 | 158667153 | 158667253 | CC |
| rs1834619 | 84 | 2 | 17901444 | 17901527 | GG |
| rs1876482 | 120 | 2 | 17362526 | 17362645 | CC |
| rs260690 | 115 | 2 | 109579681 | 109579795 | AA |
| rs3827760 | 108 | 2 | 109513546 | 109513653 | TT |
| rs6754311 | 98 | 2 | 136707920 | 136708017 | CT |
| rs798443 | 84 | 2 | 7968221 | 7968304 | AA |
| rs12498138 | 119 | 3 | 121459545 | 121459663 | GG |
| rs1919550 | 117 | 3 | 121364112 | 121364228 | AA |
| rs1229984 | 120 | 4 | 100239288 | 100239407 | GG |
| rs3811801 | 114 | 4 | 100244261 | 100244374 | CC |
| rs4833103 | 95 | 4 | 38815462 | 38815556 | AC |
| rs7657799 | 116 | 4 | 105375396 | 105375511 | TT |
| rs7722456 | 114 | 5 | 170202901 | 170203014 | TT |
| rs870347 | 119 | 5 | 6844995 | 6845113 | TT |
| rs16891982* | 108 | 5 | 33951621 | 33951728 | GG |
| rs192655 | 70 | 6 | 90518235 | 90518304 | AG |
| rs3823159 | 119 | 6 | 136482701 | 136482819 | AA |
| rs917115 | 71 | 7 | 28172543 | 28172613 | TT |
| rs1462906 | 84 | 8 | 31896545 | 31896628 | CC |
| rs1871534 | 71 | 8 | 145639652 | 145639722 | CC |
| rs2196051 | 120 | 8 | 122124216 | 122124335 | TT |
| rs6990312 | 111 | 8 | 110602270 | 110602380 | GG |
| rs3814134 | 104 | 9 | 127267664 | 127267767 | TT |
| rs4918664 | 168 | 10 | 94920962 | 94921129 | AA |
| rs1079597 | 167 | 11 | 113296227 | 113296393 | GG |
| rs174570 | 120 | 11 | 61597179 | 61597298 | CC |
| rs2238151 | 113 | 12 | 112211753 | 112211865 | CT |
| rs671 | 136 | 12 | 112241658 | 112241793 | GG |
| rs1572018 | 116 | 13 | 41715225 | 41715340 | AG |

ForenSeq DNA Signature Prep Kit with Enhanced PCR1 Buffer System Reference Guide

| Locus | Amplicon Length (bp) | Chromosome | Amplicon Start Position | Amplicon End Position | 2800M Allele Control |
|-------------|----------------------|------------|-------------------------|-----------------------|----------------------|
| rs2166624 | 71 | 13 | 42579949 | 42580019 | AG |
| rs7326934 | 96 | 13 | 49070482 | 49070577 | GG |
| rs7997709 | 85 | 13 | 34847693 | 34847777 | TT |
| rs9522149 | 119 | 13 | 111827125 | 111827243 | CC |
| rs200354 | 165 | 14 | 99375246 | 99375410 | GG |
| rs12439433 | 100 | 15 | 36219979 | 36220078 | GG |
| rs1426654 | 92 | 15 | 48426457 | 48426548 | AA |
| rs1800414 | 116 | 15 | 28196969 | 28197084 | AA |
| rs735480 | 108 | 15 | 45152321 | 45152428 | TT |
| rs12913832* | 119 | 15 | 28365523 | 28365641 | AG |
| rs459920 | 78 | 16 | 89730800 | 89730877 | TT |
| rs11652805 | 119 | 17 | 62987113 | 62987231 | TT |
| rs17642714 | 118 | 17 | 48726060 | 48726177 | AT |
| rs2593595 | 102 | 17 | 41056210 | 41056311 | TC |
| rs4411548 | 158 | 17 | 40658440 | 40658597 | GG |
| rs4471745 | 67 | 17 | 53568849 | 53568915 | GG |
| rs2042762 | 83 | 18 | 35277568 | 35277650 | AA |
| rs3916235 | 120 | 18 | 67578894 | 67579013 | AG |
| rs4891825 | 106 | 18 | 67867615 | 67867720 | AG |
| rs7226659 | 149 | 18 | 40488180 | 40488328 | GG |
| rs7251928 | 200 | 19 | 4077044 | 4077243 | AA |
| rs310644 | 89 | 20 | 62159472 | 62159560 | AA |
| rs2024566 | 88 | 22 | 41697312 | 41697399 | AA |

* Also used for phenotype prediction.

Interpreting Loci D22S1045


The following sections provide example interpretation methods to help interpret the aSTR locus D22S1045. Determine actual values and methods based on application and internal validation data.

Loci D22S1045 might indicate elevated n-1 repeat stutter, particularly with decreased marker coverage. Heterozygote imbalance might occur regardless of marker coverage. When determining the presence of a DNA mixture, consider multilocus genotype.

Data Trends for D22S1045


Elevated n-1 stutter can occur in low coverage situations, particularly for stutter in STR positions and lengths ≥ 15 . Stutter percentages increase as coverage decreases, and in extreme cases can approach or surpass the read depth of the parent allele. The following figure shows progressively increasing n-1 stutter (15 position) as locus coverage decreases.

Figure 4 Increasing stutter



Heterozygote imbalance can occur at low or high locus coverage. Imbalance increases with a larger spread between allele lengths (for example, 11,18). The following figure shows progressively decreasing intralocus balance (allele count ratio [ACR]) as the allele number spread increases.

Figure 5 Decreasing intralocus balance




Genotype Determination at D22S1045

The following flowcharts illustrate example methods of genotype determination for loci D22S1045. The read values in each flowchart are intended as conservative examples that demonstrate an interpretation method using specific read level guidelines. Base actual values and methods on laboratory application and internal validation data.

One Typed Allele

In the following figure, H is the example allele, H is a true allele, and INC is an inconclusive result. An inconclusive result is a conservative conclusion that eliminates inadvertently typing stutter position. The asterisk (*) accounts for potential drop-out due to imbalance.

Figure 6 Flowchart for one typed allele




Two Typed Alleles

In the following figure, H is the allele with highest number of reads and L is the allele with lower number of reads. An inconclusive result (INC) is a conservative conclusion that eliminates the chance of accidentally typing stutter position when two potential alleles are present with < 150 available reads.

H,L indicates a true allele and L > 50 reads outside the n-1 position = obligate sister. H,* indicates a true allele, and L < 50 reads or L in the n-1 position might be elevated stutter.

Figure 7 Flowchart for two typed alleles



Interpreting Loci DYS392


The following sections provide example interpretation methods to help interpret the Y-STR locus DYS392. Determine actual values and methods based on application and internal validation data.

Loci DYS392 might have indicate n-1 repeat stutter, particularly with decreased marker coverage. Consider multilocus genotype when determining the presence of a DNA mixture.

Data Trends for DYS392

Elevated n-1 stutter can occur when locus coverage is low. Stutter increases as coverage decreases, and in extreme cases can approach or surpass the read depth of the parent allele. The following figure shows progressively increasing n-1 stutter (12 position) as locus coverage decreases.

Figure 8 Increasing stutter




Genotype Determination at DYS392

The following flowcharts illustrate example methods of genotype determination for loci DYS392. The read values in each flowchart are intended as conservative examples that demonstrate an interpretation method using specific read level guidelines. Base actual values and methods on laboratory application and internal validation data.

One Typed Allele

In the following figure, H is the example allele, H is a true allele, and INC indicates an inconclusive result. An inconclusive result is a conservative conclusion that eliminates inadvertently typing stutter position.


Figure 9 Flowchart for one typed allele



Stutter Position Typing

In the following figure, H is the allele with the highest number of reads and L is the allele with the lowest number of reads. INC is an inconclusive result and conservative conclusion to eliminate accidentally typing stutter position. A conclusion of H or L indicates the potential for either to be elevated stutter. A conclusion of H is a true allele, even with L at a high n-1 stutter percent.

Figure 10 Flowchart for n-1 stutter position typing with parent allele



Technical Support

For technical assistance, contact Verogen Technical Support.

Contact Information

Address 11111 Flintkote Avenue
San Diego CA 92121 USA

Website www.verogen.com

Email techsupport@verogen.com

Telephone +1.833.837.6436 toll-free (North America)
+1.858.285.4101 tel
+44 (0) 208 054 8706 (United Kingdom)

Safety data sheets (SDS)—Available for download from verogen.com/product-documentation.

Product documentation—Available for download from verogen.com/product-documentation.

Meet any challenge

Verogen is a dedicated developer of human identification products for sequencing and analysis of forensic genomic samples. Working closely with the forensics community, Verogen places exceptional value on flexible, scalable solutions that deliver results when you need them most.

Verogen

+1.833.837.6436 toll-free
(North America)

+1.858.285.4101
(outside North America)

techsupport@verogen.com

www.verogen.com