# PAXgene<sup>®</sup>

## *Blood DNA Tube*

<b>Code</b>	<b>Language</b>	<b>Page</b>	<b>Code</b>	<b>Language</b>	<b>Page</b>
EN	English	1 – 6	KO	한국어	111 – 116
BG	Български	7 – 13	LT	Lietuvių k.	117 – 123
CS	čeština	14 – 20	LV	Latviski	124 – 130
DA	Dansk	21 – 27	NL	Nederlands	131 – 137
DE	Deutsch	28 – 34	NO	Norsk	138 – 144
EL	Ελληνικά	35 – 41	PL	Polski	145 – 151
ES	Español	42 – 48	PT	Português	152 – 158
ET	Eesti	49 – 55	RO	Română	159 – 165
FI	Suomi	56 – 62	RU	Русский	166 – 172
FR	Français	63 – 69	SK	Slovenský (jazyk)	173 – 179
HR	Hrvatski	70 – 76	SL	Slovenski (jezik)	180 – 186
HU	Magyar	77 – 83	SV	Svenska	187 – 193
IS	Íslenska	84 – 90	TR	Türkçe	194 – 200
IT	Italiano	91 – 97	UK	Українська	201 – 207
JA	日本語	98 – 103	ZH	中文	208 – 213
KK	Қазақша	104 – 110		<b>Symbol and Mark Key</b>	214

## Blood DNA Tube

### EN – English

#### I. Intended Use

The PAXgene Blood DNA Tube is intended to collect, anticoagulate, stabilize, transport and store a venous whole blood sample for preparation of high quality DNA for use with molecular diagnostic test methods that require DNA.

**The performance characteristics of this device have not been established for molecular diagnostic assays in general. Users must validate use of product for their specific molecular diagnostic assay.**

#### Product Features

PAXgene Blood DNA Tube • 2.5 ml • 13 x 75 mm • K<sub>2</sub>EDTA additive  
761165 • Plus • BD Hemogard™ • 100 tubes/shelf carton • 1,000 tubes/case

#### II. Summary and Explanation

The PAXgene Blood DNA Tube provides a means for the collection of whole blood for the isolation of genomic DNA in a closed, evacuated system. Blood is collected under a standard phlebotomy protocol into an evacuated tube that contains K<sub>2</sub>EDTA additive. Complete isolation of DNA is carried out using manual or automated methods such as salting out precipitation, magnetic beads, or silica membrane-based technologies.

#### III. Warnings

Contents of this tube are irritating to skin.


- After inhalation, move the exposed person to fresh air.
- After skin contact, immediately flush skin with plenty of soap and water for at least 15 minutes.
- After eye contact, immediately flush with plenty of water for at least 15 minutes.
- After swallowing, immediately give large amounts of water. Get medical attention.

#### IV. Specimen Collection and Preparation for Analysis

##### A. Required Blood Collection Accessories (Not included with the PAXgene Blood DNA Tube)

1. Blood collection needle such as the BD Vacutainer® Eclipse™ Blood Collection Needle or the BD Vacutainer Eclipse Signal™ Blood Collection Needle. See Ordering Information.
2. A BD Vacutainer Needle Holder must be used to ensure proper function. See Ordering Information.
3. Labels for positive donor identification of samples, if required
4. Alcohol swab for cleansing site
5. Dry sterile gauze
6. Tourniquet
7. Needle disposal container for used needle or needle/holder combination

##### B. Procedure for Specimen Collection

## General Instructions

WEAR GLOVES DURING VENIPUNCTURE AND WHEN HANDLING BLOOD COLLECTION TUBES TO MINIMIZE EXPOSURE HAZARD.

1. Assemble needle in holder. Be sure needle is firmly seated to ensure needle does not unthread during use.
2. Place tube into holder. Note: Do not puncture stopper.
3. Select site for venipuncture.
4. **Apply tourniquet. Prepare venipuncture site with an appropriate antiseptic. DO NOT PALPATE VENIPUNCTURE AREA AFTER CLEANSING.**
5. Place patient's arm in a downward position.
6. Remove needle shield. Perform venipuncture WITH ARM DOWNWARD AND TUBE STOPPER UPPER-MOST.
7. Center tubes in holder when penetrating the stopper to prevent sidewall penetration and resultant premature vacuum loss. Push tube onto needle, puncturing stopper diaphragm.
8. REMOVE TOURNIQUET AS SOON AS BLOOD APPEARS IN TUBE. DO NOT ALLOW CONTENTS OF TUBE TO CONTACT THE STOPPER OR END OF THE NEEDLE DURING PROCEDURE.

Note: Blood may occasionally leak from the needle sleeve. Practice Universal Precautions to minimize exposure hazard. If no blood flows into tube or if blood ceases to flow before an adequate specimen is collected, the following steps are suggested to complete satisfactory collection:

- a. Push tube forward until tube stopper has been penetrated. If necessary, hold in place to ensure complete vacuum draw.
- b. Confirm correct position of needle in vein.
9. When tube has filled to its stated volume and blood flow ceases, remove it from holder.
10. Invert 8 times.
11. Remove needle from vein, applying pressure to puncture site with dry sterile swab until bleeding stops.
12. Once clotting has occurred, apply bandage if desired.
13. After venipuncture, the top of the stopper may contain residual blood. Take proper precautions when handling tubes to avoid contact with this blood.
14. Dispose of needle and holder per your facility's policy and guidelines.
15. Store the filled PAXgene Blood DNA Tube at room temperature (18-25°C). See Performance Characteristics for DNA stability in blood samples at room temperature (18-25°C), 2-8°C, or -20°C.

See Procedure for Freezing and Thawing Specimens in the PAXgene Blood DNA Tubes for details if storage temperature of -70°C / -80°C is desired.

### The following techniques shall be used to prevent possible backflow:

- a. Place donor's arm in a downward position.
- b. Hold tube with the stopper uppermost.
- c. Release tourniquet as soon as blood starts to flow into tube.
- d. Make sure tube additives do not touch stopper or end of the needle during venipuncture.

### C. Procedure for Freezing and Thawing Specimens Collected in the PAXgene Blood DNA Tubes

1. Stand the PAXgene Blood DNA Tube upright in a wire rack. **Do not freeze tubes upright in a Styrofoam™ tray as this may cause the tubes to crack.**
2. The PAXgene Blood DNA Tubes can be stored at -20°C and below. If tubes are to be kept at temperatures below -20°C, freeze them first at -20°C for 24 hours, then transfer them to -70°C or -80°C.

(Continued)

3. Thaw the PAXgene Blood DNA Tubes in a wire rack at ambient temperature (18-25°C) for approximately two hours. **Do not thaw the PAXgene Blood DNA Tubes at temperatures above 25°C.**
4. Carefully invert the thawed PAXgene Blood DNA Tubes 10 times.

**Note:** The frozen PAXgene Blood DNA Tubes are subject to breakage upon impact. To reduce the risk of breakage during shipment, frozen tubes should be treated in the same manner as glass tubes. Users must validate their own freezing and shipping protocol for the PAXgene Blood DNA Tubes.

## D. Procedure for Specimen Preparation for Analysis

Blood samples should be processed in accordance with the instructions provided with DNA sample preparation kit. For examples, see Ordering Information.

## V. Performance Characteristics

After blood is introduced into the tube, the DNA remains stable for 14 days at 18 to 25°C (Figure 1). Based on DNA stability studies, blood can be stored in the tube for 14 days at 18 to 25°C, 28 days at 2°C to 8°C, 3 days at 35°C, or for a minimum of 6 months at -20°C with DNA concentration  $\geq 12.5$  ng/ $\mu$ l and DNA purity ( $A_{260}/A_{280}$ ) of 1.7-1.9.

For the most recent data on the stability of frozen specimens, go to [www.PreAnalytiX.com](http://www.PreAnalytiX.com).

### DNA Stability in Blood Samples at 18-25°C


Figure 1

Blood was drawn from 12 consented adult donors and stored at 18-25°C for the indicated number of days followed by total DNA purification. Blood was collected and stored in PAXgene Blood DNA Tubes, and total DNA was purified using the QIAGEN® QIA Symphony® DSP DNA Mini Kit on the QIAGEN QIA Symphony automated instrument.

## VI. Limitations

1. The quantity of blood drawn should be 2.5 ml per PAXgene Blood DNA Tube, but this volume may vary with altitude, ambient temperature, barometric pressure, tube age, venous pressure, and filling technique.

## VII. Precautions


1. Practice universal precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to bloodborne pathogens.

- Handle all biologic samples and blood collection “sharps” according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury), since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases.  
Utilize any built-in needle protector, if the blood collection set provides one. PreAnalytiX does not recommend resheilding used needles. However, the policies and procedures of your facility may differ and must always be followed.
- Discard all blood collection tubes in biohazard containers approved for their disposal.
- Do not re-use the PAXgene Blood DNA Tubes.
- Do not use the PAXgene Blood DNA Tubes after the expiration date printed on the tube label.
- Excessive centrifugation speed (over 10,000 RCF) may cause PAXgene Blood DNA Tube breakage, exposure to blood and possible injury.
- Do not transfer a sample from a syringe to a tube.

### VIII. Storage


- Store the unused PAXgene Blood DNA Tubes at 18 to 25°C. Limited excursion temperatures up to 40°C are permitted. Do not use tubes after their expiration date.

### IX. Instructions for Removal of BD Hemogard Closure

- Grasp the PAXgene Blood DNA Tube with one hand, placing the thumb under the BD Hemogard Closure. (For added stability, place arm on solid surface). With the other hand, twist the BD Hemogard Closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.
- Move thumb away before lifting closure. **DO NOT** use thumb to push closure off tube. If the tube contains blood, an exposure hazard exists. To help prevent injury during closure removal, it is important that the thumb used to push upward on the closure be removed from contact with the tube as soon as the BD Hemogard Closure is loosened.
- Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

### X. Instructions for Reinsertion of BD Hemogard Closure

- Replace closure over tube.
- Twist and push down firmly until stopper is fully resealed. Complete reinsertion of the stopper is necessary for the closure to remain securely on the tube during handling.

## Technical Assistance

If you have any questions regarding the PAXgene Blood DNA Tube, contact one of the BD Technical Services Departments listed in the BD - Customer Service section.

## Ordering Information

### PAXgene Products

#### **PAXgene Blood DNA Tube / 761165**

Content: 1,000 blood collection tubes

### BD Products\*

#### **BD Vacutainer Eclipse Blood Collection Needle / 368609 CE**

Content: 21G 1¼" inch needle, 0.8 mm x 32 mm, 48/box, 480/case

#### **BD Vacutainer Eclipse Blood Collection Needle with Pre-Attached Holder / 368650**

Content: 21G 1¼" inch needle, 0.8 mm x 32 mm, 100/case

#### **BD Vacutainer Eclipse Signal Blood Collection Needle / 368835 CE**

Content: 21G 1¼" inch needle, 0.8 mm x 32 mm, 50/box, 400/case

#### **BD Vacutainer One Use Holder / 364815**

Content: Case only for 13 mm and 16 mm diameter 1000/case

*\* These catalog numbers represent typical products that can be used with the PAXgene Blood DNA Tube.*

### To order BD Blood Collection Accessories:

Go to: <http://www.bd.com/support/contact/international.asp>

• Select country • Scroll down for BD Diagnostics - Preanalytical Systems • Dial number

### QIAGEN Products\*

#### **QIAGEN QIASymphony DSP DNA Kit / 937255 Midi (96) / 937236 Mini (192)**

Content: Reagent cartridge, enzyme rack, piercing lid, buffer ATE, reuse seal set

#### **QIAGEN QIAamp® DSP DNA Blood Mini Kit / 61104 (50)**

Content: QIAamp Mini spin columns, collection tubes, buffer AL, buffer ATL, buffer AW1, buffer AW2, buffer AE, QIAGEN protease, protease solvent

#### **QIAGEN Genra® Puregene® Blood Kit / 158422 (3 ml) / 158445 (30 ml) / 158467 (120 ml) / 158389 (1000 ml)**

Content: RBC lysis solution, cell lysis solution, protein precipitation solution, DNA hydration solution, RNase A solution

*\* These catalog numbers represent typical DNA preparation kits that can be used with the PAXgene Blood DNA Tube.*

### To order QIAGEN DNA Preparation Kits:

Go to: <http://www.qiagen.com/Products/Ordering-Information/>

## BD – Customer Service

### Argentina, Uruguay and Paraguay

Orders: 011.4551.7100  
E-mail: crc\_argentina@bd.com

### Australia

Orders: 1.800.656.100  
Fax: 1.800.656.110  
E-mail: bd\_anz@bd.com

### Austria

Orders: 43.1.7063660  
Fax: 43.1.7063660-11

### Belgium

Orders: 32.53.720.556  
Fax: 32.53.720.549

### Brazil

Orders: 0800.055.5654  
E-mail: consultoria\_vacutainer@bd.com

### Canada

Orders: 800.268.5430  
Fax: 800.565.0897  
E-mail: custsvccan@bd.com

### Denmark

Orders: 45.43.43.45.66  
Fax: 45.43.96.56.76  
E-mail: bdcs\_dk@europe.bd.com

### East Europe

Orders: 48.22.651.75.88  
Fax: 48.22.651.75.89

### Finland

Orders: 358.9.88.70.780  
Fax : 358.9.88.70.7816  
E-mail: bdorder\_finland@europe.bd.com

### France

Orders: 33.476.689496  
Fax: 33.476.683693  
E-mail: BD\_Vacutainer\_France@europe.bd.com

### Germany

Orders: 49.6221305553  
Fax: 49.6221305379

### Italy

Orders: 39.02.48.240.500  
Fax: 39.02.48.240.344  
Technical: 39.02.48.240.264  
E-mail: ordini\_medicale@europe.bd.com

### Middle East & Africa (EMA)

Orders: 971.4.3379525  
Fax: 971.4.03379551  
E-mail: EMA\_PAS@bd.com

### The Netherlands

Orders: 31.20.6545.716  
Fax: 31.20.5829.421  
E-mail:  
vacutainer\_systems\_bnl@europe.bd.com

### New Zealand

Orders: 0800.572.468  
Fax: 0800.572.469  
E-mail: bd\_anz@bd.com

### Norway

Orders: 47.73.59.12.00  
Fax: 47.73.59.12.01

### Spain

Orders: 34.902.27.17.27  
Fax: 34.91.848.81.04  
E-mail: servicio\_de\_clientes@bd.com

### Sweden

Orders: 46.8.775.51.60  
Fax 46: 8.775.51.95  
E-mail: bdorder\_sweden@europe.bd.com

### Switzerland

Orders: 41.61.4852222  
Fax: 41.61.4852200

### UK

Orders: 44.1.865.781.666  
Fax: 44.1.865.717.313  
E-mail:  
bduk\_customerservice@europe.bd.com

### USA

Orders: 888.237.2762  
Fax: 800.847.2220  
Technical: 800.631.0174  
E-mail:  
www.bd.com/vacutainer/contact/

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