

## DNeasy<sup>®</sup> Plant Maxi Kit

The DNeasy Plant Maxi Kit (cat. no. 68163) should be stored dry at 2–8°C upon arrival and are stable for up to 1 year under these conditions.

For more information, please refer to the *DNeasy Plant Handbook*, which can be found at [www.qiagen.com/handbooks](http://www.qiagen.com/handbooks).

For technical assistance, please call toll-free 00800-22-44-6000, or find regional phone numbers at [www.qiagen.com/contact](http://www.qiagen.com/contact).

### Notes before starting

- Perform all centrifugation steps at 3000–5000 x g (5000 x g is preferable) at room temperature (15–25°C) in a centrifuge with a swing-out rotor.
- If necessary, redissolve any precipitates in Buffer AP1 and Buffer AP3/E concentrates.
- Add ethanol to Buffer AW and Buffer AP3/E concentrates.
- Preheat a water bath or heating block to 65°C.
- Preheat Buffer AP1 to 65°C

1. Disrupt samples ( $\leq 1$  g wet weight or  $\leq 0.2$  g lyophilized tissue) using the TissueRuptor<sup>®</sup>, the TissueLyser, or a mortar and pestle.
2. Add 5 ml preheated Buffer AP1 and 10  $\mu$ l RNase A. Vortex vigorously until no tissue clumps are visible. Incubate for 10 min at 65°C. Invert tube 2–3 times during incubation.

**Note:** Do not mix Buffer AP1 and RNase A before use.

3. Add 1.8 ml Buffer AP2. Mix and incubate for 10 min on ice.
4. Centrifuge at 3000–5000 x g for 5 min.
5. Decant the supernatant into a QIAshredder spin column placed in a 50 ml collection tube. Centrifuge at 3000–5000 x g for 5 min.

January 2011



6. Transfer the flow-through, without disturbing the pellet if present, into a new 50 ml tube. Add 1.5 volumes Buffer AP3/E, and mix by vortexing.  
**Note:** It is important to add Buffer AP3/E to the lysate and to mix immediately.
7. Transfer the mixture from step 6 (maximum 15 ml) into a DNeasy Maxi spin column placed in a 50 ml collection tube. Centrifuge at 3000–5000 x g for 5 min. Discard the flow-through.
8. Add 12 ml Buffer AW, and centrifuge for 10 min at 3000–5000 x g. Discard the flow-through.
9. Transfer the spin column to a new 50 ml tube.
10. Add 0.75–1.0 ml Buffer AE for elution. Incubate for 5 min at room temperature (15–25°C). Centrifuge for 5 min at 3000–5000 x g.
11. Repeat step 10.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN® kit handbook or user manual.

Trademarks: QIAGEN®, DNeasy®, TissueRuptor® (QIAGEN Group). 1066958  
01/2011 © 2011 QIAGEN, all rights reserved.

