
June 2019

Epitect[®] Fast Bisulfite Conversion Handbook

For sample lysis and complete bisulfite conversion/cleanup of DNA from FFPE, blood, cultured cells, or tissue samples, optimized for methylation analysis

Contents

Kit Contents	3
Shipping and Storage	6
Intended Use.....	6
Safety Information.....	7
Quality Control.....	7
Introduction.....	8
Principle and Procedure.....	10
Description of Protocols	15
Equipment and Reagents Supplied by User.....	16
Important Notes.....	17
Yield and Size of DNA	17
Preparation of Reagents	19
Protocol: Bisulfite Conversion of Unmethylated Cytosines in DNA.....	21
Protocol: Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from FFPE Tissue Samples	26
Protocol: Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from Whole Blood, Cultured Cells, or Tissue	30
Troubleshooting Guide	34
Ordering Information	37
Document Revision History.....	39

Kit Contents

Epitect Fast DNA Bisulfite Kit	(50)	(200)
Catalog no.	59824	59826
Number of preps	50	200
Bisulfite Solution	3 x 1.5 ml	12 x 1.5 ml
DNA Protect Buffer	1.9 ml	4 x 1.9 ml
RNase-Free Water	1.9 ml	4 x 1.9 ml
MinElute® DNA Spin Columns	50	200
Collection Tubes (2 ml)	100	400
Buffer BL*	31 ml	3 x 31 ml
Buffer BW (concentrate)	2 x 13 ml	2 x 52 ml
Buffer BD (concentrate)	3 ml	4 x 3 ml
Buffer EB	1 ml	4 x 1 ml
Carrier RNA	310 µg	3 x 310 µg
Quick Start Protocol	2	2

* Contains a guanidine salt. Not compatible with disinfectants containing bleach. See "Safety Information".

Epipect Fast FFPE Bisulfite Kit	(50)
Catalog no.	59844
Number of preps	50
Epipect Fast FFPE Lysis Kit Deparaffinization Solution	8 ml
Lysis Buffer FTB	0.8 ml
Proteinase K	1.4 ml
Epipect Fast DNA Bisulfite Kit Bisulfite Solution	3 x 1.5 ml
DNA Protect Buffer	1.9 ml
RNase-Free Water	1.9 ml
MinElute DNA Spin Columns	50
Collection Tubes (2 ml)	100
Buffer BL*	31 ml
Buffer BW (concentrate)	2 x 13 ml
Buffer BD (concentrate)	3 ml
Buffer EB	1 ml
Carrier RNA	310 µg
Quick Start Protocol	3

* Contains a guanidine salt. Not compatible with disinfectants containing bleach. See "Safety Information".

Epipect Fast LyseAll Bisulfite Kit	(50)	(200)
Catalog no.	59864	59866
Number of preps	50	200
Epipect Fast LyseAll Lysis Kit Buffer EL	25 ml	4 x 25 ml
Lysis Buffer FTB	0.8 ml	4 x 0.8 ml
Proteinase K	1.4 ml	4 x 1.4 ml
Epipect Fast DNA Bisulfite Kit Bisulfite Solution	3 x 1.5 ml	12 x 1.5 ml
DNA Protect Buffer	1.9 ml	4 x 1.9 ml
RNase-Free Water	1.9 ml	4 x 1.9 ml
MinElute DNA Spin Columns	50	200
Collection Tubes (2 ml)	100	400
Buffer BL*	31 ml	3 x 31 ml
Buffer BW (concentrate)	2 x 13 ml	2 x 52 ml
Buffer BD (concentrate)	3 ml	4 x 3 ml
Buffer EB	1 ml	4 x 1 ml
Carrier RNA	310 µg	3 x 310 µg
Quick Start Protocol	3	3

* Contains a guanidine salt. Not compatible with disinfectants containing bleach. See "Safety Information".

Shipping and Storage

The EpiTect Fast DNA Bisulfite Kit, EpiTect Fast FFPE Bisulfite Kit, and EpiTect Fast LyseAll Bisulfite Kit are shipped at room temperature (15–25°C). Upon arrival, the MinElute DNA spin columns, DNA Protect Buffer, and Buffer BD should be stored at 2–8°C. However, short-term storage (up to 4 weeks) at room temperature does not affect their performance. All other buffers including the Bisulfite Solution should be stored at room temperature and are stable for at least 6 months under these conditions.

Lyophilized carrier RNA can be stored at room temperature for 1 year. Carrier RNA should only be dissolved in RNase-free water. Dissolved carrier RNA should be added immediately to Buffer BL, as described in “Preparation of Reagents” on page 19. This solution should be prepared fresh, and is stable at 2–8°C for up to 48 hours. Unused portions of carrier RNA dissolved in RNase-free water should be frozen in aliquots at –30 to –15°C and can be stored for up to 1 year.


Intended Use

The EpiTect Fast DNA Bisulfite Kit, EpiTect Fast FFPE Bisulfite Kit, and EpiTect Fast LyseAll Bisulfite Kit are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

<p>CAUTION</p> 	<p>DO NOT add bleach or acidic solutions directly to the Bisulfite (W2) Solution or to waste containing Buffer BL.</p>
---	--

Buffer BL contains a guanidine salt, which can form highly reactive compounds when combined with bleach. If liquid containing this buffer is spilled, clean with suitable laboratory detergent and water.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of EpiTect Fast DNA Bisulfite Kit, EpiTect Fast FFPE Bisulfite Kit, and EpiTect Fast LyseAll Bisulfite Kit is tested against predetermined specifications to ensure consistent product quality.

Introduction

Epigenetics describes the study of heritable changes in gene function that occur without a change in the nuclear DNA sequence. In addition to RNA-associated silencing and histone modification, a major epigenetic mechanism in higher-order eukaryotes is DNA methylation.

Methylation of DNA occurs on cytosine residues, especially on CpG dinucleotides enriched in small regions of DNA (<500 bp). These regions, with a GC content greater than 55%, are known as CpG islands. They are usually clustered around the regulatory region of genes and can affect the transcriptional regulation of these genes. Methylation of CpG islands by DNA methylases has been shown to be associated with gene inactivation and plays an important role in the development of cancer and cell aging. Reversal of DNA methylation at these sites is a potential therapeutic strategy as this reversal may restore expression of transcriptionally silenced genes. In addition to CpG, methylated cytosine residues are also found at CpNpG or CpNpN sites (N = A, T, or C) in plants.

The methylation status of a DNA sequence can best be determined using bisulfite conversion. Bisulfite treatment of the target DNA results in conversion of unmethylated cytosine residues into uracil, leaving the methylated cytosines unchanged. Therefore, bisulfite treatment gives rise to different DNA sequences for methylated and unmethylated DNA (see below).

	Original Sequence	After bisulfite sequence
Unmethylated DNA	N-C-G-N-C-G-N-C-G-N	N-U-G-N-U-G-N-U-G-N
Methylated DNA	N-C-G-N-C-G-N-C-G-N	N-C-G-N-C-G-N-C-G-N

The most critical step for the correct determination of a methylation pattern is the complete conversion of unmethylated cytosines. This is achieved by incubating the DNA in high bisulfite salt concentrations at high temperature and low pH. These harsh conditions often lead to a

high degree of DNA fragmentation and subsequent loss of DNA during purification. Purification is necessary to remove bisulfite salts and chemicals used in the conversion process that inhibit sequencing. Common bisulfite procedures usually require high amounts of input DNA to compensate for DNA degradation during conversion and DNA loss during purification that often lead to low DNA yield, highly fragmented DNA, and irreproducible conversion rates.

EpiTect Fast DNA Bisulfite Conversion Kits now provide a very fast and streamlined procedure for efficient conversion and purification of DNA prepared from FFPE, blood, cell, or tissue samples. The kits contain highly concentrated Bisulfite Solution, which reduces the time required to convert unmethylated cytosine residues into uracil from 5 hours to as little as 30 minutes, as well as preparation buffers that make it unnecessary to isolate the DNA prior to bisulfite treatment. DNA fragmentation is prevented during the bisulfite conversion by the unique DNA Protect Buffer, which contains a pH-indicator dye as a mixing control in reaction setup, allowing confirmation of the correct pH for cytosine conversion.

Furthermore, the bisulfite thermal cycling program provides an optimized series of incubation steps necessary for thermal DNA denaturation and subsequent sulfonation and cytosine deamination, enabling high cytosine conversion rates of over 99%. Desulfonation, the final step in chemical conversion of cytosines, is achieved by a convenient on-column step included in the purification procedure. The final elution volume can be as low as 10 µl though this may result in a yield reduction.


Principle and Procedure

The EpiTect Fast Bisulfite Conversion Kit procedure comprises a few simple steps: preparation of DNA from sample, bisulfite-mediated conversion of unmethylated cytosines; binding of the converted single-stranded DNA to the membrane of a MinElute DNA spin column; washing; desulfonation of membrane-bound DNA; washing of the membrane-bound DNA to remove the desulfonation agent; and elution of the pure, converted DNA from the spin column. Sample preparation is different for FFPE slices, whole blood, and cell cultures or tissues, whereas the procedure for bisulfite conversion of extracted DNA is the same for all sample types (see flowcharts on pages 11–12). The eluted, bisulfite converted DNA is suited for all techniques currently used for the analysis of DNA methylation, including PCR, real-time PCR, methylation-specific PCR, bisulfite sequencing (direct and cloning), COBRA, and Pyrosequencing®.


Bisulfite Solution

The Bisulfite Solution is conveniently provided in separate aliquots (17 conversion reactions per aliquot) that are ready to use. The bisulfite in each aliquot is supplied in a unique formulation that provides the optimal pH for complete conversion of cytosine to uracil, without the need for tedious pH adjustment. Bisulfite Solution can be stored at room temperature (15–25°C) for at least 6 months.

EpiTect Fast FFPE Bisulfite Conversion Procedure




EpiTect Fast LyseAll Bisulfite Conversion Procedure



After setting up the bisulfite reaction, continue with the procedure shown in the flowchart “EpiTect Fast DNA Bisulfite Kit Procedure”, page 11.

EpiTect Fast DNA Bisulfite Conversion Procedure



DNA Protect Buffer

DNA Protect Buffer is uniquely formulated to prevent the fragmentation usually associated with bisulfite treatment of DNA at high temperatures and low pH values. It also provides effective DNA denaturation, resulting in the single-stranded DNA necessary for complete cytosine conversion. In addition, the DNA Protect Buffer contains a pH indicator dye as a mixing control and to allow confirmation of the correct pH for cytosine conversion.

Bisulfite thermal cycling

The thermal cycling program provides an optimized series of short incubation steps necessary for thermal DNA denaturation and subsequent sulfonation and cytosine deamination, enabling the highest cytosine conversion rates.

Carrier RNA

Carrier RNA is provided to enhance the binding of small quantities of DNA to the MinElute DNA spin-column membrane. If using more than 100 ng of genomic DNA template, it is not necessary to use carrier RNA, though we strongly recommend its use when processing fragmented DNA or DNA isolated from formalin-fixed, paraffin-embedded (FFPE) tissues samples. Carrier RNA does not influence downstream applications. It should be dissolved in RNase-free water before use (see "Preparation of reagents", page **Error! Bookmark not defined.**).

Optimized Buffers

EpiTect Fast DNA Conversion Kits contain carefully optimized buffers, enabling effective lysis of samples (FFPE, cells, tissue, whole blood) in DNA isolation procedures, deparaffinization of FFPE tissues, and maximum cytosine conversion and subsequent DNA purification. Buffer BL promotes binding of the converted single-stranded DNA to the MinElute DNA spin-column membrane. Subsequently, the membrane-bound DNA is washed using Buffer BW, which efficiently removes residual Bisulfite Solution. After desulfonation using Buffer BD, the DNA is

further desalted using Buffer BW before elution from the spin-column membrane using Buffer EB.

Storage stability of converted and purified DNA

The DNA converted and purified using EpiTect Fast Bisulfite Conversion Kits can be stored at -20°C for at least 9 months without decrease in quality or conversion. Further investigations into long-term storage are ongoing. Contact QIAGEN for more information.

Description of Protocols

The choice of kit and the corresponding protocol to use is determined by the type of starting material. The protocols are interconnected and the actual bisulfite conversion reaction is the same for all samples, as illustrated in the flowcharts on pages 11–12. The standard protocol “Protocol: Bisulfite Conversion of Unmethylated Cytosines in DNA”, on page 21, can be used for conversion of 1 ng – 2 µg DNA in a volume of up to 20 µl or 1–500 ng in a maximum volume of 40 µl.

Formalin-fixed, paraffin-embedded (FFPE) tissues are processed with the EpiTect Fast FFPE Bisulfite Kit and the protocol “Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from FFPE Tissue Samples”. This kit consists of the EpiTect Fast FFPE Lysis Kit, containing specialized buffers for efficient deparaffinization and lysis of the tissue sample, and the EpiTect Fast DNA Bisulfite Kit for the bisulfite conversion of the extracted DNA. The protocol includes an optimized step to facilitate binding of DNA and can be used with single slices of FFPE tissue (10 µm in thickness).

Whole blood, cultured cells, or tissue samples are processed with the EpiTect Fast LyseAll Bisulfite Kit and the protocol “Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from Whole Blood, Cultured Cells, or Tissue”. This kit consists of the EpiTect Fast LyseAll Lysis Kit, including an innovative lysis buffer, and the EpiTect Fast DNA Bisulfite Kit for the bisulfite conversion of the extracted DNA. The protocol includes an optimized step to facilitate binding of DNA and can be used with 0.5–20 µl blood or 10–10⁵ cells (as little as 60 pg of DNA).

All protocols achieve the same cytosine conversion rates and lead to equal DNA recoveries after purification of converted DNA, independent of DNA starting amounts.

Equipment and Reagents Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- Ethanol (molecular biology grade, 96–100%)*
- Pipets and pipet tips (we recommend pipet tips with aerosol barriers for preventing cross-contamination)
- 0.2 ml reaction tubes or 8-well strips
- Thermal cycler with heated lid (since the bisulfite reaction is not overlaid with mineral oil, only thermal cyclers with heated lids are suitable for this procedure)
- 1.5 ml microcentrifuge tubes for elution steps (available from Brinkmann [Safe-Lock, cat. no. 022363204], Eppendorf® [Safe-Lock, cat. no. 0030 120.086], or Sarstedt® [Safety Cap, cat. no. 72.690])†
- Microcentrifuge
- Optional: Heating block, thermomixer, or heated orbital incubator

* Do not use denatured alcohol, which contains other substances such as methanol or methylethylketone.

† This is not a complete list of suppliers and does not include many important vendors of biological supplies.

Important Notes

Yield and Size of DNA

The yield of DNA purified after bisulfite conversion depends on the amount of DNA and source of the starting material. DNA can be prepared from FFPE tissue (10 µm in thickness), 0.5–20 µl blood, or 10–10⁵ cells.

If following the protocol “Bisulfite Conversion of Unmethylated Cytosines in DNA”, the EpiTect Fast DNA Bisulfite Kit is suited for DNA input amounts ranging from 1 ng to 2 µg, with high levels of DNA recovery throughout this range.

The size of the template DNA can vary between 500 bp (in laser microdissections) and 30 kb (fresh samples or blood). DNA purified from serum, urine, or FFPE tissue may be <500 bp in length.

If purifying bisulfite-treated DNA originating from very small sample amounts or that is very fragmented (e.g., from biopsies or FFPE tissues), we strongly recommend adding carrier RNA to Buffer BL (see “Preparation of Reagents”).

Note: The purified sample will contain considerably more carrier RNA than DNA. Carrier RNA does not influence downstream applications.

Starting material

Bisulfite conversion also depends on the nature of DNA used as starting material. Genomic DNA should be used for bisulfite treatment without any previous restriction digest step.

If working with plasmid DNA, the DNA should be linearized before starting the procedure due to the very quick reannealing of single-stranded DNA after the denaturation step.

Handling of MinElute DNA spin columns

Due to the sensitivity of nucleic acid amplification technologies, the following precautions are necessary when handling MinElute DNA spin columns to avoid cross-contamination between sample preps:

- Carefully pipet the sample or solution into the MinElute DNA spin column without wetting the rim of the column. Avoid touching the MinElute DNA spin column membrane with the pipet tip.
- Always change pipet tips between liquid transfers. We recommend the use of aerosol-barrier pipet tips.
- Open one MinElute DNA spin column at a time, and take care to avoid generating aerosols.
- Wear gloves throughout the entire procedure. In case of contact between gloves and sample, change gloves immediately.

Centrifugation

MinElute DNA spin columns will fit into most standard 1.5–2 ml microcentrifuge tubes. A set of 2 ml collection tubes is supplied for the dry centrifugation step.

All centrifugation steps should be carried out at room temperature (15–25°C).

Processing MinElute DNA spin columns in a microcentrifuge

Always close MinElute DNA spin columns before placing them in the microcentrifuge.

For efficient parallel processing of multiple samples, we recommend filling a rack with the collection tubes into which the MinElute DNA spin columns can be transferred after centrifugation. Collection tubes can be used several times.

Preparation of Reagents

Buffer BW

For 50-prep kits (cat. nos. 59824, 59844, and 59864): Add 30 ml ethanol (96–100%) to the 13 ml Buffer BW concentrate and store at room temperature (15–25°C). Invert the bottle several times before starting the procedure.

For 200-prep kits (cat. nos. 59826 and 59866): Add 120 ml ethanol (96–100%) to the 52 ml Buffer BW concentrate and store at room temperature (15–25°C). Invert the bottle several times before starting the procedure.

Buffer BD

Add 27 ml ethanol (96–100%) to Buffer BD and store at 2–8°C. Invert the bottle several times before starting the procedure and make sure to close the bottle immediately after use. White precipitates may form in the Buffer BD–ethanol mix after some storage time. These precipitates will not affect the performance of Buffer BD. However, avoid transferring precipitates to the MinElute DNA spin column.

Carrier RNA

Add 310 µl RNase-free water to the lyophilized carrier RNA (310 µg) to obtain a 1 µg/µl solution. Dissolve the carrier RNA thoroughly by vortexing. When processing 50 samples in parallel, add the complete volume of dissolved carrier RNA to the bottle of Buffer BL, and check the box on the bottle lid label. If processing fewer samples, split the dissolved carrier RNA into conveniently sized aliquots (e.g., 50 µl) and store at –20°C. Aliquots can be stored for up to 1 year. If fewer than 50 conversions will be performed in a 2-week period, only make enough Buffer BL–carrier RNA solution as required (see Table 1, for example volumes). Carrier RNA enhances binding of DNA to the spin-column membrane, especially if there are very few target molecules in the sample. Carrier RNA is not necessary if >100 ng DNA is used.

Add dissolved carrier RNA to Buffer BL. Calculate the volume of Buffer BL and dissolved carrier RNA required for the number of samples to be processed (see Table 1 for example volumes). If Buffer BL contains precipitates, dissolve by heating (maximum 70°C) with gentle agitation.

Table 1. Buffer BL and Carrier RNA Volumes

Number of samples	Volume of Buffer BL*	Volume of carrier RNA solution†
1	350 µl	3.5 µl
5	1.75 ml	17.5 µl
10	3.5 ml	35 µl
15	5.25 ml	52.5 µl
25	8.75 ml	87.5 µl
50	17.5 ml	175 µl

* The volumes given contain a 10% surplus for pipetting inaccuracies.

† Resulting in a final concentration of 10 µg/ml carrier RNA in Buffer BL.

Protocol: Bisulfite Conversion of Unmethylated Cytosines in DNA

This protocol enables bisulfite conversion of DNA amounts of 1 ng – 2 µg in a volume of up to 20 µl (high concentration), or 1–500 ng in a maximum volume of 40 µl (low concentration). Conversion of high-concentration or low-concentration DNA samples differs only in the setup of the bisulfite reactions (see Table 2). All other protocol steps are the same.

Important points before starting

- DNA Protect Buffer should turn from green to blue after addition to the DNA–Bisulfite Solution mixture, indicating sufficient mixing and correct pH for the bisulfite conversion reaction.
- Perform all centrifugation steps at room temperature (15–25°C).

Things to do before starting

- Prepare the kit reagents as described in “Preparation of Reagents”.
- Equilibrate samples and buffers to room temperature.
- Optional: Set a thermomixer, heating block, or heated orbital incubator to 60°C to dissolve the Bisulfite Solution.

Procedure

Bisulfite conversion of DNA

1. Thaw DNA to be used in the bisulfite reactions. Ensure that the Bisulfite Solution is completely dissolved.

Note: If necessary, heat the Bisulfite Solution to 60°C and vortex until all precipitates are dissolved again.

Note: Do not place dissolved Bisulfite Solution on ice.

2. Prepare the bisulfite reactions in 200 μ l PCR tubes (not provided) according to Table 2. Add each component in the order listed.

Note: The combined volume of DNA and RNase-free water must total 20 μ l for high-concentration samples, and 40 μ l for low-concentration samples.

Table 2. Bisulfite reaction components

Component	High-concentration samples (1 ng – 2 μ g)	Low-concentration samples (1–500 ng)
	Volume per reaction (μ l)	Volume per reaction (μ l)
DNA	Variable* (maximum 20 μ l)	Variable† (maximum 40 μ l)
RNase-free water	Variable*	Variable†
Bisulfite Solution	85	85
DNA Protect Buffer	35	15
Total volume	140	140

* The combined volume of DNA and RNase-free water must total 20 μ l.

† The combined volume of DNA and RNase-free water must total 40 μ l.

3. Close the PCR tubes and mix the bisulfite reactions thoroughly. Store the tubes at room temperature (15–25°C).

Note: DNA Protect Buffer should turn from green to blue after addition to the DNA–Bisulfite Solution mixture, indicating sufficient mixing and correct pH for the bisulfite conversion reaction.

4. Perform the bisulfite DNA conversion using a thermal cycler. Program the thermal cycler according to Table 3.

The complete cycle should take approximately 30 min.

Optional: In some cases, it may be necessary to extend the 60°C cycle time up to 20 min to achieve complete bisulfite DNA conversion.

Note: If using a thermal cycler that does not allow you to enter the reaction volume (140 µl), set the instrument to the largest volume setting available.

Table 3. Bisulfite conversion thermal cycler conditions

Step	Time	Temperature
Denaturation	5 min	95°C
Incubation	10 min*	60°C
Denaturation	5 min	95°C
Incubation	10 min*	60°C
Hold	Indefinite†	20°C

* In some cases, it may be necessary to extend the 60°C cycle time up to 20 min to achieve complete bisulfite DNA conversion.

† Converted DNA can be left in the thermal cycler overnight without any loss of performance.

5. Place the PCR tubes containing the bisulfite reactions into the thermal cycler. Start the thermal cycling incubation.

Important: Since the bisulfite reaction is not overlaid with mineral oil, only thermal cyclers with heated lids are suitable for this procedure. It is important to use PCR tubes that close tightly.

Converted DNA can be left in the thermal cycler overnight without any loss of performance.

Clean up of bisulfite converted DNA

1. Once the bisulfite conversion is complete, briefly centrifuge the PCR tubes containing the bisulfite reactions, and then transfer the complete bisulfite reactions to clean 1.5 ml microcentrifuge tubes.
2. Transfer of precipitates in the solution will not affect the performance or yield of the reaction.
3. Add 310 μ l freshly prepared Buffer BL containing 10 μ g/ml carrier RNA (see "Preparation of Reagents") to each sample. Mix the solutions by vortexing and then centrifuge briefly.

Note: Carrier RNA is not necessary when using >100 ng DNA.

4. Add 250 μ l ethanol (96–100%) to each sample. Mix the solutions by pulse vortexing for 15 s, and centrifuge briefly to remove the drops from inside the lid.
5. Place the necessary number of MinElute DNA spin columns and collection tubes in a suitable rack. Transfer the entire mixture from each tube from step 8 into the corresponding MinElute DNA spin column.
6. Centrifuge the spin columns at maximum speed for 1 min. Discard the flow-through, and place the spin columns back into the collection tubes.
7. Add 500 μ l Buffer BW (wash buffer) to each spin column, and centrifuge at maximum speed for 1 min. Discard the flow-through, and place the spin columns back into the collection tubes.
8. Add 500 μ l Buffer BD (desulfonation buffer) to each spin column, and incubate for 15 min at room temperature (15–25°C).
9. If there are precipitates in Buffer BD, avoid transferring them to the spin columns.
Important: The bottle containing Buffer BD should be closed immediately after use to avoid acidification from carbon dioxide in air.
Note: It is important to close the lids of the spin columns before incubation.
10. Centrifuge the spin columns at maximum speed for 1 min. Discard the flow-through, and place the spin columns back into the collection tubes.

-
11. Add 500 μ l Buffer BW to each spin column and centrifuge at maximum speed for 1 min. Discard the flow-through and place the spin columns back into the collection tubes.
 12. Repeat step 14 once.
 13. Add 250 μ l ethanol (96–100%) to each spin column and centrifuge at maximum speed for 1 min.
 14. Place the spin columns into new 2 ml collection tubes, and centrifuge the spin columns at maximum speed for 1 min to remove any residual liquid.

Optional: Place the spin columns with open lids into a clean 1.5 ml microcentrifuge tube (not provided) and incubate the columns for 5 min at 60°C in a heating block. This step ensures the evaporation of any remaining liquid.

15. Place the spin columns into clean 1.5 ml microcentrifuge tubes (not provided). Add 15 μ l Buffer EB (elution buffer) directly onto the center of each spin-column membrane and close the lids gently.
16. Incubate the spin columns at room temperature for 1 min.
17. Centrifuge for 1 min at 15,000 \times g (12,000 rpm) to elute the DNA.

Note: As little as 10 μ l Buffer EB can be used for elution if a higher DNA concentration is required, but the yield will be reduced by approximately 20%. Do not elute with less than 10 μ l Buffer EB as the spin column membrane will not be sufficiently hydrated.

Note: We recommend storing purified DNA at 2–8°C for up to 24 h. When storing purified DNA for longer than 24 h, we recommend storage at –20°C.

Protocol: Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from FFPE Tissue Samples

This protocol is designed to be used with the EpiTect Fast FFPE Bisulfite Kit for processing DNA from formalin-fixed, paraffin-embedded (FFPE) tissue samples.

Important points before starting

- If using FFPE samples on slides, scrape the FFPE slice from the slide and proceed with step 1.
- Deparaffinization Solution solidifies at temperatures below 18°C. Incubate at 30°C to resolve.
- Precipitates may form in Lysis Buffer FTB. Make sure all precipitates are dissolved at 30°C.
- DNA Protect Buffer should turn from green to blue after addition to the DNA–Bisulfite Solution mixture, indicating sufficient mixing and correct pH for the bisulfite conversion reaction.
- Perform all centrifugation steps at room temperature (15–25°C).

Things to do before starting

- Prepare the kit reagents as described in “Preparation of Reagents”.
- Equilibrate samples and buffers to room temperature.
- Optional: If performing the deparaffinization, lysis, and decrosslinking of the FFPE slice in a 1.5 ml tube (see step 1), set a heating block to 56°C.
- Optional: Set a thermomixer, heating block, or heated orbital incubator to 60°C to dissolve the Bisulfite Solution.

Procedure

Deparaffinization, lysis, and decrosslinking of FFPE slice

1. Place a FFPE slice (10 μm) in a 200 μl reaction tube or 8-well strip (not provided) and add 150 μl Deparaffinization Solution.

Optional: Deparaffinization, lysis, and decrosslinking of the FFPE slice (steps 1–5) can be performed in a 1.5 ml tube (not provided).

2. Flick or vortex the tube until all paraffin is dissolved.
3. Add 20 μl distilled water, 15 μl Lysis Buffer FTB, and 5 μl proteinase K.

Note: A master mix comprising distilled water, Lysis Buffer FTB, and Proteinase K may be prepared in advance.

4. Vortex and briefly centrifuge the samples.

Note: The Deparaffinization Solution will form a layer above the Lysis Buffer FTB with the addition of Proteinase K.

5. Perform the lysis and decrosslinking using a thermal cycler. Program the thermal cycler according to Table 4.

Optional: If using 1.5 ml tubes, perform the lysis and decrosslinking in a thermal block. Incubate the tubes in a thermal block set to 56°C for 30 min to lyse the tissues. Ensure that tissues are completely lysed (the solution will become homogeneous); If not, incubate the tubes for an additional 30 min at 56°C. Once all tissues are lysed, increase the temperature of the heating block to 95°C for 60 min for the decrosslinking step.

Table 4. Lysis thermal cycling conditions

Step	Time	Temperature
Lysis	30 min*	56°C
Decrosslinking	60 min	95°C

* Ensure that the tissue is completely lysed; if not, add an additional lysis step (30 min at 56°C).

6. Place the PCR tubes containing the lysis reactions into the thermal cycler. Start the thermal cycling incubation.

Note: Samples should be kept at room temperature (15–25°C). Proceed as soon as possible with bisulfite conversion.

Bisulfite conversion of DNA

1. Make sure the Bisulfite Solution is completely dissolved.

Note: If necessary, heat the Bisulfite Solution to 60°C and vortex until all precipitates are dissolved again.

Note: Do not place dissolved Bisulfite Solution on ice.

2. Remove the Deparaffinization Solution (approximately 130 µl) from the lysis reactions.

3. Remove as much Deparaffinization Solution as possible without disturbing the lysed sample material to make sufficient space in the reaction tube for the bisulfite reaction components (Table 5).

Note: Small amounts of remaining Deparaffinization Solution have no effect on the bisulfite reaction.

4. Prepare the bisulfite reactions by adding the reagents listed in Table 5. Add each component in the order listed.

Optional: If using 1.5 ml tubes for the deparaffinization, lysis, and decrosslinking, transfer the lysis reactions remaining from step 8 to 200 µl tubes. Add the reagents to the 200 µl tubes in the order listed in Table 5.

Table 5. Bisulfite reaction components

Component	Volume per reaction (µl)
Lysis reaction	Approx. 40
Bisulfite Solution	85
DNA Protect Buffer	15
Total volume	140

-
5. Proceed to Protocol "Bisulfite Conversion of Unmethylated Cytosines in DNA", starting with step 3, to perform the bisulfite DNA conversion.

Protocol: Sample Lysis and Bisulfite Conversion of Unmethylated Cytosines in DNA Prepared from Whole Blood, Cultured Cells, or Tissue

This protocol is designed to be used with the EpiTect Fast LyseAll Kit for processing DNA from whole blood, cultured cells, or tissue samples.

Important points before starting

- Precipitates may form in Lysis Buffer FTB. Ensure that all precipitates are dissolved at 30°C.
- DNA Protect Buffer should turn from green to blue after addition to DNA–Bisulfite Solution mixture, indicating sufficient mixing and correct pH for the bisulfite conversion reaction.
- Perform all centrifugation steps at room temperature (15–25°C).

Things to do before starting

- Prepare the kit reagents as described in “Preparation of reagents”, page **Error! Bookmark not defined.**
- Equilibrate samples and buffers to room temperature.
- Optional: Set a thermomixer, heating block, or heated orbital incubator to 60°C to dissolve the Bisulfite Solution.

Procedure

Sample lysis

Note: If working with whole blood samples, follow steps 2–11; If working with cultured cells, follow steps 12–17; If working with tissue samples, follow steps 18–22.

Whole blood samples

1. Dilute the blood sample (maximum 20 µl) 1:20 with Buffer EL (e.g., add 19 µl Buffer EL to 1 µl blood).
2. Incubate at room temperature (15–25°C) for 10–15 min.
Note: Invert tubes several times during incubation.
3. Centrifuge at maximum speed for 5 min.
4. Discard supernatant and add an additional 125 µl Buffer EL.
Important: Make sure to not disturb the blood cell pellet while removing supernatant.
5. Centrifuge at maximum speed for 1 min.
6. Resuspend pellet in 10 µl PBS and transfer into a 200 µl reaction tube or into 8-well strips (not provided).
Important: Do not place cells on ice, as this will cause Lysis Buffer FTB to precipitate.
7. Add 10 µl distilled water, 15 µl Lysis Buffer FTB, and 5 µl Proteinase K.
Note: A master mix comprising of distilled water, Lysis Buffer FTB, and Proteinase K may be prepared in advance.
8. Vortex and briefly centrifuge the samples.
9. Incubate samples for 30 min at 56°C.
10. Proceed as soon as possible with Bisulfite conversion, step 23.
Note: Samples should be kept at room temperature.

Cultured cells

11. Harvest cells according to your current protocol.
12. Resuspend cells in 10 µl PBS.
Note: Do not use more than 1×10^5 cells per 10 µl PBS.
Important: Do not place cells on ice, as this will cause Lysis Buffer FTB to precipitate.
13. Add 10 µl distilled water, 15 µl Lysis Buffer FTB, and 5 µl Proteinase K.

Note: A master mix comprising distilled water, Lysis Buffer FTB, and proteinase K may be prepared in advance.

14. Vortex and briefly centrifuge the samples.
15. Incubate samples for 30 min at 56°C.
16. Proceed as soon as possible with bisulfite conversion, step 23.

Note: Samples should be kept at room temperature (15–25°C).

Tissue samples

17. Add 20 µl distilled water to the tissue sample.

Note: Do not use more than 100 µg tissue per 20 µl distilled water.

Important: Do not place tissues on ice, as this will cause Lysis Buffer FTB to precipitate.

18. Add 15 µl Lysis Buffer FTB and 5 µl proteinase K.

Note: A master mix comprising Lysis Buffer FTB and proteinase K may be prepared in advance.

19. Vortex and briefly centrifuge the samples.
20. Incubate samples for 30 min at 56°C.
21. Proceed as soon as possible with bisulfite conversion, step 23.

Note: Samples should be kept at room temperature (15–25°C).

Bisulfite conversion

22. Ensure that the Bisulfite Solution is completely dissolved.

Note: If necessary, heat the Bisulfite Solution to 60°C and vortex until all precipitates are dissolved again.

Note: Do not place dissolved Bisulfite Solution on ice.

23. Prepare the bisulfite reactions in 200 µl PCR tubes (not provided) according to Table 6. Add each component in the order listed.

Table 6. Bisulfite reaction components

Component	Volume per reaction (µl)
Lysis reaction	Approx. 40
Bisulfite Solution	85
DNA Protect Buffer	15
Total volume	140

24. Proceed to Protocol “Bisulfite Conversion of Unmethylated Cytosines in DNA”, starting with Step 3, to perform the bisulfite DNA conversion.

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit www.qiagen.com).

Comments and suggestions

Little or no DNA recovery in purification step

- | | |
|--|--|
| a) Carrier RNA not added to Buffer BL | Prepare carrier RNA and add to Buffer BL, as described in "Preparation of Reagents". |
| b) Buffer BW or Buffer BD prepared incorrectly | Check that Buffer BW and Buffer BD concentrates were diluted with the correct volumes of ethanol (96–100%). Do not use denatured alcohol, which contains other substances such as methanol or methylethylketone. |
| c) Buffer BW or Buffer BD prepared with 70% ethanol | Check that Buffer BW and Buffer BD concentrates were diluted with 96–100% ethanol. Do not use denatured alcohol, which contains other substances such as methanol or methylethylketone. |
| d) Buffer BW and Buffer BD used in the wrong order | Ensure that Buffer BW and Buffer BD are used in the correct order in the protocol. |
| e) Sample not completely passed through the membrane | Centrifuge for 1 min at full speed or until the entire sample has passed through the membrane. |
| f) Buffer BL contains precipitates | Check Buffer BL for precipitate. Dissolve by heating (maximum 70°C) with gentle agitation. |

Comments and suggestions

Low conversion rate

- | | | |
|----|--|---|
| a) | Bisulfite reaction components not added in the correct order | Ensure that DNA, Bisulfite Solution, and DNA Protect Buffer are added in the order indicated in Table 2, Table 5, or Table 6. |
| b) | Incorrect thermal cycling conditions used | Use the thermal cycling conditions given in Table 3. |
| c) | Poor DNA quality (i.e., protein contamination) | Check that the A_{260}/A_{280} ratio of the sample DNA is between 1.7 and 1.9. |
| d) | Amount of DNA used outside recommended range | Increase or decrease the amount of starting DNA material to stay within the range of 1 ng to 2 µg for purified DNA. For FFPE samples, cells, tissue, or blood, please refer to the indicated sample amounts on page Error! Bookmark not defined.. |
| e) | DNA Protect Buffer not added | Upon addition of DNA Protect Buffer, the DNA-Bisulfite Solution mixture should turn from green to blue indicating sufficient mixing and the correct pH for DNA binding to the MinElute DNA spin column. If this color change does not occur, repeat the reaction ensuring that DNA Protect Buffer has been added. |
| f) | Bisulfite Solution contains precipitates | Heat the Bisulfite Solution to 60°C and vortex until all precipitates are dissolved again. |
| g) | Insufficient time for bisulfite DNA conversion | Extend thermal cycling step of 60°C from 10 min to 20 min. |

Incomplete deparaffinization

- | | | |
|----|--|---|
| a) | Insufficient time for deparaffinization | Vortex the sample in Deparaffinization Solution until paraffin is visibly dissolved. |
| b) | Tissue slice not exposed to Deparaffinization Solution | Vortex and invert tube to ensure entire slice is uniformly exposed to Deparaffinization Solution. |

Incomplete lysis

- | | | |
|----|---|--|
| a) | FFPE tissue slice not deparaffinized | Ensure paraffin is completely dissolved before adding Lysis Buffer FTB and proteinase K. |
| b) | Sample not exposed to Lysis Buffer FTB and proteinase K | Ensure tissue is covered by lysis reagents (e.g., tissue should not be stuck to the tube cap). It may be necessary to fold the tissue into the solution using a clean pipet tip. |
| c) | Insufficient time for lysis | Incubate the sample at 56°C for an additional 30 min. |

Comments and suggestions

-
- | | | |
|----|-------------------------------------|--|
| d) | Lysis reaction prepared incorrectly | Make sure to add all necessary components of the lysis reaction, as described on page 27 (FFPE slices) or page 31 (whole blood), page 31 (cell culture), or page 32 (tissues). |
|----|-------------------------------------|--|

Poor results in downstream PCR applications

- | | | |
|----|--|---|
| a) | Little or no PCR product, even in control reaction | If performing hot-start PCR, confirm that the initial enzyme activation step was performed.
Ensure that all PCR components were added and that suitable cycling conditions were used. |
| b) | Failure of conversion reaction | The starting DNA was not sufficiently pure. Ensure that only high-quality DNA is used for the conversion reaction.
Ensure that all steps of the modification and cleanup protocol were followed.
Sample DNA was degraded before modification reaction. Ensure that sample DNA is handled and stored correctly.
PCR primers were not appropriate or incorrectly designed. Check primer design.
Amount of template DNA used in PCR was insufficient. Increase amount of template DNA. |

Unexpected findings in buffers

- | | | |
|----|--|--|
| a) | Color of DNA Protect Buffer changes from light green to olive during storage | DNA Protect Buffer is stable at 2–8°C for one year, and a change in color within this time has no influence on performance. |
| b) | Precipitates in Buffer BD | There may be slight clouding and/or insoluble precipitates in Buffer BD during storage.
Buffer BD is stable at 2–8°C for one year, and a precipitate within this time has no influence on performance. Precipitates should not be transferred onto the membrane of the spin column. |
| c) | Bisulfite Solution contains precipitates | Heat the Bisulfite Solution to 60°C and vortex until all precipitates are dissolved again. |

Ordering Information

Product	Contents	Cat. no.
EpiTect Fast DNA Bisulfite Kit (50)	For 50 preps: Bisulfite Solution, DNA Protect Buffer, MinElute DNA Spin Columns, Carrier RNA, and Buffers	59824
EpiTect Fast FFPE Bisulfite Kit (50)	For 50 preps: Deparaffinization Solution, Lysis Buffer, Proteinase K, Bisulfite Solution, DNA Protect Buffer, MinElute DNA Spin Columns, Carrier RNA, and Buffers	59844
EpiTect Fast LyseAll Bisulfite Kit (50)	For 50 preps: Lysis Buffer, Proteinase K, Bisulfite Solution, DNA Protect Buffer, MinElute DNA Spin Columns, Carrier RNA, and Buffers	59864
Related products		
EpiTect Fast 96 DNA Bisulfite Kit	2 x EpiTect 96-well Plates, Bisulfite Solution, DNA Protect Buffer, Carrier RNA, Buffers	59720
EpiTect Fast 96 FFPE Bisulfite Kit	Deparaffinization Solution, Lysis Buffer, Proteinase K, 2 x EpiTect 96-well Plates, Bisulfite Solution, DNA Protect Buffer, Carrier RNA, Buffers	59740
EpiTect Fast 96 LyseAll Bisulfite Kit	Lysis Buffers, Proteinase K, 2 x EpiTect 96-well Plates, Bisulfite Solution, DNA Protect Buffer, Carrier RNA, Buffers	59760
EpiTect Bisulfite Kit (48)	48 EpiTect Bisulfite Spin Columns, Reaction Mix, DNA Protect Buffer, Carrier RNA, Buffers	59104
EpiTect 96 Bisulfite Kit (2)	2x EpiTect Bisulfite 96-well Plates, Reaction Mix, DNA Protect Buffer, Carrier RNA, Buffers	59110
EpiTect Whole Bisulfite Kit — for amplification of bisulfite converted DNA		
EpiTect Whole Bisulfite Kit (25)	REPLI-g® Midi DNA Polymerase, EpiTect WBA Reaction Buffer, Nuclease free water	59203
EpiTect Whole Bisulfite Kit (100)	REPLI-g Midi DNA Polymerase, EpiTect WBA Reaction Buffer, Nuclease free water	59205

Product	Contents	Cat. no.
EpiTect MSP Kit — for highly accurate methylation-specific PCR without optimization		
EpiTect MSP PCR Kit (25)	EpiTect MSP Master Mix for 25 x 50 µl reactions	59303
EpiTect MSP PCR Kit (100)	EpiTect MSP Master Mix for 100 x 50 µl reactions	59305
EpiTect MSP PCR Kit (400)	EpiTect MSP Master Mix for 400 x 50 µl reactions	59307
EpiTect MethyLight PCR Kit — for real time quantification of methylation status		
EpiTect MethyLight PCR Kit (200)	Master mix for methylation-specific real-time PCR analysis, 200 x 50 µl reactions	59436
EpiTect MethyLight PCR Kit (1000)	Master mix for methylation-specific real-time PCR analysis, 1000 x 50 µl reactions	59438
EpiTect MethyLight PCR + ROX™ Vial Kit (200)	Master mix without ROX for methylation-specific real-time PCR analysis, 200 x 50 µl reactions	59496
EpiTect MethyLight PCR + ROX Vial Kit (1000)	Master mix without ROX for methylation-specific real-time PCR analysis, 1000 x 50 µl reactions	59498

For up-to-date licensing information and product specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Document Revision History

Date	Changes
06/2019	Deleted Epitect Fast FFPE Bisulfite Kit (200) column in Kit Contents section (product discontinued); Corrected minor spelling errors in the Epitect Fast Bisulfite Conversion Procedure figure; Removed References section; Layout updates.

Limited License Agreement for Epitect Fast Bisulfite Conversion Kits

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

Trademarks: QIAGEN®, MinElute®, EpiTect®, REPLI-g®, Pyrosequencing® (QIAGEN Group); ROX™ (Applera Corporation or its subsidiaries); Eppendorf® (Eppendorf AG); Sarstedt® (Sarstedt AG & Co.) Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

06/2019 HB-1211-003 © 2019 QIAGEN, all rights reserved.

