

## RNeasy<sup>®</sup> Mini Kit, Part 2

The RNeasy Mini Kit (cat. nos. 74104 and 74106) can be stored at room temperature (15–25°C) for at least 9 months if not otherwise stated on label.

### Further information

- *RNeasy Mini Handbook*: [www.qiagen.com/HB-0435](http://www.qiagen.com/HB-0435)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

#### On-column DNase digestion

- If using the RNase-Free DNase Set for the first time, prepare DNase I stock solution by injecting 550 µl RNase-free water into the DNase I vial using an RNase-free needle and syringe. Mix gently by inverting the vial. Do not vortex.
  - For long-term storage of DNase I stock solution, divide it into single-use aliquots and store at –20°C for up to 9 months. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze aliquots after thawing.
1. Add 350 µl Buffer RW1 to RNeasy column, close lid, centrifuge for 15 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm). Discard flow-through.
  2. Add 10 µl DNase I stock solution (see above) to 70 µl Buffer RDD. Mix by gently inverting the tube. Centrifuge briefly.
  3. Add DNase I incubation mix (80 µl) directly to RNeasy column membrane, and place on benchtop (20–30°C) for 15 min.
  4. Add 350 µl Buffer RW1 to RNeasy column, close lid, centrifuge for 15 s at  $\geq 8000 \times g$ . Discard flow-through. Continue with step 5 of “RNA purification from cells/tissue samples” in *Quick-Start Protocol RNeasy Mini Kit, Part 1*, or step 4 of “RNA cleanup” (below).

## Notes before starting

### RNA cleanup

- Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
1. Adjust the sample to a volume of 100  $\mu$ l with RNase-free water. Alternatively, follow steps in “DNase digestion of RNA before RNA cleanup” in Appendix E of *RNeasy Mini Handbook*. Add 350  $\mu$ l Buffer RLT, and mix well.
  2. Add 250  $\mu$ l ethanol (96–100%) to the diluted RNA, and mix well by pipetting. Do not centrifuge. Proceed immediately to step 3.
  3. Transfer the sample (700  $\mu$ l) to an RNeasy Mini spin column placed in a 2 ml collection tube (supplied). Close the lid. Centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.  
**Optional:** If performing optional on-column DNase digestion, follow steps 1–4 of “On column DNase digestion” (above) after this step.
  4. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid. Centrifuge for 15 s at  $\geq 8000 \times g$  to wash the membrane. Discard the flow-through.
  5. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid. Centrifuge for 2 min at  $\geq 8000 \times g$  to wash the membrane.  
**Optional:** Place the RNeasy spin column in a new 2 ml collection tube (supplied). Close the lid, and centrifuge at full speed for 1 min.
  6. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50  $\mu$ l RNase-free water directly to the spin column membrane. Close the lid, and centrifuge for 1 min at  $\geq 8000 \times g$  to elute the RNA.
  7. If the expected RNA yield is  $>30 \mu\text{g}$ , repeat step 6 using another 30–50  $\mu$ l of RNase-free water. Alternatively, use the eluate from step 6 (if high RNA concentration is required). Reuse the collection tube from step 6.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, RNeasy®, Sample to Insight® (QIAGEN Group). 1102231 04/2016 HB-0571-002 © 2016 QIAGEN, all rights reserved.