



## QIAGEN Supplementary Protocol:

### Fast-forward protocol for transient transfection of NIH/3T3 cells in 96-well plates using PolyFect® Transfection Reagent

The following protocol is optimized for transient transfection of NIH/3T3 cells in 96-well plates without pre-plating of cells 24 hours prior to transfection. Cell plating and transfection are performed on the same day, making this protocol rapid and convenient. Two possibilities for transfection-complex formation (in tubes or in the wells of a 96-well plate) are provided in protocol step 2. Please read the protocol thoroughly before beginning this procedure.

**IMPORTANT:** Please consult the “General Guidelines” section in the *PolyFect Transfection Reagent Handbook* before beginning this procedure.

#### Important note before starting

- To ensure optimal results, we strongly recommend using the optimized amounts of DNA and PolyFect Reagent given in the protocol below. **The amounts given are for one well of a 96-well plate.**

#### Procedure

- 1. Dilute 0.15 µg DNA dissolved in TE buffer pH 7 to pH 8 (minimum DNA concentration: 0.1 µg/µl) with medium containing no serum or antibiotics to a total volume of 30 µl per well. Mix, then centrifuge for a few seconds to remove any liquid from the top of the tube.**  
**IMPORTANT:** Serum and antibiotics present during this step will interfere with transfection-complex formation and will significantly decrease transfection efficiency.
- 2. Dilute 1 µl PolyFect Reagent with medium containing no serum or antibiotics to a total volume of 20 µl per well. Add the diluted PolyFect Reagent to the DNA solution. Mix, then centrifuge for a few seconds to remove any liquid from the top of the tube. Alternatively, pipet the diluted DNA (step 1) and diluted PolyFect Reagent into one well of a 96-well plate. Mix by pipetting up and down 5 times.**  
**IMPORTANT:** Serum and antibiotics present during this step will interfere with transfection-complex formation and will significantly decrease transfection efficiency.  
**Note:** It is not necessary to keep PolyFect Reagent on ice at all times. 10–15 min at room temperature will not alter its stability.
- 3. Incubate the samples for 5–10 min at room temperature (15–25°C) to allow transfection-complex formation. Continue with steps 4 and 5 during this incubation.**  
**Note:** Transfection-complex formation takes a minimum of 5–10 min. The transfection complexes will remain stable during the time it takes to prepare the cells for transfection; however, avoid extending this incubation for too long.
- 4. Harvest the cells by trypsinization and suspend in growth medium (containing serum and antibiotics).**  
**Note:** The cells should be healthy and in logarithmic growth phase.

5. **Count the harvested cell suspension and adjust the cell density to 3.3–4.0 x 10<sup>5</sup> cells/ml.**
6. **If transfection-complex formation was not performed directly in a 96-well plate (step 2), pipet 50 µl of the solution containing the transfection complexes into one well of a 96-well plate.**
7. **Add 150 µl of the cell suspension (5–6 x 10<sup>4</sup> cells) to wells containing transfection complexes. Mix by pipetting up and down twice.**

At this stage, the serum and antibiotics present in the growth medium will not interfere with, but rather significantly enhance, the transfection efficiency of PolyFect Reagent.

8. **Incubate cells with the transfection complexes at 37°C and 5% CO<sub>2</sub>. Assay cells for expression of the transfected gene after an appropriate incubation time.**

For example, cells transfected with β-gal or cat reporter constructs are typically incubated for 24–48 h after transfection to obtain maximal levels of gene expression.

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor.

Selected handbooks can be downloaded from [www.qiagen.com/literature/default.aspx](http://www.qiagen.com/literature/default.aspx).

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from [www.qiagen.com/ts/msds.asp](http://www.qiagen.com/ts/msds.asp).

Trademarks: QIAGEN®, PolyFect® (QIAGEN Group).

© 2006–2007 QIAGEN, all rights reserved.