



QIAGEN Supplementary Protocol:

Fast-forward protocol for transient transfection of HeLa cells in 96-well plates using PolyFect® Transfection Reagent

The following protocol is optimized for transient transfection of HeLa cells in 96-well plates without pre-plating of cells 24 hours prior to transfection. Cell plating and transfection are performed on the same day, making this protocol rapid and convenient. Two possibilities for transfection-complex formation (in tubes or in the wells of a 96-well plate) are provided in protocol step 2. Please read the protocol thoroughly before beginning this procedure.

Please note that a separate protocol is available for HeLa-S3 cells.

IMPORTANT: Please consult the “General Guidelines” section in the *PolyFect Transfection Reagent Handbook* before beginning this procedure.

Important note before starting

- To ensure optimal results, we strongly recommend using the optimized amounts of DNA and PolyFect Reagent given in the protocol below. **The amounts given are for one well of a 96-well plate.**

Procedure

1. **Dilute 0.4 µg DNA dissolved in TE buffer pH 7 to pH 8 (minimum DNA concentration: 0.1 µg/µl) with medium containing no serum or antibiotics to a total volume of 30 µl per well. Mix, then centrifuge for a few seconds to remove any liquid from the top of the tube.**

IMPORTANT: Serum and antibiotics present during this step will interfere with transfection-complex formation and will significantly decrease transfection efficiency.

2. **Dilute 2 µl PolyFect Reagent with medium containing no serum or antibiotics to a total volume of 20 µl per well. Add the diluted PolyFect Reagent to the DNA solution. Mix, then centrifuge for a few seconds to remove any liquid from the top of the tube.**

Alternatively, pipet the diluted DNA (step 1) and diluted PolyFect Reagent into one well of a 96-well plate. Mix by pipetting up and down 5 times.

IMPORTANT: Serum and antibiotics present during this step will interfere with transfection-complex formation and will significantly decrease transfection efficiency.

Note: It is not necessary to keep PolyFect Reagent on ice at all times. 10–15 min at room temperature will not alter its stability.

3. **Incubate the samples for 5–10 min at room temperature (15–25°C) to allow transfection-complex formation. Continue with steps 4 and 5 during this incubation.**

Note: Transfection-complex formation takes a minimum of 5–10 min. The transfection complexes will remain stable during the time it takes to prepare the cells for transfection; however, avoid extending this incubation for too long.

- 4. Harvest the cells by trypsinization and suspend in growth medium (containing serum and antibiotics).**

Note: The cells should be healthy and in logarithmic growth phase.

- 5. Count the harvested cell suspension and adjust the cell density to $3.3\text{--}4.0 \times 10^5$ cells/ml.**

- 6. If transfection-complex formation was not performed directly in a 96-well plate (step 2), pipet $50 \mu\text{l}$ of the solution containing the transfection complexes into one well of a 96-well plate.**

- 7. Add $150 \mu\text{l}$ of the cell suspension ($5\text{--}6 \times 10^4$ cells) to wells containing transfection complexes. Mix by pipetting up and down twice.**

At this stage, the serum and antibiotics present in the growth medium will not interfere with, but rather significantly enhance, the transfection efficiency of PolyFect Reagent.

- 8. Incubate cells with the transfection complexes for 2–3 h at 37°C and 5% CO_2 .**

Note: Cells should have adhered to the plate by the end of this incubation.

- 9. Remove medium containing the transfection complexes from the cells by gentle aspiration, and wash cells once with $150 \mu\text{l}$ PBS.**

Note: For optimal results, we strongly recommend removing the transfection complexes when using the fast-forward transfection protocol with HeLa cells. If absolutely necessary for a particular application, it may be possible to omit this step by using less DNA and PolyFect Reagent for transfection. As a starting point, we suggest using $0.10 \mu\text{g}$ DNA and $0.75 \mu\text{l}$ PolyFect Reagent; however these amounts may need to be optimized to maximize transfection efficiency and minimize cytotoxicity.

- 10. Add fresh growth medium (containing serum and antibiotics) and incubate cells at 37°C and 5% CO_2 . Assay cells for expression of the transfected gene after an appropriate incubation time.**

For example, cells transfected with β -gal or cat reporter constructs are typically incubated for 24–48 h after transfection to obtain maximal levels of gene expression.

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor.

Selected handbooks can be downloaded from www.qiagen.com/literature/default.aspx.

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from www.qiagen.com/ts/msds.asp.

Trademarks: QIAGEN®, PolyFect® (QIAGEN Group).

© 2006–2007 QIAGEN, all rights reserved.