



How to Collect Blood Using the PAXgene® Blood RNA Tube

For Molecular Diagnostic Testing




Required Items:

- 1a.** Ensure that the PAXgene Blood RNA Tube (A) is at room temperature (18°C–25°C) prior to use and is properly labeled with patient identification.
- 1b.** If the PAXgene Blood RNA Tube is the only tube to be drawn, a small amount of blood should be drawn into a “Discard Tube” (D) prior to drawing blood into the PAXgene Blood RNA Tube. Otherwise, the PAXgene Blood RNA Tube should be the last tube drawn in the phlebotomy procedure.




Venipuncture:

- 2.** Using a blood collection set such as the BD Vacutainer UltraTouch Push Button or BD Vacutainer Safety-Lok Blood Collection Set (B) and the One-Use Holder (C), perform venipuncture using your institution’s recommended procedure for standard venipuncture.



Blood Collection:

- 3a.** Hold the PAXgene Blood RNA Tube vertically, below the blood donor’s arm, during blood collection.
- 3b.** Allow at least 10 seconds for a complete blood draw to take place. Ensure that the blood has stopped flowing into the tube before removing the tube from the holder.



After Blood Collection:

- 4a.** Gently invert the PAXgene Blood RNA Tube 8 to 10 times.
- 4b.** Store the PAXgene Blood RNA Tube upright at room temperature (18°C–25°C) for a minimum of 2 hours and a maximum of 72 hours before processing or transferring to refrigerator (2–8°C) or freezer (–20°C).*

Before Blood Collection



After Blood Collection



RNA Stabilization Reagent

RNA Stabilized Whole Blood

Ref #	Product Description
762165	PAXgene® Blood RNA Tube
367365	BD Vacutainer® UltraTouch™ Push Button Blood Collection Set
367281 (US) 367286 (CE)	BD Vacutainer Safety-Lok™ Blood Collection Set
364815	BD Vacutainer One-Use Holder
367812 (US) 368975 (CE)	BD Vacutainer Plus Plastic Serum Tube (Discard Tube)

BD Customer Service / Orders:
888.237.2762 (North America)

www.PreAnalytiX.com



*Refer to PAXgene Blood RNA Tube product circular.
 PAXgene and PreAnalytiX are registered trademarks of PreAnalytiX, Hombrechtikon, CH ©2016 PreAnalytiX GmbH
 All other trademarks are property of Becton, Dickinson and Company.
 06/16 VS9339

PreAnalytiX GmbH, Feldbachstrasse, 8634 Hombrechtikon, CH



A QIAGEN / BD Company