

RNeasy[®] Mini Kit, Part 1

The RNeasy Mini Kit (cat. nos. 74104 and 74106) can be stored at room temperature (15–25°C) for at least 9 months if not otherwise stated on label.

Further information

- *RNeasy Mini Handbook*: www.qiagen.com/HB-0435
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- If purifying RNA from cell lines rich in RNases, or tissue, add either 10 µl β-mercaptoethanol (β-ME), or 20 µl 2 M dithiothreitol (DTT), to 1 ml Buffer RLT. Buffer RLT with β-ME or DTT can be stored at room temperature for up to 1 month.
 - Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
 - Remove RNAprotect[®] stabilized tissue from the reagent using forceps.
 - For RNeasy Protect Mini Kit (cat. nos. 74124 and 74126), please start with the *Quick-Start Protocol RNAprotect Tissue Reagent, RNAprotect Tissue Tubes, and RNeasy Protect Kits*.
1. **Cells:** Harvest a maximum of 1×10^7 cells, as a cell pellet or by direct lysis in the vessel. Add the appropriate volume of Buffer RLT and select a suitable method for disruption and homogenization (see Table 1).
Tissues: Do not use more than 30 mg tissue. Disrupt the tissue and homogenize the lysate in the appropriate volume of Buffer RLT (see Table 1). Centrifuge the lysate for 3 min at maximum speed. Carefully remove the supernatant by pipetting, and use it in step 2.
 2. Add 1 volume of 70% ethanol to the lysate, and mix well by pipetting. Do not centrifuge. Proceed immediately to step 3.
 3. Transfer up to 700 µl of the sample, including any precipitate, to an RNeasy Mini spin column placed in a 2 ml collection tube (supplied). Close the lid, and centrifuge for 15 s at $\geq 8000 \times g$. Discard the flow-through.

Optional: For DNase digestion, follow steps 1–4 of “On column DNase digestion” in Quick-Start Protocol RNeasy Mini Kit, Part 2.

4. Add 700 μ l Buffer RW1 to the RNeasy spin column. Close the lid, and centrifuge for 15 s at $\geq 8000 \times g$. Discard the flow-through.
5. Add 500 μ l Buffer RPE to the RNeasy spin column. Close the lid, and centrifuge for 15 s at $\geq 8000 \times g$. Discard the flow-through.
6. Add 500 μ l Buffer RPE to the RNeasy spin column. Close the lid, and centrifuge for 2 min at $\geq 8000 \times g$.

Optional: Place the RNeasy spin column in a new 2 ml collection tube (supplied). Centrifuge at full speed for 1 min to dry the membrane.

7. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50 μ l RNase-free water directly to the spin column membrane. Close the lid, and centrifuge for 1 min at $\geq 8000 \times g$ to elute the RNA.
8. If the expected RNA yield is $>30 \mu$ g, repeat step 7 using another 30–50 μ l of RNase-free water, or using the eluate from step 7 (if high RNA concentration is required). Reuse the collection tube from step 7.

Table 1. Volumes of Buffer RLT for sample disruption and homogenization

Sample	Amount	Dish	Buffer RLT (μ l)	Disruption and homogenization
Animal cells	$<5 \times 10^6$	<6 cm	350	Add Buffer RLT, vortex ($\leq 1 \times 10^5$ cells); or use QIAshredder, TissueRuptor [®] , or needle and syringe
	$<1 \times 10^7$	6–10 cm	600	
Animal tissues	<20 mg	–	350	TissueLyser LT; TissueLyser II; TissueRuptor, or mortar and pestle followed by QIAshredder or needle and syringe
	≤ 30 mg	–	600	

Document Revision History

Date	Changes
11/2021	Removed reference to RNAlater.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN[®], Sample to Insight[®], RNAprotect[®], RNeasy[®], TissueRuptor[®] (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1126491 11/2021 HB-0570-003 © 2021 QIAGEN, all rights reserved.