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“RNAlater[®]” is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents.

Trademarks: QIAGEN[®], QIAamp[®], AllPrep[®], BioRobot[®], BioSprint[®], DNeasy[®], EpiTect[®], FastLane[®], GeneGlobe[™], HotStarTaq[®], QuantiFast[™], QuantiTect[®], RNAprotect[®], RNeasy[®], TissueRuptor[™] (QIAGEN Group); ABI PRISM[®], Applied Biosystems[®] (Applied Biosystems Corporation or its subsidiaries); PAXgene[™] (PreAnalytiX GmbH); Pyrosequencing[®] (Biotope AB); SYBR[®] (Molecular Probes, Inc.); TaqMan[®] (Roche Group).

The QuantiTect Whole Transcriptome Kit is for use only as licensed by Amersham Biosciences Corp (part of GE Healthcare Bio-Sciences) and QIAGEN GmbH. The Phi 29 DNA polymerase may not be re-sold or used except in conjunction with the other components of this kit. See U.S. Patent Nos. 5,854,033, 6,124,120, 6,143,495, 5,001,050, 5,198,543, 5,576,204, and related U.S. and foreign patents.

Purchase of this product (QuantiFast SYBR Green Kits, QuantiTect SYBR Green Kits, and FastLane Cell SYBR Green Kit) is accompanied by a limited, non-transferable immunity from suit to use it with detection by a dsDNA-binding dye as described in U.S. Patents Nos. 5,994,056 and 6,171,785 and corresponding patent claims outside the United States for the purchaser's own internal research. No real-time apparatus or system patent rights or any other patent rights, and no right to use this product for any other purpose are conveyed expressly, by implication or by estoppel.

This product (QuantiFast Probe PCR Kits, QuantiTect Probe Kits, QuantiTect Multiplex Kits, FastLane Cell Probe Kit, and FastLane Cell Multiplex Kits) is an Authorized 5' Nuclease Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under certain patents owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research when used in conjunction with Licensed Probe. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Purchase of this product (QuantiFast SYBR Green RT-PCR Kit, QuantiTect SYBR Green RT-PCR Kit, QuantiTect Probe RT-PCR Kit, QuantiTect Multiplex RT-PCR Kits, FastLane Cell SYBR Green Kit, FastLane Cell Probe Kit, and FastLane Cell Multiplex Kits) is accompanied by a limited, non-transferable license under RT and Reverse Transcription-PCR patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd to use it for the purchaser's own internal research. No real-time patent rights of any kind, no right under any other patent claims (such as apparatus or system claims), and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel.

This product (QuantiTect Primer Assays) is compatible for use in the 5' nuclease process or the dsDNA-binding dye processes covered by patents owned by Roche or owned by or licensed to Applied Biosystems. No license under these patents to practice the 5' nuclease process or the dsDNA-binding dye processes are conveyed expressly or by implication to the purchaser by the purchase of this product.

siRNA technology licensed to QIAGEN is covered by various patent applications, owned by the Massachusetts Institute of Technology, the Carnegie Institute of Washington, Alnylam Corporation, and others.

Use of methylation specific PCR (MSP) is covered by US patents 5,786,146, 6,017,704, 6,200,756, 6,265,171, and corresponding foreign patents and applications. No license under these patents to use the MSP process is conveyed to the purchaser by purchasing this product.

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弊社顧客情報の取り扱いに関しては、www.qiagen.com/goto/PrivacyPolicy をご覧ください。 © 2007 QIAGEN, all rights reserved.