

## Template File Info Note

# Investigator<sup>®</sup> Template Files for GeneMapper<sup>®</sup> ID-X v11-Standard Investigator and GeneMapper ID-X v11-CODIS

## Table of Contents

- Overview
- Template Files v11-Standard and Template Files v11-CODIS
- Template Files Content
- File Install Instructions
- Legal Disclaimer

## Overview

The Investigator Template Files are software sets for the GeneMapper ID-X Software that help simplify data analysis. Investigator Template Files may be used with ABI PRISM™ single- and multi-capillary instruments from Applied Biosystems®.

GeneMapper ID-X Software (used in conjunction with Investigator Template Files) assigns the analyzed DNA fragments, relative to their length, to the allele designation of the short tandem repeat (STR) loci. Optionally, the corresponding fragment length of peaks in base pairs (bp) or peak height in relative fluorescent units (RFU) can be indicated. These data (genotypes) can be tabulated and exported.

---

Investigator Template Files are available for all Investigator Human Identification PCR Kits for GeneMapper ID-X Software version 1.2 or higher.

**Note:** Panels and BinSets are essential for the use of Investigator Human Identification PCR Kits. The use of additional Investigator Template Files is optional.

## Template Files v11-Standard and Template Files v11-CODIS

QIAGEN offers 2 different versions of the v11 Investigator Template Files for GeneMapper ID-X: the Investigator Template Files for GeneMapper ID-X v11-Standard and v11-CODIS.

The Investigator Template Files v11-CODIS include a Panel-File, which contains a different spelling of the Amelogenin marker (i.e., "AMEL" instead of "AM" for v11-Standard). In addition, v11-CODIS contains a further analysis method (i.e., Analysis\_HID\_200rfu\_CODIS) with a threshold of 200 RFU and a different global stutter filter.

## Template Files Content

The Investigator Template Files consist of the following files:

### Size standards

- SST-BTO\_60-500bp
- SST-ROX\_50-500bp
- SST-BTO\_60-450bp

### Panels, BinSets, and Stutter

- Panel files are specific for each kit (e.g., 24plex\_Panels\_v4)
- Bin files are specific for each kit (e.g., 24plex\_Bins\_v4)
- Stutter files are specific for each kit (e.g., 24plex\_Stutter\_v4)

---

## Analysis methods

- Analysis\_HID\_310\_50rfu
- Analysis\_HID\_310\_200rfu
- Analysis\_HID\_3130\_50rfu
- Analysis\_HID\_3130\_200rfu
- Analysis\_HID\_3500\_50rfu
- Analysis\_HID\_3500\_200rfu
- Analysis\_HID\_200rfu\_CODIS

(only in Template File GeneMapper ID-X v11-CODIS)

## Plot settings

- Plots\_2dyes
- Plots\_3dyes
- Plots\_4dyes
- PlotsBT5\_4dyes
- Plots\_5dyes
- Plots\_6dyes

## File Install Instructions

Investigator Template Files for GeneMapper ID-X Software are available as a free download from the online catalog pages of all Investigator Human Identification PCR Kits ([www.qiagen.com/InvestigatorIDKits](http://www.qiagen.com/InvestigatorIDKits)).

Before using the GeneMapper ID-X Software for the first time, Investigator Template Files must be saved to the local analysis computer. Panels, BinSets, and Stutter are imported by the Panel Manager. Other template files, such as Analysis Methods, Table Settings, Plot Settings, etc., can be imported to the local PC using the GeneMapper ID-X Manager.

For detailed information, follow the installation instructions contained in the user manual of the Investigator PCR Kit you are using.

## Document Revision History

Date	Changes
10/2019	Replaced v9 with v10 for both Standard and CODIS. Replaced "GeneMapper ID-X Software version 1.0 or higher" with "GeneMapper ID-X Software version 1.2 or higher". Added SST-BTO_60-450bp to "Size standards".
05/2021	Replaced v10 with v11 for both Standard and CODIS. Updates to the Template Files are as follows: (1) Bins, Panels, and Stutter files were updated with Investigator Argus Y-28 QS information, and (2) analysis files for 3500 Genetic Analyzers are changed to include hexanucleotide STR definition for global stutter function.

### Software License Agreement and Warranty Statement

**IMPORTANT:** PLEASE READ THIS LICENSE AGREEMENT CAREFULLY. QIAGEN LICENSES THIS SOFTWARE TO YOU ("LICENSEE") ONLY UPON THE ACCEPTANCE OF ALL OF THE TERMS IN THIS LICENSE AGREEMENT. BY USING THIS SOFTWARE PACKAGE, LICENSEE CONSENTS TO BE BOUND BY THE TERMS OF THIS LICENSE.

### Software License

The software which accompanies this license (the "Software") is the property of QIAGEN and is protected by copyright law. The software under this license is licensed, not sold, to Licensee by QIAGEN upon the conditions contained in this License. QIAGEN retains full title to the Software. This license is valid for the Investigator Template Files for GeneMapper ID-X.

The Software may be used solely in accordance with the *Investigator Template Files for GeneMapper ID-X User Guide* and for use with components contained in the Software only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Software with any components not included within this Software except as described in the *Investigator Files for GeneMapper ID-X User Guide* and additional protocols available at [www.qiagen.com](http://www.qiagen.com).

Trademarks: QIAGEN®, Sample to Insight®, Investigator® (QIAGEN Group); Applied Biosystems®, ABI PRISM™, GeneMapper® (Thermo Fisher Scientific or its subsidiaries). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

05/2021 HB-1883-006 © 2021 QIAGEN, all rights reserved.