

Checklist for multiplex, real-time PCR

Before carrying out multiplex, real-time PCR with a QuantiTect® Multiplex PCR or RT-PCR Kit, please go through the checklist below. For detailed information, refer to the kit handbook or *Critical factors for success in real-time, multiplex PCR* or contact **QIAGEN Technical Service** (contact information on reverse side).

Design primers and probes

Has a primer–probe set (e.g., TaqMan® Gene Expression Assay or published in literature) already been designed for each of your target genes?

- If yes, there is no need to design new primer–probe sets (however, it may be necessary to choose new reporter dyes; see next step).
- If no, use Primer Express® or other software to design primers and probes. Use BLAST® software (www.ncbi.nlm.nih.gov/BLAST) to check that the sequences are unique.

Choose reporter dyes and quenchers

Refer to Table 1 (see reverse side) to choose the appropriate combination of reporter dyes for your probes. Check that the appropriate filters are installed on your real-time cycler. We recommend using a nonfluorescent quencher (e.g., Black Hole Quencher® dye) for all probes.

Reconstitute primers and probes

Use TE to prepare a 100 µM stock solution of each lyophilized primer or probe. Check concentrations by spectrophotometry, and store frozen in small aliquots away from light.

Set up your real-time cycler

Check whether your real-time cycler needs to be calibrated for each of the reporter dyes to be used in the multiplex assay. Be sure to activate the filter or detector for each reporter dye.

Be sure to adjust the analysis settings for each reporter dye channel in each run. Be aware that default analysis settings may not be optimal.

Evaluate the performance of the multiplex assay

Check that each primer–probe set works in single PCR. Then compare the performance of the primer–probe sets in single PCR and multiplex PCR: using serial 10-fold dilutions of template, multiplex PCR should give C_T values comparable to those obtained in single PCR.



Table 1. Combinations of Reporter Dyes for Multiplex Assays*

Cycler	Reference dye	Dye 1 [†]	Dye 2 [†]	Dye 3 [†]	Dye 4 [†]
ABI PRISM [®] 7700	ROX	FAM	HEX, JOE, VIC [®]	—	—
ABI PRISM 7000 and 7900, Applied Biosystems [®] 7300	ROX	FAM	HEX, JOE, VIC	Bodipy [®] TMR, NED	—
Applied Biosystems 7500	ROX	FAM	HEX, JOE, VIC	Bodipy TMR, NED	Alexa Fluor [®] 647, Cy [®] 5
iCycler iQ [®] and iQ5	Not required	FAM	HEX, JOE, TET, VIC	Texas Red [®] , ROX	Cy5
LightCycler [®] 2.0	Not required	FAM	HEX, JOE, VIC	Texas Red, ROX	Alexa Fluor 660, Bodipy 630/650, Pulsar [®] 650
Mx3000P™, Mx3005P™	Not required	FAM	HEX, JOE, VIC	Texas Red, ROX	Cy5
Rotor-Gene™ 3000	Not required	FAM	HEX, JOE, VIC	Texas Red, ROX	Cy5

* Visit www.qiagen.com/multiplex to view dye combinations for other cyclers.

[†] Preferably, select Dye 1 for the least abundant target, Dye 2 for the second least abundant target, and Dyes 3–4 for the most abundant targets.

Emission Maximum of Selected Reporter Dyes




Figure 1 The emission maximum (nm) of selected reporter dyes are displayed in parentheses. Emission maximum may vary depending on buffer conditions. Other dyes with similar wavelengths may not be suitable for multiplex assays due to low fluorescence and/or stability.

See multiplex PCR application data from researchers at www.qiagen.com/multiplex .

Trademarks: QIAGEN[®], QuantiTect[®] (QIAGEN Group); Alexa Fluor[®], Bodipy[®], Texas Red[®] (Molecular Probes, Inc.); ABI PRISM[®], Applied Biosystems[®], Primer Express[®], VIC[®] (Applied Biosystems Corporation or its subsidiaries); Black Hole Quencher[®], Pulsar[®] (Biosearch Technologies, Inc.); BLAST[®] (US National Library of Medicine); Cy[®] (General Electric Company); iCycler iQ[®] (Bio-Rad Laboratories, Inc.); LightCycler[®], TaqMan[®] (Roche Group); Mx3000P™, Mx3005P™ (Stratagene); Rotor-Gene™ (Corbett Research).

Each QuantiTect Multiplex Kit is an Authorized 5' Nuclease Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under certain patents owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research when used in conjunction with Licensed Probe. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA. Purchase of a QuantiTect Multiplex RT-PCR Kit is accompanied by a limited, non-transferable license under RT and Reverse Transcription-PCR patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd to use it for the purchaser's own internal research. No real-time patent rights of any kind, no right under any other patent claims (such as apparatus or system claims), and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Certain specific embodiments of the process of multiplex PCR may be covered by patents of third parties in certain countries and may require a license.

1038604 06/2006 © 2006 QIAGEN, all rights reserved.

Australia ■ 1-800-243-066
Austria ■ 0800 281011
Belgium ■ 0800 79556
Canada ■ 800-362-7737
China ■ 021-51345678
Denmark ■ 80-885942

Finland ■ 0800-914413
France ■ 01-60-920-930
Germany ■ 02103-29-12400
Ireland ■ 1800 555 061
Italy ■ 800 787980
Luxembourg ■ 8002 2067

The Netherlands ■ 0800 0229602
Norway ■ 800-18712
Sweden ■ 020-798328
Switzerland ■ 055-254-22-12
UK ■ 01293-422-999
USA ■ 800-362-7737

