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
Sample & Assay Technologies

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High-Throughput Reverse Transfection



图 1 原代细胞获得 40% 的转染效率。Primary rabbit aortic smooth muscle cells were transfected with 0.4 μg of a green fluorescent protein reporter plasmid using Effectene Transfection Reagent in 6-well plates. Cells were viewed by fluorescence microscopy 24 h after transfection. Approximately 40% of the cells were transfected as determined by FACS[®] analysis. Data kindly provided by K. Veit, 2nd Medical Clinic, Dept. Clinical Pharmacology, Mainz, Germany.


图 2 Effectene 试剂效果优于其它脂质体试剂。Murine teratocarcinoma F9 cells (5×10^5) were transfected in 6-well plates with a luciferase reporter plasmid using Effectene Reagent from QIAGEN or Reagent L or Reagent F from other suppliers, according to the manufacturers' instructions. Transfection efficiencies were determined by measuring luciferase activity 48 h after transfection and are given as relative light units (RLU). Data kindly provided by I. Clavereau, D. Petitprez, and I. Van Seuning, Unité INSERM 377, Place de Verdun, Lille Cedex, France.

Attractene Transfection Reagent NEW

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
图 1 Attractene 试剂转染效率优于其它试剂。 DNA (pGFP which expresses green fluorescent protein) was transfected into the cell types indicated using Attractene Reagent according to the conditions recommended in the handbook. Transfection was also performed with **A** Reagent F using 2 reagent volumes recommended by the manufacturer or **B** Reagent L according to the manufacturer's recommendations, with and without the enhancer provided. Transfection efficiency was estimated by counting the number of fluorescent cells by FACS[®] analysis. Transfection efficiency is expressed relative to efficiency using Attractene Reagent which was set at 100%.


图 2 Attractene 试剂毒性极低。 HepG2 cells were transfected with DNA (pGFP) using Attractene Transfection Reagent or Reagent L from another supplier. **A** Cells transfected using Attractene Reagent were healthy and viable. **B** In contrast, cells transfected using Reagent L displayed high levels of cell death.


图 3 共转染 shRNA 载体后有效沉默。 Using Attractene Reagent, HEK293 cells were **A** cotransfected with pGFP and a plasmid expressing an shRNA targeted against the green fluorescent protein gene (pGFP + pshGFP), or **B** cotransfected with pGFP and a negative control plasmid expressing a scrambled shRNA (pGFP + pNegative). After cotransfection of pGFP and the shRNA vector pshGFP, the green fluorescent protein was effectively silenced indicating effective cotransfection.

NanoFect Transfection Reagent NEW

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
图 1 纳米多聚物技术提供高效转染。Nanopolymer technology provides efficient transfection. DNA (pGFP which expresses green fluorescent protein) was transfected into the cell types indicated using NanoFect Reagent according to the conditions recommended in the handbook. Transfection efficiency was estimated by counting the number of fluorescent cells.


图 2 储存转染复合物未降低效率。Complexes were prepared using the volumes of NanoFect Reagent indicated and 0.25 μg of pCMV β, a plasmid which expresses the lacZ gene. Complexes were immediately transfected into HeLa cells in 96-well plates or stored at 4°C or -20°C. Stored complexes were transfected after 4 hours or overnight storage. 48 hours after transfection, efficiency was measured by ONPG assay. Transfection was highly efficient under all storage conditions.

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
图 1 A HiPerFect Transfection Reagent enables effective siRNA uptake and release in cells. B Suboptimal lipid reagents do not efficiently release siRNA in cells, leading to inefficient knockdown.


图 2 HeLa cells were reverse transfected in a 384-well plate with a range of concentrations of siRNA targeted against lamin A/C using HiPerFect Transfection Reagent from QIAGEN or Reagent L from another supplier. After 48 hours, gene silencing was assessed by quantitative, real-time RT-PCR.

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
图 1 MCF-7 cells were transfected with positive control MAPK1 siRNA or with nonsilencing siRNA (Control siRNA) using HiPerFect Transfection Reagent. After transfection and incubation, protein expression was determined by western blot analysis. Blots were probed using MAPK1-specific Tag • 100 antibody and anti-tubulin primary antibody. Untransfected cells were also analyzed.

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
图 1 高效共转染。 HeLa S3 cells were cotransfected with a plasmid expressing green fluorescent protein (GFP) and an siRNA targeting either GFP or lamin A/C using TransMessenger Transfection Reagent. After 20 hours, GFP expression was analyzed by **A** fluorescence microscopy and **B** fluorescence activated cell sorting (FACS, the percentage of cells with GFP fluorescence is shown on the plot). **C** Lamin A/C expression in cells was analyzed using quantitative RT-PCR. Lamin A/C expression was effectively silenced in cells that had been cotransfected with lamin A/C siRNA. GFP expression was effectively silenced in cells that had been cotransfected with GFP siRNA.

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
实验细节尽在 Transfection Cell Database

Transfection Cell Database

Our searchable database contains data from researchers like yourself who have shared their experimental results with us. All available information is shown in the detailed view. As always, QIAGEN cannot verify data supplied from outside sources.

Tell us about your [transfection results](#) and get a free gift.

Search for Experimental Data

Nucleic Acid: Transfection Reagent:
Cell Line: Type of Transfection:
Cell Line Species/Tissue:

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Transfection Cell Database - Detailed View

[Back to search results](#)

Cell Line: WI-38

Cell Line Species/Tissue:	Human / Lung fibroblast
Transfection Reagent:	HiPerFect
Nucleic Acid:	siRNA (dsRNA)
Growth Medium:	MEM with L-Glu; supplemented with Vitamins(1/100) and non essential aminoacids (1/100)
Percent Serum (%):	20 %
Reporter System:	
Plasmid Purification Method:	
Plate Format:	30 mm wells (mostly in 6-well plates)
Number of Cells:	225,000 per well
Percent Confluence(%):	Approx. 60-70 %
Amount of Nucleic Acid (µg):	2 µg per well
Amount of Enhancer (µl):	
Amount of Reagent (µl):	80 µl
Complex Incubation on Cells (hrs):	24 h
Analysis Performed Post-Transfection (hrs):	72 h
Transfection Efficiency (%):	90 %
Knockdown Efficiency (%):	80 %
Any modifications to the protocol?:	
	Adjusted the amount of HiPerFect Reagent
Notes:	
References:	

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TransFect Protocol Database

This invaluable resource for transfection experiments takes the guesswork out of transfection protocols. Rather than adapting existing protocols to fit a certain cell type or plate format, the database provides exactly the protocol needed, saving time and effort. Simply select the cell type, nucleic acid, and culture format to receive a QIAGEN transfection protocol to print out or download in convenient PDF format.


Download Transfection protocol

We've successfully generated your custom Transfection protocol for the [Attractene Transfection Reagent](#).

 [Download your protocol](#)

If you have any questions, don't hesitate to contact [Technical Service](#).

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QIAGEN Supplementary Protocol

QIAGEN

Protocol: Fast-Forward Transfection of HCT-116 cells with DNA using Attractene Transfection Reagent

This protocol uses standard parameters and might require further optimization for transfection with this cell type. The amounts given are for one well of a 24-well plate. In this protocol, cell plating and transfection are performed on the same day.

Procedure

1. Shortly before transfection, seed 8×10^4 cells per well of a 24-well plate in 500 μ l of an appropriate culture medium containing serum and antibiotics.
2. For the short time until transfection, incubate the cells under normal growth conditions (typically 37°C and 5% CO₂).

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产品	规格	货号
PolyFect Transfection Reagent — fast and easy DNA transfection of standard cell lines		
PolyFect Transfection Reagent (1 ml)*	For 25-65 transfections in 60 mm dishes or 50-100 transfections in 6-well plates	301105
PolyFect Transfection Reagent (4 x 1 ml)*	For 100-260 transfections in 60 mm dishes or 200-400 transfections in 6-well plates	301107
PolyFect Transfection Reagent (100 ml)*	For 2500-6500 transfections in 60 mm dishes or 5000-10,000 transfections in 6-well plates	301108
Effectene Transfection Reagent — the reagent of choice for DNA transfection of primary cells		
Effectene Transfection Reagent (1 ml)*	1 ml Effectene Reagent, Enhancer, Buffer EC; for 40 transfections in 60 mm dishes or 160 transfections in 12-well plates	301425
Effectene Transfection Reagent (4 x 1 ml) *	4 x 1 ml Effectene Reagent, Enhancer, Buffer EC; for 160 transfections in 60 mm dishes or 640 transfections in 12-well plates	301427
NEW Attractene Transfection Reagent — the reagent of choice for DNA transfection of adherent cells		
Attractene Transfection Reagent (0.5 ml)	Reagent for up to 330 transfections in 24-well plates	301004
Attractene Transfection Reagent (1 ml)	Reagent for up to 660 transfections in 24-well plates	301005
Attractene Transfection Reagent (4 x 1 ml)	Reagent for up to 2640 transfections in 24-well plates	301007
NEW NanoFect Transfection Reagent — DNA transfection for lipid or signal transduction research		
NanoFect Transfection Reagent (0.5 ml)	Reagent for up to 250 transfections in 24-well plates	301204
NanoFect Transfection Reagent (1 ml)	Reagent for up to 500 transfections in 24-well plates	301205
NanoFect Transfection Reagent (4 x 1 ml)	Reagent for up to 2000 transfections in 24-well plates	301207
HiPerFect Transfection Reagent — high target knockdown using siRNA with reduced off-target		
HiPerFect Transfection Reagent (0.5 ml)*	For up to 166 transfections in 24-well plates	301704
HiPerFect Transfection Reagent (1 ml)*	For up to 333 transfections in 24-well plates	301705
HiPerFect Transfection Reagent (4 x 1 ml)*	For up to 1332 transfections in 24-well plates	301707
HiPerFect Transfection Reagent (100 ml)*	For transfections in up to 1388 96-well plates	301709
RNAi Human/Mouse Starter Kit — all you need for siRNA transfection optimization and RNAi control experiments		
RNAi Human/Mouse Starter Kit	0.5 ml HiPerFect Transfection Reagent, Nonsilencing Control siRNA (Alexa Fluor 488 Labeled), Hs/Mm_MAPK1 Control siRNA	301799
TransMessenger Transfection Reagent — efficient plasmid DNA/siRNA cotransfection		
TransMessenger Transfection Reagent (0.5 ml)*	For 60 transfections in 6-well plates or 80 transfections in 12-well plates	301525

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