

HSV-2 Probe ASR

Analyte Specific Reagent. Analytical and performance characteristics are not established.

Product	Contents	Cat. no.
HSV-2 Probe ASR (24)	2 x 24 μ l HSV-2 DNA probe (2.25 μ M) in buffered solution	1066160

Product description

The HSV-2 Probe ASR contains a specific nucleic acid sequence for identification of human herpes simplex virus type 2. The HSV-2 Probe ASR is designed to detect a 154 bp region of the UL5 gene in the HSV-2 genome. The DNA is synthetically derived.

Storage conditions

Upon receipt, store the HSV-2 Probe ASR at -30°C to -15°C . Do not use after the expiration date printed on the outer box label. Repeated thawing and freezing (>2 x) should be avoided, as this may reduce reagent performance. If the reagents are to be used intermittently, they should be frozen in aliquots. Storage at $2-8^{\circ}\text{C}$ should not exceed a period of 5 hours.

Handling instructions

Store the probe away from any source of contaminating DNA, especially amplified DNA products. Use aerosol-barrier pipet tips for pipetting probe.

Quality control

This product was manufactured and released in accordance with the QIAGEN GmbH, Hilden, Quality Control and Quality Assurance procedures. The product has met QIAGEN GmbH, Hilden, quality requirements. Purity of the synthetic oligonucleotide is $>80\%$ as determined by high-performance liquid chromatography (HPLC).



Warnings and precautions

Dispose of in accordance with local regulations.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN[®] kit and kit component.

Trademarks: QIAGEN[®] (QIAGEN Group).
THE PURCHASE OF THIS PRODUCT ALONE DOES NOT IMPLY ANY
LICENSE UNDER PATENTS OWNED BY ROCHE MOLECULAR SYSTEMS,
INC. OR F. HOFFMANN-LA ROCHE LTD.

