

# Uracil DNA Glycosylase

(UDG)



1.2018

**blirt**

# Uracil DNA Glycosylase (UDG)

*E. coli* Uracil DNA Glycosylase (UDG) catalyzes the hydrolysis of the N-glycosylic bond between uracil and sugar, leaving an apyrimidinic site in uracil-containing single-stranded or double-stranded DNA. The enzyme shows no activity on RNA or oligonucleotides.

## Features

- Active over a broad pH range (optimum at pH 8.0)
- Isolated from a recombinant source (*E. coli*)

## Applications

- Helps to eliminate carry-over contamination in PCR
- As a probe for protein-DNA interaction studies
- Glycosylase mediated single nucleotide polymorphism detection (GMPD)
- For cloning of PCR products

## Usage

Treatment of 0.1 µg of uracil-containing DNA with 1 unit of UDG for 10 minutes at 37°C renders the DNA incapable of being copied by DNA polymerase.



## Heat Inactivation

The enzyme can be irreversibly inactivated by incubation at 95°C for 10 min.

## 10x UDG Reaction Buffer

250 mM Tris-HCl (pH 8.0), 1 mM EDTA, 10 mM DTT.

## Quality control

The absence of DNases has been confirmed using the relevant procedures.

## Unit definition

One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from uracil-containing dsDNA. Activity is measured by release of [<sup>3</sup>H]-uracil in a 50 µl reaction containing 0.2 µg DNA in 30 minutes at 37°C.

# Uracil DNA Glycosylase (UDG)

Component	EN19-050 500 U	EN19-250 2500 U
UDG (1 U/ $\mu$ l)	500 $\mu$ l	5 x 500 $\mu$ l
10x UDG Reaction Buffer	1 ml	5 x 1 ml

## Storage & shipping

### Storage conditions

All components should be stored at  $-20^{\circ}\text{C}$  in a freezer without a defrost cycle. When stored under optimum conditions, the reagents are stable until the expiry date.

### Shipping conditions

Shipping on dry or blue ice.

 For research use only