

March 2011

qBiomarker iPSC Expression Lentivirus and Plasmid Handbook

qBiomarker iPSC Expression Lentivirus

qBiomarker iPSC Expression Plasmid

For iPSC creation using transduction-ready
lentiviruses or transfection-ready plasmids



Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Product Use Limitations

qBiomarker iPSC Expression Lentiviruses and qBiomarker iPSC Expression Plasmids are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

CONTENTS

I.	Background and Introduction	4	
II.	Materials Provided	4	
III.	Additional Materials Required	5	
IV.	Protocol	5	
	A.	Before the Experiment	5
	B.	Transfection protocol for reprogramming plasmids	6
	C.	Transduction protocol for lentivirus-expressing reprogramming factors	7

I. Background and Introduction

The discovery that mouse fibroblasts can be reprogrammed and generate induced pluripotent stem cells (iPSCs) with qualities remarkably similar to embryonic stem cells has created a valuable new source of pluripotent cells for drug discovery, cell therapy, and basic research.

SABiosciences (a QIAGEN company) has developed expression vectors for six iPSC related proteins (Nanog, Oct4, Sox2, c-Myc, Klf4 and Lin28) to create reagents for producing induced pluripotent stem cells. Ectopic expression of these factors has been shown to create pluripotent cells which resemble embryonic stem cells. The expression vectors for the six iPSC-related proteins are available as either ready-to-transfect plasmids or ready-to-transduce lentiviral particles.

II. Materials Provided:

1. Reprogramming factors expression plasmid:

Component/Description	Amount	Concentration	Volume
Ready-to-transfect expression vector	10 µg	500 ng/µl	20 µl

2. Lentiviral particles expressing reprogramming factors:

Component/Description	Concentration	Volume
Ready-to-transduce Lentiviral particles	$\geq 0.8 \times 10^7$ TU/ml	250 µl

Storage Conditions:

Ready-to-transfect expression vectors should be stored at -20°C and ready-to-transduce lentiviral expression vectors should be stored at -80°C upon receipt.

Description

Available iPSC Induction vectors express human or mouse Nanog, Oct4, Sox2, c-Myc, Klf4 or Lin28 under the control of cytomegalovirus (CMV) promoter. For the map of parent vectors, see Appendix.

* The DYKDDDDK epitope tag is also referred to as the FLAG[®] epitope tag.

III. Additional Materials Required:

- **Mammalian cell line cultured in the appropriate growth medium**
 - **Human fibroblast medium**
 - **hESC medium**
 - **MEF**
- **Cell culture medium and standard cell culture supplies**
- **Tissue culture plates**
- **Attractene Transfection Reagent** (QIAGEN, Catalog No. 301004) for ready-to-transfect expression vectors
- **SureENTRY Transduction Reagent** (QIAGEN; Catalog No. 336921) for ready-to-transduce lentiviral expression vectors
- **Opti-MEM™ I Reduced-Serum Medium** (Invitrogen)

IV. Protocol:

A. Before the Experiment:

Optimization of transfection conditions (For Plasmids): Optimizing transfection conditions for each cell type is important for the success of an experiment. Variables to consider when optimizing the transfection conditions include cell density, cell viability, amount of DNA, ratio of DNA to transfection reagent, transfection complex formation time, and transfection incubation time (see the detailed protocols for recommendations). The constitutive expressing CMV-mGFP vector (QIAGEN; Cat # 336881) can be used to determine the optimal transfection conditions.

Optimization of transduction conditions (For Lentiviral Particles): Optimization of transduction conditions for each cell type is important for successful lentivirus particle use. Variables to consider, when optimizing the transduction conditions include Multiplicity of Infection (MOI), concentration of SureENTRY Transduction Reagent used, time of assay development and the cell density. The Signal Lenti Reporter Control, Positive Control (GFP) (QIAGEN; Cat # 336891) can be used for determining the optimal transduction conditions.

Multiplicity of Infection (MOI) (For Lentiviral Particles): The transduction efficiency of iPSC Induction expression vectors varies significantly for different cell types. Users should determine the Multiplicity of Infection (MOI), which is the number of transducing lentiviral

qBiomarker iPSC Expression Lentiviruses and Plasmids

particles per cell, required for desired transduction efficiency of a new cell type. The MOI is typically adjusted by increasing or decreasing the amount of virus added per well to a series of wells containing the same number of cells. We recommend testing the Cignal Lenti Reporter Control, Positive Control (Cat # 336891) at MOIs of 5, 10, and 50 (each MOI in triplicate), in order to establish the optimal MOI for each cell type to be studied.

To calculate:

$$\text{Multiplicity of Infection (MOI)} = \frac{\text{Number of transducing units (TU) deposited in a well}}{\text{Number of target cells present in that well}}$$

Total transducing units needed per well (TU) = (Total number of cells per well) x (Desired MOI)

$$\text{Total mL of lentiviral particles to add to each well} = \frac{\text{Total TU needed per well}}{\text{reported on Certificate of Analysis}} \text{ TU/mL}$$

SABiosciences has found that some commonly used cell lines (like HT1080, HEK293 and HepG2 etc.) can be effectively transduced using an MOI between 10 and 25, however, some cell types (like primary cells) are more resistant to transduction and efficient transduction of these cell types may require a higher MOI (~ 50).

B. Transfection protocol for reprogramming plasmids:

The following protocol is designed to transfect the adherent cell line using Attractene Transfection Reagent (Cat. # 301004) in a 96-well plate format. Transfection reagent from other vendors can be used with the iPSC Induction plasmids. *If you are using a transfection reagent other than Attractene Transfection Reagent, follow the manufacturer's protocol.* Moreover, the use of plates or wells of different size require the user to adjust the components in proportion to the surface area of the new plate or well. **The protocol below is just a general guideline; the optimal conditions/amounts should be adjusted according to the cell type and the study requirements.**

Read the protocol completely before starting the experiment.

1. One day before transfection, seed (1×10^5 cells) in human fibroblast medium to each well of a 6-well plate.
2. On the day of transfection, add 5 μg of expression vector into separate 250 μl aliquots of Opti-MEM I Reduced-Serum Medium. Mix gently and incubate mixture for 5 minutes at room temperature (15–25°C).
3. For each well, add 4.5 μl of Attractene Transfection Reagent into 250 μl of Opti-MEM separately. Mix gently and incubate mixture for 5 minutes at room temperature.

IMPORTANT: The amount of transfection reagent and DNA required for optimal performance may vary, depending on the cell line.

4. Add 250 μ l of Attractene mix to 250 μ l expression vector mix. Mix gently and incubate for 20 minutes at room temperature.
 5. Add 500 μ l expression vector-Attractene complexes in medium to the appropriate well containing cells and 1 ml of normal growth medium. This gives a final volume of 1.5 ml. Mix gently by rocking the plate back and forth.
 6. Incubate cells at 37°C in a 5% CO₂ incubator for 16-24 hours.
 7. One day after transduction, remove the medium containing lentiviral particles from wells. Wash the cells with 3 ml of 1 X PBS three times and add 3 ml of human fibroblast medium.
 8. Four days after transduction, plate 5.5×10^5 MEFs in MEF medium on a 10-cm dish coated with 0.1% gelatin. Incubate until the next day.
 9. Five days after transduction, incubate the infected fibroblasts with 1 ml 0.05% trypsin/EDTA for 3 min at 37°C. Stop the trypsinization with 11 ml human fibroblast medium and place all cells into one 10-cm dish plated with MEFs (prepared at step 8).
 10. Seven days after infection remove human fibroblast medium and add 10 ml hESC medium (containing KOSR and 10 ng/ml of β FGF).
 11. Replenish the cells with 10 ml hESC medium every day and observe the cells for any sign of colony formation.
- Typically, sometime around three weeks after infection, various types of colonies will appear.
12. One day before the colonies are big enough to pick, plate 1.9×10^4 MEFs per well of a 24-well plate that have been coated with 0.1% gelatin.
 13. Mark the colonies that show morphology similar to hES cells. Pick the marked colonies under a dissection microscope using a 20- μ l pipette. Put each colony into one well of the 24-well plate prepared at step 12.
 14. Follow standard human ES cells culture procedure to expand and maintain the colonies.

C. Transduction protocol for lentivirus-expressing reprogramming factors:

The following protocol is designed to transduce fibroblast using lentiviral particles expressing reprogramming factors in a 6-well plate format. If you are using plates or wells of different size, adjust the components in proportion to the surface area. **This is just a general guideline; the optimal transduction conditions and cell culture medium for the parental cells should be optimized according to the cell type and the study requirements.**

Read the protocol completely before starting the experiment.

1. One day before transduction, seed (1×10^5 cells) in human fibroblast medium to each well of a 6-well plate.

qBiomarker iPSC Expression Lentiviruses and Plasmids

2. On the day of transduction, replace the old medium from with 2 ml of fresh human fibroblast medium.
3. Add 125 μl of each lentiviral particles (at 0.8×10^7 TU/ml) expressing the reprogramming factors together. Adjust to a final volume of 2 ml using human fibroblast medium.

The amount of lentiviral particles should be adjusted accordingly if the concentration is not 0.8×10^7 TU/ml.

4. (Optional) Add SureENTRY Transduction Reagent to a final concentration of 8 $\mu\text{g/ml}$ in each well. Gently swirl the plate to mix. Incubate at 37°C in a humidified incubator in an atmosphere of 5% CO_2 .
5. One day after transduction, remove the medium containing lentiviral particles from wells. Wash the cells with 3 ml of 1 X PBS three times and add 3 ml of human fibroblast medium.
6. Four days after transduction, plate 5.5×10^5 MEFs in MEF medium on a 10-cm dish coated with 0.1% gelatin. Incubate until the next day.
7. Five days after transduction, incubate the infected fibroblasts with 1 ml 0.05% trypsin/EDTA for 3 min at 37°C. Stop the trypsinization with 11 ml human fibroblast medium and place all cells into one 10-cm dish plated with MEFs (prepared at step 6).
8. Seven days after infection remove human fibroblast medium and add 10 ml hESC medium (containing KOSR and 10 ng/ml of βFGF).
9. Replenish the cells with 10 ml hESC medium every day and observe the cells for any sign of colony formation.

Typically, sometime around three weeks after infection, various types of colonies will appear.

10. Under a dissection microscope and using a 20- μl pipette, pick a colony that shows morphology similar to hES cells and put the colony into one well of a 24-well plate that has been preplated with MEFs at a density of 1×10^4 cells/ cm^2 on 0.1% gelatin-coated wells.
11. Follow standard human ES cells culture procedure to expand the colonies.

Note:

SureENTRY Transduction Reagent enhances transduction of most cells, however, some cells like primary neurons are sensitive to the SureENTRY Transduction Reagent. Do not add SureENTRY Transduction Reagent to these types of cells. If working with a cell type for the first time, a SureENTRY Transduction Reagent control only well should be used to determine cell sensitivity.

V. Frequently Asked Questions:

1. Which reprogramming system should I chose, plasmids or lentiviral particles?

Both systems work for iPSC induction. Plasmids would not integrate into the host genome therefore is the better choice if the genome integrity is your first concern. Lentiviral particles generally have higher expression of the genes and is the better choice if the iPSC induction rate if your first concern.

2. I see DYKDDDDK on the vector maps. What is it?

DYKDDDDK is the epitope tag that is also referred to as the FLAG[®] epitope tag.

3. What is MOI?

MOI is an abbreviation for Multiplicity of Infection or the number of viral particles exposed to a cell.


4. Where I can get more information about the plasmids or lentiviral particles?

More information can be found at:

http://www.sabiosciences.com/support_faq.php?target=reporter


Appendix: Vector Maps

1. Vector map of reprogramming factor expression plasmid:



MCS: Multiple cloning site; ORF: open reading frame

2. Vector map of lentiviral plasmid expressing reprogramming factors:



MCS: Multiple cloning site; ORF: open reading frame

Ordering Information

Product	Contents	Cat. no.
qBiomarker iPSC Expression Lentivirus	$\geq 0.8 \times 10^7$ TU/ml ready-to-transduce lentiviral particles	337402
qBiomarker iPSC Expression Plasmid	10 μ g ready-to-transfect expression vector	337401

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

qBiomarker iPSC Expression Lentiviruses and Plasmids

Notes

Notes

qBiomarker iPSC Expression Lentiviruses and Plasmids

Notes

Trademarks: QIAGEN® (QIAGEN Group); Opti-MEM™ (Life Technologies, Inc.); FLAG® (Sigma-Aldrich Co.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the qBiomarker iPSC Expression Lentivirus and qBiomarker iPSC Expression Plasmid to the following terms:

1. The qBiomarker iPSC Expression Lentivirus and qBiomarker iPSC Expression Plasmid may be used solely in accordance with the *qBiomarker iPSC Expression Lentivirus and Plasmid Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *qBiomarker iPSC Expression Lentivirus and Plasmid Handbook* and additional protocols available at www.qiagen.com.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see www.qiagen.com.

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replace this product in the event the product fails to perform due to any manufacturing defect. SABiosciences Corporation makes no other warranties of any kind, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose. SABiosciences Corporation shall not be liable for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use this product.

NOTICE TO PURCHASER

The purchase of a qBiomarker iPSC Expression product includes a limited, nonexclusive license to use the kit components for research use only. This license does not grant rights to use the kit components to manufacture commercial products without written approval of SABiosciences Corporation. No other license, expressed, implied or by estoppel, is granted.

Signal™ Lenti Reporter Assay System Limited License.

The purchaser of this product agrees to only use the lentiviral particles in this kit in cell-based reporter assays for in vivo and in vitro internal research.

Use of this product for Commercial Purposes requires a license from Sigma-Aldrich Corporation. The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party, or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration, but excludes not-for-profit core facilities providing services within their own research institutions at cost.

This product is licensed under U.S. Pat. Nos. 5,817,491; 5,591,624; 5,716,832; 6,312,682; 6,669,936; 6,235,522; 6,924,123 and foreign equivalents from Oxford BioMedica (UK) Ltd., Oxford, UK, and is provided for use in academic and commercial in vitro and in vivo research for elucidating gene function, and for validating potential gene products and pathways for drug discovery and development, but excludes any use of LentiVector® technology for: creating transgenic birds for the purpose of producing useful or valuable proteins in the eggs of such transgenic birds, the delivery of gene therapies, and for commercial production of therapeutic, diagnostic or other commercial products not intended for research use where such products do not consist of or incorporate a lentiviral vector. Information about licenses for commercial uses excluded under this license is available from Oxford BioMedica (UK), Ltd., Medawar Center, Oxford Science Park, Oxford OX4 4GA UK enquiries@oxfordbiomedica.co.uk or BioMedica Inc., 11622 El Camino Real #100, San Diego CA 92130-2049 USA. LentiVector is a registered US and European Community trademark of Oxford BioMedica plc.

www.qiagen.com

Australia ■ Orders 1-800-243-800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

Austria ■ Orders 0800-28-10-10 ■ Fax 0800-28-10-19 ■ Technical 0800-28-10-11

Belgium ■ Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

Brazil ■ Orders 0800-557779 ■ Fax 55-11-5079-4001 ■ Technical 0800-557779

Canada ■ Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

China ■ Orders 86-21-3865-3865 ■ Fax 86-21-3865-3965 ■ Technical 800-988-0325

Denmark ■ Orders 80-885945 ■ Fax 80-885944 ■ Technical 80-885942

Finland ■ Orders 0800-914416 ■ Fax 0800-914415 ■ Technical 0800-914413

France ■ Orders 01-60-920-926 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930 ■ Offers 01-60-920-928

Germany ■ Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

Hong Kong ■ Orders 800 933 965 ■ Fax 800 930 439 ■ Technical 800 930 425

Ireland ■ Orders 1800 555 049 ■ Fax 1800 555 048 ■ Technical 1800 555 061

Italy ■ Orders 800-789-544 ■ Fax 02-334304-826 ■ Technical 800-787980

Japan ■ Telephone 03-6890-7300 ■ Fax 03-5547-0818 ■ Technical 03-6890-7300

Korea (South) ■ Orders 080-000-7146 ■ Fax 02-2626-5703 ■ Technical 080-000-7145

Luxembourg ■ Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

Mexico ■ Orders 01-800-7742-639 ■ Fax 01-800-1122-330 ■ Technical 01-800-7742-436

The Netherlands ■ Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

Norway ■ Orders 800-18859 ■ Fax 800-18817 ■ Technical 800-18712

Singapore ■ Orders 1800-742-4362 ■ Fax 65-6854-8184 ■ Technical 1800-742-4368

Spain ■ Orders 91-630-7050 ■ Fax 91-630-5145 ■ Technical 91-630-7050

Sweden ■ Orders 020-790282 ■ Fax 020-790582 ■ Technical 020-798328

Switzerland ■ Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

UK ■ Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

USA ■ Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

