

Novel solutions for molecular food safety testing



Sarah Fakh¹, Petra Dehring¹, Corinna Küppers¹, Sandra Luley¹, Wolfgang Hauser², Wolfgang Weber², and Holger Engel¹

¹ QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, Germany

² ifp Institut für Produktqualität GmbH, Teltowkanalstrasse 2, 12247 Berlin, Germany


Introduction

Food safety concerns have steadily increased in recent years and are now one of the most important challenges for food authorities. There is a requirement for highly sensitive pathogen detection and most food testing laboratories rely on time-consuming culturing methods, with associated high storage and production costs.

QIAGEN has recently introduced the *mericon*[™] Food Safety Testing product line, covering reagent kits and instruments for all major food applications – from DNA extraction to highly specific real-time PCR detection.

The data shown here illustrate the sensitivity of the new *mericon* PCR assays, as well as tolerance of the Multiplex PCR Master Mix to PCR inhibitors. Representative pathogen DNA extractions from food enrichment cultures and detection with the *mericon* assays are shown for *Salmonella enterica* and *Listeria monocytogenes*. Efficiency of 2 new DNA extraction systems for Gram-positive and -negative bacteria are compared for cultured samples with and without food matrices.

We also present the QIASymphony[®] Rotor-Gene[®] Q, a fully automated system for DNA extraction directly from enrichment cultures, followed by automated reaction setup, and detection on the Rotor-Gene Q real-time PCR instrument. The protocol setup and results of QIASymphony extraction detection runs for food enrichment cultures are described.



The QIAGEN food safety testing pathogen workflow

Sample preparation

mericon DNA Bacteria Kit

- Thermal lysis for Gram-negative bacteria.

mericon DNA Bacteria Plus Kit

- Mechanical lysis for Gram-positive and Gram-negative bacteria.

The TissueLyser LT, for simultaneous mechanical disruption of 12 samples.

The QIASymphony Rotor-Gene Q system

- Fully automated instrument for pathogen DNA or RNA extraction, assay setup, and detection.

The QIASymphony Rotor-Gene Q.

Real-time PCR assays

mericon PCR detection kits

- 11 real-time PCR assays for all major food-borne pathogens.

The *mericon* portfolio includes:

- Salmonella* spp.
- L. monocytogenes*
- Listeria* spp.
- C. sakazakii*
- Campylobacter* spp.
- C. jejuni, laridis, coli*
- E. coli* VTEC stx 1/2
- S. aureus*
- Y. enterocolitica*
- Shigella* spp.
- Legionella* spp.
- L. pneumophila*

Detection

Rotor-Gene Q real-time PCR cycler


- All *mericon* assays are optimized for use with the Rotor-Gene Q, but are also compatible with other real-time (block) cyclers.

The Rotor-Gene Q is a precise and versatile real-time PCR cycler with a unique centrifugal rotary design.



mericon assays.

Real-time PCR pathogen detection in complex food matrices: *mericon* Assay systems

- Assay portfolio includes all major food-borne pathogens.
- Highly sensitive target detection down to 10 copies.
- Internal Control to indicate possible PCR inhibition.
- PCR Multiplex Master Mix highly tolerant to PCR inhibitors.
- Combine several assays in one run; all *mericon* PCR assays use the same cycling conditions.



Comparison of resistance to residual PCR inhibitors of the *mericon* Multiplex PCR Master Mix and a master mix provided by supplier AB. A plasmid carrying an artificial DNA sequence (DNA model TaqMan[®] system, 10⁶ copies) was amplified in reactions spiked with inhibitors potentially carried over from sample preparations. The change in threshold cycles of each inhibitor compared to that of a water negative control is shown.



The *mericon* Salmonella spp Kit detection sensitivity with purified *Salmonella* DNA on the Rotor-Gene Q. Titration was carried out in 1:10 dilutions from 10,000 to 10 copies per PCR reaction. Above: Template DNA amplified in duplicates, detection in the green channel (R²=0.993). Below: Internal Control, MAX-labeled, detection in the yellow channel.

Fully automated pathogen DNA extraction and detection: QIASymphony Rotor-Gene Q

QIASymphony Rotor-Gene Q System

Sample preparation

- Loading of untreated enrichment culture
- Thermal lysis at 90°C for 5 min
- DNA sequestration via magnetic, silica-coated beads
- Automated transfer of elution plate to assay setup module

Assay setup

- Automated data transfer for continuous sample tracking
- Highly accurate pipetting into Rotor-Gene Q ring*

Protocol time: 25 min

Real-time PCR detection



- PCR target detection via integrated Rotor-Gene Q instrument

Protocol time: 74 min

Processing of 48 samples:
1 ml protocol → 140 min
200 µl protocol → 96 min

Automated pathogen DNA extraction and detection with the QIASymphony Rotor-Gene Q. 750 g chocolate with coffee filling was homogenized in 6.75 l buffered peptone water and inoculated with ~10 cfu *Salmonella*. Detection was performing using the *mericon* Salmonella spp Kit.

Time-efficient manual pathogen DNA extraction and detection: *mericon* DNA Bacteria Kits




Comparison of DNA extraction efficiency for Gram-negative and Gram-positive bacteria. **A** 25 g of milk or hazelnut spread was homogenized in 225 ml buffered peptone, inoculation with ~10 cfu *Salmonella* and enriched for 20 h at 37°C. Undiluted and diluted enrichment culture samples were processed using both *mericon* DNA Bacteria Kits and tested with the *mericon* Salmonella spp Kit on the Rotor-Gene Q. Overall, the efficiency of DNA extraction from Gram-negative bacteria is similar with both preparation protocols. **B** Serial dilutions of enrichment cultures of *S. enterica* (buffered peptone water) and *L. monocytogenes* (1/2 Fraser medium) were processed using both *mericon* DNA Bacteria Kits. For *L. monocytogenes* (Gram-positive), mechanical lysis delivered a more effective extraction, most evident at low bacterial concentrations where the average C_t was 4.9 cycles lower with mechanical lysis.

- QIAGEN introduces 2 new manual systems for pathogen DNA extraction from enrichment cultures.
- The *mericon* DNA Bacteria Plus Kit is based on a mechanical lysis using glass beads. The method was developed for Gram-positive or difficult to lyse bacteria, but is universally applicable.
- The *mericon* DNA Bacteria Kit is based on a thermal pathogen lysis for Gram-negative or easy to lyse bacteria.
- Pathogen DNA extraction is reduced to a minimum of simple, time-efficient steps; 30 samples can be processed in ~30 minutes.

Conclusions

- QIAGEN introduces a comprehensive product portfolio for all major applications in molecular food safety testing.
- The new *mericon* products cover reagent kits and instruments from DNA extraction to real-time PCR detection.
- The portfolio contains sample preparation and detection solutions for all major food pathogens, as well as complete solutions for GMO, animal, plant, and allergen detection.



Trademarks: QIAGEN[®], QIASymphony[®], *mericon*[™], Rotor-Gene[®] (QIAGEN Group); FAM[®] (Life Technologies); MAX[™] (Integrated DNA Technologies, Inc.); TaqMan[®] (Roche Group).