

September 2024

Important Note

Use of USB to serial adapter on QIAxcel® Advanced

Dear valued QIAxcel Advanced customer,

When using the QIAxcel Advanced instrument with USB to serial adapter, QXA (cat. no. 9027105) (Roline® USB to Serial Converter, 1.8m, no. 12.02.1163) download of a driver and installation on the operating computer is required.

Operating system specific driver can be downloaded from the **Resources** tab of this page: <https://prolificusa.com/product/pl2303gt-usb-rs232-serial-bridge-controller-built-rs232-xcvr/>

In case of any question, please contact our Technical Support Center at www.qiagen.com/support/technical-support

Best regards,

QIAGEN