


pQE-16 Vector

Positions of elements in bases

Vector size (bp)	3996
Start of numbering at <i>Xho</i> I (CTCGAG)	1–6
T5 promoter/lac operator element	7–87
T5 transcription start	61
Dehydrofolate reductase coding sequence	133–690
6xHis-tag coding sequence	703–720
Lambda <i>t</i> ₀ transcriptional termination region	743–837
<i>rrnB</i> T1 transcriptional termination region	1599–1697
ColE1 origin of replication	2173
β-lactamase coding sequence	3791–2931



pQE-16

