

QIAGEN® Proteinase K

Version 1

Proteolytic Enzyme

For Laboratory Use

Active ingredient: Proteinase K solution (>600 mAU/ml)

QIAGEN Proteinase K enables efficient lysis of sample material.

REF	19134
------------	-------

MAT	1094085
------------	---------



QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY

Specifications

Volume: 10 ml

Activity: >600 mAU/ml (1 mAU is the activity that releases folin-positive amino-acids and peptides from denatured hemoglobin in 1 minute corresponding to 1 μ mol tyrosine)

Storage

QIAGEN Proteinase K can be stored dry at room temperature (15–25°C) until the expiration date.

Symbols

	Catalog number
	Material number
	Lot number
	Global Trade Item Number
	To be used by
	Temperature limitations
	Legal manufacturer



Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of QIAGEN Proteinase K is tested against predetermined specifications to ensure consistent product quality.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN® kit and kit component.

Ordering Information

Product	Contents	Cat. no.
QIAGEN Proteinase K (10 ml)	10 ml (>600 mAU/ml, solution)	19134

Trademarks: QIAGEN® (QIAGEN Group).

1094085 03/2015 © 2015 QIAGEN, all rights reserved.

www.qiagen.com

Canada = 800-572-9613

China = 8621-3865-3865

Denmark = 80-885945

Australia = 1-800-243-800

Finland = 0800-914416

Austria = 0800/281010

France = 01-60-920-930

Belgium = 0800-79612

Germany = 02103-29-12000

Brazil = 0800-557779

Hong Kong = 800 933 965

Ireland = 1800 555 049

Italy = 800-787980

Japan = 03-6890-7300

Korea (South) = 080-000-7145

Luxembourg = 8002 2076

Mexico = 01-800-7742-436

The Netherlands = 8000 0229592

Norway = 800-18859

Singapore = 1800-742-4368

Spain = 91-630-7050

Sweden = 020-790282

Switzerland = 055-254-22-11

UK = 01293-422-911

USA = 800-426-8157



Sample & Assay Technologies