

Discover a new level of reliability

in your manual lab procedure



QIAGEN[®] digitalized protocols and Gilson[®] digital pipetting system

Refine pipetting accuracy. Promote traceability. Achieve consistently precise and reproducible manual procedures with Gilson pipettes and QIAGEN kits.

- **Skip the paper-based protocols.**
Get real-time guidance on your QIAGEN protocol through Gilson's TRACKMAN[®] Connected tablet-based system.
- **Say goodbye to tedious manual pipette adjustments.**
Gilson's PIPETMAN[®] M Connected pipettes follow the QIAGEN protocol running on TRACKMAN Connected, configuring your pipette settings automatically so you don't have to do it yourself.
- **Report and trace experiment runs automatically.**
Let TRACKMAN Connected generate your run reports for you. Analyze your team's output and processes with full data on hand.

Four simple steps

To run a QIAGEN protocol using Gilson's TRACKMAN Connected system



1

Choose the QIAGEN kit you need



2

Answer questions to define your experiment*



3

Follow real-time guidance on the procedure†



4

Download a comprehensive run report to accelerate reporting‡

* Based on your answers, TRACKMAN Connected automatically selects the protocol you need.

† TRACKMAN Connected guides you through the entire procedure while automatically adjusting the PIPETMAN M Connected pipettes.

‡ The run report can be uploaded to ELN systems or stored locally.

➔ Find all QIAGEN kits that you can use with the TRACKMAN Connected system at www.qiagen.com/trackman-connected. For more information on the TRACKMAN Connected system and how it can link with QIAGEN nucleic acid isolation kits, visit www.gilson.com/gilson-connect/qiagen.

Ordering Information

Product	Cat. no.	Product	Cat. no.
DNeasy® Blood & Tissue Kit (50)*	69504	QIAamp ccfDNA/RNA Kit (50)	55184
QIAamp® DNA Mini Kit (50)*	51304	RNeasy® Plus Mini Kit (50)*	74134
QIAprep® Spin Miniprep Kit (50)*	27104	QIAamp Viral RNA Mini Kit (50)*	52904
QIAquick® PCR Purification Kit (50)*	28104	QIAamp 96 Viral RNA Kit (10)	52962

* Also available in other sizes.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAprep®, QIAquick®, DNeasy®, RNeasy® (QIAGEN Group); Gilson®, PIPETMAN®, TRACKMAN® (Gilson Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2021 QIAGEN, all rights reserved. PROM-17952-001

In collaboration with



Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com