


NeXtal Stock Kit Salt

For recreating Screening Suite conditions



The NeXtal Stock Kit Salt provides:

- Stock solutions ideal for developing grids around original hit conditions
- Compatibility with all EasyXtal Screening Suites
- Maximized reproducibility through online access to production reports

No.	Salt	Concentration
1	Ammonium acetate	1 M
2	Ammonium chloride	5 M
3	Ammonium dihydrogen phosphate	2.5 M
4	Ammonium fluoride	10 M
5	Ammonium formate	10 M
6	Ammonium iodide	1 M
7	Ammonium nitrate	10 M
8	Ammonium sulfate	3.5 M
9	Calcium acetate hydrate	1 M
10	Calcium chloride dihydrate	2 M
11	di-Ammonium hydrogen citrate	2.5 M
12	di-Ammonium hydrogen phosphate	3.5 M
13	di-ammonium tartrate	2 M
14	di-Potassium hydrogen phosphate	3 M
15	di-Sodium hydrogen phosphate dihydrate	0.5 M
16	di-Sodium tartrate dihydrate	1.5 M
17	Lithium acetate dihydrate	5 M
18	Lithium chloride anhydrous	10 M
19	Lithium nitrate	8 M
20	Lithium sulfate monohydrate	1 M
21	Magnesium acetate tetrahydrate	1 M
22	Magnesium chloride hexahydrate	2 M
23	Magnesium formate dihydrate	1 M
24	Magnesium nitrate hexahydrate	3 M

No.	Buffer	Concentration
25	Magnesium sulfate heptahydrate	2.5 M
26	Potassium acetate	5 M
27	Potassium chloride	4 M
28	Potassium dihydrogen phosphate	1.5 M
29	Potassium fluoride	6 M
30	Potassium formate	14 M
31	Potassium iodide	1 M
32	Potassium nitrate	3 M
33	Potassium sodium tartrate tetrahydrate	1.5 M
34	Potassium sulfate	0.5 M
35	Potassium thiocyanate	8 M
36	Sodium acetate trihydrate	1.6 M
37	Sodium chloride	5 M
38	Sodium dihydrogen phosphate monohydrate	5 M
39	Sodium fluoride	0.8 M
40	Sodium formate	7 M
41	Sodium iodide	1 M
42	Sodium nitrate	7 M
43	Sodium sulfate decahydrate	1 M
44	Sodium thiocyanate	8 M
45	tri-Lithium citrate tetrahydrate	1.5 M
46	tri-Potassium citrate monohydrate	2.5 M
47	tri-Sodium citrate dihydrate	1.6 M
48	Zinc acetate dihydrate	1 M

Find out more and order EasyXtal and NeXtal products online at www.qiagen.com/crystallization

Trademarks: QIAGEN® (QIAGEN Group) 10/2006 © 2006 QIAGEN, all rights reserved

Australia ■ 03-9840-9800

Austria ■ 0800/28-10-10

Belgium ■ 0800-79612

Canada ■ 800-572-9613

China ■ 021-51345678

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-920

Germany ■ 02103-29-12000

Ireland ■ 1800 555 049

Italy ■ 02-33430411

Japan ■ 03-5547-0811

Luxembourg ■ 8002-2076

The Netherlands ■ 0800-0229592

Norway ■ 800-18859

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

