



July 2024

Product Sheet

Buffer ASL

Contents

Buffer ASL (560 mL)

Catalog no.

19082

Buffer ASL

4 x 140 mL

For Laboratory Use

REF

19082



QIAGEN GmbH
QIAGEN Strasse 1, 40724 Hilden, GERMANY

R1

MAT

1135486

Shipping and Storage

The Buffer ASL 4 × 140 mL (cat. no. 19082) is shipped at ambient temperature. It should be stored immediately upon receipt at ambient temperature. Keep container tightly closed in a dry and well-ventilated place.










When stored correctly at 15–25°C, the Buffer ASL 4 × 140 mL is stable until the expiration date printed on the kit label as well as on the bottle labels.

Do not use if stored outside the specifications, if the packaging has been damaged, or if other signs of deterioration or malfunction are visible.

Disposal

Dispose used and unused product in compliance with local and national regulations. Follow the recommendations in the Safety Data Sheet (SDS).

Symbols

	Catalog number
	Consult instructions for use
	Material number
	Lot number
	Global Trade Item Number
	Unique Device Identifier
	To be used by
	Temperature limitations
	Manufacturer

Intended Use

The Buffer ASL 4 × 140 mL enables pre-treatment of human stool samples to support lysis and solubilization of the samples. Buffer ASL can be used before a subsequent manual or automated extraction of nucleic acids.

The Buffer ASL is intended for laboratory use.

Description and Principle

Buffer ASL (Stool Lysis Buffer) enables pre-treatment of human stool samples to support lysis and solubilization of the samples.

Buffer ASL (Stool Lysis Buffer) can be used before a subsequent manual or automated extraction of nucleic acids.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety, where you can find, view, and print the SDS for each QIAGEN® kit and kit component.

Specimens and samples are potentially infectious. Discard sample and assay waste according to your local safety procedures.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of Buffer ASL 4 × 140 mL is tested against predetermined specifications to ensure consistent product quality.

Procedure

Refer to the instructions for use (handbook) of the respective product requiring use of Buffer ASL as part of a product workflow.

Ordering Information

Product	Contents	Cat. no.
Buffer ASL 4 × 140 mL	Buffer for pre-treatment of stool samples	19082

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN® products to follow any applicable local regulation and we also recommend to follow any applicable standards and guidelines.

Document Revision History

Date	Changes
07/2024	Initial release

Trademarks: QIAGEN®, Sample to Insight® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1135486 07/2024 HB-3602-001 © 2024 QIAGEN, all rights reserved.