


NeXtal Stock Kit Buffer

For recreating Screening Suite conditions



The NeXtal Stock Kit Buffer provides:

- Stock solutions ideal for developing grids around original hit conditions
- Compatibility with all EasyXtal Screening Suites
- Maximized reproducibility through online access to production reports

No.	Buffer	Concentration	Final pH
1	Citric Acid	1.0 M	2.2
2	Citric Acid	1.0 M	2.4
3	Citric Acid	1.0 M	2.6
4	Citric Acid	1.0 M	2.8
5	Citric Acid	1.0 M	3
6	Citric Acid	1.0 M	3.2
7	Citric Acid	1.0 M	3.4
8	Citric Acid	1.0 M	3.6
9	Citric Acid	1.0 M	3.8
10	Sodium Acetate	1.0 M	4
11	Sodium Acetate	1.0 M	4.2
12	Sodium Acetate	1.0 M	4.4
13	Sodium Acetate	1.0 M	4.6
14	Sodium Acetate	1.0 M	4.8
15	tri-Sodium Citrate	1.0 M	5
16	tri-Sodium Citrate	1.0 M	5.2
17	tri-Sodium Citrate	1.0 M	5.4
18	tri-Sodium Citrate	1.0 M	5.6
19	tri-Sodium Citrate	1.0 M	5.8
20	Sodium Cacodylate	1.0 M	6
21	Sodium Cacodylate	1.0 M	6.2
22	Sodium Cacodylate	1.0 M	6.4
23	Sodium Cacodylate	1.0 M	6.6
24	Sodium Cacodylate	1.0 M	6.8

No.	Buffer	Concentration	Final pH
25	Sodium HEPES	1.0 M	7
26	Sodium HEPES	1.0 M	7.2
27	Sodium HEPES	1.0 M	7.4
28	Sodium HEPES	1.0 M	7.6
29	Sodium HEPES	1.0 M	7.8
30	Tris Hydrochloride	1.0 M	8
31	Tris Hydrochloride	1.0 M	8.2
32	Tris Hydrochloride	1.0 M	8.4
33	Tris Hydrochloride	1.0 M	8.6
34	Tris Hydrochloride	1.0 M	8.8
35	CAPSO	0.5 M	9
36	CAPSO	0.5 M	9.2
37	CAPSO	0.5 M	9.4
38	CAPSO	0.5 M	9.6
39	CAPSO	0.5 M	9.8
40	CAPS	0.5 M	10
41	CAPS	0.5 M	10.2
42	CAPS	0.5 M	10.4
43	CAPS	0.5 M	10.6
44	CAPS	0.5 M	10.8
45	CAPS	0.5 M	11
46	CAPS	0.5 M	11.2
47	CAPS	0.5 M	11.4
48	CAPS	0.5 M	11.6

Find out more and order EasyXtal and NeXtal products online at
www.qiagen.com/crystallization

Trademarks: QIAGEN® (QIAGEN Group) 10/2006 © 2006 QIAGEN, all rights reserved

Australia ■ 03-9840-9800

Austria ■ 0800/28-10-10

Belgium ■ 0800-79612

Canada ■ 800-572-9613

China ■ 021-51345678

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-920

Germany ■ 02103-29-12000

Ireland ■ 1800 555 049

Italy ■ 02-33430411

Japan ■ 03-5547-0811

Luxembourg ■ 8002-2076

The Netherlands ■ 0800-0229592

Norway ■ 800-18859

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



WWW.QIAGEN.COM