


NeXtal Stock Kit Cryo

For recreating Screening Suite conditions



The NeXtal Stock Kit Cryo provides:

- Stock solutions ideal for developing grids around original hit conditions
- Compatibility with all EasyXtal Screening Suites
- Maximized reproducibility through online access to production reports

No.	Cryoprotectant	Concentration
1	1,6-Hexanediol	6 M
2	D-(+)-Glucose	70% (w/v)
3	D-(+)-Sucrose	70% (w/v)
4	Glycerol	100% (v/v)
5	Isopropanol	100% (v/v)
6	Jeffamine M-600 pH 7.0	40% (v/v)
7	Lithium acetate	5 M
8	Lithium chloride	10 M
9	Lithium nitrate	8 M
10	Lithium sulfate	2 M
11	Magnesium acetate	3.4 M
12	MPD	100% (v/v)
13	NVH Oil	100% (v/v)
14	Paratone-N	100% (v/v)
15	PEG 200	100% (v/v)
16	PEG 400	100% (v/v)
17	PEG 4000	50% (w/v)
18	Pentaerythritol ethoxylate (15/4 EO/OH)	50% (v/v)
19	Pentaerythritol propoxylate (5/4 PO/OH)	50% (v/v)
20	Polyvinylpyrrolidone K15	50% (w/v)
21	Propylene glycol	100% (v/v)
22	Sodium formate	7 M
23	Sodium nitrate	7 M
24	Xylitol	70% (w/v)

Protein Crystallization Suites and Formats

	EasyXtal Microplate	NeXtal Deep-Well Block	EasyXtal Tool X-Seal	EasyXtal DG Tool X-Seal	NeXtal Tubes
The Classics Suite	■	■	■	■	■
The Classics Lite Suite	■	■	■	■	■
The Cryos Suite	■	■	■	■	■
The PEGs Suite	■	■	■	■	■
The AmSO ₄ Suite	■	■	■	■	■
The MPD Suite	■	■		■	■
The Anions Suite	■	■	■	■	■
The Classics Suite	■	■	■	■	■
The pHClear Suite	■	■	■	■	■
The pHClear II Suite	■	■	■	■	■
The MbClass Suite	■	■		■	■
The MbClass II Suite	■	■		■	■
The Protein Complex Suite	■	■	■	■	■
The PEGs II Suite	■	■	■	■	■
The CompPAS Suite	■	■	■	■	■
The PACT Suite	■	■	■	■	■
The Nucleix Suite	■	■	■	■	■
The JCSG+ Suite	■	■	■	■	■
The Opti-Salts Suite	■		■	■	
Pre-Screen Assay				■	

Find out more and order EasyXtal and NeXtal products online at www.qiagen.com/crystallization

Trademarks: QIAGEN® (QIAGEN Group) 10/2006 © 2006 QIAGEN, all rights reserved

Australia ■ 03-9840-9800

Austria ■ 0800/28-10-10

Belgium ■ 0800-79612

Canada ■ 800-572-9613

China ■ 021-51345678

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-920

Germany ■ 02103-29-12000

Ireland ■ 1800 555 049

Italy ■ 02-33430411

Japan ■ 03-5547-0811

Luxembourg ■ 8002-2076

The Netherlands ■ 0800-0229592

Norway ■ 800-18859

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

