

Single Tube Magnet

for separation of BioMag® particles in a 1.5 ml or 2 ml tube

The Single Tube Magnet is a convenient separation device for cell separation, magnetocapture, or immunoassays using BioMag® magnetic particles in a 1.5 ml or 2 ml tube. The beads are pulled to the side of the tube by the magnetic field, which has a strength of greater than 20 megaoersted, and held in place while solutions are exchanged or removed. The Single Tube Magnet can separate volumes between 40 µl and 1 ml that contain <5 mg BioMag particles. Removal of the magnetic field allows easy resuspension of the magnetic particles.

Important notes about the Single Tube Magnet

- Do not immerse the Single Tube Magnet
- Wipe clean only using a mild detergent
- Harsh chemicals and organic solvents may damage the magnet
- **WARNING:** Strong magnetic field. Store away from magnetic media, e.g., computers, videotapes, or credit cards

Ordering Information

Product	Contents	Cat. No.
Single Tube Magnet	Magnet for separating magnetic particles in a 1.5 ml or 2 ml tube	36910
Related products		
12-Tube Magnet	Magnet for separating magnetic particles in 12 x 1.5 ml or 2 ml tubes	36912
15 ml/50 ml Tube Magnet	Magnet for separating magnetic particles in 5 x 15 ml and 3 x 50 ml tubes	36935
Flask Magnet	Magnet for separating magnetic particles in a cell culture flask	36937

Patented or patent-pending and/or registered or registration-pending trademark of the QIAGEN Group: QIAGEN®, BioMag®.
Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2002 QIAGEN, all rights reserved.

www.qiagen.com

QIAGEN:

Australia Tel. 03-9489-3666 Fax 03-9489-3888	Canada Tel. 800-572-9613 Fax 800-713-5951	France Tel. 01-60-920-930 Fax 01-60-920-925
Germany Tel. 02103-29-12400 Fax 02103-29-22000	Italy Tel. 02-33430411 Fax 02-33430426	Japan Tel. 03-5547-0811 Fax 03-5547-0818
Switzerland Tel. 061-319-30-31 Fax 061-319-30-33	UK and Ireland Tel. 01293-422-999 Fax 01293-422-922	USA Tel. 800-426-8157 Fax 800-718-2056

Distributors:

Argentina Tecnolab S.A. [011] 4555 0010 **Austria/Slovenia** Merck Eurolab GmbH [01] 576 00 0 **Belgium/Luxemburg** Westburg b.v. 0800-1-9815 **Brazil** Unisciencia do Brasil 011 3622 2320 **China** Gene Company Limited [852]2896-6283 **Cyprus** Scientronics Ltd [02] 765 416 **Czech Republic** BIOCONSULT spol. s.r.o. [420] 2 417 29 792 **Denmark** Merck Eurolab A/S 43 86 87 88 **Egypt** Clinilab 52 57 212 **Finland** Merck Eurolab Oy [09]804 551 **Greece** BioAnalytica S.A. [01]640 03 18 **India** Genetix [011]542 1714 or [011]515 9346 **Israel** Westburg [Israel] Ltd. 08 6650813/4 or 1-800 20 22 20 **Korea** LRS Laboratories, Inc. [02] 924-86 97 **Malaysia** RESEARCH BIOLABS SDN. BHD. (603)8070 3101 **Mexico** Quimica Valaner S.A. de C.V. [55] 55 25 57 25 **The Netherlands** Westburg b.v. [033]4950094 **New Zealand** Biolab Scientific Ltd. [09] 980 6700 or 0800 933 966 **Norway** Merck Eurolab AS 22 90 00 00 **Poland** Syngen Biotech Sp.z.o.o. [071] 351 41 06 or 0601 70 60 07 **Portugal** IZASA PORTUGAL LDA [21] 424 7312 **Singapore** Research Biolabs Pte Ltd 2731066 **Slovak Republic** BIOCONSULT Slovakia spol. s.r.o. [02] 5022 1336 **South Africa** Southern Cross Biotechnology [Pty] Ltd [021] 671 5166 **Spain** IZASA, S.A. [93] 902.20.30.90 **Sweden** Merck Eurolab AB (08) 621 34 00 **Taiwan** TAIGEN Bioscience Corporation [02] 2880 2913 **Thailand** Theera Trading Co. Ltd. [02] 412-5672 **In other countries contact:** QIAGEN, Germany



1019511 01/2002