



QIAGEN GmbH
QIAGEN Str. 1
40724 Hilden, Germany

December 2022

Release Note: EZ2[®] Connect, EZ2 Connect Fx, and EZ2 Connect MDx Protocol Package v7

Dear valued customer,

At QIAGEN[®], we are continuously making improvements to provide the highest level of performance with our instruments and protocols. Consequently, we are pleased to inform you about the opportunity to upgrade your EZ2 Connect, EZ2 Connect Fx, or EZ2 Connect MDx instrument to Protocol Package v7. This protocol package includes the improvements and changes to the protocols listed below.

| Protocol name | New protocol version | Changes |
|---------------------------------|----------------------|--|
| EZ1&2 [™] DNA FFPE Kit | 2 | Bug fix: improved elution in 100 µL |
| EZ1&2 ccfDNA Kit | 4 | Implemented additional measures to reduce the risk of magnet contamination |
| AdnaTest CTC Select Kit | 3 | Implemented additional measures to reduce the risk of magnet contamination |

If you have further questions, please contact your local QIAGEN representative or contact our Technical Support Center at www.qiagen.com/support/technical-support.

Best regards,

Your EZ2 Connect Team