

July 2011

---

## Important Note

Dear QIAquick Nucleotide Removal Kit user,

We wish to inform you that, due to regulatory guidelines, we have changed the composition of Buffer PN. The new buffer supplied with the QIAquick® Nucleotide Removal Kit — Buffer PNI — shows the same performance as Buffer PN. The quality and yield achieved and overall kit performance are not affected by this change.

Before using Buffer PNI, please add the amount of 100% isopropanol stated on the buffer bottle label. All other protocol steps remain unchanged.

We apologize for any inconvenience that this may cause.

If you have any questions, please call one of the QIAGEN Technical Service Departments or distributors (see back cover of handbooks or visit [www.qiagen.com](http://www.qiagen.com)).

Best regards,

QIAGEN

Trademarks: QIAGEN®, QIAquick® (QIAGEN Group).

