



October 2022

Important Note

QIAxcel® Advanced Power Requirements

Dear valued QIAxcel Advanced customer,

This note is to share with you an important piece of information that you must consider when working with your QIAxcel Advanced:

In addition to the "Power requirements" information listed on page 16 of the *QIAxcel Advanced User Manual*, please be aware that main supply voltage fluctuations are not to exceed $\pm 10\%$ of the nominal voltages.

In case you have any questions, please reach out to QIAGEN Technical Services or your local QIAGEN representative.

Best regards,

QIAGEN

Trademarks: QIAGEN®, Sample to Insight®, QIAxcel® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

10/2022 HB-3148-001© 2022 QIAGEN, all rights reserved.