

QIAGEN Product List

EZ1[®] Advanced Instruments for DNA extraction from human identity and forensic samples

QIAGEN is a leading provider of integrated solutions for human identity and forensic testing. We have an extensive portfolio of technology for extracting DNA from forensic samples suitable for a broad range of throughput and application needs.

Automated sample preparation enables efficient, standardized purification of nucleic acids. Manual sample handling is minimized, ensuring safety to the user and reliable processing of samples.

The EZ1 Advanced Instruments use proven EZ1 Kits for purification of DNA from a wide range of starting materials. Purified nucleic acids are ready for immediate use in downstream applications like quantification or genetic fingerprinting of forensic samples. The combination of easy-to-use instruments and prefilled, sealed reagent cartridges makes efficient nucleic acid purification easy and reliable.



QIAGEN Product List

Application	Product name	Cat. no.
Sample preparation		
For fully automated sample preparation from 1-6 samples	EZ1 Advanced	9001410
For fully automated sample preparation from 1-14 samples	EZ1 Advanced XL	9001492
For unmatched sensitivity and performance for the automated extraction of forensic sample types on the EZ1 system	EZ1 DNA Investigator Kit (48)	952034
Related Products		
For complete tracking of samples and reagents throughout the entire purification process	Bar Code Scanner	9241226
For reproducible high performance by using easy-to-run protocols	Investigator Card	9016387
For easy, trouble-free installation	Installation Service	9240098

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, EZ1®, Investigator® (QIAGEN Group)

© 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia = 1-800-243-800

Austria = 0800/281010

Belgium = 0800-79612

Canada = 800-572-9613

China = 021-51345678

Denmark = 80-885945

Finland = 0800-914416

France = 01-60-920-930

Germany = 02103-29-12000

Hong Kong = 800 933 965

Ireland = 1800 555 049

Italy = 800 787980

Japan = 03-5547-0811

Korea (South) = 1544 7145

Luxembourg = 8002 2076

The Netherlands = 0800 0229592

Norway = 800-18859

Singapore = 65-67775366

Spain = 91-630-7050

Sweden = 020-790282

Switzerland = 055-254-22-11

UK = 01293-422-911

USA = 800-426-8157

