



QIAGEN GmbH  
QIAGEN Str. 1  
40724 Hilden, Germany

August 2022

## Release of EZ2<sup>®</sup> Connect Fx Protocol Package v6

Dear valued customer,

At QIAGEN<sup>®</sup>, we are continuously making improvements to provide the highest level of performance with our instruments and protocols. Consequently, we are pleased to inform you about the opportunity to upgrade your EZ2 Connect Fx instrument to Protocol Package v6. This protocol package includes the improvements and changes to the protocols below.

Protocol name	New protocol version	Changes
DNA Investigator LV Heated	2	<ul style="list-style-type: none"><li>• GUI text update including clarity of which tube to load in the heater well 11</li><li>• Improved aspiration and dispensation to avoid dripping at the heater positions</li></ul>
DNA Investigator LV RT	2	<ul style="list-style-type: none"><li>• Protocol cleanup</li><li>• Relocation of micro step: Collecting Buffer</li><li>• Delete duplicated micro step: Resuspension of beads</li></ul>

If you have further questions, please contact your local QIAGEN representative or contact our Technical Support Center at [www.qiagen.com/support/technical-support](http://www.qiagen.com/support/technical-support).

Best regards,

Your EZ2 Connect Team